


AD-A100 799

AIR FORCE INST OF TECH WRIGHT-PATTERSON AFB OH SCH00--ETC F/6 6/5
ANALYSIS AND PERFORMANCE EVALUATION OF ELECTROCARDIOGRAM DATA C--ETC(U)
DEC 80 M D TOWNSEND
AFIT/GE/EE/80D-46

UNCLASSIFIED

NL

1 OF 3
ADA
100799



AFIT/GE/EE/80D-46

Accession For	
NTIS GRA&I	X
DTIC TAB	
Unannounced	
Justification	
By	
Distribution/	
Availability	
Annals	
Dist	Special
A	

ANALYSIS AND PERFORMANCE EVALUATION

OF ELECTROCARDIOGRAM

DATA COMPRESSION TECHNIQUES

THESIS

AFIT/GE/EE/80D-46

Melvin D. Townsend
Captain USAF

Approved for public release; distribution unlimited

AFIT/GE/EE/80D-46

ANALYSIS AND PERFORMANCE EVALUATION
OF ELECTROCARDIOGRAM
DATA COMPRESSION TECHNIQUES

THESIS

Presented to the Faculty of the School of Engineering
of the Air Force Institute of Technology
Air University
in partial Fulfillment of the
Requirements for the Degree of
Master of Science

by

Melvin D. Townsend, B.S.E.E.

Captain USAF

Graduate Electrical Engineering

December 1980

Approved of public release; distribution unlimited.

Acknowledgements

This thesis is the result of a great deal of effort and could not have been accomplished without the assistance and understanding of many people. I first want to thank my fellow students, especially Capt. ''Chip'' Lutz and Capt. Lee Baker, for their assistance in microcomputer interface problems. Next I would like to thank my thesis committee, and in particular Dr. Rustan my thesis advisor, for their assistance and counsel during this thesis. Finally, I wish to extend my sincerest thanks to my wife, Gail, without whose love and support this thesis would never have come to completion.

Contents

	Page
Acknowledgements	i
List of Figures	iv
Abstract	vi
I. Introduction	1
Background	1
Problem	3
Scope	3
General Approach	4
Sequence of Presentation	5
II. Data Compression and the EKG	8
Entropy Reduction	9
Sampling	10
Orthogonal Expansions and Filtering	11
EKG Filtering Compression	13
Transforms	14
Fast Fourier Transform	15
Discrete Karhunen-Loeve Transform	17
Redundancy Reduction	23
Predictors/Interpolators	24
Difference Reductions	26
Time Compression	27
Entropy Encoding	29
Summary	31
III. Tolan and Dower Compression Techniques	33
Tolan EKG Data Compressor	39
Tolan Decorrelator	39
Tolan Entropy Encoder	40
Dower EKG Data Compressor	47
Dower Decorrelator	48
Dower Entropy Encoder	51
Tolan versus Dower	57
Other Microcomputer EKG Compressors	59
Turning Point Algorithm	59
Marquette Algorithm	61
Chapter III Summary	63
IV. EKG-Data Acquisition and Analysis System	66

EKG-DAAS Hardware	66
EKG-DAAS Software	69
EKG-EXEC	74
NOCPRS	76
DISPLAY	78
TOLAN	80
Summary	83
V. Experimental Procedure, Data Analysis, and Results	84
Experimental Procedure	84
EKG Compression Measurement Parameters	87
Data Analysis and Results	92
TOLAN and DOWER Performance Comparison	98
Summary	99
VI. Summary, Conclusions and Recommendations	101
Bibliography	104
Appendix A: The Electrocardiogram	108
The Physiology and Electrical Characteristics of the Heart	108
The Cardiac Cycle	111
Heart Disease and the EKG	113
Appendix B: Fundamentals of Information Theory	117
Information Source	118
Transmission Channel	121
Encoder	125
Decoder	127
Appendix C: Software Program Listings	129
EKG-EXEC	130
DISPLAY	143
PRSTAT	151
NOCPRS	166
TOLAN-A	176
DECPRS	189
ENTROPY	199
Appendix D: EKG Experimental Data Listing	206
Appendix E: Equipment Specifications	231
Vita	252

List of Figures

Figure	Page
1 Typical EKG Waveform with Arrhythmia	13
2 Typical EKG Redundancy Reduction Data Compressor.	23
3 Second Order Difference Operator.	28
4 Relative Frequency of 2nd Order Difference Operator	28
5 Nonstationary Waveform With Random Level and Slope	34
6 EKG Waveform With Moving Average Slope	35
7 Simulated Variance Reduction Via Second Difference	37
8 Tolan Collection and Decorrelation Algorithm	41
9 Tolan Code and Code Tree	42
10 Tolan Variable Length Encoder Algorithm	44
11 Dower Collection and Decorrelation Algorithm	45
12 Dower EKG Data Compression System	47
13 Three Hypothetical Sample Sequences	50
14 Example of State Space Partition for 4 Symbols	52
15 Dower Encoder Accumulator-Memory Buffer Interface	53
16 Dower VLC Algorithm	56
17 Turning Point Patterns	60
18 EXORCISER Component Module Layout	68
19 EKG-DAAS Hardware Configuration	69
20 EKG-DAAS Software Control Flowgraph	71
21 EKG-DAAS Overlay Structure and Memory Map	72
22 EKG-EXEC Functional Flowchart	75
23 NOCPRS Functional Flowchart	77

24	DISPLAY Functional Flowchart	79
25	TOLAN Functional Flowchart	81
26	Experimental Data Collection Setup	86
27	Original and Reconstructed EKG With Best Compression Ratio	89
28	Original and Reconstructed EKG With Average Compression Ratio	90
29	Original and Reconstructed EKG With Worst Compression Ratio	91
30	Tolan Compression Ratio Breakdown	95
31	Second Difference Distributions	97
A1	The Heart Cross-Section	109
A2	Wilson EKG Electrode System	111
A3	Typical EKG Waveform	113
B1	The Communication System Model	117
B2	A Discrete Memoryless Channel	121
B3	A Typical Rate Distortion Function	124

Abstract

EKG data compression techniques were investigated for potential real time implementation on an 8 bit Motorola 6800 microprocessor. Research indicated entropy reduction transform techniques such as the Fast Fourier Transform and the discrete Karhunen-Loeve Transform were not feasible for implementation on the 6800. Two redundancy reduction (RR) techniques (TOLAN and DOWER) utilizing 2nd order difference operations in conjunction with variable length encoding were studied in detail. One such RR technique (TOLAN) was fully implemented and tested with ''in vivo'' EKG data. Analysis revealed compression ratios ranging from 1.25:1 to 2.26:1. Investigation of the poor performance of the compression algorithm showed significant degradation of the 2nd order difference ''decorrelator'' due to a noisy collection environment. It was concluded that real time EKG data compression is feasible on the 6800 but that time compression techniques which store a zero value sequence counter versus the value of zero are not efficient in a high noise environment.

ANALYSIS AND PERFORMANCE EVALUATION
OF ELECTROCARDIOGRAM DATA COMPRESSION TECHNIQUES

I. Introduction

Background

Currently the USAF School of Aerospace Medicine (USAFSAM) is receiving more than thirty thousand electrocardiogram (EKG) data records a year from Air Force flight personnel worldwide. This data is presently recorded on paper strip charts for immediate medical analysis and long term storage (via microfilm).

Advances in computerized biomedical analysis has generated a need for digitized storage of the EKG waveform to allow computerized interpretation and comparison of present and past cardiographic data. As computerized diagnosis becomes more accurate, a long term digital record of historical EKG data will enable the cardiology staff at USAFSAM to identify developing heart disease before it becomes a danger to a flight crew member or his fellow crewman.

Data compression of the sampled EKG has been an area of active research since the late nineteen sixties. References (1), (7), (12), (26), (28), (29), (32), (33), and (35) are representative of the research efforts performed in the last ten years. The reasons for compressing EKG data are twofold: 1) digital storage costs are rapidly approaching

analog storage costs and ; 2) increased use of computer aided diagnosis requires large digital data bases.

In the past, EKG data compression has generally been performed at a central computer facility. The data is normally collected in a physician's office or clinic and transmitted in an analog format over a standard telephone link to the central computer. The data is then digitized and input to a computer for diagnosis and storage (in compressed form).

With the current advances in microprocessor technology, sampling and compression at the collection site is now a viable alternative. State of the art digital communication systems can operate at 9600 bits per second over the standard 3 kilohertz (KHz) bandwidth telephone channel. With error detection and correction protocols (Ref 27), high fidelity digital transmission of the compressed EKG to the central processing center appears to be the wave of the future (Ref 21:253-254).

The Department of Defense (DOD) is currently installing a computerized EKG interpretation system at centralized US military medical centers worldwide. To collect and transmit this EKG data, remote medical clinics will use a commercial EKG ''cart'' containing a Motorola 6800 microprocessor. This EKG cart performs internal data compression, record formatting, and error protective ''channel'' encoding for digital transmission to the central medical center.

As an independent study, USAFSAM has sponsored this

thesis to investigate the general field of microprocessor based EKG data compression. The results of this study are to be used for comparison against the data compression achieved in the DOD system and to create a measurement baseline to evaluate future microcomputer based EKG data compression systems.

Problem

The problem addressed by this thesis was the analysis of currently available EKG data compression algorithms for implementation in a microprocessor based computational environment. Additionally, a performance measure was to be developed by which differing EKG data compression techniques could be compared.

Scope

In this thesis, an EKG Data Acquisition and Analysis System (EKG-DAAS) was constructed. The EKG-DAAS collects 3 leads of an amplified EKG (i.e., Analog/Digital (A/D) dynamic range $\approx \pm 5$ volts), and samples the data at a operator controlled rate between 300-700 hertz. Data is digitized at 12 bit precision and subsequently rounded to 8 bits for uniform truncation error performance (Ref 24:413-418).

The EKG-DAAS hardware was developed around a Motorola Exorciser microcomputer (6800 microprocessor) with associated A/D converter, disk memory, computer terminal, and

hard copy printer. The EKG-DAAS hardware is controlled via a software program called EKG-EXEC. EKG-EXEC is an assembly language program which performs terminal, printer, and disk input/output (I/O) operations as well as providing a supporting structure for the data compression and analysis software. Because of the lack of a high order language (e.g. FORTRAN, PASCAL, etc.) and the desire for maximum program execution speed, all programming done in this thesis is in assembly language. This has proven to be a major limitation.

Only one EKG data compression algorithm was completed and implemented on the EKG-DAAS. "In vivo" EKG data was taken, however, and EKG data compression performance analyzed.

General Approach

To accomplish the objectives in the problem statement, the literature was first searched for EKG compression algorithms whose implementation and execution on a 6800 microprocessor was considered feasible. As a result of this literature search and private correspondence (Ref 11 and 31), two compression routines were found. These two EKG data compression algorithms will hereafter be referred to as the Tolan and Dower methods (Ref 31,12).

With the Tolan and Dower algorithms identified, construction of the EKG-DAAS was begun by assembling the data acquisition subsystem. This data acquisition system

consists of an Analog/Digital converter and a interrupt sampling clock. The A/D was calibrated and coding of the EKG-EXEC software started. Several months of effort resulted in the EKG-EXEC software and programming of the Tolan compression algorithm initiated. Upon completion of the Tolan compression routines, the Dower algorithm was analyzed but not implemented.

Based on the results of the Tolan compression algorithm, a performance measure was formulated which compares achieved compression against an approximate maximum compression computed from the data statistics. The Dower compression algorithm was analyzed for similarity with the Tolan technique and an estimated performance figure calculated with respect to the compression measure or ''metric''.

Sequence of Presentation

Chapter 2 begins the thesis development with a general survey of the field of data compression. Data compression is shown to be divided into two subclasses (Entropy Reduction and Redundancy Reduction) and each subclass is defined. Several EKG compression techniques found in the literature were included in the ER category. These ER techniques were described and their performance advantages and limitations analyzed. Chapter 2 proceeds by describing redundancy reduction and several RR electrocardiogram compression algorithms are also analyzed. Chapter 2

concludes with the determination that implementation of the identified ER compression routines would require programming efforts beyond the scope of this thesis. Hence this thesis is limited to redundancy reduction compression.

Chapter 3 describes, in detail, the Tolan and Dower redundancy reduction algorithms whose implementation on the Exorciser was considered feasible. The Tolan algorithm is described first followed by a discussion of the Dower method. Two other EKG compression techniques, discovered during the research of this thesis, are also summarized.

Chapter 4 reviews the hardware and software configuration of the EKG-DAAS. The hardware system is described first and the specifications of the Exorciser microcomputer, A/D converter, and I/O peripherals are presented. Following the hardware description, the EKG-EXEC program is documented and the software design philosophy examined. Finally the Tolan compression module is discussed.

Chapter 5 outlines the results obtained using the EKG-DAAS with the Tolan algorithm. A detailed description of the experimental setup is presented along with a discussion of the performance parameters measured by EKG-EXEC. Finally a performance metric is described and the Tolan compression results compared against this metric. Chapter five ends with an analysis of how well the Dower compression technique would have performed against the performance metric. Chapter 6 concludes the thesis and

presents recommendations for future study.

This thesis contains five appendices. Appendix A surveys basic electrocardiology and is recommended to readers unfamiliar with this subject. Appendix B is a tutorial on Information and Coding Theory and is likewise recommended to those readers unfamiliar with this field. Appendix C contains a listing of the EKG-EXEC assembly language software as well as one BASIC program used on the Exorciser for data analysis. Appendix D contains data printouts of the data taken in the data collection experiment. Appendix E presents photocopies of the specification sheets for the equipment used in this thesis.

With the sequence of presentation outlined, attention now turns to the theoretical section of this thesis. The first subject is data compression theory.

II. Data Compression and the EKG

Data compression is an operation in which data from an information source is ''simplified'' or ''filtered'' in a manner that produces an approximation of the original with at most some predefined amount of distortion. Some form of data compression is usually necessary when storage limitations, bandwidth requirements, or transmission channel capacity prohibit operation on the original data.

In general, data compression can be divided into two types of operations (Ref 10:4). The first operation, called entropy reduction (ER), is an irreversible transformation which reduces or compresses the data by mapping a source into an approximation of itself with a lower entropy rate. Sampling is an example of such a transformation. The second operation, known as redundancy reduction (RR) compresses the sampled data train by reducing, or eliminating, the redundancy existing in digital sequence. Since the redundant components of the data train contain no information about the source, RR is an ''exact'' data compression operation.

This chapter discusses data compression as it has been applied to the electrocardiogram (EKG). Two types of data compression strategies prevail. The first strategy involves two (or more) ER operations on the EKG data and attempts to compress the data by filtering selected components of discrete transforms. Because ER operations are irreversible, this type of data compression is sometimes

referred to as "inexact" compression. The second strategy performs data compression by reducing the redundancy present in the sampled EKG data train. This redundancy arises from two causes : 1) neighboring signal samples are not statistically independent and ; 2) the quantized signals amplitudes do not occur with equal probability.

Because of the speed limitations of the Motorola Exorciser microcomputer used in this thesis, only redundancy reduction algorithms were tested. This chapter, therefore, is intended as a review of the work that has been done in EKG data compression to allow comparison with the results obtained by this author.

Entropy Reduction

As is known from information theory (Ref 19), entropy $\{H(X)\}$ is defined as a measure of the "randomness" or "uncertainty" of an information source. If the symbols emitted by the source are statistically independent, then the source is said to be "memoryless" and entropy is given by the equation:

$$H(X) = - \sum_{i=1}^N p_i \log_2 p_i \quad (\text{bits}) \quad (1)$$

where p_i represents the probability of occurrence of the i th symbol and N is the number of distinct symbols output by the source. For those "symbols" in a data set which occur

with zero probability, the term $p_i \log_2 p_i$ is defined equal to zero.

Immediately obvious from Eq.(1) is the observation that $H(X)$ is defined on a discrete probability distribution. In fact, Shannon (Ref 30) defined the entropy of a continuous source as equal to positive infinity. Thus an ER transformation must, in some manner, ''discretize'' a continuous waveform into a countable set of components and attach some probability to the elements of the set.

Sampling. The clearest example of this ''discretization'' is sampling, obviously the most important operation necessary for digital signal processing. In sampling, an electrical circuit (such as an Analog-to-Digital (A/D) converter) periodically measures the value of a signal $x(t)$ and records the data as a numeric (usually binary) number. If the signal is sampled at least twice the highest frequency component of $x(t)$, and the duration of the sampling operation is long enough that ''aliasing'' effects (Ref 36:68-72) of ''windowing'' are negligible, then all of the ''frequency domain'' components of interest in $x(t)$ will be preserved. Physical constraints, however, limit the accuracy of the amplitude measurement to some finite precision. Amplitude information residing below the sensitivity of the A/D converter is irretrievably lost.

As an example, let the A/D converter digitize to 8 bits. With 8 levels there are 256 possible outputs, each with a certain probability of occurrence. If the continuous

waveform fed to the A/D is the output of a stationary stochastic process with uniform statistics, then each numeric value will occur with probability of $1/256$ (assuming the source max/min deviation \geq the A/D dynamic range). Assuming sample-to-sample independence, the source entropy is then calculated as

$$H(X) = \sum_{256} (1/256) \log_2 256 = 8. \quad (2)$$

Entropy has been ''reduced'' from $+\infty$ to 8. If the stochastic source has ''less random'' statistics (like gaussian), entropy would be reduced even further.

All of the EKG compression routines discussed in this thesis are implemented on digital computers, hence the ER operation of sampling is always performed. In ''inexact'' EKG data compression, the next operation is another ER mapping in which the sample sequence is transformed into an alternate domain and filtered.

Orthogonal Expansions and Filtering. A bandlimited waveform $x(t)$ can be expanded (over a given interval T) as a linear combination of orthonormal basis functions $\phi_n(t)$ (Ref 36:20-21). That is

$$x(t) = \sum_{n=0}^{\infty} a_n \phi_n \quad (3)$$

The functions $\phi_n(t)$, as $n=0,1,\dots,\infty$, are defined as orthonormal if

$$\int_T \phi_i(t) \phi_k(t) dt = 1, \quad i=k$$

$$=0 \quad \text{otherwise} \quad (4)$$

The orthonormal set $\{\phi_n(t)\}$ is called complete (C.O.N.) if, for any given $\varepsilon > 0$, there exists an N and a finite expansion

$$\tilde{x}(t) = \sum_{n=0}^{N-1} a_n \phi_n(t) \quad (5)$$

such that

$$\int_T |x(t) - \tilde{x}(t)|^2 dt < \varepsilon \quad (6)$$

Given a complete, orthonormal set $\{\phi_n(t)\}$, then representation of any physical, noise-limited signal (over a prescribed interval T) is possible with a finite set of weighting coefficients $\{a_0, a_1, \dots, a_{N-1}\}$.

If k ($k < N$) terms in Eq.(5) are selectively suppressed (i.e., $a_0, a_1, \dots, a_k = 0$), then the signal $x(t)$ is said to be filtered of those components and a compressed representation of $x(t)$ is produced.

EKG Filtering Compression. A great deal of work (1),(2),(26),(33),(35),(21) has been done in attempting to "compress" EKG data by storing, or transmitting, selective components of an orthonormal expansion. This type of compression strategy is complicated by the complex characteristics of the EKG waveform. Foremost are the limitations imposed due to the "nonstationarity" or variability of the EKG source, both within any given waveform and the general population as a whole.




Fig 1. Typical EKG waveform with arrhythmia.

During an EKG collection, the heart rate of the person undergoing testing can vary. This variation, or arrhythmia, requires that a suitable reference point be established on

the P-QRS-T complex (Fig 1) around which the orthogonal expansion can be performed. Without this reference, the data falling within an expansion interval T (based purely on a constant Δt) will appear to be ''almost'' random. If the ''filtering'' used to effect compression is based on an expected waveform within the transform window, improper registration will cause severe degradation in compression efficiency. The above situation is known as the ''epoch'' problem and is discussed by Womble (Ref 35:703).

Assuming that the expansion interval is aligned properly on some feature of the EKG (e.g. R wave), then the problem of variable P, QRS, and T complexes must be considered.

If the individual under test has a heart disease, then in many cases the P, QRS, and T segments of that person's EKG may vary significantly from the population ''norm''. Again, if the compression filtering is based on an assumed, or ''normal'' waveform, then significant reduction in compression efficiency can be expected.

These problems, among others, require that the EKG ''filtering'' compressor be robust enough to handle the variations possible within the waveform. With these limitations in mind, discussion will now proceed to compression transforms.

Transforms. Two major types of discrete filtering compression strategies, as applied to EKGs, will now be discussed. These are the Fast Fourier Transform (FFT) and

the discrete Karhunen-Loeve Transform (DKLT). These two techniques are representative of current research into EEG filtering compression and will serve as good basepoints against which to compare the redundancy reduction compression techniques outlined in chapter 3.

Fast Fourier Transform. The fast Fourier transform is a computationally efficient algorithm for calculating the discrete Fourier transform (DFT). The following discussion involves the DFT, but calculation by the FFT is implied.

Data compression using the DFT is achieved by zonal filtering in which "zones" of the transform sequence are selectively discarded. Usually these zones are defined by a "cutoff frequency" f_c and only those components less than or equal to f_c are saved. Data reconstruction is performed by computing the inverse DFT (FFT) on the filtered sequence and replacing the filtered components with zeroes. Using the symmetry relationships described in Oppenheim and Shaffer (Ref 24:103-105), then the transform sequence for $N=8$ is:

$$X(k) = \begin{bmatrix} X(0) \\ X(1) \\ X(2) \\ X(3) \\ X(4) \\ X(5) \\ X(6) \\ X(7) \end{bmatrix} = \begin{bmatrix} R(0) + jIm(0) \\ R(1) + jIm(1) \\ R(2) + jIm(2) \\ R(3) + jIm(3) \\ R(4) + jIm(4) \\ R(3) - jIm(3) \\ R(2) - jIm(2) \\ R(1) - jIm(1) \end{bmatrix} \quad (7)$$

where R is the real part, Im the imaginary part of the complex number X , and $j=\sqrt{-1}$.

Zonal filtering takes advantage of the symmetry in Eq.

(7) by saving those positive frequency components less than or equal to a certain cutoff frequency f_c . If f_c were chosen as equal to $X(2)$, then only $X(0)$, $X(1)$, and $X(2)$ would be saved. On reconstruction, $X(3)$, $X(4)$, and $X(5)$ would be set equal to zero (i.e. they have been filtered out) while $X(6)$ and $X(7)$ would be recreated from $X(1)$ and $X(2)$ using symmetry. The inverse DFT (FFT) now produces a "filtered" approximation of $x(n)$.

EKG data compression via filtered FFT spectra has been studied by Womble (Ref 35) and the TRW corporation (Ref 33). In both studies the distortion criterion used in establishing f_c was visible reproducibility with no detectable distortion.

In the TRW study, compression ratios (defined as bits in : bits out) from 2:1 to as high as 17:1 were obtained using zonal filtering. The high compression ratios were measured using highly rhythmic EKGs taken from an individual (an astronaut) in a low noise environment. This extraordinary compression ratio is far from "normal", however. Womble (Ref 35) has shown that, on the average, 40-80 terms of a 512 sample FFT (sampled at 500 Hz) are necessary to reproduce the EKG with acceptable visual distortion. Since the FFT requires storage of 2 numbers per term (real and imaginary components) then, in general, FFT zonal filtering compression ratios of 5:1 to 3:1 are more common.

The FFT is often used in digital signal processing

because of its speed and representation in the frequency domain. As will be seen from the data in chapter 5, compression ratios of 5 : 1 are usually higher than those obtainable by redundancy reduction techniques. Even so, the FFT is not the "optimal" transform for representing a sample sequence.

Discrete Karhunen-Loeve Transform. Because all of the filtering compression techniques are ER transformations, distortion upon reconstruction is inevitable. One would like a performance measure against which this distortion could be evaluated.

One common measure of performance is the mean square error defined by the equation:

$$e(M) = E\{(\underline{X} - \tilde{\underline{X}})^2\} \quad (8)$$

where E is the statistical expectation operator, \underline{X} is a discrete signal composed of N sample values, and $\tilde{\underline{X}}$ is the estimate of \underline{X} via some orthogonal coordinate system. Given that a representation of \underline{X} is desired with less than N components, Ahmed and Rao (Ref 2:200-203) have shown that the DKLT is the optimum transform to minimize mean square error.

The DKLT can be described as follows: Let the orthonormal transform matrix [T] be defined as

$$[T]' = [\phi_1, \phi_2, \dots, \phi_N] \quad \{[T]' = [T] \text{ transpose}\} \quad (9)$$

and ϕ_i are N dimensional, real valued basis vectors. Hence a transform vector \underline{Y} can be created by

$$\underline{Y} = [T] \underline{X} \quad (10)$$

where $\underline{X} = [x_1, x_2, \dots, x_N]$; $\underline{Y} = [y_1, y_2, \dots, y_N]$. Since the ϕ_i are orthonormal, then

$$\underline{X} = y_1 \phi_1 + y_2 \phi_2 + \dots + y_N \phi_N = \sum_{i=1}^N y_i \phi_i \quad (11)$$

The goal is to optimally represent \underline{X} by a subset $\{y_1, y_2, \dots, y_M\}$ where $M < N$. If the remaining $N-M$ terms are represented by the constants b_i , then an estimate $\tilde{\underline{X}}$ is defined such that

$$\tilde{\underline{X}} = \sum_{i=1}^M y_i \phi_i + \sum_{i=M+1}^N b_i \phi_i \quad (12)$$

An error vector, $\Delta \underline{X}$ is now created where

$$\Delta \underline{X} = (\underline{X} - \tilde{\underline{X}}) = \sum_{i=M+1}^N (y_i - b_i) \phi_i \quad (13)$$

Now the mean square error defined in Eq.(8) can be redefined

as

$$e(\underline{X}) = E\{(\Delta \underline{X})'(\Delta \underline{X})\}$$

$$= E\left\{ \sum_{i=M+1}^N \sum_{j=M+1}^N (y_i - b_i)(y_j - b_j)(\phi'_i \phi_j) \right\}$$

$$= \sum_{i=M+1}^N E\{(y_i - b_i)^2\} \quad (14)$$

Ahmed and Rao (Ref 2:202) show that by minimizing Eq.(14)

$$e(M) = \sum_{i=M+1}^N [\phi'_i E\{(\underline{X} - \bar{\underline{X}})(\underline{X} - \bar{\underline{X}})'\} \phi_i] \quad (15)$$

where $\bar{\underline{X}}$ is the mean of \underline{X} .

The expectation in Eq. (15) is recognized as the covariance of \underline{X} , hence Eq. (15) becomes

$$e(M) = \sum_{i=M+1}^N (\phi'_i [K_{\underline{X}}] \phi_i) \quad (16)$$

where $[K_{\underline{X}}]$ is the covariance matrix of \underline{X} .

By minimizing Eq.(16), Ahmed and Rao (Ref 2:200-205)

show that the orthonormal functions ϕ_i become the eigenfunctions of the covariance matrix $[K_x]$ and the mean square error becomes

$$e(M) = \sum_{i=M+1}^N \lambda_i \quad (17)$$

where the λ_i are the eigenvalues of $[K_x]$. By expanding \underline{X} with the eigenfunctions corresponding to the M largest eigenvalues, then \underline{X} is filtered of $N-M$ components with minimum mean square error. This is the discrete Karhunen-Loeve transform.

The DKLT has been studied by Ahmed (Ref 1) and Womble (Ref 35). Both studies compared the performance of EKG data compression using the DKLT against other orthogonal transforms.

Ahmed's tests utilized canine EKGs and clearly showed the optimality (in the 'mean square sense') of the DKLT. Ahmed, however, deemed the DKLT to be too complex for practical implementation, and proceeded with the development of EKG compression strategies using suboptimal transforms. The details of these suboptimal expansions will not be given here, but the reader is referenced to Ahmed's paper (Ref 1).

Womble, on the otherhand, demonstrated that the DKLT is not too difficult to implement. In his experiments, Womble

took 3 lead vectorcardiogram (VCG) data (see appendix A), hereafter referred to as X,Y,Z, and first transformed the data into a different, orthogonal coordinate system (U,V,W). The transformation was determined by solving for the eigenvalues of the 3 X 3 matrix:

$$S = 1/N \sum_{i=1}^N \begin{bmatrix} X(i) \\ Y(i) \\ Z(i) \end{bmatrix} \begin{bmatrix} X(i), Y(i), Z(i) \end{bmatrix} \quad (18)$$

where N was chosen equal to 200 and X(i),Y(i),Z(i) are the X,Y,Z components of the VCG in the Frank (Ref 14) coordinate system.

Next an "average", or mean, heartbeat was calculated using 900 patients and the data (in the eigenvector coordinate system) was subtracted from the mean forming a "patient" vector ρ . The vector ρ is defined as

$$\rho = \begin{bmatrix} u(1) - u_m(1) \\ u(2) - u_m(2) \\ \vdots \\ u(N) - u_m(N) \\ v(1) - v_m(1) \\ v(2) - v_m(2) \\ \vdots \\ v(N) - v_m(N) \\ z(1) - z_m(1) \\ z(2) - z_m(2) \\ \vdots \\ z(N) - z_m(N) \end{bmatrix} \quad (19)$$

This patient vector ρ is expanded using the eigenfunctions

of the matrix

$$L = 1/N \sum_{i=1}^N (\rho_i)(\rho_i)' \quad (20)$$

where N is now a large number of patients being "averaged" (up to 300 reported in Ref 33). The eigenvectors of this tremendous matrix are calculated (once for all time), and then M ($M \ll 200$) are used to expand the data vector ρ .

The above approach has allowed accurate representation of EKG sample sequences, properly aligned on the heartbeat, with an $M=20$ DKLT coefficients per lead per second (Ref 35). At Womble's sample rate of 250 Hz, this represents a compression ratio of over 12 : 1. With the most significant eigenfunctions stored from the solution of Eq. (20), then the compression of digitized EKG data by filtering the DKLT expansion on ρ is approaching feasibility on a microcomputer.

The preceeding discussion was intended as a basic tutorial on some of the entropy reducing (ER) transformations currently under test for EKG data compression. As stated earlier, ER operations are "inexact" because some "information" is always discarded. By discarding data efficiently, zonal filtering can achieve compression ratios larger than those obtained using the redundancy reduction techniques implemented in this thesis.

Redundancy Reduction

Redundancy reduction is an operation which performs data compression by indentifying and removing the redundant components of a digital sequence. This redundancy exists for two reasons : 1) the sample points are not statistically independent and ; 2) the quantized amplitude values of the sample train do not occur with equal probability. The redundant components of the data sequence carry no ''information'' (see appendix B), hence their removal does not affect the ''message'' content of the data.

This section of chapter two discusses redundancy reduction (RR) operations which have been applied to the EKG. This is done because the compression algorithms examined in this thesis are of the RR type, and the background given here will aid in the descriptions given in chapter three.




Fig 2. Typical EKG Redundancy Reduction Data Compressor.

The sampled EKG contains redundancy resulting from both

sample-to-sample correlation and unequal amplitude probability. As is illustrated by Figure 2, each of these redundancy components can be reduced by a different operator.

The decorrelator, as inferred by its title, transforms the data in a manner which decorrelates the data stream. The decorrelated residual sequence is then efficiently compressed by means of the entropy encoder.

To reduce the redundancy arising from intersample correlation, two basic approaches are used. The first approach utilizes linear predictor/interpolators (Ref 28,33) to produce an information bearing "residual" sequence. The second method generates the residual sequence by taking successive differences on the data stream. As will be shown, these residual sequences carry all of the source information with minimum intersample redundancy. The predictor/interpolator method is discussed first.

Predictors and Interpolators. A predictor estimates the next sample value of a sequence (i.e. x_n) based on a linear combination of k past samples. That is

$$\tilde{x}_n = \sum_{i=1}^k a_i x_{n-i} \quad (21)$$

where a_i are coefficients chosen to minimize the mean square error $\sigma_e^2 = E\{(X - \tilde{X})^2\}$ between the sample sequence X_n and the

predicted sequence \tilde{X}_n . An error , or residual, sequence e_n can then be formed where

$$e_n = X_n - \tilde{X}_n \quad (22)$$

thus

$$X_n = \tilde{X}_n + e_n \quad (23)$$

Conceptually, this technique decomposes a sample value into a part which is correlated with k past sample values (i.e. the redundant component) and a part which is uncorrelated with them. Linear mean square estimation theory (Ref 24:385-430) shows that the uncorrelated part (e_n) may be retained alone with no loss of information.

In their research in EKG data compression, Ruttimann and Pipberger (Ref 28:616) prove that since $E\{e_n\}=0$, then the mean square error σ_e^2 is equal to the variance of e_n . With a second order predictor ($k=2$), Ruttimann and Pipberger have demonstrated a variance reduction (σ_x^2/σ_e^2) of over 25 : 1. This reduction in variance implies that the EKG residual sequence e_n is much more tightly clustered around a given mean than is the original sequence X_n . This clustering, as will be discussed later, enhances data compression by means of entropy , or source encoding.

In an analogous way, interpolators estimate a value of x_n . In the case of the interpolator, however, the estimate

of x_n consists of a linear combination of past and future samples. That is

$$\tilde{x}_n = \sum_{i=1}^k a_i x_{n-i} + \sum_{i=1}^m b_i x_{n+i} \quad (24)$$

where k past and m future values are used. As with the predictor, the coefficients a_i and b_i are again chosen to minimize the mean square error. As was the case with the second order predictor, Ruttiman and Pipberger (Ref 28:617) have found that a second order interpolator ($k=1, m=1$) yields EKG residual sequences e_n with the most significant variance reductions (greater than 31 : 1).

EKG data compression utilizing predictors/interpolators have been studied by the TRW Corporation (Ref 33), as well as Ruttimann and Pipberger (Ref 28). In both cases, predictors/interpolators of order 2 seem to prevail. As is shown by the two groups above, second order systems have the smallest residual sequence variance hence are most amenable to entropy encoding. An alternate approach used to reduce intersample redundancy is by means of difference operations.

Difference Reduction. In difference reduction, the residual sequence is formed by taking successive differences of the sample data train. Because there is no multiplication (as in Eq.(24)), the difference operation is inherently a

simpler procedure than is prediction or interpolation. Difference sequences, however, are not optimized with respect to mean square error and less ''efficient'' redundancy reduction (decorrelation) occurs.

An example of an EKG difference reduction compressor is illustrated in Figure 3. This technique was implemented by Cox and Ripley (Ref 7) and tested against a data base of 45 patients. Figure 4 shows that Cox and Ripley's second difference decorrelator produced a sharply peaked relative frequency distribution. This ''peakedness'' or clustering of the residual sequence permits efficient entropy encoding. A variant of the difference reduction above is time compression.

Time Compression. EKG time compression, as implemented by Dower, Berghofer, and Stewart (Ref 12,29), uses a second order difference reduction but instead of keeping the value of zero (the most common), a run length counter (Δt) is kept. This run length counter measures the number of repetitive $\Delta^2 X$ sequence terms equal to zero. By saving (then encoding) only those $\Delta^2 X \neq 0$, and the Δt between the nonzero second differences, data compression can be realized.

Time compression is the method that has been implemented and studied extensively by this author. Chapter 3 discusses two different RR compression algorithms which both use second order time compression for the decorrelator. For now, discussion is focused on the second

operation identified in Figure 2, the entropy encoder.




Fig 3. Second Order Difference Operator (from Ref 7:336).




Fig 4. Relative Frequency of 2nd Order Difference Operator (Ref 7:336)

Entropy Encoding. Entropy, or source, encoding operates on the remaining redundancy in the residual sequence e_n resulting from unequal source symbol probability. How well e_n can be compressed is given by Shannon's noisless source coding theorem (Ref 30) which states that the average number of binary symbols per source output can be made to approach the entropy of the source and no less. Unfortunately, the EKG "'source'" has memory, hence the entropy of this source will be less than that which would be calculated from Eq.(1). In general, the true value of a nonstationary, memory source is difficult (or impossible) to calculate.

To circumvent this "'memory'" problem, the assumption is made in the literature (Ref 29,12, 7,35) that the "'uncorrelated'" nature of the residual sequence e_n approaches independence, hence the entropy of e_n as calculated by Eq.(1) is a good bound on the possible compression. Since the variance of e_n is sharply peaked around zero, entropy encoding via variable length coding (VLC) appears attractive.

In variable length coding, those values of e_n which occur most often are assigned the shortest code words (i.e. fewest code symbols per source symbol). As example, assume that both the source symbol alphabet and the code alphabet are binary. The values in the sequence e_n are the result of algebraic operations on fixed length numeric sample values and are hence fixed length binary numbers (c.g. $L=8$). The variable length coder maps these fixed length numeric values

(binary) into uniquely decodable (UD) binary code words whose bit length is a function of the frequency of occurrence of the values in e_n .

Coding theory (Ref 19:237-248) shows that, in principle, it is possible to have variable length codes with an average bit length (\bar{l}) equal to the entropy of source. Hence with the right VLC, a compression ratio of L/\bar{l} is achievable.

A common VLC used in EKG data compression research (Ref 7,28) is the Huffman code. Huffman, in 1952, developed an algorithm (Ref 16) for generating the optimal UD code (assuming stationary, memoryless source). This code is in a class of UD codes called prefix codes in which no code word is the prefix of any other. Unfortunately, a codeword must be assigned to every possible symbol which occurs, regardless of how infrequently. This means that although the average code word bit length approaches the entropy of the source, the longest code word can be substantially larger. In a straight Huffman code, these "'long'" code words can induce severe problems due to "'buffer overflow'" (Ref 17).

For practical implementation in EKG data compression, Ruttimann and Pipberger (Ref 28) and Cox and Ripley (Ref 7) constructed a modified Huffman code. In this modified code, the residual sequence source words are partitioned into a frequent set and an infrequent set (known as "'else'"). A Huffman code was then formed with all of the residual words

in the first set plus a special code word used as a prefix for any source word from the infrequent set. The prefix, when it occurs, is followed by a fixed length suffix which contains the value of the infrequent source word. The probability that "'else'" will occur must be kept small enough to maintain the efficiency of the truncated Huffman code.

Using the modified Huffman code in conjunction with a second order interpolator, Ruttimann and Pipberger (Ref 28) have attained compression ratios as high as 9 : 1. To get this ratio, however, significant digital signal processing was performed on the "'raw'" 8 bit data. This preprocessing involved digital filtering for noise reduction and Lagrange interpolation/decimation to produce an effective 200 Hz sample rate. Nonetheless, this represents a significant approach to the compression ratios attainable by the discrete Karhunen-Loeve transform mentioned previously.

Although optimal, Huffman codes are not the only VLC used in EKG data compression. Two different codes will be described in chapter three, one of which (Dower code) may possibly be "'more optimal'" than Huffman in the compression of quantized EKGs.

Summary

This chapter has reviewed the theory of data compression and how this theory has been applied to the electrocardiogram. It was shown that data compression could

be partitioned into two types of operations: 1) entropy reducing transformations which map a data source into an approximation of itself with a lower entropy rate and ; 2) redundancy reduction operations which compress by removing redundancy resulting from sample-to-sample dependence and unequal source symbol probabilities.

The first operation was applied to the EKG via "filtering" of orthogonal expansions of the digitized EKG. Filtering transforms, especially the discrete Karhunen-Loeve transform, were shown to be very efficient "compressors" if the loss of the filtered components were tolerable. If this loss was not acceptable, then redundancy reduction operations are used.

In redundancy reduction techniques, the EKG is first processed by passing the digitized data through a decorrelator which reduced the redundant component caused by source symbol dependence. This decorrelator could be implemented with predictor/interpolators or difference operations. Next the data is encoded via entropy encoding operations, usually utilizing variable length codes. It was finally shown that a second order interpolator, followed by an "optimal" Huffman encoder could achieve compression ratios which approach those obtained by the discrete Karhunen-Loeve transform. Both of these last two techniques, however, require substantial computational overhead.

III. Tolan and Dower EKG Compression Techniques

Two EKG data compression algorithms are studied in detail in this chapter. The first compression technique, hereafter referred to as the Tolan method, was conceived by Dr. Gil Tolan, MC from the USAF School of Aerospace Medicine (USAFSAM), Brooks AFB, Texas. The second compression approach, referred to as the Dower method, was developed by Mr. Roger Dower and Mr. Dave Berghofer from Shaughnessy Hospital, Vancouver, B.C., Canada. Both of the EKG compression procedures are redundancy reduction operations using time compression for decorrelation and variable length codes for entropy encoding.

This author had originally intended to implement and test both the Tolan algorithm and the Dower algorithm on the Motorola microcomputer (See Chapter 4). Unfortunately, difficulties with hardware failures and insufficient software tools (i.e., a high order language) prevented implementation of the Dower algorithm. Nonetheless, this chapter compares the design of both the Tolan and the Dower algorithms and illustrates their differences as well as their similarities. Based on the results of the Tolan compression approach (chapter 5), this comparison will illustrate the potential performance of the Dower compression strategy in a real time microcomputer environment. Before continuing with the Dower and Tolan algorithm descriptions, a short digression will be made.

In both the Tolan and The Dower compressors, a second

difference operation is used to produce an information bearing sequence with reduced component-to-component dependence. Why the second difference operation results in the lowest correlation is a question for closer scrutiny.




Fig 5. Nonstationary waveforms with random level and slope.
(From ref 6:91)




Fig 6. EKG waveform with moving average slope.

A nonstationary process (like the EKG) exhibits a wide amplitude distribution due to the random ''wander'' of the waveform as shown in Figure 6. This wander, or ''moving average'', can be induced by patient variation, EKG apparatus drift, or both. Time series analysis (Ref 6) shows that a nonstationary process which has sample functions which are ''locally stationary'' or homogenous (Figure 5) can be represented by a process model which calls for the d'th difference of the process to be stationary. The proof of this assertion is given by Box and Jenkins (Ref 6:85-125).

If the nonstationary process sample sequence (time series) exhibits a random level as illustrated in Figure 5a, then a first difference operation will remove this ''moving average'' and force the resulting difference sequence to be centered around zero. If the waveform exhibits a random slope, as shown in Figure 5b, then a second difference operation will remove this quadratic ''bias'' with a corresponding reduction in amplitude distribution variance (Figure 7.) Comparison of Figures 5 and 6 shows that an EKG trace can ''look'' similar to the example in Figure 5b.

From the theory of stochastic processes (Ref 9:330-331), it is known that the expected value of a sample mean obtained by sampling a wide-sense stationary random process along a sample function in time is equal to the constant mean value of that random process.




Fig 7. Simulated variance reduction via second differencing.

In addition, the variance of the sample mean is inversely proportional to the number of samples taken when the samples are pairwise uncorrelated. Another fact about wide-sense stationary processes is that their autocorrelation functions are dependent only on the time difference between observation of the sample sequence and that $R(t_1 - t_2) \leq R(0)$ for $t_1 - t_2 \neq 0$. Since the second difference is the first "difference" which can be modelled as coming from a "stationary" process, and the maximum correlation occurs for zero time difference, then the correlation must be reduced for adjacent sequence values in the second difference operation.

This heuristic argument is far from complete and does not explain why the third difference exhibits worse behavior than the second difference. The EKG waveform is a complex function from a complex source and higher order effects could begin to dominate with the third difference operation. Further analysis of this anomaly is left for future study.

The remainder of this chapter is organized in the following manner. First the Tolan algorithm is described, followed by a description of the Dower procedure. Next a short synopsis will be made of three other EKG data compression algorithms uncovered during the research of this thesis.

Tolan EKG Data Compressor

The Tolan EKG data compression algorithm is a redundancy reduction procedure which processes a three lead EKG (VCG) and produces a compressed, digital output. This digital output sequence could subsequently be channel encoded for ''errorless'' transmission (Ref 27 and Appendix B) or stored for later retrieval and reconstruction.

As was the case for the RR techniques described in chapter 2, the Tolan compressor is subdivided into a data decorrelator and an entropy encoder. The decorrelator is discussed first.

Tolan Decorrelator. The Tolan decorrelator is a second order difference reduction operation which utilizes a time compression approach to form a decorrelated output sequence. A second order system was chosen based on the experimental results of Dower and Berghofer (Ref 12) and Cox and Ripley (Ref 7) which showed maximum compression gain with a second order difference operation. The Tolan decorrelator works on a three lead EKG (VCG) signal set, assumed to be sampled at a constant rate. The data compression is achieved in real time between successive samples.

The Tolan decorrelator algorithm is defined in Figure 8. Close examination of the algorithm in Figure 8 reveals that the second difference data is stored only if : 1) any of the three Δ^2 values are nonzero or ; 2) if the Δt counter

records more than 127 repetitive cases where $\Delta^2 x - \Delta^2 y - \Delta^2 z = 0$. Rule 1 was chosen due to the high degree of correlation between the leads of an FKG (i.e. when one lead is changing, so are the others). The capacity of the time counter in rule 2 was arbitrarily chosen to be sufficient to record the long quiescent periods which occur in the FKG (see Figure 1) and short enough to be efficiently stored.

Step 12 of the Tolan algorithm in Figure 8 calls the variable length encoding subroutine (Figure 10) which encodes the Δ^2 terms calculated by the second difference decorrelator. This encoding procedure is the next subject to be discussed.

Tolan Entropy Encoder. The Tolan code is an uniquely decodable variable length code which stores the Δ^2 values as a contiguous sequence of binary 1's. The length of this "run" of binary 1's is equal to the magnitude of the Δ^2 term. To delineate between the "runs", binary 0's are used as codeword delimiters. Since the second difference has both negative and positive values, a sign bit (0 for positive, 1 for negative) immediately follows the 0 bit delimiter. The three values $\Delta^2 x, \Delta^2 y, \Delta^2 z$ are encoded and stored sequentially followed by a delimited, 7 bit, uncoded Δt value. A Δ^2 value of zero is indicated by 3 successive 0 bits.

```

1:  $\Delta x=0, \Delta y=0, \Delta z=0, \Delta t=1$ 

2:  $x(i)=(a/d \text{ ch } 0), y(i)=(a/d \text{ ch } 1), z(i)=(a/d \text{ ch } 2)$ 

3:  $\text{mem}(0)=x(i), \text{mem}(1)=y(i), \text{mem}(2)=z(i)$ 

4:  $\tilde{x}(i+1)=x(i)+\Delta x, \tilde{y}(i+1)=y(i)+\Delta y, \tilde{z}(i+1)=z(i)+\Delta z$ 

5: if ready for next sample then GOTO 6 else GOTO 5

6:  $i=i+1, x(i)=(a/d \text{ ch } 0), y(i)=(a/d \text{ ch } 1), z(i)=(a/d \text{ ch } 2)$ 

7:  $\Delta^2 x=x(i)-\tilde{x}(i), \Delta^2 y=y(i)-\tilde{y}(i), \Delta^2 z=z(i)-\tilde{z}(i)$ 

8: if  $\Delta^2 x \neq 0$  or  $\Delta^2 y \neq 0$  or  $\Delta^2 z \neq 0$  then GOTO 11 else GOTO 9

9:  $\Delta t=\Delta t + 1$ 

10: if  $\Delta t \leq 127$  then GOTO 4 else GOTO 12

11:  $\Delta x=\Delta x + \Delta^2 x, \Delta y=\Delta y + \Delta^2 y, \Delta z=\Delta z + \Delta^2 z$ 

12: go subroutine coder {  $\Delta^2 x, \Delta^2 y, \Delta^2 z, \Delta t$  }

13: if memory is full then STOP else GOTO 14

14:  $\Delta t=1$ 

15: GOTO 4

```

where

$x(i), y(i), z(i)$ = sampled, 8 bit precision, EKG data

$\Delta x, \Delta y, \Delta z$ = first difference { $\Delta x(n)=x(n)-x(n-1)$ }

$\Delta^2 x, \Delta^2 y, \Delta^2 z$ = second difference { $\Delta^2 x(n)=\Delta x(n)-\Delta x(n-1)$ }

Δt = time difference between nonzero Δ^2 values

$\tilde{x}(i+1), \tilde{y}(i+1), \tilde{z}(i+1)$ = next predicted data points

Fig 8. Tolan Collection and Decorrelation Algorithm.

codewords are prefix codes as shown in the tree diagram of Figure 9. To decode the sequence correctly, however, it is imperative that codeword synchronization be maintained in order to determine when the coded Δ^2 values end and the uncoded, 7 bit Δt variable starts. This loss of synchronization is detected by the "error" state shown in the code tree of Figure 8. The decoder which implements the code tree in Figure 9 will not be discussed here but is listed in Appendix C.

In contrast with the Huffman code (see chapter 2 and Ref 16) and the Dower code to be discussed, the Tolan code is not constructed using the apriori knowledge of the source word (i.e. $\Delta^2 x, \Delta^2 y, \Delta^2 z$) relative frequency of occurrence. This means that the Tolan code will only produce bit compression (bits out/bits in < 1) if the Δ^2 values are sharply peaked around a mean of zero such that few code words exceed the 8 bit, fixed length value of the Δ^2 terms. As will be shown in chapter 5, the second order difference decorrelator does produce such a sharply peaked relative distribution.

The algorithm which implements the Tolan variable length coder is shown in Figure 10. This algorithm works in conjunction with the Tolan decorrelator in Figure 8.


```

1: entry ( $\Delta^2x, \Delta^2y, \Delta^2z, \Delta t$ )
2: cnt  $\leftarrow$  0
3: reset current memory bit, increment bit pointers
4: if end of memory, set eom flag and RETURN else GOTO 5
5: cnt  $\leftarrow$  cnt + 1
6: if cnt=1 then tvar  $\leftarrow$   $\Delta^2x$ 
7: if cnt=2 then tvar  $\leftarrow$   $\Delta^2y$ 
8: if cnt=3 then tvar  $\leftarrow$   $\Delta^2z$ 
9: if cnt  $\geq$  4 then GOTO 18 else GOTO 10
10: if tvar  $\geq$  0 then GOTO 11 else GOTO 12
11: reset memory bit, increment bit counters and GOTO 13
12: set memory bit, increment bit counters
13: if end of memory, set eom flag and RETURN else GOTO 14
14: if tvar=0 then GOTO 5 else GOTO 15
15: set memory bit, increment bit counters
16: if end of memory, set eom flag and RETURN else GOTO 17
17: tvar  $\leftarrow$  tvar - 1 and GOTO 14
17: store 7 bit  $\Delta t$  counter to memory, update bit pointers
18: if end of memory, set eom flag and RETURN else RETURN

```

Fig 10. Tolan Variable Length Encoder Algorithm.

```

0: START OF DOWER DATA COLLECTION AND DECORRELATION

1: FLAG=TIM=DT=DX1=DY1=DZ1=DX2=DY2=DZ2=DDX=DDY=DDZ=CNT=0

2: If sample interrupt detected then GOTO 3 else GOTO 2

3: X(n) ← A/D CH 0, Y(n) ← A/D CH 1, Z(n) ← A/D CH 2

4: CNT=CNT + 1

5: If CNT =1 then GOTO 6 else GOTO 7

6: X(n-1) ← X(n), Y(n-1) ← Y(n), Z(n-1) ← Z(n) and GOTO 2

7: If CNT =2 then GOTO 8 else GOTO 13

8: DX1 ← X(n)-X(n-1), DY1 ← Y(n)-Y(n-1), DZ1 ← Z(n)-Z(n-1), DT ← DT+1

9: If DX1≠0 OR DY1≠0 OR DZ1≠0 OR DT > 127 then GOTO 11 else GOTO 10

10: CNT ← 1 and GOTO 2

11: TIM ← DT-63, X(n-2) ← X(n-1), Y(n-2) ← Y(n-1), Z(n-2) ← Z(n-1)

12: X(n-1) ← X(n), Y(n-1) ← Y(n), Z(n-1) ← Z(n) DT ←0 and GOTO 2

13: DX2 ← X(n)-X(n-1), DY2 ← Y(n)-Y(n-1), DZ2 ← Z(n)-Z(n-1), DT ← DT+1

14: If DX2≠0 OR DY2≠0 OR DZ2≠0 OR DT > 127 then GOTO 16 else GOTO 15

15: CNT ← 2 and GOTO 2

16: DDX ← DX2-DX1, DDY ← DY2-DY1, DDZ ← DZ2-DZ1

17: If |DDX| > 63 then GOTO 18 else GOTO 21

18: DUMX(n)=.25X(n-1)+.75X(n)

    DUMX(n-1)=.25X(n-2)+.5X(n-1)+.25X(n)

    DUMX(n-2)=.75X(n-2)+.25X(n-1)

19: X(n) ← DUMX(n), X(n-1) ← DUMX(n-1), X(n-2) ← DUMX(n-2)

20: DX1=X(n-1)-X(n-2), DX2=X(n)-X(n-1) and GOTO 16

```

Fig 11-a. Dower Collection and Decorrelation Algorithm.

```

21: If |DDY| > 63 then GOTO 22 else GOTO 25
22: DUMY(n)=.25Y(n-1)+.75Y(n)
    DUMY(n-1)=.25Y(n-2)+.5Y(n-1)+.25Y(n)
    DUMY(n-2)=.75Y(n-2)+.25Y(n-1)
23: Y(n) ← DUMY(n), Y(n-1) ← DUMY(n-1), Y(n-2) ← DUMY(n-2)
24: DY1=Y(n-1)-Y(n-2), DY2=Y(n)-Y(n-1) and GOTO 16
25: If |DDZ| > 63 then GOTO 26 else GOTO 29
26: DUMZ(n)=.25Z(n-1)+.75Z(n)
    DUMZ(n-1)=.25Z(n-2)+.5Z(n-1)+.25Z(n)
    DUMZ(n-2)=.75Z(n-2)+.25Z(n-1)
27: Z(n) ← DUMZ(n), Z(n-1) ← DUMZ(n-1), Z(n-2) ← DUMZ(n-2)
28: DZ1=Z(n-1)-Z(n-2), DZ2=Z(n)-Z(n-1) and GOTO 16
29: If FLAG=0 then GOTO 30 else GOTO 31
30: VLC ← (0:DX1,DY1,DZ1)
31: VLC ← (TDT:DDX,DDY,DDZ)
32: If MEMORY FULL then STOP else GOTO 33
33: X(n-2) ← X(n-1), Y(n-2) ← Y(n-1), Z(n-2) ← Z(n-1)
34: X(n-1) ← X(n), Y(n-1) ← Y(n), Z(n-1) ← Z(n)
35: TDT ← DT-63, FLAG=1 and GOTO 2

```

Fig 11-b. Dower Collection and Decorrelation Algorithm.

Dower EKG Data Compressor

The Dower EKG compression technique, like the Tolan method, is a redundancy reduction procedure. The Dower compressor combines a zero order time compression operation with a second order difference reduction transformation to produce a decorrelated residual frame sequence. This residual frame sequence is then compressed by a specially tailored variable length code whose average codeword bit length approaches the entropy bound without the buffer overflow problem encountered with the ''optimal'' Huffman code.




Fig 12. Dower EKG Data Compression System.

Dower Decorrelator. The Dower data collection and decorrelation operation (Figure 11) begins with a zero order time compression process where only those EKG sample values (x,y,z) which differ from the previous sample (Δx or Δy or $\Delta z \neq 0$) are saved along with a run length counter Δt . This operation produces a sequence of data frames ($\Delta t: x,y,z$). Following the creation of the data frames, a second difference with respect to frame number is performed generating a sequence of second order difference frames ($\Delta t: \Delta_f^2 x, \Delta_f^2 y, \Delta_f^2 z$).

The Dower algorithm operates on 8 bit data, hence there are 256 potential source symbols for each lead. When second differences are taken, the range increases to 1024 potential Δ^2 symbols. By experimental evidence, Dower and Berghofer (Ref 12) have found that $0 \pm 63 \Delta^2$ values are sufficient to reproduce all but the fastest EKG artifacts (e.g., pacemaker spikes). The variable length encoder, therefore, is designed to expect 127 source symbols (values) and no more. To insure the Δ^2 dynamic range of 0 ± 63 is not exceeded, the Δ^2 values are limited by a preprocessor shown in steps 18-26 of Figure 11. This preprocessor is iterative, and irreversably modifies the three sample points which produced the $|\Delta^2| > 63$ until the Δ^2 value falls within the encoder range. With care taken to record the necessary initial conditions (first x,y,z and first $\Delta x, \Delta y, \Delta z$), this second difference frame sequence now contains all the significant information in the original sample sequence with

reduced interframe correlation.

The mechanics of the lower decorrelator are best illustrated by an example. Figure 13 shows three hypothetical sample sequences. By applying the rules in the preceeding paragraph, a set of zero order data frames ($\Delta t: x, y, z$) is formed. That is:

(0:0,-1,2), (0:1,1,1), (0:3,2,3), (0:4,3,1), (3:3,1,-3), (0:1,0-1),
(0:-1,-1,-2), (0:-2,-2,-1), (4:-3,1,-1)

To complete the correlation reduction process, the second difference with respect to frame number must now be performed. This results in the set of second difference frames ($\Delta t: \Delta_f^2 x, \Delta_f^2 y, \Delta_f^2 z$) shown below.

(0:0,-1,2), (0:1,2,-1) | (0:1,-1,3), (0:-1,0,0), (3:-2,-3,-2),
(0:-1,1,6), (0:0,0,-3), (0:1,0,2), (4:0,4,-1)

The data sets preceeding the verticle bar are the initia
conditions necessary for reconstruction. The first data set is the first zero order time compression sample frame. The second frame is the first difference between the first two zero order data frames.




Fig 13. Three hypothetical sample sequences.

Dower Entropy Encoder. The Dower entropy encoder also uses variable length codes for the second stage of the FKG data compressor. The Dower code, however, is significantly different than the prefix code used in the Tolan compression method.

As was mentioned previously, the Dower VLC is configured for 127 source symbols (8 bit Δ^2 values). Unlike the Tolan code, however, the Dower VLC does not perform a 1:1 mapping between a single source symbol and a single code symbol.

In the Dower coder, the source symbols (Δ^2_i) are mapped onto the state space of a 14 bit accumulator (A). This state space is partitioned into symbol regions R_i which are assigned according to the probability of occurrence of the source symbols. The size ($||R_i||$) of the symbol region R_i is given by the relation

$$||R_i|| = ds_i = s_{i+1} - s_i = \text{Pr}(\Delta^2_i) 16384 \quad (24)$$

where s_i is the initial state of R_i and s_{i+1} is the initial state of R_{i+1} . The initial state, s_i , is also determined by symbol probability. The symbol regions R_i corresponding to those symbols which occur least often are assigned to the low end of the state space range such that if:

$$\text{Pr}(\Delta^2_k) > \text{Pr}(\Delta^2_n) > \text{Pr}(\Delta^2_m) > \text{Pr}(\Delta^2_j)$$

then

$$||R_k|| > ||R_n|| > ||R_m|| > ||R_j||$$

and

$$s_k > s_n > s_m > s_j$$

This relationship is illustrated in Figure 14. For those potential symbols Δ_i^2 which do not occur, 1 state is always assigned.




Fig 14. Example of State Space Partition for 4 Symbols.

The physical configuration of the Dower encoder is shown in Figure 15. As can be seen from Figure 15, the memory storage buffer for the encoded data is contiguous with accumulator and all accumulator bit shifts (right or left) also shift the entire memory buffer. Right shifting operations always shift a binary 0 into the most significant bit (MSB) of A.




Fig 15. Dower Encoder Accumulator-Memory Buffer Interface.

A source symbol Δ_i^2 is encoded by adding the value s_i to the current contents of the accumulator. If the value in the accumulator is such that $s_i + A \geq 2^{14}$, then an arithmetic overflow condition would occur with a corresponding loss of data. To insure this overflow does not occur, the accumulator is first right shifted k_i times into memory. The value k_i is the maximum number of right shifts which would ever be necessary and can be found from the inequality:

$$2^{14} \leq 2^{k_i}(s_{i+1} - s_i) < 2^{15} . \quad (25)$$

The k_i values are tabulated in Table I for various values of $ds = s_{i+1} - s_i$. This right shifting diminishes the value in the accumulator (i.e., divides A by 2 per right shift) until the value in A can be safely mapped into a region R_i . The initial state, s_i can now be added to A without overflow occurring.

In some instances, the maximum shift k_i would not be

necessary (i.e., $k_i - 1$ right shifts would still prevent overflow). To determine when this situation occurs, a third variable y_i is defined where

$$y_i = ds - 2^{14 - k_i} \quad (26)$$

represents the excess states by which ds exceeds the next lowest integral power of 2.

Table I

Dower Code Variables Per Accumulator Partition Size

Partition Size	k	y
$ds=1$	14	0
$1 < ds < 4$	13	$ds-2$
$4 < ds < 8$	12	$ds-4$
$8 < ds < 12$	11	$ds-8$
$16 < ds < 32$	10	$ds-16$
$32 < ds < 64$	9	$ds-32$
$64 < ds < 128$	8	$ds-64$
$128 < ds < 256$	7	$ds-128$
$256 < ds < 512$	6	$ds-256$
$512 < ds < 1024$	5	$ds-512$
$1024 < ds < 2048$	4	$ds-1024$
$2048 < ds < 4096$	3	$ds-2048$
$4096 < ds < 8192$	2	$ds-4096$
$8192 < ds < 16384$	1	$ds-8192$

If the value in $A < y_i$ after the initial k_i shifts, then one too many right shifts occurred. A is then left shifted once and s_i added to A without fear of overflow. If $A \geq y_i$ then a full k_i right shifts were required. To insure maximum efficiency in the next encoding, however, the position of the accumulator in state space map should be

right justified within the current symbol region. This right justification is obtained by adding $y_i + s_i$ to A whenever $A \geq y_i$.

By use of the variable y_i maximum efficiency in the number of right shifts (i.e. information bits stored to memory) can be obtained. Dower and Berghofer (Ref 12) assert that, on the average, the number of right shifts per source symbol (Δ_i^2) will approach quite closely to the value $-\log_2 \Pr(\Delta_i^2)$ which is the self information 'content' of the Δ_i^2 value (see appendix B).

The algorithm which implements the Dower encoder is shown in Figure 15. This algorithm assumes that the variable tables for $\Delta_{fx}^2, \Delta_{fy}^2, \Delta_{fz}^2$ and Δt have been calculated from the second difference histograms.

The Dower decoding operation is just the inverse of the encoding process. For decoding, the current state of the accumulator is mapped into the state space table of the current frame variable ($\Delta t, \Delta_{fx}^2, \Delta_{fy}^2, \Delta_{fz}^2$). Once the correct region R_i is determined, then s_i is subtracted from A . If $A < 2y_i$ then one right shift is performed followed by k_i left shifts. If $A \geq 2y_i$, then y_i is subtracted from A followed by k_i left shifts. For a more detailed example of the Dower encoding/decoding operation, the reader is referred to Dower and Berghofer (Ref 12).

Although conceptually more difficult than the Tolan VLC, the algorithmic structure of the Dower entropy encoder is only slightly more complex. As was shown above, the

Dower VLC can be implemented solely with shifting and table look up operations. The Dower method does, however, require that code tables (s_i, k_i, y_i) be constructed prior to the encoding process.

```

1: Entry  $(\Delta t, \Delta_f^2 X, \Delta_f^2 Y, \Delta_f^2 Z)$ 
2: IF this is the first frame stored then  $\text{Acc} \leftarrow 0$  else GOTO 3
3:  $\text{DDVAR} \leftarrow \Delta t, k \leftarrow k(\Delta t), s \leftarrow s(\Delta t), y \leftarrow (\Delta t)$ 
4: RIGHT SHIFT Acc k times
5: IF  $\text{Acc} < y$  then GOTO 8 else GOTO 6
6: LEFT SHIFT Acc once
7:  $\text{Acc} \leftarrow \text{Acc} + s$  and GOTO 9
8:  $\text{Acc} \leftarrow \text{Acc} + s + y$ 
9: IF memory is full, then SET eom flag and RETURN else GOTO 10
10: IF  $\text{DDVAR} = \Delta t$  then GOTO 11 else GOTO 12
11:  $\text{DDVAR} \leftarrow \Delta_f^2 X, k \leftarrow k(\Delta_f^2 X), s \leftarrow s(\Delta_f^2 X), y \leftarrow y(\Delta_f^2 X)$  and GOTO 4
12: IF  $\text{DDVAR} = \Delta_f^2 X$  then GOTO 13 else GOTO 14
13:  $\text{DDVAR} \leftarrow \Delta_f^2 Y, k \leftarrow k(\Delta_f^2 Y), s \leftarrow s(\Delta_f^2 Y), y \leftarrow y(\Delta_f^2 Y)$  and GOTO 4
14: IF  $\text{DDVAR} = \Delta_f^2 Y$  then GOTO 15 else RETURN
15:  $\text{DDVAR} \leftarrow \Delta_f^2 Z, k \leftarrow k(\Delta_f^2 Z), s \leftarrow s(\Delta_f^2 Z), y \leftarrow y(\Delta_f^2 Z)$  and GOTO 4

```

Fig 16. Dower VLC algorithm.

Tolan versus Dower

The similarities and differences between the Tolan and Dower EKG compression techniques are now discussed. The decorrelators are compared first.

Both the Dower and the Tolan decorrelators use a second difference technique where time compression is used to eliminate the storing of the most common Δ^2 value of zero. The Tolan algorithm, however, continuously calculates the second difference and applies these differences to the Tolan VLC. If a long "run" of $\Delta^2 = 0$ values occur, the Tolan encoder overflows at $\Delta t = 128$. This run counter overflow forces a storage "dump" to the VLC with a resultant loss in efficiency.

The Dower decorrelator approaches time compression slightly differently. In the Dower algorithm, only those sample points where the EKG data was changing (i.e., Δx .or. Δy .or. $\Delta z \neq 0$) are saved forming data frames. A second difference with respect to frame number is performed and the Δ_f^2 's and Δt are fed to the VLC. Although this frame methodology appears more efficient than the Tolan technique, the Dower decorrelator is still constrained by a maximum Δt of 127 (i.e., maximum Dower VLC code range). On the basis of this analysis, it appears that both the Tolan and the Dower decorrelators have similar performance.

Since the Tolan and Dower decorrelators appear about equally efficient, the real compression payoff is in the variable length encoders. In the Tolan VLC, code words as

long as 1024+2 bits are feasible. These long codewords would only occur if an extremely large and fast "spike" (e.g. a pacemaker pulse) appeared in the data. In the population as a whole, spikes of this magnitude occur very infrequently. Nevertheless, the Tolan compression efficiency will degrade seriously in an environment where "impulsive" artifacts appear.

The Dower VLC is designed to minimize sensitivity to impulses. In the Dower system, the VLC is configured such that the longest codeword (i.e., the number of right shifts to memory) is 14 bits. Extremely large signal spikes are numerically filtered to reduce their second difference within the 0 ± 63 range. Such drastic limiting action would occur infrequently, however, and not seriously affect the reproduction of the EKG.

As an example of their performance, let a second difference of 45 be encountered by both the Tolan and Dower encoders. The Tolan VLC would require 47 bits to encode this data. The exact codeword size of the Dower routine is dependent on the $\Delta^2 = 45$ probability of occurrence. Nonetheless, the "codeword" is always ≤ 14 bits long; an obvious increase in efficiency over 47 bits.

A calculation of the performance differential between the Tolan and Dower routines will be estimated in chapter 5. For now it is sufficient to say that the Dower EKG compressor outperforms the Tolan method. A short synopsis will now be made of two other EKG compression techniques.

6.1 Microcomputer FKG Compressors.

Two other FKG compression techniques were discovered by this author. The first is referred to as the Turning Point Method and is discussed in reference 32. The second method is currently in use by Marquette Electronics for data compression in a commercial EKG 'cart'. The Turning Point technique is presented first.

Turning Point Algorithm. The turning point algorithm is, by definition, a 2:1 data compressor where one of two consecutive sample points is discarded. The algorithm which determines which sample point is discarded is as follows. The first sample point is stored and assigned as the reference point (X_0). The next two consecutive points become X_1 and X_2 . With 3 sample points there are 8 'patterns' or combinations which reflect the 'trends' in data (see Figure 17). The Turning Point algorithm stores the circled point (Figure 17) which becomes the new reference point X_0 . The point not circled (X_1 or X_2) is discarded. The next two points are sampled, their values are assigned to X_1 and X_2 , and the process repeated.




























PATTERN	x_0	x_1	x_2
1			
2			
3			
4			
5			
6			
7			
8			
9			

Fig 17. Turning Point Patterns (from Ref 32:6.61).

It can be shown (Ref 32:6.59-6.65) that :

$$\text{if } (X_2 - X_1) * (X_1 - X_0) < 0 \quad X_0 = X_1$$

$$\text{if } (X_2 - X_1) * (X_1 - X_0) \geq 0 \quad X_0 = X_2$$

hence the sign of the product of consecutive first differences determines the ''significant'' point to be saved.

The advantage of the Turning Point compressor is speed of execution. Execution of this algorithm on even the slowest microprocessor would pose no problem at all. The disadvantages of the turning point routine are numerous, however. First the routine discards data so that reconstruction of the original sample sequence is not possible. In this regard, the Turning Point algorithm falls in the class of Entropy Reducing (ER) techniques as was discussed in chapter 2. Second no attempt is made to use the probabilistic distribution of the EKG to enhance data compression as is done with the Tolan and Dower routines. This obviously leads to inefficiency. Finally, this routine only produces a 2:1 compression ratio; incredibly poor in relation to the other techniques already discussed in this thesis. Only where the simplest technique is necessary, would the Turning Point technique be beneficial.

Marquette Algorithm. The Marquette algorithm was developed by Marquette Electronics, Milwaukee, Wisconsin. This is similar to the Dower and Tolan algorithms in that variable length encoding is used to compress the output of a

"difference" decorrelator.

In contrast to the Dower and Tolan techniques, the Marquette compressor only calculates the first difference for input into the VLC. The Marquette variable length encoder stores the first difference data as 1,3,5, or 7 nibbles where a nibble is defined as 4 bits. These nibble codes are arranged as follows:

<u>Range of Difference</u>	<u>Code Length in Nibbles</u>
(-7,7)	1
(-127,127)	3
(-2047,2047)	5
(-32767,32767)	7

To encode the first differences, and delineate between code words, the following rules apply:

- 1) Differences must be coded on the smallest possible range, and attempts to encode a difference of +5, for example, using more than 1 nibble will result in a decoding error.
- 2) Single nibble codes are difference plus 8. A nibble value of zero does not occur.
- 3) Three nibble codes start with a single zero nibble. The remaining two nibbles are obtained as follows:
 - (1) Positive Differences +8
 - (2) Negative Differences +7
- 4) Five nibble codes start with two zero nibbles. The remaining three nibbles are obtained as follows:
 - (1) Positive Differences +128
 - (2) Negative Differences +127
- 5) Seven nibble codes start with three zero nibbles. The remaining four nibbles are obtained as follows:

- (1) Positive Differences +2048
- (2) Negative Differences +2047

6) The sequence of nibbles in a code starts with the zero flag nibbles (if any) followed by the most significant through least significant nibble.

The above set of encoding rules were obtained from Mr. Tom Divers, Marquette Electronics project engineer (Ref 11).

From their own analysis, Marquette has shown (Ref 11) that at a 250 Hz sampling rate, 89.1 percent of the first differences fall within the ± 7 range with 99.8 percent falling within a ± 127 range. Marquette reports that at an A/D precision of 10 bits, an average of 4.89 bits/sample (across the total EKG population) is obtained with their compression routine.

The Marquette EKG compression appears to work well, even with a first order difference correlation reducer. The Marquette variable length encoder, however, is tailored to the 8 bit ASCII data communications environment and is not "optimum" in any sense. The Marquett VLC does perform "exact" redundancy reduction entropy compression with sufficient efficiency to make this encoding scheme commercially viable.

Chapter III Summary.

This chapter has looked in detail at two EKG data compression techniques which perform "exact" redundancy reduction. The Tolan routine, which was implemented by this author (see chapters 4 and 5) decorrelated the sample

sequence data by performing a second order time compression operation. The residual sequence resulting from the decorrelator was ''compressed'' by a uniquely decodable variable length code which was shown to be ''suboptimal'' with respect to the ''optimal'' Huffman code.

The Dower compression technique also performed a second difference operation, but preceeded the second order ''differencer'' by a zero order time compressor. The Dower zero order time compressor produced a sequence of ''data frames'' which in turn were converted to a sequence of second order difference frames with respect to frame number. Since the Dower VLC limited the Δt time compression counter to a maximum of 127 (as did the Tolan VLC), input symbols, the Dower and Tolan decorrelators were considered to perform equally well. The Dower entropy encoder was also shown to be a variable length coding operation but not a ''prefix'' code as was the Tolan VLC. The Dower VLC maps decorrelator ''symbols'' into the ''state space'' of 14 bit accumulator which encodes data by adding the initial address of a symbol's state space region to the accumulator. To prevent accumulator overflow, the accumulator data is shifted out to a memory storage buffer. The ''number of shifts'' necessary to prevent accumulator overflow represent the codeword size and it was shown that the Dower VLC approached the ''entropy'' bound of the decorrelated input sequence.

The chapter was concluded by a synopsis of two other

EKG compression techniques. The first of these two, the Turning Point algorithm, was shown to be of marginal use because it is not an 'exact' technique and produces compression of only 2:1. The second technique, however, was the Marquette compression system and it was shown to produce acceptable compression worthy of commercial application.

From the available algorithms, the Dower compression technique has the capacity to produce the best 'exact' compression of any of the techniques studied in this chapter. An interesting experiment would be the combining of the second order interpolator (discussed in chapter 2) used by Ruttiman and Pipberger (Ref 28) and the variable length encoder used by Dower (Ref 12). This combination should prove to be very powerful and effective and is left for further study.

The next chapter in this thesis discusses the configuration of the EKG Data Acquisition and Analysis System assembled by this author to test the Tolan EKG compression algorithm.

IV. EKG-Data Acquisition and Analysis System

The EKG-Data Acquisition and Analysis System (EKG-DAAS) was assembled for this thesis as the testbed on which experimental EKG data could be acquired, compressed, analyzed, stored, and reconstructed. The EKG-DAAS design can be separated into the categories of hardware and software. The hardware is discussed first.

EKG-DAAS Hardware.

The EKG-DAAS was constructed around the Motorola Exorciser microcomputer (appendix E) which uses a 6800 microprocessor for its central processing unit (CPU). In the EKG-DAAS, the Exorciser is configured with 32 kilobytes (K) of read/write (RAM) memory and 16 K of read only memory (ROM). In addition, the Exorciser was equipped with the EXBUG debugging module which allowed interactive program debugging with preselectable software breakpoints, execution tracing modes, and CPU register display.

To provide extended memory, a Midwest Scientific Instruments (MSI) FD-8 Disk Memory unit was interfaced to the Exorciser and provides approximately 290 K of online user memory. The FD-8 is accessed by a MSI Disk Operating System (DOS) and communicates with the CPU via a MEX6820 Input/Output Module installed in the Exorciser chassis.

Terminal input/output (I/O) is accomplished by means of a standard RS-232 serial interface which has switch selectable baud rate from 110 to 9600 baud. In the EKG-DAAS configuration, the Exorciser serial ''port'' was connected in parallel with a Heathkit H-14 dot matrix line printer for program listings and data printouts.

Data I/O was accomplished by means of a Sinetrac ST-6800 Analog/Digital-Digital/Analog (A/D-D-A) converter module which samples and digitizes analog data to 12 bit precision. The ST-6800 has the capability of sampling 32 distinct analog channels (A/D) as well as output 2 channels (D/A) with simple memory addressed LDA (load) and STA (store) instructions. For the EKG-DAAS, the ST-6800 was setup for ± 5 volt, 2's compliment data and was addressed (A/D ch 0) at E400 Hexadecimal (Hex). The internal configuration of the Exorciser as used in the EKG-DAAS is illustrated in Figure 18.

To uniformly sample the EKG input data via the ST-6800 required the use of an external interrupt timer as is shown in Figure 19. This timer allowed data sampling rates between 300 and 700 Hertz but for the duration of this research was set, and calibrated, at 500 Hertz. The interrupt was interfaced to the Exorciser via an interrupt line on the ST-6800.

Considerable problems arose in this thesis due to hardware problems associated with the FD-8 Disk Memory. The original configuration of the EKG-DAAS used two FD-8 systems




Fig 18. EXORCISER Component Module Layout (From Ref 22).

but one failed about midway through the software development. This failure caused a major rewrite of the thesis software and destroyed several weeks of work. Although the EKG-DAAS can now operate in a one disk environment, considerable "manhandling" of the data diskettes is necessary.




Fig 19. EKG-DAAS Hardware Configuration

Further information concerning the Exorciser's hardware capability can be obtained in appendix E and reference 22. Attention is now turned to the EKG-DAAS software.

EKG-DAAS Software.

The EKG-DAAS software was written in 6800 assembly language and controls all aspects of terminal, disk memory, and A/D-D/A operation. The EKG-DAAS programs consist of approximately 4300 lines of assembly language and are listed in appendix C.

The software used in the EKG-DAAS was written in 6800

assembly language for two reasons: 1) no high order language (e.g., PASCAL, FORTRAN) is available on the Exerciser and 2) speed limitations imposed by the thesis requirement for online, real time EKG data compression made it imperative that the compression programs run as fast as possible. The laborious task of writing and testing assembly language slowed software development to the point where only one EKG data compression/reconstruction routine (Tolan) was completed.

The EKG-DAAS software is integrally tied to the MSI DOS (Ref 23). All disk I/O operations initiated by the EKG-DAAS routines flow through the MSI-DOS and hence the DOS must be "live" somewhere in memory. To insure that the DOS routines are always available, the DOS was disassembled and relocated in high memory ROM (C400 Hex).

The basic flow of EKG-DAAS program control is illustrated in Figure 20. The EKG-DAAS software is broken into overlaid modules which are called into memory and executed by EKG-EXEC and DISPLAY. The basic memory map and overlay structure is illustrated in Figure 21.

As can be seen in Figure 21, extensive memory management was required in order to allow a sufficiently large memory buffer for the EKG data. As configured in Figure 21, the EKG-DAAS could collect 11.6 seconds of uncompressed (3 leads, 8 bits/lead, 500 samples/sec) EKG data. For compression with the Tolan algorithm, a maximum of 26.2 seconds of data (TA1359PA, appendix D) was




Fig 20 EKG-DAAS Software Control Flowgraph.




Fig 21. EKG-DAAS Overlay Structure and Memory Map.

collected.

With the basic structure of the EKG-DAAS software defined, this chapter will now examine the components of the EKG-DAAS software in more detail. In the next section of this chapter, a software module will be described along with a simplified flowchart of that module's operation.

EKG-EXEC. EKG-EXEC is the executive command module which controls the execution flow of the EKG-LAB system. Upon input of a command number, EKG-EXEC loads the appropriate overlay routine into the program work buffer (0100-1000 Hex) and then passes control to that overlay. The above command and control operation is illustrated in Figure 22.

In addition to the command 'handler', EKG-EXEC contains the utility subroutines FILHDR, SAVFIL, HXASC, OVLAY, and PDFPRT. These subroutines are described as follows:

FILHDR. FILHDR clears the memory data buffer, initializes the statistics buffer variables, and queries the console for data such as FILENAME, SUBJECT, DATE, etc.

SAVFIL. SAVFIL reads the filename in the memory buffer header and then writes the memory file to disk. Disk I/O is passed through EOS subroutines.

HXASC. This subroutine converts hexadecimal data to ASCII for display on the terminal and printer devices.

OVLAY. OVLAY is the routine which actually performs the overlay function. After an overlay is loaded into memory, OVLAY jumps program control to the overlay program.

PDFPRT. PDFPRT prints the amplitude distribution to the terminal device (printer). Although resident in EKG-EXEC, PDFPRT is called only by the DISPLAY module.




Fig 22. EKG-EXEC Functional Flowchart

NOCPRS. NOCPRS is a data acquisition module in which the EKG waveform is sampled and stored without data compression. The data is rounded to 8 bits from 12 bits. The reason for this rounding is explained in chapter 5.

This module was constructed for two reasons : 1) uncompressed data was considered useful for doing experimental studies on potential data compression techniques implemented after the original data collection session and ; 2) this module was the structure around which the Tolan (and potentially other) compression routines were built. The basic operation of NOCPRS is illustrated in Figure 23.

As is seen in Figure 23, NOCPRS does more than just sample the EKG. Statistical parameters are collected and updated during an EKG data collection. These parameters are used to measure the real time performance of the compression (no compression) software. A detailed description of these measurement parameters is described in chapter 5.




Fig 23. NOCPRS Functional Flowchart.

DISPLAY. DISPLAY is a command module similar in construction to EKG-EXEC. DISPLAY controls the data output modules which display the sampled EKG data and data statistics to two output devices. The first device is the terminal (printer) and the second device is an oscilloscope. The basic command structure of DISPLAY is shown in Figure 24. DISPLAY is broken into 5 working submodules. These modules are described as follows:

PRSTAT. This submodule reads the memory file header (3C00-3D00 Hex), formats the statistical data found there, and prints this data to the terminal (printer). The statistical data in appendix D was generated by PRSTAT.

PDFPRT. This submodule prints the memory file histogram tables to the terminal (printer). PDFPRT output is also listed in appendix D.

PDFOUT. PDFOUT scans a user selected lead histogram (X,Y,Z) and formats the data for display to an oscilloscope. The data is output via D/A ch 0. An example of PDFOUT output is found in Figure 31 in chapter 5.

DECPRS. The DECPRS module scans the memory file header and identifies the compression technique which was used to encode the data in memory. The appropriate decompression algorithm is then selected and the data decoded and output on D/A channel 0.

LOAD. LOAD initiates a data file load from disk memory to RAM. This load is performed by EOS routines called by load.




Fig 24. DISPLAY Functional Flowchart.

TOLAN. The TOLAN module is the ''heart'' of the EKG-DAAS and performs the actual data compression on the sampled EKG. As was described in chapter 3, the TOLAN algorithm first performs a second difference time compression operation followed by a variable length encoder. The operation of the TOLAN compression module is shown in Figure 25.

To detect sample clock (or CPU clock) drift, time calibration operations are performed prior to and after the data collection run. Other statistical data parameters are also collected allowing post collection measurement of compression performance.




Fig 25a. TOLAN Functional Flowchart.




Fig 25b. TOLAN Functional Flowchart

Summary

This chapter has presented an overview of the complex software which makes up the EKG-DAAS. The reader is referred to the program listings in appendix C for a more thorough description of the program operation.

The software structure in the EKG-DAAS was written in a "Top Down" manner and all attempts have been made to document the operation of each routine. Since assembly language is difficult to read, this chapter was written to assist the reader in understanding the basic structure of the EKG-DAAS. The next chapter presents the results of the EKG experiment where data was collected, analyzed, and compressed by the EKG-DAAS.

V. Experimental Procedure, Data Analysis, and Results

This chapter presents the results of an EKG collection experiment where ''in vivo'' EKG data was taken from test subjects in real time. Data was taken and stored in both compressed and uncompressed formats for later analysis and reconstruction.

Chapter 5 is organized as follows. First the ''experimental'' setup is described along with a description of the EKG equipment, collection environment, and subject personnel. Next the parameters used to determine compression performance are defined followed by the analysis and results of the experimental data. Finally, the chapter concludes with a comparison between the results obtained using the Tolan compression algorithm and the estimated performance of the Dower technique. Discussion now turns to the experimental procedure.

Experimental Procedure

The data was taken from a set of nine fellow students during a laboratory course on electrocardiograms. The equipment used to produce the EKG was the model DR-12 research recorder built by Electronics for Medicine, Inc. (see appendix E). The DR-12 is a vintage medical recording system built in the late 1950's and is constructed with vacuum tube amplification circuitry.

The personnel used for test subject were all Air Force

officers in good physical health. No test subjects with obvious heart disease were used although significant ECG variations between subjects was noted (Fig 27-29). Output of the DR-12 was limited to one EKG signal which could be switched to any of the 6 "limb leads" (Ref 13:29-34) by controls on the DR-12. The electrodes of the EKG were applied to the wrists of the test subjects and in some cases, not all, an electrode jelly was applied to reduce skin-electrode resistance.

Since the EKG-Data Acquisition and Analysis System (EKG-DAAS) was configured for a 3 lead system, the X,Y,Z inputs were connected in common and the single signal available from the DR-11 applied to this connection. A Brush Instruments Mark II recorder and a Tektronics Model 465M Oscilloscope were used as analog output devices (Appendix E). The display instruments were connected in common with the A/D inputs as is shown in Figure 26.

Prior to the data recording session, the A/D was calibrated in accordance with the operating manual (Ref 8). The A/D was configured for a dynamic range of of -5.000 volts to +4.9976 volts with 2's complement binary representation.

The actual data collection proceeded as follows. First the subject was connected to the DR-12 and the amplitude of the resulting EKG signal adjusted to fall within the A/D dynamic range. A test run of all six limb leads (I,II,III,AVL,AVR,AVF) was then taken (without storage by




Fig 26. Experimental Data Collection Setup.

EKG-DAAS) and the lead with the least 'noisy' signal was selected. EKG-DAAS was then executed and a data collection made without compression. This uncompressed data collection was followed as soon as possible with another 'run' in which the Tolan compression algorithm was enabled. Following both data collections (uncompressed and compressed), the raw data traces from the Brush recorder were annotated with the time and subject and filed for later data comparison with the reconstructed waveform.

The data was stored on 'floppy' diskettes and processed post collection for the entropy and maximum compression statistics. Before the results of these data

collections are presented, the compression measurement parameters will be defined.

FKG Compression Measurement Parameters

To permit determination of compression performance, a set of statistical parameters was calculated and saved during each data collection. These statistical parameters are described as follows:

Number of Samples. This statistic was saved to determine the total number of bits that were input to the compression stages. The total number of bits were calculated by $(8 \text{ bits/sample}) * 3 \text{ leads} * (\text{num of samples/lead})$.

2nd Difference Frequency of Occurrence. Four frequencies of occurrence tables were kept with double precision binary counters. Following the data collections, these $\Delta^2x, \Delta^2y, \Delta^2z, \Delta t$ histogram tables were input to a BASIC program (ENTROPY) where the entropy of second difference "source" was calculated.

Total Waiting Loop Counts During Collection. A counting loop was established in the TOLAN compression module which allowed determination of the percent of the sampling period used for the compression and statistics calculations. One circuit of this counting loop takes 46 machine cycles of the 6800 microprocessor. A count of the total number of loop cycles completed following the sampling/compression interrupts is kept in the collection statistics buffer.

Maximum Loop Count Per Interrupt. To offset the inaccuracy which would develop if the interrupt clock period changed between runs (or if the master clock in the Exerciser drifted), a loop count calibration was performed immediately before and after each collection run. This calibration was accomplished by performing 256 sequential interrupts with no interrupt processing except return-from-interrupt. The before and after calibration counts were then averaged and the maximum loop counts per interrupt calculated. Time Efficiency of the compression operation was

AD-A100 799


AIR FORCE INST OF TECH WRIGHT-PATTERSON AFB OH SCHOO--ETC F/G 6/5
ANALYSIS AND PERFORMANCE EVALUATION OF ELECTROCARDIOGRAM DATA C--ETC(U)
DEC 80 M D TOWNSEND
AFIT/GE/EE/80D-46

UNCLASSIFIED

NL

2 OF 3

AD A
100799



then calculated from the equation:

$$T.E. = (1 - (\text{Total Loop Count} / ((\text{Num of Samples}) * (\text{Max Cnt})))) * 100$$

Channel Maximums and Minimums. The channel maximums and minimums were retained to allow determination if the analog inputs exceeded the A/D dynamic range.

Number of Memory Bits Available. This number was constant and was determined by the amount of Read/Write (RAM) memory available for data storage. For the current configuration of the Exorciser and the EKG-EXEC program this was 139248 bits (17406 bytes).

Number of Bits Available to Variable Length Coder. This counter measured the number of bits out of the data decorrelator and allowed calculation of the decorrelator's compression ratio (bits out/bits in).

Number of Bits Used to Store Channel X,Y,Z. These counters measured the number of bits used to store the data from the three input leads. This count is the number of code bits out of the variable length encoder.

Number of Bits Used to Store Time. This counter was identical to the channel counters above but measured the number of code bits used to store the Δt run counts.

Total Compression Ratio Achieved. This figure was calculated post collection by dividing the total code bits stored by the total R/W memory bits available.

The statistical data defined above was compiled by the EKG-EXEC program and is printed by the DISPLAY software module as illustrated in appendix D.




Fig 27. Original and Reconstructed EKG with
Best Compression Ratio (2.260 : 1)




Fig 28. Original and Renconstructed EKG with
Average Compression Ratio (1.60 : 1)




Fig 29 . Original and Reconstructed EKG with
Worst Compression Ratio (1.252 : 1)

Data Analysis and Results

To determine the compression efficiency of a redundancy reduction EKG compression algorithm is, in general, a difficult job to accomplish. The reasons for this difficulty are: 1) the sampled EKG data sequence contains significant correlation (memory) between sample values making calculation of the absolute bound of the entropy extremely difficult (Ref 34:479-489) and ; 2) the techniques used for both the data decorrelator and entropy encoder vary significantly from algorithm to algorithm.

In this thesis, as is done in the literature (Ref 7,12,28), the assumption is made that the output of the decorrelator is ''almost decorrelated''. Decorrelated is not ''independent'' (unless the source was statistically gaussian), but it is assumed that true entropy of the 2nd difference sequence approaches the value which would be calculated by Eq.(1) (reproduced below).

$$\Delta^2 \text{ entropy} \approx - \sum -p_i \log_2 p_i \quad (27)$$

(p_i is the probability of the i 'th second difference). This second difference entropy can then be used as an upper bound on the potential entropy encoding compression of the Tolan redundancy reduction technique.

The Δ^2 entropy values tabulated in Table II were calculated by Eq.(27). To calculate the approximate limit on the entropy encoder's compression ratio, the uncompressed data word length of 8 bits was divided by the ''lowest''

average code word (i.e., the 2nd difference entropy). A look at Table II shows that the entropy encoder compression ratio varied between 58% and 71% of this entropy bound. Since the Tolan variable length code is suboptimal, a lower efficiency is expected. Nevertheless, the Tolan entropy encoder performed more-or-less consistently across the data set. This last observation would imply that the data decorrelator influences the overall compression ratio more

Table II
Experimental Data Summary

Subject Id	2'nd Difference Entropy	Maximum Compression Possible (Approx)	Achieved 2nd Difference Compression	Percent of Max Encoder Compr	Percent of Sample Interval	Achieved Total Compr
TA1545T	3.2801	2.42 : 1	1.55 : 1	64.0%	93.9%	1.36 : 1
TA1548T	3.0171	2.63 : 1	1.61 : 1	61.2%	91.4%	1.50 : 1
TA1559B	3.016	2.64 : 1	1.72 : 1	65.1%	93.0%	1.53 : 1
TA1511S	3.323	2.39 : 1	1.65 : 1	69.0%	94.0%	1.43 : 1
TA1520B	2.601	3.06 : 1	1.81 : 1	59.1%	90.9%	1.72 : 1
TA1448L	2.930	2.72 : 1	1.73 : 1	63.6%	90.6%	1.60 : 1
TA1439S	3.267	2.43 : 1	1.63 : 1	67.1%	93.4%	1.43 : 1
TA1359P	2.487	3.22 : 1	1.88 : 1	58.3%	77.8%	2.26 : 1
TA1413L	3.783	2.10 : 1	1.50 : 1	71.4%	95.4%	1.25 : 1

than the entropy encoder.

As was described in chapter 3, the Tolan data decorrelator uses time compression in conjunction with a 2nd

difference operation. The Tolan time compression technique makes the assumption that the second difference value of zero occurs so frequently that the encoding of 0 would be less efficient than the storing of a zero value run counter. In the experimental situation in this thesis, signal noise was quite evident in the EKG traces (Fig 29). The sharp ''spikes'' induced by noise are accentuated by the second difference operation, hence a $\Delta^2 \neq 0$ was a common occurrence. This forced the storage of a lot of 7 bit time counters.

In 8 out of the 9 compression runs made, the frequent storage of time counts actually caused the Tolan data decorrelator to produce negative compression (i.e. more bits out than went in). Since the entropy encoder was producing a larger positive compression ratio, the overall compression figure remained positive. The effect of this operation is graphically illustrated in Figure 30.

A look at the original and reconstructed EKG traces (Fig 27-29) in conjunction with the data in Table II, shows that as the ''noise'' level increased on the signal the compression became progressively worse. Since it was concluded that the entropy encoder performed approximately the same across the data set, the degradation in total compression must be due, in large degree, to the degradation in the Tolan time compression data decorrelator performance.

A clear effect of the noise is illustrated in Figure




Fig 30 . Tolan Compression Ratio Breakdown.

31. From Figure 31, it can be seen that as the "noise" level increased the variance of the 2nd difference distribution became increasingly larger. As the distribution became less peaked, the efficiency of the variable length encoder decreased with a subsequent loss in achieved compression ratio.

The nonzero component at $\Delta^2 = 0$ in Figure 31 can be explained as follows. The decision was made early in the design of the EKG-DAAS that 8 bit data would be used versus 9 or 10 as recommended by the American Heart Association (Ref 3). This decision was made to simplify the software (i.e. single precision could be used). Since the A/D converter has 12 bit resolution, the sample was rounded to 8 bits for uniform error distribution (Ref 24:424-432). Unfortunately because of this rounding action, small differences in the least significant bits of the A/D converter may have affected the rounding operation. Since the A/D cannot sample all three channels simultaneously, the probability that a "noisy" signal will change the least significant bit (or bits) is high. With a rounding operation, these changes may ripple to affect the least significant bit of the 8 bit data.

At most this effect would only cause a change in the least significant bit of the 8 bit rounded values. This would cause, however, the second differences between channels X,Y,Z to be different even though they were connected in common. As was described in chapter 3, any of




Figure 1. The distribution of the number of segments in the top plot is the same as the distribution of the number of segments in the middle plot. The distribution of the number of segments in the bottom plot is the same as the distribution of the number of segments in the top plot.

the three leads with a non zero second difference forced the storage of all three data points. Hence many zeros data points were stored because of 2nd difference asymmetry which would not have occurred had the 12 bit A/D data been truncated instead of rounded.

This fact undoubtedly affected the overall compression efficiency of the Tolan algorithm. Nonetheless, the distribution in Figure 31 would not have changed (except the zero value count) significantly and noise would still have broadened the 2nd difference distribution.

Tolan and Dower Performance Comparison

As has been reiterated several times in the text of this thesis, the author's original intention was to implement both the Dower and the Tolan algorithms for experimental test and comparison. Since time did not permit the Dower implementation, an experimental comparison was not possible.

As an attempt to compare the results of the Dower and Tolan compression routines, the results quoted from the papers by Dower, Berghofer, and Stewart (Ref 12,29) will be used. Dower states that his state space variable length encoder approaches the entropy bound of the second difference source ''with about 1.65% wastage'' (Ref 12:3). This value is significantly higher than the approximately 30% ''wastage'' observed with the Tolan variable length encoder.

Assuming similar operations of the second difference time compression (in this author's 'noisy' environment), the Dower decorrelator is expected to have negative compression. With the Tolan and Dower decorrelators assumed 'equal', then the real gain of the Dower technique over the Tolan approach is in the VLC.

Assuming the Dower entropy encoder approached closely to the entropy bound, it can be extrapolated that the Dower algorithm would have achieved a maximum compression of approximately $(1.2) * (3.22) = 3.86:1$ for the 'best' EKG in Table II and Figure 30. The worst compression ratio would have been $(.84) * (2.1) = 1.76:1$ for the worst (noisest) EKG.

Chapter 5 Summary

This chapter began with a description of the EKG collection experiment. Although the EKG apparatus was limited to one channel, successful collection and compression of EKG data was performed. Analysis of the data revealed that time compression is inefficient in a 'noisy' environment and that 'rounding' of the 12 bit A/D samples in conjunction with placing all three sample lead in common accentuated the degradation caused by the time compression data decorrelator. Nonetheless, the Tolan algorithm did achieve an overall positive data compression figure but significantly lower than the 9:1 ratio achieved by Ruttiman and Pipberger (Ref 28) or the average value of 7.3:1 reported by Stewart, Berghofer, and Dower (Ref 29). Finally

it was heuristically shown that if the Dower compression algorithm lived up to the statements by Dower, then a compression ratio gain of 3.86:1 to 2.26:1 could have been achieved with the Dower EKG data compression technique. This thesis will now proceed to provide conclusions and recommendations

VI. Summary, Conclusions, and Recommendations

Summary and Conclusions

This thesis has investigated the field of electrocardiogram data compression with the objective of evaluating compression algorithms on a 6800 microprocessor based computer system. Accomplishment of this goal required the construction of the EKG-Data Acquisition and Analysis System utilizing the Motorola Exorciser microcomputer.

To determine those EKG compression algorithms which had potential for Exorciser implementation, a literature search was made to locate EKG data compression techniques. In addition to the literature search, personal correspondence (Ref 11,31) yielded several EKG compression algorithms. The results of this research is presented in chapter 2.

Since thesis requirements dictated the need for an online, real time data compression algorithm, only the fastest EKG compressors could be considered. Two routines were selected for detailed analysis and implementation. These two compression algorithms (Tolan and Dower) were discussed at length in chapter 3.

It was deduced from the methodology of the two data compression techniques that the Dower algorithm would perform better than the Tolan procedure. To test this

hypothesis, the EKG-DAAS was constructed and the Tolan algorithm implemented. The implementation of the EKG-DAAS is documented in chapter 4.

Time constraints prohibited completion of the Dower compression algorithm, but data was successfully compressed, analyzed, and decompressed with the Tolan algorithm. The results of this analysis are presented in chapter 5.

The conclusions of this research effort are as follows. First, EKG data compression can be accomplished in real time by a microprocessor based computer system. The Exerciser is a slow microcomputer (1 MHz cycle time) yet it was still possible to implement the Tolan algorithm with a 500 Hz sample rate. The second conclusion is that signal noise can dramatically affect the efficiency of the EKG routines in the same class as the Tolan algorithm. The expense and implementation difficulties of prefiltering, low electromagnetic noise environment, and proper EKG lead attachment are well worth the gain in compression achieved. Finally, a software project of this magnitude should not be attempted totally in assembly language. Although assembly language offers the greatest flexibility and speed, algorithm implementation and debugging efforts are enormous. A high order language would have allowed this author to complete his original thesis objectives.

Recommendations

Several recommendations are offered for further study in microcomputer based EKG data compression. First, implementation of the EKG compression algorithms using a 16 bit microprocessor (e.g. 6 MHz Intel 8086) would permit an order of magnitude improvement in speed of execution. Hardware multiply and divide along with 16 bit arithmetic registers would permit easy implementation of the Dower algorithm and would even make use of the Transform compressors (i.e. FFT) feasible. Second, further study is needed on determination of decorrelator inefficiency on the overall data compression. A large study of different decorrelators such as 1, 2, 3 difference operations with and without time compression is needed. Next, an EKG compression algorithm implemented using the Ruttiman and Pipberger 2nd order interpolator (Ref 28) for the decorrelator along with the Dower variable length encoder (Ref 12) should be built. This combination should prove to be very efficient. Finally, programming and experimental testing should be done on a full scale microcomputer development system, complete with a high order language, A/D-D/A capability, and flexible disk file manipulation software. Such a system is the Zilog MCZ1/25 microcomputer resident here at A.F.I.T.

Bibliography

1. Ahmed, Nasir et al. ''Electrocardiographic Data Compression Via Orthogonal Transforms,'' IEEE Transactions on Biomedical Engineering, BME-22, (6) : 484-487 (November 1975).
2. Ahmed, N and Rao, K.R. Orthogonal Transforms for Digital Signal Processing. New York: Springer-Verlag, 1975.
3. American Heart Association, Rep. Comm. ECG. ''Recommendations for Standardization of Leads and Specifications for Instruments in Electrocardiography and Vectorcardiography,'' Circulation, 54: 11-31 (1975).
4. Berger, Toby. Rate Distortion Theory: A Mathematical Basis For Data Compression. Englewood Cliffs, N.J.: Prentice-Hall, Inc., 1971.
5. Blasbalg, H. and Van Blerkom, R. ''Message Compression,'' IRE Transaction on Space Electronics and Telemetry, 8 : 228-238 (September 1962).
6. Box, E and Jenkins, G.W. Time Series Analysis: Forecasting and Control. San Francisco, California : Holden-Day Inc, 1976.
7. Cox, J.R. and Ripley, K.L. ''Compact Digital Coding of Electrocardiographic Data,'' Proceedings Sixth Hawaii International Conference on System Sciences. 333-336. Honolulu, Hawaii : University of Hawaii, January 9-11, 1973.
8. DATEL SYSTEMS, Inc. Sinetrac Series Model ST-6800 A/D-D/A Peripheral Systems Instruction Manual, Part No. 58-12140-25, Mansfield, Mass. January, 1979.
9. Davenport, W.B. Probability and Random Processes: An Introduction For Applied Scientists and Engineers. New York : McGraw-Hill Book Company, 1970.
10. Davisson, Lee D. and Gray, Robert M., editors. Benchmark Papers In Electrical Engineering and Computer Science, Volume 14 : Data Compression. Stroudsburg, Pennsylvania : Dowden, Hutchinson, and Rose, Inc., 1976.
11. Diver, Tom. Project Engineer. Personal Correspondence. Marquette Electronics, Inc., Milwaukee, Wisconsin, 15 May 1980.

12. Dower, R.G. with appendix by Berghofer, D. "'Optimal ECG Data Compression,'" working paper submitted for publication in Proceedings of the Engineering Foundation Conferences: Computerized Interpretation of the ECG V. Asilomar, California. April 27-March 2, 1980.
13. Dubin, Dale, M.D. Rapid Interpretation of EKG's (Second Edition). Tampa, Florida: COVER Publishing Company, 1970.
14. Frank, Ernest. "'An Accurate, Clinical Practical System For Spatial Vectorcardiography,'" Circulation, 13 : 737-749 (May 1956).
15. Gray, R.M. and Davisson, L.D. "'A Mathematical Theory of Data Compression?'" , Benchmark Papers In Electrical Engineering and Computer Science, Volume 14 : Data Compression, edited by L.D. Davisson and R.M. Gray. 21-25. Stroudsburg, Pennsylvania : Dowden, Hutchinson, and Rose, Inc., 1976.
16. Huffman, David A. "'A Method For The Construction of Minimum Redundancy Codes,'" Proceedings of the IRE, 40 : 1098-1101 (September 1952).
17. Jelinek, Frederick. "'Buffer Overflow in Variable Length Coding of Fixed Rate Sources,'" IEEE Transactions on Information Theory, IT-14 (3) : 490-501 (May 1968).
18. Makhoul, John. "'Linear Prediction: A Tutorial Review,'" Proceedings of the IEEE, 63 (4) : 561-580 (April 1975).
19. McEliece, Robert J. The Theory of Information and Coding : Encyclopedia of Mathematics and its Applications, Volume 3. Gian-Carlo Rota, editor. Reading, Massachusetts : Addison-Wesley Publishing Company, 1977.
20. McFec, R. and Baule, G.M. "'Research in Electrocardiography and Magnetocardiography,'" Proceedings of the IEEE , 60 (3) : 290-303 (March 1972).
21. McGibbon, C.I. "'Data Communication Networks of Interest to Medical Information Processing,'" Proceedings of the IFIP TC 4 Working Conference on Optimization of Computer ECG Processing, edited by H.K. Wolf and P.W. MacFarlane. 253-254. Halifax, Nova Scotia, Canada. June 5-7, 1979.

22. Motorola Inc., M6800 UNO-Peiser Users Guide, Phoenix, Arizona: Second Edition, 1975.
23. Midwest Scientific Instruments, Operating and Assembly Manual, The FD-8 Floppy Disk Memory System, Olathe, Kansas.
24. Oppenheim, Alan V. and Shafer, Ronald W. Digital Signal Processing. Englewood Cliffs, N.J. : Prentice-Hall, Inc. 1975.
25. Papoulis, Athanasios. Probability, Random Variables, and Stochastic Processes. New York: McGraw Hill Book Company, 1965.
26. Rautaharju, P.M., Warren, J., and Wolf, H. "'Waveform Vector Analysis of Orthogonal Electrocardiograms: Quantification and Data Reduction, Journal of Electrocardiography, 6: 103-111 (March 1973).
27. Ristenbatt, Marlin P. "'Alternatives in Digital Communications,'" Proceedings of the IEEE, 61 (6):703-721 (June 1973).
28. Ruttimann, Urs E. and Pipberger, Hubert V. "'Compression of the ECG by Prediction or Interpolation and Entropy Encoding,'" IEEE Transactions on Biomedical Engineering, BME-26, 11 : 613-622 (November 1979).
29. Stewart, D., Berghofer, D. with appendix by Dower, R.G. "'Data Compression of ECG Signals,'" Proceedings of the Engineering Foundation Conferences, Computerized Interpretation of ECG IV. Asilomar, California. January, 1979.
30. Shannon, Claude E. "'A Mathematical Theory of Communication,'" Bell Systems Technical Journal, 27 :379-423 (July 1948).
31. Tolan, Gil. (MD), Lt. Col. USAFMC, Personal Correspondence, USAF School of Aerospace Medicine, Brooks AFB, Texas. April, 1980.
32. Tompkins, W. and Webster, J.G., editors. Design of Microcomputer-Based Medical Instrumentation. Englewood Cliffs, N.J.: Prentice-Hall, Inc., 1980.
33. TRW Systems Group, Houston, Texas. "'Investigation of Data Compression Techniques,'" Final Report, NASA-CR-115177. September, 1971, N71-35323.
34. Viterbi, Andrew and Omura, Jim. Principles of Digital Communication and Coding. New York : McGraw Hill Book Company, 1979.

35. Womble, Edward M. et al. "'Data Compression For Storing and Transmitting ECG/VCG's,'" Proceedings of the IEEE, 65 (5) : 702-706 (May 1977).
36. Ziemer, R.E. and Tranter, W.H. Principles of Communications: Systems, Modulation, and Noise. Boston, Massachusetts : Houghton Mifflin Company, 1976.

Appendix A

The Electrocardiogram

Introduction

The electrocardiogram (or EKG) is a record of the electrical activity of the heart as measured from the body surface. The magnitude, shape, and timing of the electrical potentials generated by the heart reveal a great deal of information concerning the health of the cardiac system. This appendix will describe how the electrical signals from the heart are generated, how the cardiac cycle is coordinated and controlled, and finally how the EKG can be used as a diagnostic tool.

The Physiology and Electrical Characteristics of the Heart

The heart (Fig. A1) is an organ about the size of a fist with four main pumping chambers and a specialized electrical conduction system. At the top are two thin walled pumps called the atrium which prime the main pumps of the heart, the ventricles. The ventricles are separated by a thick wall of muscle tissue called the septum. Blood from the right ventricle goes to the lungs and blood from the left ventricle goes to the rest of the body.




Fig A1. The Heart Cross-Section (From Ref 20:291).

During the resting state between contractions, the cells of the heart are ''polarized'' with an electrical potential existing between the inside (-) and the outside (+) of the cell. This potential is generated by an ionic gradient across the cellular membrane and is normally maintained for approximately .2 to .4 of a second before spontaneous ''depolarization'' occurs. Depolarization (caused by an intrush of sodiun ions into the cell) induces the cell to contract for approximately 1/4 of a second. Because all heart cells contain specialized conducting fibers, the depolarization of one cell initiates the depolarization of neighboring cells and a ''wave of excitation'' sweeps across the myocardium (heart) at a rate

of about 1 meter per second.

The voltages measured at the body surface are the superposition of thousands of heart cells depolarizing (or polarizing) as the wave of excitement flows through the myocardium. Early work by Wilson (Ref 20:292) showed that the heart could be represented by an equivalent electrical dipole whose vector orientation sweeps through a closed loop during one cardiac cycle. In simple terms, as the wave of excitement flows toward a positive skin electrode, a positive slope is generated on the EKG record representing the projection of the heart vector onto the axis of the EKG lead.

To observe this sweeping dipole vector, Wilson developed the 12 lead EKG system in almost universal use today. This system (Fig. A2) attempts to measure the heart vector from a variety of vantage points in hope of determining the actual direction of propagation of the wave of excitement. Unfortunately, the leads of the Wilson system are not orthogonal and reconstruction of the actual heart vector orientation and amplitude is difficult. To overcome this problem, Frank (Ref 14:737-749) developed the vector cardiogram which combines 7 leads in a summing network to produce three orthogonal components of the heart vector. Frank VCG systems are popular in heart diagnosis and research because all of the information is contained in three leads of data versus 12.




Fig A2. Wilson EKG Electrode System (From Ref 20)

The Cardiac Cycle

For effective pumping action, the heart muscle must contract in a controlled, coordinated way. This means that the stimulating wave of excitation must follow a well defined conduction path to allow the heart to contract in the most efficient manner. In "normal" hearts, the contraction sequence begins in the right atrium where

special cells in an area known as the Sino-Atrium (SA) node spontaneously depolarize faster than the rest of the heart tissue. The SA node hence initiates a wave of excitation which covers the atrium in about 80 milliseconds. This atrial contraction produces the electrical signal called the P wave (Fig. A3) on a typical EKG record.

The wave of excitation started by the SA node then reaches another specialized receptor known as the Atrial-Ventricular (AV) node. Here connecting fibers delay the excitation impulse for about 50 milliseconds to allow the ventricles to fill with blood. After the 50 millisecond delay, the excitation signal is relayed to special conducting fibers in the septum known as the Bundle of His. These conducting fibers rapidly (30 ms) transmit the excitation wave to the interior (endocardium) wall of the ventricles where the wavefront propagates radially to the outer wall (epicardium) in another 30 milliseconds. The ventricular wave of excitation produces the QRS waveform complex seen on the EKG.

Following the ventricular contraction, the muscle cells of the ventricles repolarize over a period of 100 milliseconds. No muscular action is occurring but the EKG responds to this electrical activity and the T wave is noted on the EKG record. Atrial repolarization occurs during the QRS hence it is generally invisible on the EKG.




Fig A3. Typical EKG Waveform (From Ref 20)

Finally the heart rests for approximately .2 of a second and the cycle starts again. The above cardiac cycle is typical of a healthy heart. Disease, however, can effect this sequence dramatically.

Heart Disease and the EKG

The variety of ailments which plague the human heart are numerous and no attempt will made to describe the spectrum of diseases possible. There are, however, several common heart defects which routine EKG analysis usually detects. These include premature contractions (Atrial and Ventricle), bundle branch blocks, hypertrophy, and infarction.

Premature Contractions. Premature contractions are caused by the spontaneous depolarization of heart tissue outside of the SA node. This depolarization initiates a wave of excitement causing the atrium, or more notably, the

ventricles to contract out of rhythm with their normal cycle. Points where this "unscheduled" depolarization occurs are called ectopic foci and can arise because of infarct damage (to be discussed later), coronary heart disease producing oxygen starvation and a number of other causes. Premature Ventricular Contractions (PVC's) are highly visible on the EKG record. Normally the ventricles contract simultaneously and the voltage vectors generated tend to cancel keeping the QRS amplitude relatively small. A PVC, however, causes depolarization of one ventricle before the other generating an unbalanced, hence larger, voltage output.

Bundle Branch Blocks. A bundle branch block is caused by a block of the impulse of the right or left Bundle Branch. This causes a delay in the transmission of the stimulation impulse to ventricle blocked and forces the two ventricles to contract at slightly different times. This difference in ventricular contraction time shows up on the EKG as a double humped QRS.

Hypertrophy. Hypertrophy is an enlargement of one section of the heart muscle tissue. This enlargement affects the duration and strength of the wave of excitement and is visible on the EKG record as a diphasic trace if atrial hypertrophy is present. If ventricular hypertrophy exists, the QRS amplitudes are much larger than normal due to the fact that more tissue is depolarizing.

Infarction. Myocardial infarction is an injury to the heart tissue caused by an occlusion of a coronary artery. An area of the heart is then without a blood supply and often permanent damage occurs. This heart disease is often the one most commonly called a "heart attack" and, as is well known, is many times fatal. If a person survives the original "attack", then this permanent damage shows up in the EKG as a change in the QRS and T waves. This change occurs because the infarcted tissue no longer responds to the excitation wave and the wave moves around, not through, the damaged tissue.

Summary.

This appendix has briefly examined the physiology of the heart and discussed how the myocardial tissue generates the electrical fields measured by the electrocardiogram. Though EKG analysis has been practiced for fifty years, intense research continues in improving EKG diagnosis. Dramatic improvement in computer aided EKG analysis and better understanding of the electro-physiology of the heart is leading to improved cardiac health care worldwide.

The heart diseases discussed above are only a small subset of the problems which can afflict the human heart. Should the reader desire a more thorough background on heart physiology and cardiac disease, Dubin's book (Ref 13) is highly recommended. This programmed text carefully leads the reader through EKG analysis and is easily read. For a

more firm background on the electro-physics of the heart,
the article by McFee and Baule (Ref. 20) provides a good
tutorial review on EKG history and research.

Appendix B

Fundamentals of Information Theory

Introduction

In 1948, Claude Shannon published a classic paper (Ref 30:379-423) titled "'A Mathematical Theory of Communication'" in which he laid the foundation of modern information theory. Shannon used his "'information theory'" to describe, mathematically, the interrelationships between the components of a "'typical'" communication system as illustrated in Figure B-1.




Fig B-1. The communication system model (From Ref 4).

This appendix is written as a basic tutorial on information theory, and is intended to acquaint the reader

with the terminology used in describing the data compression techniques in chapters 2 and 3. The appendix begins with a review of the information source, proceeds to a discussion of the transmission channel and Shannon's rate distortion theory, and concludes with a description of the system encoder/decoder.

Information Source

For his initial analysis, Shannon proposed modelling the information source as a discrete stochastic process whose output is governed by known statistics. As pointed out by Davisson and Gray (Ref 10: 2-4), the discrete-time model was commonly used for any or all of the following reasons: (1) digital communication links have become common place; (2) a continuous time process can be modelled as discrete by sampling, orthogonal function expansion, or waveform segmentation; (3) greater simplicity.

The source is characterized by a finite set of possible outcomes known as its alphabet A . The occurrence of a particular alphabet symbol is governed by probabilistic descriptors (i.e. probability density functions) and it is assumed that the source produces only one symbol from the alphabet every T_s seconds. Hence the information source has a symbol rate of $R_s = 1/T_s$ symbols per second.

The next question of interest is how much information is conveyed by the occurrence of a given source symbol? If X is a discrete random variable occurring at time t , and x is

an element of A, then the random variable self information can be defined. That is

$$I(x) = -\log p\{X=x\} \quad (B.1)$$

where $p\{X=x\}$ is the probability that $X=x$.

According to this description, the less probable an event is, the more information is conveyed when it occurs. The base of the logarithm is unspecified, but in this thesis it is assumed to be base 2. Hence the occurrence of symbol x reveals $I(x)$ bits of information.

The amount of information received per observation is of interest, but one would like a measure of the "uncertainty" or "randomness" of the source. If the stochastic process defining the source is considered stationary (a pretentious assumption but one generally made) then the output of the source is a sequence of random variables with identical probabilistic descriptors. The probabilistic descriptor will be defined as $\{u\}$ and could represent the moments of the random variable or its probability density function (PDF). If, in addition, the source is considered ergodic, then statistical averages equal time averages and calculation of the set $\{u\}$ is greatly simplified.

Given that the discrete process is ergodic, or at least stationary, then the source output at any time is described by the random variable X with range $\Lambda = \{x(1), x(2), \dots, x(n)\}$. The measure of the "uncertainty" or "randomness" is

defined as the entropy of the source and is given by the equation

$$H(X) = - \sum_{i=1}^n p_i \log_2 p_i \quad (B.2)$$

where p_i is the probability of occurrence of the discrete value x_i . If $p_i = 0$, then the term $\log_2 1/p_i$ is defined equal to 0 (i.e. no contribution to the entropy). Should the range A be of infinite extent (i.e. a continuous source), then the above series may not converge nor would $p(x)$ necessarily be defined. In this case, $H(X)$ is defined as positive infinity.

As example (Ref 19:15), let X represent the outcome of a single roll of a fair die. Then $A = \{1, 2, 3, 4, 5, 6\}$ and $p_i = 1/6$ for each i . Here $H(X) = \sum_6 1/6 \log_2 6 = 2.58$ bits.

In the above example, the statistics governing the occurrence of a given outcome were uniform. This represents the "most random" case with a resultant maximum of the entropy function. Should the die be "loaded", then the predicted outcome is "less random" and hence the value of $H(X)$ would be reduced.

The next important component in Shannon's communication model is the transmission channel.

Transmission Channel

The transmission channel is also assumed to be a discrete time, finite alphabet device which accepts and transmits to the receiver one symbol in a finite alphabet B each T_c seconds. The alphabet B is often binary, hence the dimension of the symbol space (defined as $\{|B|\}$) is 2 and $B=\{0,1\}$. The transmission rate of the channel is defined as $R_c=1/T_c$ channel symbols (bits if binary) per second. If the output of the source process is defined as $\{X\}$, and the channel is 'noisy', then the received symbol $\{Y\}$ might not equal the transmitted symbol $\{X\}$. A noiseless channel transmits symbols with no error in which case $\{X\}=\{Y\}$.




Fig B-2. A discrete memoryless channel (From Ref 19:19)

For simplicity, let the channel be modelled as a

discrete memoryless channel (DMC) as depicted in Figure B-2. The DMC is described mathematically by a conditional probability which relates the chance that a given output y was the result of a given input x . If the input to the channel is a random variable X , and the output is a random variable Y , then a quantity called conditional entropy can be defined. The equation defining conditional entropy is:

$$H(X|Y) = \sum_{x=1}^n \sum_{y=1}^n p(x,y) \log_2 \frac{1}{p(x|y)} \quad (B.3)$$

For a given pair X, Y of random variables, $H(X|Y)$ represents the amount of uncertainty remaining about X after Y has been observed.

Now since $H(X)$ represents the uncertainty about X before X is known and $H(X|Y)$ represents the uncertainty after, the difference $H(X) - H(X|Y)$ must represent the amount of information provided about X by Y . This quantity is called the mutual information between X and Y , and is denoted by :

$$I(X;Y) = H(X) - H(X|Y) \quad (B.4)$$

With the above definitions in hand, the most important quantity of a communications channel can be described; that quantity is the channel capacity. Channel capacity is

defined as the maximum amount of information, per unit of time, which can be 'reliably' transmitted over the channel. That is :

$$C = \max\{I(X;Y)\} \quad (B.6)$$

Channel capacity is closely related to another important parameter of a communication system known as the rate-distortion function.

Rate-Distortion. Shannon postulated the existence of a mathematical distortion measure, $d(X,Y)$, to measure the distortion or loss resulting if a source symbol X is reproduced as Y . Unfortunately this abstract distortion measure is difficult to quantify. As stated by Berger (Ref 4:6), 'the unavailability of a distortion measure that is both physically meaningful and analytically tractable constitutes one of the major obstacles to progress in (communication) system design.'

Assuming one has such a distortion measure, then associated with most source-user pairs is a function $R(D)$ called the rate distortion function. The rate distortion is important in that it gives the designer a mathematical tool to measure the amount of distortion that can be expected for a given transmission rate. A communication system can achieve a given fidelity D if and only if the capacity C exceeds $R(D)$. Hence $R(D)$ is the effective rate at which the

source produces information subject to the constraint that the user can tolerate an average distortion of D .




Fig B-3. A typical rate distortion function (From Ref 4:7).

The simplest source-user pair is a discrete memoryless source (DMS) and a single letter fidelity criteria. The DMS produces statistically independent, identically distributed, discrete random variables. Assuming a single letter fidelity criterion, every time the system presents the letter y to the user when the source output was actually x , a nonnegative penalty $p(x, y)$ is determined. The rate-distortion function for the above case is plotted in Figure B-3. As can be seen from the Figure, $R(0)$ is equal to $H(X)$. This last result leads to Shannon's famous channel coding theorem which states that if the symbol rate is less than the channel capacity, it is possible to transmit with perfect fidelity. That is, if the entropy of the source,

$H(X)$, is less than or equal to the channel capacity C , then transmission with zero error can be obtained. This amazing result is not without cost, however. To reduce the source entropy (as seen by the channel) requires coding. Encoder

In most circumstances, the output of the information source is not suitable for direct input into the channel. To match the source to the channel, much like matching impedances in circuit theory, is the job of the encoder.

The encoder incorporates all of the functions which process the source data for transmission over the communication channel. This includes coding, analog-to-digital conversion, and modulation. In order to transmit without loss, integers K and L must exist such that $KT_s = LT_c$. This guarantees that the received sequence $\{Y\}$ has the same symbol rate as the transmitted sequence $\{X\}$.

As illustrated in Figure B-1, the encoder is divided into two functional subunits. The first of these subunits is the source encoder.

Source Encoder. Source encoding is the operation by which the source output is mapped into an alternate symbol set with the goal of reducing source sequence dependence (i.e. reduce redundancy). Source coding is the transformation in which data compression occurs.

The first source codes, as conceived by Shannon, were block codes where blocks of source symbols were mapped into a single representative channel symbol. This type of

encoding is simple, but efficient only if the source repeats a given, constant length, sequence of symbols on a frequent basis. Another type of source coding is variable length coding.

In variable length codes, those source symbols which occur most frequently are assigned the shortest channel codewords.

One of the most popular variable length codes in use was originally proposed by Huffman (Ref 16:31-34) in the early fifties. This code is in a class known as uniquely decodable (UD) codes. UD codes imply that the codewords, regardless of length, are unique sequences. Hence decoding is instantaneous upon codeword reception (i.e. the decoding does not depend upon reception of the next codeword). McEliece has shown (Ref 19:244-245) that the Huffman code is optimal in the class of UD source codes.

Many techniques for source coding are available. Chapter 2 of this thesis discusses some approaches used in "compressing" electrocardiogram waveforms.

The output of the source coder is generally a sequence of discrete symbols (i.e. binary 1|0) suitable for further processing by the channel coder.

Channel Encoder. Shannon proved that if the rate of the source, as seen by the channel, is less than the channel capacity C , then "noiseless" or error free transmission is theoretically possible. The goal of channel encoder is to combat channel noise to achieve a probability of error (P_e)

which approaches zero.

Channel encoders allow a reduction of P_e by selectively reinserting redundancy which has been removed by the source coder. This redundancy allows error detection, and with the proper codes, error correction by the decoder. How close the P_e approaches zero is purely a function of the effort (and money!) spent on channel coding.

A common channel coding technique is the use of 'parity' bits in digital communication. With parity checking, a bit (or bits) is added to the source word which represents the number (odd or even) of 'ones' in the word. For example, if odd parity is defined and the source word is 1100101 then a single parity bit of 0 implies an even number of ones (1100101|0).

Parity checking in the example above will detect a single bit error in transmission. A more powerful technique is Hamming codes.

The Hamming code is in a large class of codes known as linear codes. Hamming not only detects errors in transmission, but will correct errors to a certain level. This correction is accomplished by multiplying the received codeword by a matrix known as the syndrome. In the case of Hamming codes, the output from this transformation is the bit position in error.

Decoder

As would be expected, the decoder is the inverse

operation which outputs an estimate of the data input to the encoder. The decoder is composed of the channel decoder and the source decoder. The channel decoder uses the redundancy added by the channel encoder to perform error checking and/or correction. The source decoder takes the "correct" data output from the channel decoder and "decompresses" the data to produce an estimate of the information source. Information theory has shown that if enough time, complexity, and money is spent on channel encoding, and the information rate is below the channel capacity, then "error free" transmission is possible.

Conclusion

This appendix has been a very brief summary of a very large field of study. The key words underlined throughout this text are terminology which appear in the theory chapters (2 and 3) of the thesis. For a more thorough, mathematical treatment of the fields of Information Theory and Coding, the reader is encouraged to refer to the textbooks by McEliece (Ref 19) and Berger (Ref 4). The IEEE Press and Benchmark book (Ref 10) is an excellent source for a survey of the key papers in the field of Data Compression.

Appendix C

This appendix contains the 6800 assembly language source programs of the EKG-DAAS.

```

00030      *:
00040      *
00050      * PROGRAM NAME: EKG-EXEC
00060      * AUTHOR : CAPT. TEL. TOLAN
00070      * VERSION : 1.8
00080      * VERSION DATE 2 OCT 80
00090      *
00100      *
00110      * PROGRAM DESCRIPTION
00120      *
00130      * THIS PROGRAM IS THE EXECUTIVE ROUTINE WHICH
00140      * CONTROLS THE EXERCISER EKG DATA ACQUISITION
00150      * SYSTEM. THIS ROUTINE CALLS OVERLAYED PROGRAMS
00160      * WHICH PERFORM DATA COLLECTION, COMPRESSION, STORAGE
00170      * AND RECONSTRUCTION.
00180      * THIS SOFTWARE IS IN SUPPORT OF THESIS RESEARCH
00190      * TO IDENTIFY THE MOST EFFICIENT EKG DATA COMPRESS-
00200      * ION ALGORITHM IN THE TEST SET.
00210      *
00220      * COMMAND OPTIONS
00230      *
00240      * 0=STORAGE WITHOUT COMPRESSION (10 BIT)
00250      *
00260      * 1=COMPRESS WITH ALGOR TOLAN-A
00270      *
00280      * 2=COMPRESS WITH ALGOR TOLAN-B
00290      *
00300      * 3=COMPRESS WITH ALGOR DOWER
00310      *
00320      * 4=COMPRESS WITH ALGOR 2ND ORDER INTERPOL
00330      *
00340      * 5=COMPRESS WITH ALGOR TURNPT
00350      *
00360      * 6=DISPLAY COLLECTED DATA & STATS
00370      *
00380      * 7=JUMP PROGRAM CONTROL TO DOS
00390      *
00400      * 8=JUMP PROGRAM CONTROL TO EXBUG
00410      *
00420      * 9=LOAD & EXECUTE OTHER OVERLAYS
00430      *
00440      * S=SAVE CUR MEM FILE TO DISK
00450      *
00460      * START OF PROGRAM
00470      *
00480      *
00490 1D00      ORG      $1D00      PROGRAM START LOCATION
00500      *
00510      OPT      0          ASSB OPT TO CREATE CUR FILE
00520      OPT      NOG       ASSB OPT TO SUP FCC LIST
00530      *
00540      *
00550      * LABEL DECLARATIONS
00560      *

```



```

00570      * SUPPORT SUBROUTINE ADDRESSES
00580      *
00590      C27F  QWTRD0 EQU  $C27F  EOS. QWTRD0 INPUT ROUTINE
00610      CA36  KEYRDO EQU  $CA36  EOS. KEYRD INPUT ROUTINE
00620      F000  EXBUG  EQU  $F000  EXBUG. EXBUG ENTRY PT
00630      2800  DOS    EQU  $2800  DOS. DOS ENTRY PT
00640      C75B  WRTE0  EQU  $C75B  EOS. CLR PASSWD BUFFER
00650      CC87  CLRPAS EQU  $CC87  EOS. CLR PASSWD BUFFER
00660      C7C8  DRIVE  EQU  $C7C8  EOSIO. DSK DRIVE SELECT
00670      C803  RLIB   EQU  $C803  EOSIO.
00680      C807  WLIB   EQU  $C807  EOSIO.
00690      CE52  LOAD2  EQU  $CE52  EOS. LOAD PRGM ROUTINE
00700      CBDA  DOSTR2 EQU  $CBDA  EOS. ALT EOS ENTRY LOC
00710      E055  BYTE  EQU  $E055  MIREUG. GET TO HEX DIG FRM TE
00720      *
00730      * DATA BUFFERS
00740      *
00750      3066  NAME   EQU  $3066  NAME BUFFER FOR FILE I/O
00760      3060  TEMPX1 EQU  $3060  TEMP 2 BYTE STOR BUFFER
00770      3062  TEMPX2 EQU  $3062  TEMP 2 BYTE STOR BUFFER
00780      3058  STARTX EQU  $3058  LOWEST ADDR USED IN PRGM
00790      305A  ENDX   EQU  $305A  HIGHEST ADDR USED IN PRGM
00800      305C  GOX    EQU  $305C  START EXECUTE ADDR
00810      001B  PROGX  EQU  $001B  DISK I/O ERR VEC ADDR
00820      0008  RMEMH  EQU  $0008
00830      3002  ENDBUF EQU  $3002  BUF WITH ADDR OF LAST CHAR IN
00840      3400  HDRSTR EQU  $3400  START OF MEM FILE HDR SECTOR
00850      3004  BUFFER EQU  $3004  KEYED INPUT BUFFER START
00860      3057  TYPE   EQU  $3057  TYPE OF FILE FOR I/O
00870      1C96  STKSAV EQU  $1C96  TEMP STACK SAVE BUFFER
00880      1C98  CPRTYP EQU  $1C98  COMPRESSION TYPE BUFFER
00890      1C9A  COUNT  EQU  $1C9A  GENERAL 8 BIT COUNTER
00900      1C9B  NAMPTR EQU  $1C9B  NAME POINTER FOR CONSOLE I/O
00910      1C9D  VECsav EQU  $1C9D  I/O VEC SAVE BUFFER
00920      1C9F  OLAYGO EQU  $1C9F  OVERLAY EXECUTE FLAG
00930      1CA0  LGOFIG EQU  $1CA0  PRSTAT LOAD VS EXECUTE FLAG
00940      1CA1  FILHLC EQU  $1CA1  FILHDR SUB ADDR PASS BUF
00950      1CA3  SAVEDL  EQU  $1CA3  SAVEFL SUB ADDR PASS BUF
00960      1CA5  HXASLC EQU  $1CA5  HXASC SUB ADDR PASS BUF
00970      1CA7  HXBUFF EQU  $1CA7  HXASC PARAMETER BUFFER
00980      1CA9  PDFPLC EQU  $1CA9  PDFPRT SUB ADDR PASS BUF
00990      1CAB  OVRLLC EQU  $1CAB  OVERLAY SUB ADDR PASS BUF
01000      1CAD  OVPRBF EQU  $1CAD  OVERLAY PARAMETER PASS BUFFER
01010      3490  LOCPCT EQU  $3490  TOTL CAL LOOPS EXEC IN WIT
01020      3494  SAMPNO EQU  $3494  NUM OF SAMPLES TAKEN
01030      3496  LPCAL  EQU  $3496  NUM OF CAL LOOPS PER 1 INTR
01040      3497  MAXZ    EQU  $3497  MAX VLU IN CH Z
01050      3498  MAXZLO EQU  $3498  LOC OF MAX VLU IN CH Z
01060      349A  MINZ    EQU  $349A  MIN VLU IN CH Z
01070      349B  MINZLO EQU  $349B  LOC OF MIN VLU IN CH Z
01080      349D  MAXY    EQU  $349D  MAX VLU IN CH Y
01090      349E  MAXYLO EQU  $349E  LOC OF MAX VLU IN CH Y
01100      34A0  MINY    EQU  $34A0  MIN VLU IN CH Y

```

01110	34A1	MINLO EQU	\$34A1	LOC OF MIN VLU IN CH Y
01120	34A3	MAXX EQU	\$34A3	MAX VLU IN CH X
01130	34A4	MAXLO EQU	\$34A4	LOC OF MAX VLU IN CH X
01140	34A6	MINX EQU	\$34A6	MIN VLU IN CH X
01150	34A7	MINLO EQU	\$34A7	LOC OF MIN VLU IN CH X
01160	34A9	MEMBIT EQU	\$34A9	NUM OF BITS AVAIL FOR STO
01170	34AC	DTABIT EQU	\$34AC	NUM OF BITS USED TO STO DTA
01180	34B0	XBITS EQU	\$34B0	NUM OF BITS USED TO STO X
01190	34B3	YBITS EQU	\$34B3	NUM OF BITS USED TO STO Y
01200	34B6	ZBITS EQU	\$34B6	NUM OF BITS USED TO STO Z
01210	34B9	TBITS EQU	\$34B9	NUM OF BITS USED TO STO TIME
01220	34BC	ACELCT EQU	\$34BC	# BITS FED TO VAR LEN CODER
01230	34C2	BASSAV EQU	\$34C2	SAVE LOC FOR VLU \$3620 FROM
01240	3460	MAXMIN EQU	\$3460	START OF MAX,MIN ASCII BUFFER
01250	34C3	ENTRPY EQU	\$34C3	START OF ENTROPY TABLE IN ASC
01260	3C00	SECZRO EQU	\$3C00	SEC 0 OF MEM FILE
01270	0019	SAVEX EQU	\$0019	TEMP LOC TO SAVE INDEX REG
01280	3500	XPDF EQU	\$3500	LOC OF X PDF BUFFER
01290	3700	YPDF EQU	\$3700	LOC OF Y PDF BUFFER
01300	3900	ZPDF EQU	\$3900	LOC OF Z PDF BUFFER
01310	3B00	TPDF EQU	\$3B00	LOC OF TIME PDF BUFFER
01320	8000	BUFEND EQU	\$8000	END OF MEM BUF
01330		*		
01340		* HARDWARE ADDRESSES		
01350		*		
01360	E500	DACZRO EQU	\$E500	DAC 0 ADDRESS
01370	E502	DACONE EQU	\$E502	D :THIS ROUTINE IS THE EXECUTIVE
01460		* CONTROLLER OF THE EKG DATA ACQ SYS		
01470		*		
01480		*		
01490	1D00 0F	START SET		STOP POSSIBLE INTR ON RESET
01500	1D01 CE 4000	LDX	#\$4000	
01510	1D04 FF E500	STX	DACZRO	CLR DACS & SET SEL 1 HIGH
01520	1D07 FF E502	STX	DACONE	
01530	1D0A CE 214A	LDX	#FILHDR	
01540	1D0D FF 1CA1	STX	FILHLC	PUT FILHDR ADDR IN PASS BUF
01550	1D10 B6 3620	LDA A	\$3620	
01560	1D13 B7 34C2	STA A	BASSAV	SAVE CURRENT \$3620 IN BUF FO
01570	1D16 CE 208E	LDX	#SAVFIL	
01580	1D19 FF 1CA3	STX	SAVFILC	PUT SAVFIL ADDR IN PASS BUF
01590	1D1C CE 23DB	LDX	#HXASC	
01600	1D1F FF 1CA5	STX	HXASLC	PUT HXASC ADDR IN PASS BUF
01610	1D22 CE 22C2	LDX	#PDFPRT	
01620	1D25 FF 1CA9	STX	PDFPLC	PUT PDFPRT ADDR IN PASS BUF
01630	1D28 CE 2067	LDX	#OVLAY	
01640	1D2B FF 1CAB	STX	OVLILC	PUT OVLAY ADDR IN PASS BUF
01650	1D2E CE 1D00	LDX	#START	
01660	1D31 FF 2F01	STX	DOS+1	
01670	1D34 7F 1C9F	CLR	OLAYGO	CLR OVLAY HCN FLAG
01680	1D37 7F 1CA0	CLR	LGOFILG	CLR FRSTAT LOAD FIG
01690	1D3A 4F	CLR A		
01700	1D3B BD C7C8	JSR	DRIVE	INSURE EXEC TLE TO DRIVE 0
01710	1D3E BD 1F6F	JSR	CHDIR	PROMPT & GET IN CHDIR
01720	1D41 C1 30	CMP B	#'0	IS CHDIR ASCII 0?

01730	1D43	26	09	BNE	EXCM1	NO. CHECK OTHER CLEDS
01740	1D45	CE	1DC7	LDX	#NOCPR	LOAD OVERLAY NAME PTR
01750	1D47	FF	1CND	STX	OVRBUF	
01760	1D4E	7E	2067	JMP	OVRLAY	LOAD & EXECUTE OVRLAY
01770	1D4E	C1	31	EXCM1	CMP B	#'1
01780	1D50	26	09	BNE	EXCM2	
01790	1D52	CE	1DCF	LDX	#TOLAN1	
01800	1D55	FF	1CAD	STX	OVRBUF	
01810	1D58	7E	2067	JMP	OVRLAY	
01820	1D5B	C1	32	EXCM2	CMP B	#'2
01830	1D5D	26	09	BNE	EXCM3	
01840	1D5F	CE	1DD7	LDX	#TOLAN2	
01850	1D62	FF	1CAD	STX	OVRBUF	
01860	1D65	7E	2067	JMP	OVRLAY	
01870	1D68	C1	33	EXCM3	CMP B	#'3
01880	1D6A	26	09	BNE	EXCM4	
01890	1D6C	CE	1DDF	LDX	#DOVER	
01900	1D6F	FF	1CAD	STX	OVRBUF	
01910	1D72	7E	2067	JMP	OVRLAY	
01920	1D75	C1	34	EXCM4	CMP B	#'4
01930	1D77	26	09	BNE	EXCM5	
01940	1D79	CE	1DE6	LDX	#INTER	
01950	1D7C	FF	1CAD	STX	OVRBUF	
01960	1D7F	7E	2067	JMP	OVRLAY	
01970	1D82	C1	35	EXCM5	CMP B	#'5
01980	1D84	26	09	BNE	EXCM6	
01990	1D86	CE	1DEE	LDX	#TURNPT	
02000	1D89	FF	1CAD	STX	OVRBUF	
02010	1D8C	7E	2067	JMP	OVRLAY	
02020	1D8F	C1	36	EXCM6	CMP B	#'6
02030	1D91	26	09	BNE	EXCM7	
02040	1D93	CE	1DF6	LDX	#DISPI	
02050	1D96	FF	1CAD	STX	OVRBUF	
02060	1D99	7E	2067	JMP	OVRLAY	
02070	1D9C	C1	37	EXCM7	CMP B	#'7
02080	1D9E	26	09	BNE	EXCM8	
02090	1DA0	CE	2PCD	LDX	#S2BCD	
02100	1DA3	FF	2801	STX	DOS+1	
02110	1DA6	7E	2800	JMP	DOS	
02120	1DA9	C1	38	EXCM8	CMP B	#'8
02130	1DAB	26	03	BNE	EXCM9	
02140	1DAD	7E	F000	JMP	EXBUG	
02150	1DB0	C1	39	EXCM9	CMP B	#'9
02160	1DB2	26	09	BNE	EXCMS	
02170	1DB4	CE	1DDE	LDX	#OVLMSG	
02180	1DB7	ED	CA87	JSR	OUTPUT	
02190	1DBA	7E	CFDA	JMP	DOSTR2	
02200	1DBD	C1	53	EXCMS	CMP B	#'S
02210	1DBF	26	03	BNE	EXCEND	
02220	1DC1	ED	208E	JSR	SAVEFL	
02230	1DC4	7E	1D00	EXCEND	JMP	START
02240				*		
02250	1DC7	4E		NOCPR	FCC	/NOCPRS* /
02260	1DCF	54		TOLAN1	FCC	/TOLAN-A* /

```

02297 1DD7 54      TOLAN2 FCC      /TOLAN-B*/
02298 1DDF 44      DOWER FCC      /DOWER* /
02299 1DD7 49      INTER FCC     /INTERPL*/
02300 1DD7 54      TURNPT FCC    /TURNPT* /
02310 1DF6 44      DISPL FCC     /DISPLAY*/
02320 1DFE 1A      OVLMSG FCB     $1A
02330 1DFF 45      FCC          /ENTER "RUN FILENAME" TO LOAD & EXECUT
02340 1E25 4F      FCC          /OVERLAY "FILENAME"./
02350 1E38 0D0A    FDB          $0D0A
02360 1E3A 28      FCC          /(WARNING: OVERLAY MUST FIT BETWEEN /
02370 1E5D 30      FCC          /0100-1400)/
02380 1E67 0D0A    FDB          $0D0A,$0D0A,$003F,$2004
02390              *
02400              * END EKG-EXEC
02410              *
02420              *FUNCTION :CMDIN
02430              *INPUTS (REG) :NONE
02440              *OUTPUTS (REG):B
02450              *CALLS :OUTPUT,KEYBD0
02460              *DESTROYS :ALL REGS
02470              *PURPOSE :TO PRINT PROMPT MESSAGE AND READ COMMAND
02480              * FROM CONSOLE.
02490              *
02500 1E6F CE 1E7F CMDIN LDX      #CMDMSG  LOAD OUT PTR
02510 1E72 BD CA8F   JSR        OUT1CR
02520 1E75 BD CA36   JSR        KEYBD0
02530 1E78 FE 3002   LDX        ENDBUF   SET PRT TO LAST CHAR ENTERED
02540 1E7B 09       DEX
02550 1E7C E6 00     LDA B  0,X
02560 1E7E 39       RTS
02570              *
02580 1E7F 1A07      CMDMSG FDB     $1A07,$0D0A,$0A0A,$0A0A,$0A0A
02590 1E89 45        FCC          /EKG DATA ACQUISITION SYSTEM/
02600 1EA4 0D0A      FDB          $0D0A,$0D0A
02610 1EA8 43        FCC          /COMMAND OPTIONS:/
02620 1EB8 0D0A      FDB          $0D0A,$0D0A
02630 1EBC 44        FCC          /DATA ACQUISITION/
02640 1ECC 0D0A      FDB          $0D0A,$0D0A
02650 1ED0 20        FCC          / 0=STORAGE WITH NO COMPRESSION/
02660 1EEE 0D0A      FDB          $0D0A
02670 1EF0 20        FCC          / 1=COMPRESSION VIA TOLAN ALGORITHM A
02680 1F14 0D0A      FDB          $0D0A
02690 1F16 20        FCC          / 2=COMPRESSION VIA TOLAN ALGORITHM B
02700 1F3A 0D0A      FDB          $0D0A
02710 1F3C 20        FCC          / 3=COMPRESSION VIA DOWER ALGORITHM/
02720 1F5E 0D0A      FDB          $0D0A
02730 1F60 20        FCC          / 4=COMPRESSION VIA 2ND ORDER INTERPL
02740 1F89 0D0A      FDB          $0D0A
02750 1F8B 20        FCC          / 5=COMPRESSION VIA TURNING POINT ALG
02760 1FB1 0D0A      FDB          $0D0A,$0D0A
02770 1FB5 44        FCC          /DATA DISPLAY/
02780 1FC1 0D0A      FDB          $0D0A,$0D0A
02790 1FC5 20        FCC          / 6=PRINT DATA STATISTICS OR /
02800 1FE1 44        FCC          /DISPLAY DATA/

```

```

02810 1FED CD0A      FDB      $0D0A,$0D0A
02820 1FF1 50        FCC      /PROGRAM CONTROL/
02830 2000 CD0A      FDB      $CD0A,$CD0A
02840 2004 20        FCC      / 7=RETURN TO DOS/
02850 2014 0D0A      FDB      $0D0A
02860 2016 20        FCC      / 8=RETURN TO EXBUG/
02870 2028 0D0A      FDB      $0D0A
02880 202A 20        FCC      / 9=LOAD & EXECUTE USER ENTERED OVERL
02890 2050 0D0A      FDB      $0D0A,$0D0A
02900 2054 45        FCC      /ENTER COMMAND NOW=/
02910 2066 04        FCB      4
02920
02930                *
02940                * END CMDIN
02950                *
02960                *FUNCTION :OVLAY
02970                *INPUTS (REG) :X
02980                *OUTPUTS (REG) :NONE
02990                *CALLS :CLRNAM,PUTNAM,OVLAY AT 0,X
03000                *DESTROYS :ALL REGISTERS
03010                *PURPOSE :TO LOAD AND THEN EXECUTE OVERLAYS
03020                * ASSOCIATED WITH PKG-EXEC &
03030                * SPECIFIED BY THE INPUT COMMAND.
03040                *
03050                *
03060                *
03070 2067 BD 20E1 OVLAY JSR      CLRNAM      CLR NAME BUFFER
03080 206A BD CC87      JSR      CLRPAS      CLEAR PASSWORD BUFFER
03090 206D CE 3C66      LDX      #NAME      PUT DOS FILENAME BUF IN NAMPT
03100 2070 FF 1C9B      STX      NAMPTR
03110 2073 FE 1CAD      LDX      OVRBUF      POINT X TO OVLAY NAME
03120 2076 BD 20EF      JSR      PUTNAM      PUT OVLAY NAME IN NAME
03130 2079 86 22        LDA A      #$22      DEFINE FILE TYPE
03140 207B B7 3057      STA A      TYPE
03150 207E BD CE52      JSR      LOAD2      LOAD FILE POINTED TO BY STX
03160 2081 86 AA        LDA A      #$AA      SET UP COMPARE FOR OVLAY FCN
03170 2083 B1 1C9F      CMP A      OVLGO      IS FLAG TRUE?
03180 2086 27 05        BEQ      OVRRTS     YES. TREAT OVLAY AS SUBR
03190 2088 8E A049      LDS      #$A049     NO. TREAT OVLAY AS ABS JUMP
03200 208B 6E 00        JMP      0,X      JUMP TO OVLAY
03210 208D 39          OVRRTS RTS
03220                *
03230                * END OVLAY
03240                *
03250                *
03260                *FUNCTION :SAVEFL
03270                *INPUTS (REG) :NONE
03280                *OUTPUTS (REG) :NONE
03290                *CALLS :CLRPAS,DRIVE,WLIB,WRITE0,OUTPUT
03300                *DESTROYS :ALL REGISTERS
03310                *PURPOSE :TO SAVE PKG DATA FILES RESIDING IN
03320                *IN MEMORY LOCATIONS 3A00-7FFF ON DISK. FILE NAME IS
03330                * PICKED OFF OF FILE HEADER AT LOC 3A02-3A03.
03340                *

```

```

03350
03360 2084 8D 51      *      SAVEFL NSR      CLRNAM      CLR NAME BUFFER
03370 208D C7 2141      LDX      #NAME      GET NAME COUNT
03380 2093 8D 4F      JSR      CLRNAM      CLEAR NAME
03390 2095 B6 3620      LDA A      S3620      GET NEW LOC FOR SAVE FROM PAS
03400 2098 B7 34C2      STA A      BASSAV      SAVE IN BUFFER
03410 209B CE 2141      LDX      #NAMEMSG      SET UP PTR TO FRT NAME
03420 209E FF 1C9B      STX      NAMEPTR
03430 20A1 CE 3402      LDX      #HDRSTR+2      POINT X TO FILENM IN DATA
03440 20A4 8D 49      JSR      PUTNAM      PUT FILENM IN "NAME" BUF
03450 20A6 86 22      LDA A      #S22      SET TYPE OF FILE
03460 20A8 B7 3057      STA A      TYPE
03470 20AB CE 3400      LDX      #HDRSTR      LOAD BUFFER START
03480 20AE FF 3058      STX      STARTX
03490 20B1 CE 7FFF      LDX      #BUFEND-1      LOAD BUFFER END
03500 20B4 FF 305A      STX      ENDX
03510 20B7 CE 2800      LDX      #S2800      LOAD LGO ADDRESS
03520 20BA FF 305C      STX      GOX
03530 20BD DF 1B      STX      PROGX
03540 20BF BD CC87      JSR      CLRPAS      CLR PASSWRD BUFFER
03550 20C2 CE 2119      LDX      #DRVMSG
03560 20C5 BD CA8F      JSR      OUTNCR
03570 20C8 BD E055      JSR      BYTE
03580 20CB BD C7C8      JSR      DRIVE      SELECT DRIVE 01 FOR STORE
03590 20CE BD C807      JSR      FLIB
03600 20D1 DE 08      LDX      ENEMH
03610 20D3 BD C75B      JSR      WRITE0      WRITE BUFFER TO DISK
03620 20D6 4F      CLR A
03630 20D7 BD C7C8      JSR      DRIVE      RESET DRIVE BACK TO 0
03640 20DA CE 2123      LDX      #CATMSG      LOAD OUTPUT POINTER
03650 20DD BD CA87      JSR      OUTPUT
03660 20E0 39      RTS
03670
03680 20E1 CE 3066      *      CLRNAM LDX      #NAME      POINT TO NAME BUFFER
03690 20E4 C6 08      CLRNA0 LDA B      #S8      SET UP CHAR COUNT
03700 20E6 86 20      LDA A      #S20      ASCII SPACE
03710 20E8 A7 00      CLRNA1 STA A      0,X      FILL NAME BUF WITH SPACES
03720 20EA 08      INX
03730 20EB 5A      DEC B
03740 20EC 26 FA      BNE      CLRNA1
03750 20EE 39      RTS
03760
03770 20EF BF 1C96      *      PUTNAM STS      STKSAV      SAVE STACK POINTER
03780 20F2 7F 1C9A      CLR      COUNT      CLEAR COUNTER
03790 20F5 35      TXS      TRANSFER X TO COUNT
03800 20F6 CE 3066      LDX      #NAME      PICK UP POS FROM NAME BUFFER
03810 20F9 C6 08      LDA B      #8      SET UP NAME BUFFER COUNTER
03820 20FB 32      PUTNAM LUL A      FILL CHAR FROM COUNT
03830 20FC 81 04      CIP A      #S04      IS IT 4?
03840 20FE 27 15      BQJ      PUTNA2      YES, NAME ENDED
03850 2100 A7 00      STA A      0,X      NO, STORE CHAR IN "NAME"
03860 2102 08      INX
03870 2103 FF 3060      STX      TEMPX1      SAVE INDEX PTR
03880 2106 FE 1C9B      LDX      NAMEPTR      LOAD CURRENT NAME PTR

```

```

03890 2109 A7 00      STA A 0,X      STORE IN NAMEMS
03900 210B 08          INX            INCR POINTER
03910 210C FF 1C93     STX          NAMEPTR
03920 210F FF 3060     LDX          TEMPX1  RESTORE INDEX REG
03930 2112 5A          DEC B
03940 2113 26 F6       RNE          PUTNA1  8 CHARS YET?
03950 2115 FE 1C96     PUTNA2 LDS      STKSAV  YES. RESTORE STACK
03960 2118 39          RTS
03970                  *
03980 2119 0D0A        DRVMSG FDB      $0D0A
03990 211B 44          FCC          /DRIVE? /
04000 2122 04          FCB          4
04010 2123 1A07        CATMSG FDB      $1A07,$0D0A,$0D0A
04020 2129 44          FCC          /DATA SAVED ON DISK FILE /
04030 2141 0008        NAMEMSG RMB      8
04040 2149 04          FCB          4
04050                  *
04060                  * END SAVEFL
04070                  *
04080                  *
04090                  *FUNCTION :FILHDR
04100                  *INPUTS (REG) :NONE
04110                  *OUTPUTS (REG) :NONE
04120                  *CALLS :OUTPUT,KEYBD0
04130                  *DESTROYS :ALL REGISTERS
04140                  *PURPOSE :TO INIT DATA BUFFER WITH ACCOUNT
04150                  * DATA (NAME,DATE,TIME,ETC) AND CLEAR FOR NEW
04160                  * DATA STORAGE.
04170                  *
04180 214A BD 21D7     FILHDR JSR      CLRBUF  FILL DATA BUF WITH NULS
04190 214D BD 21C9     JSR      CLRSEC  FILL MEM FILE SEC 3400 WITH A
04200 2150 BD 2288     JSR      PDFCLR  INITIALIZE COLL PARAM
04210 2153 FE 1C98     LDX      CPRTYP  PICK UP COLL TYPE
04220 2156 FF 3400     STX      HDRSTR  STR IN HEADER SECTOR
04230 2159 CE 21E5     LDX      #NAMEMS OUT PTR
04240 215C BD CA8F     JSR      OUTNCR  "NAME FILE PLEASE?"
04250 215F BD CA36     JSR      KEYBD0
04260 2162 C6 02       LDA B  $S02    HEADER STOR OFFSET
04270 2164 8D 42       RSR      TXTSTR  STORE FILNAME IN HEADER SEC
04280 2166 37          PSN B          SAVE NEXT AVAIL HEAD LOC
04290 2167 CE 2201     LDX      #SUBJMS
04300 216A BD CA8F     JSR      OUTNCR  "TEST SUBJECT ID?"
04310 216D BD CA36     JSR      KEYBD0
04320 2170 33          PUL B          RETRIEVE STOR OUTPUT
04330 2171 8D 35       RSR      TXTSTR
04340 2173 37          PSN B
04350 2174 CE 221F     LDX      #DATEMS
04360 2177 BD CA8F     JSR      OUTNCR  "COLLECTION SUBJECT DATE?"
04370 217A BD CA36     JSR      KEYBD0
04380 217D 33          PUL B
04390 217E 8D 28       RSR      TXTSTR
04400 2180 37          PSN B
04410 2181 CE 2247     LDX      #DATMSG
04420 2184 BD CA8F     JSR      OUTNCR  "DATE (05 JUL 80)?"

```

04470	2187	ED	CA36	JSR	KEYBD0	
04480	218A	33		PUL	B	
04490	218D	BD	1B	INR	TXTSTR	
04500	218F	37		PSH	B	
04470	218E	CE	225E	LDX	#TIMFMS	
04480	2191	BD	CA8F	JSR	OUTNCR	"TIME (1420)?"
04490	2194	BD	CA36	JSR	KEYBD0	
04500	2197	33		PUL	B	
04510	2198	8D	0E	BSR	TXTSTR	
04520	219A	37		PSH	B	
04530	219B	CE	2270	LDX	#CONTMS	
04540	219E	BD	CA8F	JSR	OUTNCR	"COMMENTS (MAX 80)?"
04550	21A1	ED	CA36	JSR	KEYBD0	
04560	21A4	33		PUL	B	
04570	21A5	8D	01	BSR	TXTSTR	
04580	21A7	39		RTS		
04590				*		
04600	21A8	BF	1C96	TXTSTR	STS	STKSAV SAVE STACK
04610	21AB	8E	3003	LDS	#BUFFER-1	POINT TO INPUT BUFFER
04620	21AF	CE	3400	LDX	#HDRSTR	POINT TO HEADER
04630	21B1	4F		CLR	A	SET A TO ZERO
04640	21B2	4C		TXTST1	INC	A
04650	21B3	08		INX		
04660	21B4	11		CBA		IS A = OFFSET?
04670	21B5	2D	FB	PLT	TXTST1	NO. KEEP INCREASING
04680	21B7	32		TXTST2	PUL	A
04690	21B8	81	04	CHP	A	#S04 IS IT 4
04700	21BA	27	06	BEQ	TXTSTR	YES. END TEXT, RETURN
04710	21BC	A7	00	STA	A	0,X NO. STORE CHAR IN HEADER
04720	21BE	5C		INC	B	
04730	21BF	08		INX		
04740	21C0	20	F5	BRA	TXTST2	
04750	21C2	A7	00	TXTSTR	STA	A
04760	21C4	5C		INC	B	
04770	21C5	BE	1C96	LDS	STKSAV	RESTORE STACK
04780	21C8	39		RTS		
04790	21C9	CE	3400	CLRSEC	LDX	#HDRSTR POINT TO DATA BUF START
04800	21CC	86	20	LDA	A	#S20 ASCII SPACE
04810	21CE	A7	00	CLRSEL	STA	A
04820	21D0	08		INX		
04830	21D1	8C	3500	CPX	#XPDF	IS SECTOR SPACED?
04840	21D4	26	F8	BNE	CLRSEL	COUNT > ZERO?
04850	21D6	39		RTS		NO. RTN
04860				*		
04870	21D7	CE	3400	CLRBUF	LDX	#HDRSTR POINT TO DATA BUF
04880	21DA	86	00	LDA	A	#C0 ASCII NULL
04890	21DC	A7	00	CLRBU1	STA	A
04900	21DE	08		INX		
04910	21DF	8C	3000	CPX	#BUFIND	IS INDEX REG AT END?
04920	21E2	26	F8	BNE	CLRBU1	NO. KEEP GOING
04930	21E4	39		RTS		YES. RETURN
04940				*		
04950	21E5	0D0A		NAMEMS	FDB	\$0D0A,\$0D0A
04960	21E9	4E		FCC		/NAME DATA FILE PLEASE? /


```

04970 2200 04          FCB      4
05000 2201 0D0A      SURIMS FDB      $0D0A,$0D0A
05010 221E 04          FCB      4
05010 221F 0D0A      RATESMS FDB      $0D0A,$0D0A
05020 2223 43          FCC      /COLLECTION SAMPLING RATE (500 HZ)? /
05030 2246 04          FCB      4
05040 2247 0D0A      DATEMS FDB      $0D0A,$0D0A
05050 224B 44          FCC      /DATE (05 JUL 80)? /
05060 225D 04          FCB      4
05070 225E 0D0A      TIMEMS FDB      $0D0A,$0D0A
05080 2262 54          FCC      /TIME (1423)? /
05090 226F 04          FCB      4
05100 2270 0D0A      CONTIMS FDB      $0D0A,$0D0A
05110 2274 43          FCC      /COMMENTS (MAX=80)? /
05120 2287 04          FCB      4
05130                  *
05140                  * END FILHDR
05150                  *
05160                  *FUNCTION :PDFCLR
05170                  *INPUTS (REG) :NONE
05180                  *OUTPUTS (REG):NONE
05190                  *CALLS :NOTHING
05200                  *DESTROYS :A,B,X,CC
05210                  *PURPOSE :TO INITIALIZE STATISTIC VAR & CLEAR
05220                  * PDF BUFFERS
05230                  *
05240 2288 4F          PDFCLR CLR A          FILL PARAM BUFF WITH 0
05250 2289 CE 3490      LDX      #LOOPCT
05260 228C A7 00      PDFCLO STA A      0,X
05270 228E 08          INX
05280 228F 8C 34C3      CPX      #ENTRPY  FILL WITH 0 UP TO ENTROPY BUF
05290 2292 26 F8          BNE      PDFCLO
05300 2294 CE 021F      LDX      #$021F  STORE MAX POSSIBLE BITS IN ME
05310 2297 86 F0          LDA A      #$F0
05320 2299 FF 34A9      STX      MEMBIT
05330 229C B7 34AB      STA A      MEMBIT+2
05340 229F 4F          CLR A
05350 22A0 CE 3500      LDX      #XPDI  NOW CLR PDF BUF AREA
05360 22A3 A7 00      PDFCLI STA A      0,X
05370 22A5 08          INX
05380 22A6 8C 3C00      CPX      #SECZRO
05390 22A9 26 F8          BNE      PDFCLI
05400 22AB 86 F0          LDA A      #$80  NOW SET UP INITIAL MAX & MINS
05410 22AD C6 7F          LDA B      #$7F
05420 22AF B7 3497      STA A      MAXZ
05430 22B2 B7 349D      STA A      MAXY
05440 22B5 B7 34A3      STA A      MAXX
05450 22B8 F7 349A      STA B      MINZ
05460 22BB F7 34A0      STA B      MINY
05470 22BE F7 34A6      STA B      MINX
05480 22C1 39          RTS
05490                  *
05500                  * END PDFCLR

```

```

06050 2327 DE 19      LDX      SNVLA      GET CUR PDF ADDR
06060 2329 A6 00      LDA A      0,X      LDA MSB BYTE OF PDF VLU
06070 232B FA 01      LDX B      1,X      LD LCA
06080 232D CA 23CD     LDX      #PDFMS
06090 2330 FF 1CA7     STX      HXBUF
06100 2333 BD 23DB     JSR      HXASC      CONV & STO MSB IN ASC STR
06110 2336 17         TBA
06120 2337 CE 23CF     LDX      #PDFMS5+2
06130 233A FF 1CA7     STX      HXBUF
06140 233D BD 23DB     JSR      HXASC      CONV & STORE LSB IN ASC STR
06150 2340 CE 23BB     LDX      #PDFMS3
06160 2343 BD CA87     JSR      OUTPUT     PRINT STRING
06170 2346 32         PUL A      RETRIEV VLU INDEX
06180 2347 DE 19      LDX      SAVEX     RETRIEC CUR PDF ADDR
06190 2349 39         RTS
06200                *
06210 234A 1A0C     PDFMS1 FDB      $1A0C,$0D0A,$0A0A,$0A0A,$0A0A
06220 2354 3A         FCC      /:FILE /
06230 235A 0008     PDFNAM RMB      8
06240 2362 0D0A     FDB      $0D0A,$0D0A
06250 2366 43         FCC      /CHANNEL /
06260 236E 0001     PDFMS2 RMB      1
06270 236F 20         FCC      / AMPLITUDE DISTRIBUTION/
06280 2386 0D0A     FDB      $0D0A,$0D0A
06290 238A 44         FCC      /DATA VALUE/
06300 2394 20         FCC      /          NUMBER OF OCCURENCES/
06310 23F6 0D0A     FDB      $0D0A,$0D0A
06320 23FA 04         FCB      4
06330 23BB 20         PDFMS3 FCC      / /
06340 23BF 0002     PDFMS4 RMB      2
06350 23C1 2E         FCC      /...../
06360 23CD 0004     PDFMS5 RMB      4
06370 23D1 20         FCC      / (HEX)/
06380 23D7 04         FCB      4
06390 23D8 0D23     PDFMS6 FDB      $0D23
06400 23DA 04         FCB      4
06410                *
06420                *
06430                * END PDFPRT
06440                *
06450                *FUNCTION :HXASC
06460                *INPUTS (REG) :A,X,DTA FOR CONV,LOC OF ASC STOR
06470                *OUTPUTS(REG) :NONE
06480                *CALLS :NOTHING
06490                *DESTROYS :A,CC
06500                *PURPOSE :TO CONVERT 1 BYTE OF HEX DATA
06510                * INTO 2 BYTES OF ASCII
06520                *
06530                *
06540 23DB FE 1CA7     HXASC     LDX      HXBUF
06550 23DE 36         PSH A      SAVE DTA ON STK
06560 23DF 44         LSR A      SHIFT TOP 4 BITS TO 1-4
06570 23E0 44         LSR A
06580 23E1 44         LSR A

```

```

05510      *
05520      *
05530      *FUNCTION :PRINT*
05540      *INPUTS (REG):NONE
05550      *OUTPUTS (REG):NONE
05560      *CALLS :OUTPUT
05570      *DESTROYS(REG):A,B,X,CC
05580      *PURPOSE :TO PRINT AMPLITUDE DISTRIBUTION
05590      * DATA COMPILED DURING DATA COLLECTION.
05600      *
05610 22C2 CE 235A PDFPRT LDX #PDFNAM PUT CUR MEM FILE NAME IN MSG
05620 22C5 FF 1C9B STX NAMPTR
05630 22C8 CE 3402 LDX #HDESTR+2
05640 22CB BD 20EF JSR PUTNAM
05650 22CE CE 234A LDX #PDFMS1 LOAD MSG POINTER
05660 22D1 85 58 LDA A #'X
05670 22D3 B7 236E STA A PDFMS2 STOR ASCII X TO MSG
05680 22D6 BD CA87 JSR OUTPUT "CHANNEL X AMPLITUDE .."
05690 22D9 CE 3500 LDX #XPDF POINT TO X PDF DTA
05700 22DC 86 7F LDA A #$7F INIT COUNTER
05710 22DE 8D 25 BSR PDFPRT1 PRINT X PDF DATA
05720 22E0 CE 234A LDX #PDFMS1 STOR ASCII Y TO MSG
05730 22E3 86 59 LDA A #'Y
05740 22E5 B7 236E STA A PDFMS2
05750 22E8 BD CA87 JSR OUTPUT
05760 22EB CE 3700 LDX #YPDF POINT TO Y PDF DTA
05770 22EE 86 7F LDA A #$7F
05780 22F0 8D 13 BSR PDFPRT1
05790 22F2 CE 234A LDX #PDFMS1
05800 22F5 86 5A LDA A #'Z
05810 22F7 B7 236E STA A PDFMS2
05820 22FA BD CA87 JSR OUTPUT
05830 22FD CE 3900 LDX #ZPDF
05840 2300 86 7F LDA A #$7F
05850 2302 8D 01 BSR PDFPRT1
05860 2304 39 RTS
05870      *
05880 2305 4C PDFPRT1 INC A
05890 2306 81 7F CMP A #$7F 256 VALUES PRINT'D YET?
05900 2308 27 07 BEQ PDFPRT2 YES. EXIT
05910 230A BD 231B JSR PRTVLU NO. PRINT VALUE TO CONSOLE
05920 230D 08 INX
05930 230E 08 INX
05940 230F 20 F4 BRA PDFPRT1 PRINT NEXT VLU
05950 2311 BD 231B PDFPRT2 JSR PRTVLU
05960 2314 CE 23D8 LDX #PDFMS6
05970 2317 BD CA8F JSR OUTINCR
05980 231A 39 RTS
05990      *
06000 231B 36 PRTVLU PSH A SAVE VALUE INDEX IN A
06010 231C DF 19 STX SAVEX SAVE LOC POINTED TO BY X
06020 231E CE 23BF LDX #PDFMS4 PICK UP STRING PTR
06030 2321 FF 1CA7 STX HXBUF
06040 2324 BD 23DB JSR HXASC CONV VLU INDEX TO ASCII & STR

```

EKG-EX-8

06590	23E2	44	LSR A		
06600	23E3	8D 08	BSR	HXTOAS	CONV SHFT'D BITS TO ASC
06610	23E5	A7 00	STA A	0,X	STR IN ALPH STRING
06620	23E7	32	PUL A		
06630	23E8	8D 03	BSR	HXTOAS	CONV LOW 4 TO ASC
06640	23EA	A7 01	STA A	1,X	STR IN ALPH STRING
06650	23EC	39	RTS		
06660			*		
06670	23ED	84 0F	HXTOAS	AND A	#\$0F CLR TOP 4 BITS
06680	23EF	81 0A	CMP A	#\$0A	IS NUM GE HEX A?
06690	23F1	2C 03	BGE	HXTOA1	YES. BRANCH & ADD BIAS OF 37
06700	23F3	8B 30	ADD A	#\$30	NO. ADD BIAS OF 30 & RET
06710	23F5	39	RTS		
06720	23F6	8B 37	HXTOA1	ADD A	#\$37
06730	23F8	39	RTS		
06740			*		
06750			*	END HXASC	
06760			*		
06770			*	END OF EKG-EXEC ROUTINES	
06780			*		
06790				END	

```

00070
00080
00090
00100
00110
00120
00130
00140
00150
00160
00170
00180
00190
00200
00210
00220
00230
00240
00250
00260
00270
00280
00290
00300
00310
00320
00330
00340
00350
00360
00370
00380
00390
00400
00410
00420
00430
00440
00450
00460
00470
00480
00490
00500
00510 0100
00520
00530
00540
00550
00560

```

```

*****
*
*   OVERLAY NAME      :DISPLAY
*   ALIAS             :CH1 CH2 EXTRACTED
*   VERSION            :1.5
*   VERSION DATE       :4 NOV 80
*
*****
*
*           PROGRAM DESCRIPTION
*
*   THIS PROGRAM DISPLAYS THE STATISTICAL DATA
*   CALCULATED DURING THE COLLECTION OF AN EKG
*   RECORDING. THE USER HAS SEVERAL MODES OF DATA
*   OUTPUT AND THEY ARE LISTED IN THE COMMAND STRING
*   BELOW. THIS PROGRAM ALSO PRODUCES PROBABILITY
*   DENSITY FUNCTION (PDF) PLOTS OF THE DATA COLLECTED
*   ON CHANNEL X,Y,OR Z VIA D/A CONVERTERS. THE DATA
*   (PDF) IS DISPLAYED ON AND OSCILLISCOPE.
*   IN ADDITION, THE COMMAND STRING ALLOWS RECONSTRUCT
*   AND DISPLAY OF PREVIOUSLY COLLECTED DATA ON THE
*   OSCILLISCOPE. THE RECONSTRUCTION ROUTINES ARE
*   OVERLAYED OVER THE PRINT STATISTICS (PRSTAT)
*   ROUTINES.
*
*   COMMAND OPTIONS
*
*       0= RETURN TO EKG-EXEC
*       1= PRINT CURRENT MEMORY FILE STATISTICS
*       2= PRINT PDF AMPLITUDE TABLES
*       3= DISPLAY XPDF ON THE OSCILLISCOPE
*       4= DISPLAY YPDF ON THE OSCILLISCOPE
*       5= DISPLAY ZPDF ON THE OSCILLISCOPE
*       6= RECONSTRUCT & DISPLAY CH X ON
*           THE OSCILLISCOPE
*       7= RECONSTRUCT & DISPLAY CH Y ON
*           THE OSCILLISCOPE
*       8= RECONSTRUCT & DISPLAY CH Z ON
*           THE OSCILLISCOPE
*       9= LOAD MEMORY FILE FROM DISK
*
*****
*
*           START OF DISPL
*
*****
*
*       ORG      $0100      PROGRAM ORIGIN
*
*       OPT      O          ASSB OPT. LIST ASSEMBLY
*       OPT      NOG        ASSB OPT. SUPPRESS FCC LIST
*
*****

```

```

00570
00580
00590
00600
00610
00620      CA8F      OUTNCR EQU      SCA8F      EPRODOS. OUTPUT STRING WITH
00630      CA36      KEYPDO EQU      SCA36      EPRODOS. INPUT ALPH STRING
00640      1D00      START EQU      $1D00      EKG-NEW C. START ADDRESS
00650      CCA0      MAKNAE EQU      SCCA0      EPRODOS. GET FILE NAME FROM
00660      C803      RLIB EQU      $C803      EPRODOS. SET UP READ I/O
00670      2800      DOS EQU      $2800      DOS. RUSTART LOC
00680      C43C      READ0 EQU      $C43C      EPRODOS. READ DISK ROUTINE
00690      E055      BYTE EQU      $E055      MIKBUG. INPUT 2 HEX DIG FROM T
00700      C7C8      DRIVE EQU      $C7C8      EPRODOSIO. READ DISK ROUTINE
00710      0500      PRSTAT EQU      $0500      PRINT COLLEC STAT ROUTINE
00720
00730      * DATA BUFFERS
00740      *
00741      0020      CHNLEF EQU      $0020      CHANNEL DESIGNATOR FOR DECPRE
00750      3400      HDRSTR EQU      $3400      MEM FILE HDR ADDR
00760      3500      XPDF EQU      $3500      LOC OF CH X PDF BUFFER
00770      3700      YPDF EQU      $3700      LOC OF CH Y PDF BUFFER
00780      3900      ZPDF EQU      $3900      LOC OF CH Z PDF BUFFER
00790      FFF8      IRQVEC EQU      $FFF8      LOC OF INT VECTOR ADDRESS
00800      1C9D      VFC SAV EQU      $1C9D      TEMP IRQVEC SAVE LOC
00810      3002      ENDBUF EQU      $3002      BUF WITH ADDR OF LAST CHAR IN
00820      1C96      STR SAV EQU      $1C96      TEMP STACK STORAGE BUFFER
00830      0019      SAVEX EQU      $0019      TEMP INDPX REG STORAGE BUFFER
00840      000A      EMETH EQU      $000A      END MEM ADDR FOR DISK I/O
00850      001B      PROGX EQU      $001B      ERR VFC ADDR
00860      1C9F      OLAYCO EQU      $1C9F      OVERLAY MODE FLAG
00870      3057      TYPE EQU      $3057      FILE TYPE BUFFER
00880      1CA9      PDFPLC EQU      $1CA9      PDFPRT ADDR PASS BUFFER
00890      1CAB      OVRLLC EQU      $1CAB      OVRLAY ADDR PASS BUFFER
00900      1CAD      OVRBUF EQU      $1CAD      OVRLAY PARAMETER BUFFER
00910      *
00920      * HARDWARE ADDRESSES
00930      *
00940      E500      DACZRO EQU      $E500      D/A CONVERTER ZERO
00950      E400      ADCZRO EQU      $E400      A/D CONVERTER ZERO
00960      *
00970      *
00980      *FUNCTION      :DISPL
00990      *INPUTS (REG) :NONE
01000      *OUTPUTS (REG) :NONE
01010      *CALLS      :OUTNCR,KEYPDO
01020      *DESTROYS      :A,B,CC,X
01030      *PURPOSE      :DISPL IS THE COMMAND EXEC FOR PROMPT
01040      *              :COMMAND INPUT FROM THE TERMINAL AND T
01050      *              :DIRECTING EXECUTION OF THE FOLLOWING O
01060      *              :IONS.
01070      *
01080      *
01090 0100 CE 0100 DISPL ID# #DISPL

```

DISPLA-5

01100	0103	FF	2801	STX	DOS+1	
01110	0106	CE	025B	LDX	#DISPM5	
01120	0109	2D	0A0F	JSR	PDFPRT	"GET LAST CHAR INPUT FROM TERM"
01130	010C	ED	CA36	JSR	KEYBD0	GET LAST CHAR INPUT
01140	010F	FE	3002	LDX	ENDBUF	SET UP IN DATA BUF POINTER
01150	0112	09		DEX		
01160	0113	E6	00	LDA B	0,X	GET LAST CHAR INPUT FROM TERM
01170	0115	C1	30	CMP B	#'0	WAS IT 0?
01180	0117	26	03	BNE	DISPL1	NO. CHECK IF 1
01190	0119	7E	1D00	JMP	START	YES. RET TO EKG-EXEC
01200	011C	C1	31	DISPL1	CMP B	#'1
01210	011E	26	0C	BNE	DISPL2	NO. CHECK IF 2
01220	0120	CE	01E4	LDX	#PRSNAM	GET "PRSTAT*" TO SEND TO OVRL
01230	0123	FF	1CAD	STX	OVREBUF	STORE "PRSTAT" TO OVRLAY BUFF
01240	0126	7F	1C9F	CLR	OLAYGO	RESTE OLAYGO FLAG
01250	0129	7E	04D6	JMP	OVRLAY	
01260	012C	C1	32	DISPL2	CMP B	#'2
01270	012E	26	0F	BNE	DISPL3	NO. CHECK IF 3
01280	0130	ED	04DB	JSR	PDFPRT	PRINT PDF AMPLITUDE TABLES
01290	0133	CE	0463	LDX	#DISMS1	
01300	0136	ED	CA8F	JSR	OUTNCR	"PRESS RETURN"
01310	0139	BD	CA36	JSR	KEYBD0	
01320	013C	7E	0100	JMP	DISPL	
01330	013F	C1	33	DISPL3	CMP B	#'3
01340	0141	26	0F	BNE	DISPL4	NO. CHECK 4
01350	0143	CE	3500	LDX	#XPDF	GET ADDRESS OF X PDF DATA
01360	0146	BD	0470	JSR	MAXSCN	SCAN PDF & FIND SHIFT NECESSA
01370			*			TO ALLOW OUTPUT OF PDF VIA 12
01380	0149	BD	04A0	JSR	PDFTRN	TRANSFER X PDF TO WORK BUFFER
01390	014C	BD	04BA	JSR	SCALE	SCALE WORK BUFFER BY SHIFT FR
01400	014F	7E	01F4	JMP	PDFOUT	NOW DISPLAY X PDF TO OSCOPE.
01410	0152	C1	34	DISPL4	CMP B	#'4
01420	0154	26	0F	BNE	DISPL5	NO. CHECK 5
01430	0156	CE	3700	LDX	#YPDF	GET ADDRESS OF Y PDF DATA
01440	0159	BD	0470	JSR	MAXSCN	
01450	015C	BD	04A0	JSR	PDFTRN	
01460	015F	BD	04BA	JSR	SCALE	
01470	0162	7E	01F4	JMP	PDFOUT	
01480	0165	C1	35	DISPL5	CMP B	#'5
01490	0167	26	0F	BNE	DISPL6	NO. RET & PRINT COMMAND PROMP
01500	0169	CE	3900	LDX	#ZPDF	
01510	016C	BD	0470	JSR	MAXSCN	
01520	016F	BD	04A0	JSR	PDFTRN	
01530	0172	BD	04BA	JSR	SCALE	
01540	0175	7E	01F4	JMP	PDFOUT	
01550	0178	C1	36	DISPL6	CMP B	#'6
01560	017A	26	0F	BNE	DISPL7	
01561	017C	7F	0020	CLR	CHNLRF	SET CHNL DESIGNATOR FOR CH X
01570	017F	7F	1C9F	CLR	OLAYGO	
01580	0182	CE	01EC	LDX	#DECPRS	LOAD & RUN DECPRS ROUTINE
01590	0185	FF	1CAD	STX	OVREBUF	STORE "DECPRS" TO OVRLAY BUFF
01600	0188	7E	04D6	JMP	OVRLAY	
01610	018B	C1	37	DISPL7	CMP B	#'7
01620	018D	26	10	BNE	DISPL8	

DISPLA-5

01621	018F	86	01	LDA A	#1	SET CHNL DESIGNATOR FOR CH Y
01622	0191	97	20	STA A	CHNLBF	
01623	0193	7F	1C9F	CLR	OLAYGO	
01624	0195	01	019C	INC		INCR INTR. ADDR. FOR CH Y
01650	0199	FF	1CAD	STX	OVREBUF	STORE "DECPRS" TO OVERLAY BUFF
01660	019C	7E	04D6	JMP	OVRLAY	
01670	019F	C1	38	DISPL8	CMP B	#'8
01680	01A1	26	10	BNE	DISPL9	
01681	01A3	86	02	LDA A	#2	SET CHNL DESIGNATOR FOR CH Z
01682	01A5	97	20	STA A	CHNLBF	
01690	01A7	7F	1C9F	CLR	OLAYGO	
01700	01AA	CE	01EC	LDX	#DECPRS	LOAD & RUN DECPRS ROUTINE
01710	01AD	FF	1CAD	STX	OVREBUF	STORE "DECPRS" TO OVERLAY BUFF
01720	01D0	7E	04D6	JMP	OVRLAY	
01730	01B3	C1	39	DISPL9	CMP B	#'9
01740	01B5	26	2A	BNE	CMDERR	
01750	01E7	86	22	LDA A	#\$22	SET UP FILE TYPE FOR LOAD
01760	01B9	B7	3057	STA A	TYPE	
01770	01BC	CE	024C	LDX	#NATMSG	"ENTER FILENAME?"
01780	01BF	BD	CA8F	JSR	OUTNCR	
01790	01C2	BD	CCA0	JSR	MAKHA	
01800	01C5	ED	C803	JSR	RLIB	
01810	01C8	CE	FFFF	LDX	#\$FFFF	
01820	01CB	DF	0A	STX	EMEMH	
01830	01CD	CE	2800	LDX	#DOS	
01840	01D0	DF	1B	STX	PROGX	
01850	01D2	CE	3400	LDX	#HDRSTR	PUT ADDR TO LOAD IN INDEX
01860	01D5	ED	C43C	JSR	READ0	
01870	01D8	CE	0456	LDX	#FILEMS	
01880	01DB	BD	CA8F	JSR	OUTNCR	
01890	01DE	BD	CA36	JSR	KEYED0	
01900	01E1	7E	0100	CMDERR	JMP	DISPL
01910				*		
01920	01E4	50		PRSNAM	FCC	/PRSTAT* /
01930	01EC	44		DECPRS	FCC	/DECPRS* /
01960				*		
01970	01F4	86	AA	PDFOUT	LDA A	#SAA
01980	01F6	B7	0248	STA A	DONTST	SET INTERRUPT DONE TEST FLAG
01990	01F9	0E		PDFOU1	CLI	CLEAR INTERRUPT MASK FOR STOP
02000	01FA	7F	024B	CLR	KOUNT	CLR BUFFER PRT COUNTER
02010	01FD	FE	FFF8	LDX	IRQVLC	GET CUR IRQ VECTOR
02020	0200	FF	1C9D	STX	VECSAV	SAVE IN TRIP BUFFER
02030	0203	CE	0240	LDX	#STOP	LOAD ADDRESS OF STOPPING INTR
02040	0206	FF	FFF8	STX	IRQVEC	
02050	0209	CE	07FF	LDX	#\$07FF	SEND MAX POS PULSE FOR SCOPE
02060	020C	FF	E500	STX	DACZRO	
02070	020F	CE	3200	LDX	#\$3200	GET START ADDR OF BUFFER
02080	0212	7C	024B	PDFCU2	INC	KOUNT
02090	0215	27	0E	BEQ	PDFOU3	BUMP UP PTR COUNTER
02100	0217	A6	00	LDA A	0,X	HAS WHOLE BUFFER BEEN USED
02110	0219	E6	01	LDA B	1,X	NO. GET MSB BYTE OF INDEXED D
02120	021B	B7	E500	STA A	DACZRO	GET LSB
02130	021E	F7	E501	STA B	DACZRO+1	SEND TO D/A CONV
02140	0221	08		INX		SEND TO D/A CONV
						INC INDEXED ADDRESS


```

02150 0222 08          INX
02160 0223 20 ED      BRA    PDFOU2    LOOP AGAIN
02170 0224 20 02A8     PDFOU3    LOOP AGAIN
02180          *
02190 0228 27 02      BEQ    PDFOU4    YES. EXIT DISPLAY LOOP
02200 022A 20 CD      BRA    PDFOU1    NO. KEEP SENDING DISPLAY
02210 022C 0F          PDFOU4 SEI
02220 022D CE 4000    LDX    #S4000    RESET INTR GATE.
02230 0230 FF E5C0    STX    DACZRO
02240 0233 B7 E400    STA    A    ADCZRO    RESET ST6800 INT FLIP FLOP
02250 0236 01          NOP
02260 0237 FE 1C9D    LDX    VECSAV    GET ORIG INT VEC ADDRESS
02270 023A FF FFF8    STX    IRQVEC    PUT BACK IN INT VEC BUFFER
02280 023D 7E 0100    JMP    DISPL    RET & PRINT COMMAND PROMPT
02290          *
02300 0240 7F 0248    STOP    CLR    DONTST    CLEAR DONE TEST FLAG
02310 0243 B7 E400    STA    A    ADCZRO    CLEAR ST6800 IN FLIP FLOP
02320 0246 01          NOP
02330 0247 3B          RTI
02340 0248 C001      DONTST RMB 1    RETURN FROM INTERRUPT AND STO
02350 0249 C001      SHFCNT RMB 1    INTR DONE TEST FLAG
02360 024A C001      MAXCNT RMB 1    BUF WITH NUM OF SHP FOR 12 BI
02370 024B C001      COUNT RMB 1    TEMP BUF TO COUNT 256 THRU 10
02380          *
02390 024C 45          NAMEMSG FCC    /ENTER /
02400 0252 04          FCB    4
02410 0253 44          DRVMSG FCC    /DRIVE? /
02420 025A 04          FCB    4
02430          *
02440 025B 1A07      DISPM5 FDB    S1A07,$0C0D,$0A0A,$0A0A,$0A0A,$0A0A
02450 0267 0D0A      FDB    $0D0A
02460 0269 44          FCC    /DATA DISPLAY/
02470 0275 0D0A      FDB    $0D0A,$0D0A
02480 0279 43          FCC    /COMMAND OPTIONS:/
02490 0289 0D0A      FDB    $0D0A,$0D0A
02500 028D 20          FCC    /    0=RETURN TO ERG-EXEC/
02510 02A6 0D0A      FDB    $0D0A,$0D0A
02520 02AA 20          FCC    /    1=PRINT CURRENT MEMORY FILE STA
02530 02D5 0D0A      FDB    $0D0A,$0D0A
02540 02D9 20          FCC    /    2=PRINT PDF AMPLITUDE TABLES/
02550 02FA 0D0A      FDB    $0D0A,$0D0A
02560 02FE 20          FCC    /    3=DISPLAY CH X PDF ON OSCILLISC
02570 0325 0D0A      FDB    $0D0A,$0D0A
02580 0329 20          FCC    /    4=DISPLAY CH Y PDF ON OSCILLISC
02590 0350 0D0A      FDB    $0D0A,$0D0A
02600 0354 20          FCC    /    5=DISPLAY CH Z PDF ON OSCILLISC
02610 037B 0D0A      FDB    $0D0A,$0D0A
02620 037F 20          FCC    /    6=RECONSTRUCT & DISPLAY CH X ON
02630 03A4 4F          FCC    /OSCILLOSCOPE/
02640 03P0 0D0A      FDB    $0D0A,$0D0A
02650 03F4 20          FCC    /    7=RECONSTRUCT & DISPLAY CH Y ON
02660 03D9 4F          FCC    /OSCILLOSCOPE/
02670 03E5 0D0A      FDB    $0D0A,$0D0A
02680 03E9 20          FCC    /    8=RECONSTRUCT & DISPLAY CH Z ON

```

```

02690 040E 4F      FCC      /OSCILLOSCOPE/
02700 041A 0D0A     FDB      $0D0A,$0D0A
02710 041F 20      DEC      /SELECT MEMORY FROM DISK/
02720 0427 000A     DEC      /COUNT/
02730 0443 45      FCC      /ENTER COMMAND NO.=/
02740 0455 04      FCB      4
02750 0456 46      FILEMS FCC /FILE LOADED. /
02760 0463 50      DISMS1 FCC /PRESS RETURN/
02770 046F 04      FCB      4
02780              *
02790              *
02800      *FUNCTION      :MAXSCN
02810      *INPUTS (REG)   :X
02820      *OUTPUTS        :SHIFT REQUIRED IN SHFCNT
02830      *CALLS          :NOTHING
02840      *DESTROYS        :A,B,CC
02850      *PRUPOSE          :THIS ROUTINE SCANS THE PDF BUFFER PO
02860      *                :TO BY X AND RETURNS THE MINIMUM
02870      *                :SHIFT RIGHTS NECESSARY TO BRING THE
02880      *                :AMPLITUDE COUNTS DOWN TO 11 BITS MAG
02890      *                :FOR OUTPUT FROM THE ST6800 D/A CONVE
02900      *
02910      *
02920 0470 DF 19    MAXSCN STX    SAVEX    SAVE PDF ADDR IN X
02930 0472 09      DEX          ADJUST X FOR PROPER INDEXED L
02940 0473 09      DEX
02950 0474 7F 0249  CLR          SHFCNT
02960 0477 7F 024A  CLR          MAXCNT
02970 047A 7C 024A  MAXSC1 INC     MAXCNT    256 BYTES SCANNED YET?
02980 047D 27 1E    BEQ          MAXSC4    YES. EXIT WITH CORRECT SHIFT
02990 047F 08      INX          NO. INC X & LOOK AT NEXT 2 BY
03000 0480 08      INX
03010 0481 A6 00    LDA A 0,X    GET MSB BYTE OF VALUE
03020 0483 81 08    CMP A #08    IS IT LESS THAN 8?
03030 0485 2D F3    BLT          MAXSC1    YES. NO SHIFT NECESS, GO CHEC
03040 0487 5F      CLR B        NO. SET UP COUNTER
03050 0488 0C      CLC          INSURE CARRY CLEAR
03060 0489 5C      MAXSC2 INC B   NOW SHIFT UNTIL MSB 1 BIT HIT
03070 048A 48      ASL A
03080 048B 25 02    BCS          MAXSC3    CARRY SET?
03090 048D 20 FA    BRA          MAXSC2    NO. KEEP SHIFTING LEFT
03100 048F 17      MAXSC3 TRA     YES. PUT SHIFT COUNT IN A
03110 0490 86 06    LDA A #6     GET MAX SHIFT
03120 0492 10      SRA          AND GET DIFFERENCE
03130 0493 B1 0249  CMP A SHFCNT IS SHIFT LE CURRENT SHIFT COU
03140 0496 2F E2    BLE          MAXSC1    YES. IGNORE THIS VALUE & CHEC
03150 0498 B7 0249  STA A SHFCNT NO. SAVE THIS SHIFT COUNT
03160 049B 20 DD    BRA          MAXSC1    NOW GO CHECK NEXT PDF WORD
03170 049E DE 19    MAXSC4 LDX    SAVEX    RETRIEVE INPUT X
03180 049F 39      RTS
03190      *
03200      *
03210      *FUNCTION      :PDFTRN
03220      *INPUTS        :X

```

DISPLA-5

```

03230      *OUTPUTS      :NONE
03240      *CALLS        :NOTHING
03250      *DESTROYS      :A,CC
03260      *PURPOSE        :THIS ROUTINE TRANSFERS THE PDF DATA
03270      *                :POINTED TO BY X INTO A WORK BUFFER
03280      *                :AREA AT 3200-33FF.
03290      *
03300      *
03310 04A0 DF 19 PDFTRN STX SAVEX  SAVE INDEX IN TEMP BUFFER
03320 04A2 BF 1C96      STS  STKSAV  SAVE STACK POINTER
03330 04A5 35          TXS          SET UP STACK FOR BUFF TRANSF
03340 04A6 CE 3200      LDX  #S3200  POINT TO WORK BUFF START
03350 04A9 8C 3400 PDFTR1 CPX #S3400  DONE TRANSFERING YET?
03360 04AC 27 06      BEQ  PDFTR2    YES. RET
03370 04AE 32          PUL A        NO. GET NEXT BYTE
03380 04AF A7 00      STA A 0,X     STORE IN WORK BUFFER
03390 04B1 08          INX
03400 04B2 20 F5      BRA  PDFTR1    GO GET NEXT BYTE TO TRANSFER
03410 04B4 DE 19 PDFTR2 LDX SAVEX  RETRIEVE ENTRY X
03420 04B6 BE 1C96      LDS  STKSAV  RETRIEVE STACK PTR
03430 04B9 39          RTS
03440      *
03450      *
03460      *FUNCTION      :SCALE
03470      *INPUTS        :DATA IN WORK BUFFER
03480      *OUTPUTS        :DATA IN WORK BUFFER SCALED BY SHFCNT
03490      *CALLS          :NOTHING
03500      *DESTROYS      :A,B,CC
03510      *PURPOSE        :THIS ROUTINE SCALES THE PDF DATA IN
03520      *                :THE WORK BUFFER (3200-33FF) BY THE
03530      *                :AMOUNT CALCULATED BY PARSON TO ENABL
03540      *                :THE ST6800 D/A CONVERTER TO OUTPUT
03550      *                :A POSITIVE VOLTAGE PROPORTIONAL TO
03560      *                :THE FREQUENCY OF OCCURENCE OF EACH
03570      *                :OF THE POSSIBLE 256 AMPLITUDE VALUES
03580      *
03590      *
03600 04BA DF 19 SCALE STX SAVEX  SAVE ENTRY X
03610 04BC CE 3200      LDX  #S3200  LOAD WORK BUFFER START ADDRES
03620 04BF 5F          CLR B
03630 04C0 BC 0249 SCALE1 IDA A SHFCNT  PICK UP SCALE SHIFT COUNT
03640 04C3 27 0E      BEQ  SCALE4    IF SHFCNT = 0, NO SCALING NEC
03650 04C5 27 07 SCALE2 BEQ  SCALE3    IS SHFCNT = 0?
03660 04C7 67 00      ASR  0,X       NO. SCALE DATA POINTED TO BY
03670 04C9 66 01      ROR  1,X
03680 04CB 4A          DEC A
03690 04CC 20 F7      BRA  SCALE2    GO SHIFT DOWN AGAIN
03700 04CE 08 SCALE3 INX          INCREMENT WORK BUFFER DATA TO
03710 04CF 08          INX
03720 04D0 5C          INC B
03730 04D1 26 ED      BNE  SCALE1     256 2 BYTE WORDS SHIFTED YET?
03740 04D3 DE 19 SCALE4 LDX SAVEX  NO. KEEP LOADING & SHIFTING
03750 04D5 39          RTS  RETRIEVE ENTRY X
03760      *

```

DL 11-11

```
03770 *FUNCTION :OVRPLAY
03780 *INPUTS :X
03790 *OUTPUTS :NONE
03800 *CALLS :ROUTINE POINTED TO BY X
03810 *DESTROYS :ALL REG
03820 *PURPOSE :THIS IS A VECTOR ROUTINE TO CALL OVRPLAY
03830 * EKG-EXEC TO LOAD ROUTINES FOR EXECUTION.
03840 *
03850 04D6 FE 1CAB OVRPLAY LDX OVRLLC
03860 04D9 6E 00 JMP 0,X
03870 *
03880 04DB FE 1CA9 PDFPRT LDX PDFPLC
03890 04DE 6E 00 JMP 0,X
03900 *
03910 *****
03920 *
03930 * END OF OVERLAY DISPLAY ROUTINES
03940 *
03950 *****
03960 *
03970 END
```

```

00000
00040
00070
00080
00090
00100
00110
00120
00130
00140
00150
00160
00170
00180
00190
00200
00210
00220
00230
00240
00250
00260
00270
00280
00290 0500
00300
00310
00320
00330
00340
00350
00360
00370
00380
00390
00400
00410
00420
00430
00440
00450
00460
00470
00480
00490
00500
00510
00520
00530
00540
00550
00560

*
*
* OVERLAY NAME : PRSTAT
* ZONE : 00000001 TO 00000001
* VERSION : 1:9
* VERSION DATE : 2 OCT 80
*
*
* PROGRAM DESCRIPTION
*
* THIS ROUTINE, WHEN LOADED AND EXECUTED, TAKES
* THE STATUS DATA CURRENTLY PRESENT IN THE
* MEMORY FILE HEADER SECTOR AND PRINTS THIS DATA
* TO THE CONSOLE. THIS ROUTINE CAN BE CALLED
* IMMEDIATELY FOLLOWING A DATA COLLECTION OR
* AFTER A PREVIOUSLY COLLECTED DATA FILE HAS BEEN
* LOADED INTO MEMORY FROM DISK.
*
*
*
* START OF PRSTAT
*
*
*
*      ORG      $0500      OVERLAY START ADDRESS
*
*      OPT      0          ASSB OPT-GIN ORJ FILE
*      OPT      NOG        ASSB OPT-SUPPRESS FCC LIST
*
*
*
* LABEL DESCRIPTIONS
*
* SUPPORT SUBROUTINE ADDRESSES
*
*
*      CA8F      OUTOCR EQU      $CA8F      EPROC.DOS. ALPH STRING TO CONS
*      CA36      KEYFDO EQU      $CA36      EPROC.DOS. ALPH STRING IN FRM
*      1D00      START EQU      $1D00      EKG-EXEC. SOFT EKG-EXEC START
*      2800      DOS EQU      $2800      DOS ERROR VECTOR, USED FOR JM
*
*
* DATA BUFFERS
*
*
*      1CA0      LGOFLG EQU      $1CA0      LOAD VERSUS EXECUTE FLAG
*      1CA5      HXASC EQU      $1CA5      HXASC ADDR PASS BUFFER
*      1CA7      HXPBF EQU      $1CA7      HXASC PARAMETER BUFFER
*      1C96      STKSAV EQU      $1C96      TEMP STACK SAVE BUFFER
*      3400      HDRSTR EQU      $3400      MEM FILE HEADER START
*      3490      LOOPCT EQU      $3490      TOTL CAL LOOPS EXEC IN WAIT
*      3494      SAMPLN EQU      $3494      NUM OF SAMPLES TAKEN
*      3496      LPCAL EQU      $3496      NUM OF CAL LOOPS PER 1 INTR
*      3497      MAXZ EQU      $3497      MAX VLU IN CH Z
*      3498      MAXZLO EQU      $3498      LOCATION OF MAX VLU IN CH Z

```

00570	349A	MINZ EQU	\$349A	MIN VLU IN CH Z
00580	349B	MINZLO EQU	\$349B	LOCATION OF MIN VLU IN CH Z
00590	349C	MAX EQU	\$349C	MAX VLU IN CH Y
00600	349D	MAXLO EQU	\$349D	LOCATION OF MAX VLU IN CH Y
00610	34A0	MINY EQU	\$34A0	MIN VLU IN CH Y
00620	34A1	MINYLO EQU	\$34A1	LOCATION OF MIN VLU IN CH Y
00630	34A3	MAXX EQU	\$34A3	MAX VLU IN CH X
00640	34A4	MAXXLO EQU	\$34A4	LOCATION OF MAX VLU IN CH X
00650	34A6	MINX EQU	\$34A6	MIN VLU IN CH X
00660	34A7	MINXLO EQU	\$34A7	LOCATION OF MIN VLU IN CH X
00670	34A9	MEMBIT EQU	\$34A9	NUM OF BITS AVAILABLE FOR STO
00680	34AC	DTABIT EQU	\$34AC	NUM OF BITS USED TO STOR DIA
00690	34B0	XBITS EQU	\$34B0	NUM OF BITS USED TO STOR X DT
00700	34B3	YBITS EQU	\$34B3	NUM OF BITS USED TO STOR Y DT
00710	34B6	ZBITS EQU	\$34B6	NUM OF BITS USED TO STOR Z DT
00720	34B9	TBITS EQU	\$34B9	NUM OF BITS USED TO STOR TIME
00730	34BC	ACELECT EQU	\$34BC	# BITS FED TO VAR LEN CODER
00740	34C2	BASEAV EQU	\$34C2	SAVE MEM LOC \$3620 FROM BASIC
00750	34C3	ENTREPY EQU	\$34C3	START OF ENTROPY BUFFER
00760	3460	MAXMIN EQU	\$3460	START OF MAX,MIN ASCII DATA
00770		*		
00780		*		
00790		*FUNCTION :PRSTAT		
00800		*INPUTS (REG) :NONE		
00810		*OUTPUTS (REG) :NONE		
00820		*CALLS :OUTPUT,FXASC		
00830		*DESTROYS (REG) :A,B,CC,X		
00840		*PURPOSE : TO LIST STATUS INFORMATION FROM DATA		
00850		* COLLECTION FILE TO TERMINAL.		
00860		*		
00870		*		
00880	0500 86 08	PRSTAT LDA A	#8	CLR MSG FIELDS IN STAT OUTPUT
00890	0502 CE 0A66	LDX	#NAME	
00900	0505 BD 0930	JSR	MSGCLR	
00910	0508 86 35	LDA A	#53	
00920	050A CE 0A8C	LDX	#SUBJ	
00930	050D BD 0930	JSR	MSGCLR	
00940	0510 86 0A	LDA A	#10	
00950	0512 CE 0ADF	LDX	#RATE	
00960	0515 BD 0930	JSR	MSGCLR	
00970	0518 86 14	LDA A	#20	
00980	051A CE 0B07	LDX	#DATE	
00990	051D BD 0930	JSR	MSGCLR	
01000	0520 86 08	LDA A	#8	
01010	0522 CE 0B39	LDX	#TIME	
01020	0525 BD 0930	JSR	MSGCLR	
01030	0528 86 19	LDA A	#25	
01040	052A CE 0B5F	LDX	#TYPE	
01050	052D BD 0930	JSR	MSGCLR	
01060	0530 86 35	LDA A	#53	
01070	0532 CE 0F29	LDX	#COMMENT	
01080	0535 BD 0930	JSR	MSGCLR	
01090	0538 86 0E	LDA A	#8	
01100	053A CE 0B96	LDX	#XINT	

01110 053D BD 0930	JSR	MSGCLR
01120 053E 86 08	LDA A	#8
01130 053F 86 08	LDA A	#8
01140 0545 BD 0930	JSR	MSGCLR
01150 0548 86 08	LDA A	#8
01160 054A CE 00FA	LDX	#ZINT
01170 054D BD 0930	JSR	MSGCLR
01180 0550 86 08	LDA A	#8
01190 0552 CE 002C	LDX	#TINT
01200 0555 BD 0930	JSR	MSGCLR
01210 0558 86 08	LDA A	#8
01220 055A CE 00CE	LDX	#MANCPR
01230 055D BD 0930	JSR	MSGCLR
01240 0560 86 08	LDA A	#8
01250 0562 CE 0D11	LDX	#CTACH
01260 0565 BD 0930	JSR	MSGCLR
01270 0568 86 05	LDA A	#5
01280 056A CE 0D77	LDX	#CTINTF
01290 056D BD 0930	JSR	MSGCLR
01300 0570 86 05	LDA A	#5
01310 0572 CE 0D99	LDX	#TINTF
01320 0575 BD 0930	JSR	MSGCLR
01330 0578 86 05	LDA A	#5
01340 057A CE 0D8B	LDX	#COMPR
01350 057D BD 0930	JSR	MSGCLR
01360 0580 86 08	LDA A	#8
01370 0582 CE 0F2F	LDX	#AMFN
01380 0585 BD 0930	JSR	MSGCLR
01390 0588 86 08	LDA A	#8
01400 058A CE 0D7D	LDX	#AMAXX
01410 058D BD 0930	JSR	MSGCLR
01420 0590 86 08	LDA A	#8
01430 0592 CE 0E93	LDX	#AMINX
01440 0595 BD 0930	JSR	MSGCLR
01450 0598 86 08	LDA A	#8
01460 059A CE 0F61	LDX	#AMANY
01470 059D BD 0930	JSR	MSGCLR
01480 05A0 86 08	LDA A	#8
01490 05A2 CE 0E93	LDX	#AMINX
01500 05A5 BD 0930	JSR	MSGCLR
01510 05A8 86 08	LDA A	#8
01520 05AA CE 0FC5	LDX	#AMAXX
01530 05AD BD 0930	JSR	MSGCLR
01540 05B0 86 08	LDA A	#8
01550 05B2 CE 0FF7	LDX	#AMINX
01560 05B5 BD 0930	JSR	MSGCLR
01570 05B8 CE 3402	LDX	#DIRECTR42
01580 05BB FF 003D	STX	#FUTAC
01590 05BE CE 0A66	LDX	#AMINX
01600 05C1 BD 0930	JSR	MSGCLR
01610 05C4 CE 0A66	LDX	#AMINX
01620 05C7 BD 0930	JSR	MSGCLR
01630 05CA CE 0A66	LDX	#AMINX
01640 05CD BD 0930	JSR	MSGCLR

01600 05D0 CE 0F07	LDX	#DATE	
01660 05D3 BD 093C	JSR	MSGPUT	
01680 05D5 CE 0F39	LDX	#TIME	
016A0 05D7 BD 093C	JSR	MSGPUT	
016E0 05DC CE 0F29	LDX	#COMMENT	
01700 05DF BD 093C	JSR	MSGPUT	
01710 05E2 FE 34C3	LDX	#ENTROPY	GET LOC OF ENTROPY BUFFER
01720 05E5 8C 2020	CPX	#\$2020	CHK IF ENTROPY CALC MADE?
01730 05E8 27 03	BEQ	PRST00	NO. STORE NOT CALC MSG
01740 05EA 7E 06A4	JMP	PRST00	YES. GET ASCII STRING FROM HD
01750 05ED CE 09FD PRST00	LDX	#ENTCLR	
01760 05F0 FF 095D	STX	FRMLC	
01770 05F3 CE 0E96	LDX	#XENT	
01780 05F6 BD 093C	JSR	MSGPUT	
01790 05F9 CE 09FD	LDX	#ENTCLR	
01800 05FC FF 095D	STX	FRMLC	
01810 05FF CE 0FC8	LDX	#YENT	
01820 0602 BD 093C	JSR	MSGPUT	
01830 0605 CE 09FD	LDX	#ENTCLR	
01840 0608 FF 095D	STX	FRMLC	
01850 060B CE 0EFA	LDX	#ZENT	
01860 060E BD 093C	JSR	MSGPUT	
01870 0611 CE 09FD	LDX	#ENTCLR	
01880 0614 FF 095D	STX	FRMLC	
01890 0617 CE 0C2C	LDX	#TEXT	
01900 061A BD 093C	JSR	MSGPUT	
01910 061D CE 09FD	LDX	#ENTCLR	
01920 0620 FF 095D	STX	FRMLC	
01930 0623 CE 0CCE	LDX	#MAXCPR	
01940 0626 BD 093C	JSR	MSGPUT	
01950 0629 CE 09FD	LDX	#ENTCLR	
01960 062C FF 095D	STX	FRMLC	
01970 062F CE 0D11	LDX	#CPRACH	
01980 0632 BD 093C	JSR	MSGPUT	
01990 0635 CE 09FD	LDX	#ENTCLR	
02000 0638 FF 095D	STX	FRMLC	
02010 063B CE 0D57	LDX	#CPREFE	
02020 063E BD 093C	JSR	MSGPUT	
02030 0641 CE 09FD	LDX	#ENTCLR	
02040 0644 FF 095D	STX	FRMLC	
02050 0647 CE 0D99	LDX	#TIMEFF	
02060 064A BD 093C	JSR	MSGPUT	
02070 064D CE 09FD	LDX	#ENTCLR	
02080 0650 FF 095D	STX	FRMLC	
02090 0653 CE 0DCB	LDX	#COLNUR	
02100 0656 BD 093C	JSR	MSGPUT	
02110 0659 CE 09FD	LDX	#ENTCLR	
02120 065C FF 095D	STX	FRMLC	
02130 065F CE 0D8D	LDX	#AMAXX	
02140 0662 BD 093C	JSR	MSGPUT	
02150 0665 CE 09FD	LDX	#ENTCLR	
02160 0668 FF 095D	STX	FRMLC	
02170 066B CE 0E2F	LDX	#AMINX	
02180 066E BD 093C	JSR	MSGPUT	

PRSTAT-9

02190	0671	CE	09FD	LDX	#ENTCLR	
02200	067A	FF	095D	STX	FRMLOC	
02210	0680	FF	095D	STX	FRMLOC	
02220	067A	BD	093C	JSR	MSGPUT	
02230	067D	CE	09FD	LDX	#ENTCLR	
02240	0680	FF	095D	STX	FRMLOC	
02250	0683	CE	0E93	LDX	#AMINX	
02260	0686	BD	093C	JSR	MSGPUT	
02270	0689	CE	091D	LDX	#ENTCLR	
02280	068C	FF	095D	STX	FRMLOC	
02290	068F	CE	0EC5	LDX	#AMAXZ	
02300	0692	BD	093C	JSR	MSGPUT	
02310	0695	CE	09FD	LDX	#ENTCLR	
02320	0698	FF	095D	STX	FRMLOC	
02330	069B	CE	0EF7	LDX	#AMINZ	
02340	069E	BD	093C	JSR	MSGPUT	
02350	06A1	7E	070A	JMP	PRSTAL	
02360	06A4	CE	34C3	LDX	#ENTRYP	GET ADDRESS OF ENTRYP BUFFER
02370	06A7	FF	095D	STX	FRMLOC	
02380	06AA	CE	0E96	LDX	#XENT	NOW PRINT APPROPRIATE DAT TO
02390	06AD	BD	093C	JSR	MSGPUT	
02400	06B0	CE	0FC8	LDX	#YENT	
02410	06B3	BD	093C	JSR	MSGPUT	
02420	06B6	CE	0EFA	LDX	#ZENT	
02430	06B9	BD	093C	JSR	MSGPUT	
02440	06BC	CE	0C2C	LDX	#TENT	
02450	06BF	BD	093C	JSR	MSGPUT	
02460	06C2	CE	0CCE	LDX	#MAXCPR	
02470	06C5	BD	093C	JSR	MSGPUT	
02480	06C8	CE	0D11	LDX	#CPRACH	
02490	06CB	BD	093C	JSR	MSGPUT	
02500	06CE	CE	0D57	LDX	#CPRFTF	
02510	06D1	BD	093C	JSR	MSGPUT	
02520	06D4	CE	0D99	LDX	#TIMEFF	
02530	06D7	BD	093C	JSR	MSGPUT	
02540	06DA	CE	0DCB	LDX	#COLDUR	
02550	06DD	BD	093C	JSR	MSGPUT	
02560	06E0	CE	3460	LDX	#MAXMIN	NOW GET ADDR WHERE ASCII MAX,
02570	06E3	FF	095D	STX	FRMLOC	
02580	06E6	CE	0DFD	LDX	#AMAXX	
02590	06E9	BD	093C	JSR	MSGPUT	
02600	06EC	CE	0E2F	LDX	#AMINX	
02610	06EF	BD	093C	JSR	MSGPUT	
02620	06F2	CE	0E61	LDX	#AMAXY	
02630	06F5	BD	093C	JSR	MSGPUT	
02640	06F8	CE	0E93	LDX	#AMINZ	
02650	06FB	BD	093C	JSR	MSGPUT	
02660	06FE	CE	0FC5	LDX	#AMAXZ	
02670	0701	BD	093C	JSR	MSGPUT	
02680	0704	CE	0EF7	LDX	#AMINZ	
02690	0707	BD	093C	JSR	MSGPUT	
02700	070A	FE	3400	LDX	HDRSTR	NOW DECODE CONTRS TYPE & PUT
02710	070D	8C	4E43	CPX	#S4E43	(HC) NOCPRS?
02720	0710	26	08	ENE	PRSTA2	NO

02740 0712 CE 09F9	LDX	#NOCCOUP	YES. POINT TO NOCCOUP MSG
02740 0715 FF 095D	STX	FRMLLOC	
02750 0718 20 42	BRA	PRSTA7	
02770 071D 26 08	BNE	PRSTA3	NO.
02780 071F CE 09C2	LDX	#TOLA	YES. POINT TO TOLA-A MSG
02790 0722 FF 095D	STX	FRMLLOC	
02800 0725 20 35	BRA	PRSTA7	
02810 0727 8C 5442 PRSTA3	CPX	#S5442	(TB) TOLAN-B?
02820 072A 26 08	BNE	PRSTA4	NO
02830 072C CE 09CA	LDX	#TOLB	YES. POINT TO TOLAN-B MSG
02840 072F FF 095D	STX	FRMLLOC	
02850 0732 20 28	BRA	PRSTA7	
02860 0734 8C 444F PRSTA4	CPX	#S444F	(DO) DOWER?
02870 0737 26 08	BNE	PRSTA5	NO
02880 0739 CE 09D2	LDX	#DOWR	YES. POINT TO DOWER MSG
02890 073C FF 095D	STX	FRMLLOC	
02900 073F 20 1B	BRA	PRSTA7	
02910 0741 8C 5450 PRSTA5	CPX	#S5450	(TP) TURNING POINT?
02920 0744 26 08	BNE	PRSTA6	NO.
02930 0746 CE 09D8	LDX	#TURNPT	YES. POINT TO TURNING POINT M
02940 0749 FF 095D	STX	FRMLLOC	
02950 074C 20 0E	BRA	PRSTA7	
02960 074E 8C 494E PRSTA6	CPX	#S494E	(IN) 2ND ORDER INTERPOLATOR?
02970 0751 27 03	BEQ	PRST61	YES. TYPE IS RECOGNIZED. CONT
02980 0753 7E 095F	JMP	ERROR	NO. TYPE NOT RECOGNIZED. EPR
02990 0756 CE 09E6 PRST61	LDX	#INTER	
03000 0759 FF 095D	STX	FRMLLOC	
03010 075C CE 0B5F PRSTA7	LDX	#TYPE	PUT COMERS TYPE IN MSG
03020 075F BD 093C	JSR	MSGPUT	
03030 0762 CE 0A12	LDX	#HDRMSG	NOW PT TO STATMSG PAGE 1
03040 0765 BD CA8F	JSR	OUTNCR	AND PRINT PAGE 1
03050 0768 BD CA36	JSR	KEYBD0	
03060 076B CE 0C6B	LDX	#HDRMS2	NOW PT TO STATMSG PAGE 2
03070 076E BD CA8F	JSR	OUTNCR	AND PRINT PAGE 2
03080 0771 BD CA36	JSR	KEYBD0	
03090 0774 CE 0FE4	LDX	#NUMEMP+2	PICK UP NUM OF SEPLS,
03100 0777 F6 3495	LDA A	SAMR10+1	CONVERT TO ASCII, & STORE
03110 077A BD 0955	JSR	HXASC	TO OUTPUT STRING.
03120 077D CE 0FE2	LDX	#NUMSMP	
03130 0780 B6 3494	LDA A	SAMR10	
03140 0783 BD 0955	JSR	HXASC	
03150 0786 CE 1019	LDX	#NUMLPS	PICK UP VLU OF CAL TST
03160 0789 B6 3496	LDA A	LPCAL	& STORE TO OUTPUT STRING
03170 078C BD 0955	JSR	HXASC	
03180 078F CE 1060	LDX	#TOTLUP+6	PICK UP TOTL LOOP COUNT
03190 0792 B6 3493	LDA A	LOOPCT+3	,CONVERT TO ASCII & STORE
03200 0795 BD 0955	JSR	HXASC	TO OUTPUT STRING
03210 0798 CE 105E	LDX	#TOTLUP+4	
03220 079B B6 3492	LDA A	LOOPCT+2	
03230 079E BD 0955	JSR	HXASC	
03240 07A1 CE 105C	LDX	#TOTLUP+2	
03250 07A4 B6 3491	LDA A	LOOPCT+1	
03260 07A7 BD 0955	JSR	HXASC	

03270 07AA CE 105A	LDX	#TOTLUP	
03280 07AB B6 3490	LDA A	LOOPCT	
03290 07AC BD 095E	JSR	HXASC	
03300 07AD CE 10C5	LDX	#XMAXNM	PICK UP MAX VLU ON CH X
03310 07AE B6 34A3	LDA A	MAXX	CONV TO ASCII & STR IN
03320 07AF BD 0955	JSR	HXASC	TO OUTPUT STRING
03330 07B0 CE 10E1	LDX	#XMAXLO+2	
03340 07B1 B6 34A5	LDA A	MAXXLO+1	
03350 07B2 BD 0955	JSR	HXASC	
03360 07C5 CE 10DF	LDX	#XMINLO	
03370 07C8 B6 34A4	LDA A	MAXXLO	
03380 07CB BD 0955	JSR	HXASC	
03390 07CE CE 10F2	LDX	#XMINNM	PICK UP MIN VLU ON CH X
03400 07D1 B6 34A6	LDA A	MINX	CONV TO ASCII & STORE IN
03410 07D4 BD 0955	JSR	HXASC	OUTPUT STRING
03420 07D7 CE 110E	LDX	#XMINLO+2	
03430 07DA B6 34A8	LDA A	MINXLO+1	
03440 07DD BD 0955	JSR	HXASC	
03450 07E0 CE 110C	LDX	#XMINLO	
03460 07E3 B6 34A7	LDA A	MINXLO	
03470 07E6 BD 0955	JSR	HXASC	
03480 07E9 CE 111F	LDX	#YMAXNM	PICK UP MAX VLU ON CH Y
03490 07EC B6 349D	LDA A	MAXY	CONV TO ASCII & STORE IN
03500 07EF BD 0955	JSR	HXASC	OUTPUT STRING
03510 07F2 CE 113B	LDX	#YMAXLO+2	
03520 07F5 B6 349F	LDA A	MAXYLO+1	
03530 07F8 BD 0955	JSR	HXASC	
03540 07FB CE 1139	LDX	#YMAXLO	
03550 07FE B6 349E	LDA A	MAXYLO	
03560 0801 BD 0955	JSR	HXASC	
03570 0804 CE 114C	LDX	#YMINNM	PICK UP MIN VLU ON CH Y
03580 0807 B6 34A0	LDA A	MINY	CONV TO ASCII & STORE IN
03590 080A BD 0955	JSR	HXASC	OUTPUT STRING
03600 080D CE 1168	LDX	#YMINLO+2	
03610 0810 B6 34A2	LDA A	MINYLO+1	
03620 0813 BD 0955	JSR	HXASC	
03630 0816 CE 1166	LDX	#YMINLO	
03640 0819 B6 34A1	LDA A	MINYLO	
03650 081C BD 0955	JSR	HXASC	
03660 081F CE 1179	LDX	#ZMAXNM	PICK UP MAX VLU ON CH Z
03670 0822 B6 3497	LDA A	MAXZ	CONV TO ASCII & STORE IN
03680 0825 BD 0955	JSR	HXASC	OUTPUT STRING
03690 0828 CE 1195	LDX	#ZMAXLO+2	
03700 082B B6 3499	LDA A	MAXZLO+1	
03710 082E BD 0955	JSR	HXASC	
03720 0831 CE 1193	LDX	#ZMAXLO	
03730 0834 B6 3498	LDA A	MAXZLO	
03740 0837 BD 0955	JSR	HXASC	
03750 083A CE 11A6	LDX	#ZMINNM	PICK UP MIN VLU ON CH Z
03760 083D B6 349A	LDA A	MINZ	CONV TO ASCII & STORE IN
03770 0840 BD 0955	JSR	HXASC	OUTPUT STRING
03780 0843 CE 11C2	LDX	#ZMINLO+2	
03790 0846 B6 349C	LDA A	MINZLO+1	
03800 0849 BD 0955	JSR	HXASC	

03810 084C CE 11C0	LDX #MINZLO
03820 084F B6 349B	LDA A MINZLO
03830 0851 BD 0955	JSR HXASC
03840 0853 CE 120E	LDX #MEMBIT+4 PICK UP NUM OF BITS
03850 0858 B6 34AB	LDA A MEMBIT+2 AVAILABLE FOR DATA STOR
03860 085B BD 0955	JSR HXASC , CONV TO ASCII & STOR
03870 085E CE 120C	LDX #AVABIT+2 IN OUTPUT STRING
03880 0861 B6 34AA	LDA A MEMBIT+1
03890 0864 BD 0955	JSR HXASC
03900 0867 CE 120A	LDX #AVABIT
03910 086A B6 34A9	LDA A MEMBIT
03920 086D BD 0955	JSR HXASC
03930 0870 CE 124B	LDX #TCPBIT+6
03940 0873 B6 34BF	LDA A ACELCT+3
03950 0876 BD 0955	JSR HXASC
03960 0879 CE 1249	LDX #TCPBIT+4
03970 087C B6 34BE	LDA A ACELCT+2
03980 087F BD 0955	JSR HXASC
03990 0882 CE 1247	LDX #TCPBIT+2
04000 0885 B6 34ED	LDA A ACELCT+1
04010 0888 BD 0955	JSR HXASC
04020 088B CE 1245	LDX #TCPBIT
04030 088E B6 34EC	LDA A ACELCT
04040 0891 BD 0955	JSR HXASC
04050 0894 CE 127F	LDX #TOTBIT+6 PICK UP NUM OF BITS
04060 0897 B6 34AF	LDA A DTABIT+3 ACTUALL STORID
04070 089A BD 0955	JSR HXASC BECAUSE OF DATA
04080 089D CE 127D	LDX #TOTBIT+4 COMPRESSION, CONV TO
04090 08A0 B6 34AE	LDA A DTABIT+2 ASCII & STORE IN
04100 08A3 BD 0955	JSR HXASC OUTPUT STRING
04110 08A6 CE 127B	LDX #TOTBIT+2
04120 08A9 B6 34AD	LDA A DTABIT+1
04130 08AC BD 0955	JSR HXASC
04140 08AF CE 1279	LDX #TOTBIT
04150 08B2 B6 34AC	LDA A DTABIT
04160 08B5 BD 0955	JSR HXASC
04170 08B8 CE 12B7	LDX #XBITS+4 PICK UP NUM OF BITS
04180 08BB B6 34B2	LDA A XBITS+2 ACTUALLY STORED FOR
04190 08BE BD 0955	JSR HXASC DATA IN CHANNEL X,
04200 08C1 CE 12B5	LDX #XBITS+2 CONV TO ASCII &
04210 08C4 B6 34B1	LDA A XBITS+1 STORE TO OUTPUT STRING
04220 08C7 BD 0955	JSR HXASC
04230 08CA CE 12B3	LDX #XBITS
04240 08CD B6 34B0	LDA A XBITS
04250 08D0 BD 0955	JSR HXASC
04260 08D3 CE 12EF	LDX #YBITS+4 PICK UP NUM OF BITS
04270 08D6 B6 34B5	LDA A YBITS+2 ACTUALLY STORED FOR
04280 08D9 BD 0955	JSR HXASC DATA IN CHANNEL Y,
04290 08DC CE 12DD	LDX #YBITS+2 CONV TO ASCII &
04300 08DF B6 34B4	LDA A YBITS+1 STORE IN OUTPUT STRING
04310 08E2 BD 0955	JSR HXASC
04320 08E5 CE 12EB	LDX #YBITS
04330 08E8 B6 34B3	LDA A YBITS
04340 08EB BD 0955	JSR HXASC

```

04350 08EE CE 1327      LDX      #ZBITNM+4  PICK UP NUM OF BITS
04360 08F1 B6 34B8      LDA A    ZBITS+2  ACTUALLY STORED FOR
04370 08F4 FD 0955      JSR      HXASC    DATA IN CHANNEL 1,
04380 08F7 C1 1325      LDX      #ZBITNM+2  CONV TO ASCII &
04390 08FA B6 34B7      LDA A    ZBITS+1  STORE IN OUTPUT STRING
04400 08FD BD 0955      JSR      HXASC
04410 0900 CE 1323      LDX      #ZBITNM
04420 0903 B6 34B6      LDA A    ZBITS
04430 0906 BD 0955      JSR      HXASC
04440 0909 CE 136D      LDX      #TBITNM+4  PICK UP NUM OF BITS
04450 090C B6 34EB      LDA A    TBITS+2  ACTUALLY STORED FOR
04460 090F BD 0955      JSR      HXASC    TIME (OR OTHER PARAMETER)
04470 0912 CE 136B      LDX      #TBITNM+2  CONV TO ASCII &
04480 0915 B6 34BA      LDA A    TBITS+1  STORE IN OUTPUT STRING
04490 0918 BD 0955      JSR      HXASC
04500 091B CE 1369      LDX      #TBITNM
04510 091E B6 34B9      LDA A    TBITS
04520 0921 BD 0955      JSR      HXASC
04530 0924 CE 0F89      LDX      #COSTAT  GET OUTPUT STRING ADDR
04540 0927 BD CA8F      JSR      OUTNCR  & SEND TO CONSOLE
04550 092A BD CA36      JSR      KEYBD0  WAIT FOR TERM INPUT
04560 092D 7E 2800      JMP      DOS      JUMP BACK TO CALLING ROUTINE
04570
04580      *
04590      *FUNC: MSGCLR
04600      *INPUTS: A (# OF SPACES TO CLR),X (LOC TO PUT SPACES)
04610      *OUTPUTS: ASCII $20 TO MEM AT X
04620      *CALLS: NOTHING
04630      *DESTROYS: A,B,CC
04640      *PURPOSE: THIS ROUTINE CLRS THE MSG BUFFER EACH TIME
04650      * PRSTAT IS CALLED FOR NEW INFO OFF OF THE MEM DATA
04660      *
04660 0930 4C      MSGCLR INC A      INC COUNTER
04670 0931 C6 20      LDA B    #$20    ASCII SPACE
04680 0933 4A      MSGCLR DEC A
04690 0934 27 05      BEQ      MSGCLR2  IS BUFFER CLRD?
04700 0936 E7 00      STA B    0,X      NO. KEEP LOOPING
04710 0938 08      INX
04720 0939 20 F8      BRA      MSGCLR1
04730 093B 39      MSGCLR2 RTS
04740      *
04750      * END OF MSGCLR
04760      *
04770      *****
04780      *
04790      *FUNC: MSGPUT
04800      *INPUTS: X (ADDR WHERE DATA IN FRMLLOC GOING)
04810      *OUTPUTS: ASCII DATA TO ADDR IN X
04820      *CALLS: NOTHING
04830      *DESTROYS: A,X,CC
04840      *PURPOSE: THIS ROUTINE TRANSFER ASCII TEXT
04850      * FROM MEM FILE HDR TO STAT MSG BUFFER
04860      *
04870 093C BF 1C96 MSGPUT STS      STKSAV  SAVE CURRENT STACK PTR
04880 093F BE 095D      LDS      FRMLLOC  GET LOC OF DATA FOR TRANSFER

```

```

04880 0942 34          DES
04900 0943 32      MSGPU1 PUL A
04910 0944 31 04      CMP A  #4      END OF MSG?
04920 0945 27 05      BQ      MSGPU2  YES. EXIT
04930 0948 A7 00      STA A  0,X      NO. PUT CHAR IN MSG
04940 094A 08          INX
04950 094B 20 F6      BRA      MSGPU1
04960 094D 31      MSGPU2 INS      UPDATE STK WITH NEXT TRANSFER L
04970 094E BF 095D      STS      FRMLOC  STR BACK IN FRMLOC
04980 0951 BE 1C96      LDS      STKSAV  PICK UP ENTRY STACK PTR
04990 0954 39      RTS
05000          *
05010          * END MSGPUT
05020          *
05030          *****
05040          *
05050          *FUNC: HXASC
05060          *INPUTS: A,X (DATA TO BE CONV'D, ADDR TO STORE ASCII O
05070          *CALLS: HXASC ROUTINE IN EKG-EXEC
05080          *DSTROYS: A,X,CC
05090          *PURPOSE: THIS IS A RELOC PASS ROUTINE TO PICK UP
05100          * ADDR OF HXASC IN EKG-EXEC AND JUMP TO IT
05110          *
05120 0955 FF 1CA7  HXASC STX      HXBUF      SAVE X IN PARAMETER BUFFER
05130 0958 FE 1CA5          LDX      HXASLC  GET ADDR OF CONV ROUTINE IN E
05140 095B 6E 00          JMP      0,X      JUMP TO IT
05150          *
05160          * END OF HXASC
05170          *
05180          *****
05190          *
05200 095D 0002  FRMLOC RMB      2      DATA FROM LOC FOR MSGPUT
05210          *
05220 095F CE 096B  ERROR LDX      #ERRMSG
05230 0962 BD CA8F          JSR      OUTNCR
05240 0965 BD CA36          JSR      KEYED0
05250 0968 7E 1D00          JMP      START
05260          *
05270 096B 0707  ERRMSG FDB      $0707,$0D0A
05280 096F 4D          FCC      /MEMORY FILE COMPRESSION TYPE /
05290 098C 4E          FCC      /NOT RECOGNIZED. PRESS RETURN/
05300 09A8 04          FCB      4
05310          *
05320          *****
05330          * END PRSTAT
05340          *****
05350          *
05360          *****
05370          *
05380          * OUTPUT STRING TO LIST TO CONSOLE
05390          *
05400          *****
05410          *
05420 09A9 4E      NOCOMP FCC      /NO COMPRESSION PERFORMED/

```

PRSTAT-9

05430 09C1 04		FCB	4
05440 09C2 54	TOLA	FCC	/TOLAN-A/
05450 09C9 04		FCB	4
05460 09CA 54	TOLB	FCC	/TOLAN-B/
05470 09D1 04		FCB	4
05480 09D2 44	DOWR	FCC	/DOWR/
05490 09D7 04		FCB	4
05500 09D8 54	TURNPT	FCC	/TURNING POINT/
05510 09E5 04		FCB	4
05520 09E6 32	INTER	FCC	/2ND ORDER INTERPOLATOR/
05530 09FC 04		FCB	4
05540 09FD 43	ENTCLR	FCC	/CALCULATION NOT MADE/
05550 0A11 04		FCB	4
05560 0A12 1A07	HDRMSG	FDB	\$1A07,\$0CCD,\$0A0A,\$0ACA,\$0ACA,\$0ACA
05570 0A1E 3A		FCC	/:EKG SAMPLE COLLECTION STATISTICS/
05580 0A3F 20		FCC	/ : PAGE 1/
05590 0A48 0D0A		FDB	\$0D0A,\$0D0A
05600 0A4C 46		FCC	/FILENAME. /
05610 0A66 0008	NAME	RMB	8
05620 0A6E 0D0A		FDB	\$0D0A,\$0D0A
05630 0A72 53		FCC	/SUBJECT /
05640 0A8C 0035	SUBJ	RMB	53
05650 0AC1 0D0A		FDB	\$0D0A,\$0D0A
05660 0AC5 53		FCC	/SAMPLING RATE /
05670 0ADF 000A	RATE	RMB	10
05680 0AE9 0D0A		FDB	\$0D0A,\$0D0A
05690 0AED 44		FCC	/DATE OF COLLECTION. . . . /
05700 0B07 0014	DATE	RMB	20
05710 0B1B 0D0A		FDB	\$0D0A,\$0D0A
05720 0B1F 54		FCC	/TIME OF COLLECTION. . . . /
05730 0B39 0008	TIME	RMB	8
05740 0B41 0D0A		FDB	\$0D0A,\$0D0A
05750 0B45 43		FCC	/COMPRESSION USED. /
05760 0B5F 0019	TYPE	RMB	25
05770 0B78 0D0A		FDB	\$0D0A,\$0D0A
05780 0B7C 43		FCC	/CHANNEL X ENTROPY /
05790 0B96 0008	XENT	RMB	8
05800 0B9E 20		FCC	/ BITS /
05810 0BAA 0D0A		FDB	\$0D0A,\$0D0A
05820 0BAE 43		FCC	/CHANNEL Y ENTROPY /
05830 0BC8 0008	YENT	RMB	8
05840 0BD0 20		FCC	/ BITS /
05850 0BDC 0D0A		FDB	\$0D0A,\$0D0A
05860 0BE0 43		FCC	/CHANNEL Z ENTROPY /
05870 0BFA 0008	ZENT	RMB	8
05880 0C02 20		FCC	/ BITS /
05890 0C0E 0D0A		FDB	\$0D0A,\$0D0A
05900 0C12 54		FCC	/TOTAL SOURCE ENTROPY. . . . /
05910 0C2C 0008	TENT	RMB	8
05920 0C34 20		FCC	/ BITS /
05930 0C40 0D0A		FDB	\$0D0A,\$0D0A
05940 0C44 50		FCC	/PRESS RETURN FOR PAGE 2 OF STATISTIC
05950 0C6A 04		FCB	4
05960 0C6B 1A07	HDRMS2	FDB	\$1A07,\$0DCA

PRSTAT-9

05970	0C6F	45	FCC	/EKG SAMPLE COLLECTION STATISTICS/
05980	0C8F	20	FCC	/ : PAGE 2/
06000	0C9C	41	FCC	/RATIO POSSIBLE. /
06010	0CB2	0D0A	FDB	\$0D0A
06020	0CB4	52	FCC	/RATIO POSSIBLE. /
06030	0CCE	0008	MAXCPR RMB	8
06040	0CD6	20	FCC	/ : 1 /
06050	0CE2	0D0A	FDB	\$0D0A
06060	0CE4	43	FCC	/COMPRESSION RATIO/
06070	0CF5	0D0A	FDB	\$0D0A
06080	0CF7	41	FCC	/ACHIEVED. /
06090	0D11	0C08	CPRACH RMB	8
06100	0D19	20	FCC	/ : 1 /
06110	0D25	0D0A	FDB	\$0D0A
06120	0D27	41	FCC	/ACHIEVED COMPRESSION/
06130	0D3B	0D0A	FDB	\$0D0A
06140	0D3D	45	FCC	/EFFICIENCY. /
06150	0D57	0005	CPREFF RMB	5
06160	0D5C	20	FCC	/ % OF MAXIMUM /
06170	0D6B	0D0A	FDB	\$0D0A
06180	0D6D	43	FCC	/COMPRESSION TIME/
06190	0D7D	0D0A	FDB	\$0D0A
06200	0D7F	45	FCC	/EFFICIENCY OBTAINED . . . /
06210	0D99	0005	TIMEFF RMB	5
06220	0D9E	20	FCC	/ % SMP INTERVAL/
06230	0DDD	0D0A	FDB	\$0D0A, \$0D0A
06240	0DE1	43	FCC	/COLLECTION DURATION . . . /
06250	0DCB	0005	COLDUR RMB	5
06260	0DD0	20	FCC	/ SECONDS /
06270	0DDF	0D0A	FDB	\$0D0A, \$0D0A
06280	0DE3	43	FCC	/CHANNEL X MAXIMUM /
06290	0DED	0008	AMAXX RMB	8
06300	0E05	20	FCC	/ VOLTS /
06310	0E11	0D0A	FDB	\$0D0A, \$0D0A
06320	0E15	43	FCC	/CHANNEL X MINIMUM /
06330	0E2F	0008	AMINX RMB	8
06340	0E37	20	FCC	/ VOLTS /
06350	0E43	0D0A	FDB	\$0D0A, \$0D0A
06360	0E47	43	FCC	/CHANNEL Y MAXIMUM /
06370	0E61	0008	AMAXY RMB	8
06380	0E69	20	FCC	/ VOLTS /
06390	0E75	0D0A	FDB	\$0D0A, \$0D0A
06400	0E79	43	FCC	/CHANNEL Y MINIMUM /
06410	0E93	0008	AMINY RMB	8
06420	0E9B	20	FCC	/ VOLTS /
06430	0EA7	0D0A	FDB	\$0D0A, \$0D0A
06440	0E7B	43	FCC	/CHANNEL Z MAXIMUM /
06450	0FC5	0008	AMAXZ RMB	8
06460	0FCD	20	FCC	/ VOLTS /
06470	0FD9	0D0A	FDB	\$0D0A, \$0D0A
06480	0FDD	43	FCC	/CHANNEL Z MINIMUM /
06490	0FE7	0008	AMINZ RMB	8
06500	0FEF	20	FCC	/ VOLTS /

06510	0F0B	CD0A	FDB	\$OD0A,\$OD0A
06520	0F0F	43	FCC	/COMMENTS. /
06530	0F11	CD0A	FDB	\$OD0A,\$OD0A
06550	0F62	50	FCC	/PRESS RETURN FOR PAGE 3 OF STATISTIC
06560	0F88	04	FCB	4
06570	0F89	1A07	COSTAT FDB	\$1A07,\$OCCD,\$OAOA,\$OAOA,\$OAOA,\$OAOA
06580	0F95	45	FCC	/ENG SAMPLE COLLECTION STATISTICS: PA
06590	0FED	OD0A	FDB	\$OD0A CR/LF
06600	0FDF	20	FCC	/ NUMBER OF SAMPLES TAKEN /
06610	0FD8	28	FCC	/ (SAMPLE) = /
06620	0FE2	0004	NUMSMP RFB	4
06630	0FE6	20	FCC	/ (HEX)/
06640	0FEC	OD0A	FDB	\$OD0A CP/LF
06650	0FEE	20	FCC	/ MAXIMUM LOOP COUNT PER/
06660	1005	20	FCC	/ INTERRUPT (LPCAL) = /
06670	1019	0002	NUMIPS RFB	2
06680	101B	20	FCC	/ (HEX)/
06690	1021	OD0A	FDB	\$OD0A CR/LF
06700	1023	20	FCC	/ TOTAL WAITING LOOP COUNTS DURING/
06710	1044	20	FCC	/ COLLECTION (LOOPCT) = /
06720	105A	0008	TOTLUP RFB	8
06730	1062	20	FCC	/ (HEX)/
06740	1068	OD0A	FDB	\$OD0A
06750	106A	54	FCC	/TIME EFFICIENCY = (1-(LOOPCT (SAMPLE
06760	1097	2A	FCC	/*100/
06770	109B	OD0A	FDB	\$OD0A,\$OD0A
06780	109F	43	FCC	/CHANNEL MAXIMUMS AND MINIMUMS/
06790	10BC	OD0A	FDB	\$OD0A
06800	10FE	20	FCC	/ XMAX= /
06810	10C5	0002	XMAXTM RFB	2
06820	10C7	20	FCC	/ (HEX) AT SAMPLE NUMBER /
06830	10DF	0004	XMAXLO RFB	4
06840	10E3	20	FCC	/ (HEX)/
06850	10E9	OD0A	FDB	\$OD0A
06860	10EB	20	FCC	/ XMIN= /
06870	10F2	0002	XMINTM RFB	2
06880	10F4	20	FCC	/ (HEX) AT SAMPLE NUMBER /
06890	110C	0004	XMINLO RFB	4
06900	1110	20	FCC	/ (HEX)/
06910	1116	OD0A	FDB	\$OD0A
06920	1118	20	FCC	/ YMAX= /
06930	111F	0002	YMAXTM RFB	2
06940	1121	20	FCC	/ (HEX) AT SAMPLE NUMBER /
06950	1139	0004	YMAXLO RFB	4
06960	113D	20	FCC	/ (HEX)/
06970	1143	OD0A	FDB	\$OD0A
06980	1145	20	FCC	/ YMIN= /
06990	114C	0002	YMINTM RFB	2
07000	114F	20	FCC	/ (HEX) AT SAMPLE NUMBER /
07010	1166	0004	YMINLO RFB	4
07020	116A	20	FCC	/ (HEX)/
07030	1170	OD0A	FDB	\$OD0A
07040	1172	20	FCC	/ ZMAX= /

PRSTAT-9

07050 1179 0002	ZMAXNM RMB	2	
07060 117B 20	FCC	/ (HEX) AT SAMPLE NUMBER /	
07070 117C 20	ZMINNM RMB	1	
07080 1197 20	FCC	/ (HEX)/	
07090 119D 0D0A	FDB	\$0D0A	
07100 119F 20	FCC	/ ZMIN= /	
07110 11A6 0002	ZMINNM RMB	2	
07120 11A8 20	FCC	/ (HEX) AT SAMPLE NUMBER /	
07130 11C0 0004	ZMINLO RMB	4	
07140 11C4 20	FCC	/ (HEX)/	
07150 11CA 0D0A	FDB	\$0D0A,\$0D0A	
07160 11CE 43	FCC	/COMPRESSION STATISTICS:/	
07170 11E5 0D0A	FDB	\$0D0A	
07180 11E7 20	FCC	/ NUMBER OF MEMORY BITS /	
07190 11FE 41	FCC	/AVAILABLE = /	
07200 120A 0006	AVABIT RMB	6	
07210 1210 20	FCC	/ (HEX)/	
07220 1216 0D0A	FDB	\$0D0A	
07230 1218 20	FCC	/ NUMBER OF BITS AVAILABLE/	
07240 1231 20	FCC	/ TO VAR LEN CODER = /	
07250 1245 0008	TCPBIT RMB	8	
07260 124D 20	FCC	/ (HEX)/	
07270 1253 0D0A	FDB	\$0D0A	
07280 1255 20	FCC	/ TOTAL NUMBER OF DATA BITS/	
07290 126F 20	FCC	/ STORED = /	
07300 1279 0008	TOTBIT RMB	8	
07310 1281 20	FCC	/ (HEX)/	
07320 1287 0D0A	FDB	\$0D0A	
07330 1289 20	FCC	/ NUMBER OF BITS USED TO/	
07340 12A0 20	FCC	/ STORE CHANNEL X = /	
07350 12E3 0006	XBITNM RMB	6	
07360 12B9 20	FCC	/ (HEX)/	
07370 12BF 0D0A	FDB	\$0D0A	
07380 12C1 20	FCC	/ NUMBER OF BITS USED TO/	
07390 12D8 20	FCC	/ STORE CHANNEL Y = /	
07400 12EB 0006	YBITNM RMB	6	
07410 12F1 20	FCC	/ (HEX)/	
07420 12F7 0D0A	FDB	\$0D0A	
07430 12F9 20	FCC	/ NUMBER OF BITS USED TO/	
07440 1310 20	FCC	/ STORE CHANNEL Z = /	
07450 1323 0006	ZBITNM RMB	6	
07460 1329 20	FCC	/ (HEX)/	
07470 132F 0D0A	FDB	\$0D0A	
07480 1331 20	FCC	/ NUMBER OF BITS USED TO/	
07490 1348 20	FCC	/ STORE TIME OR OTHER LAGARY	
07500 1362 45	FCC	/ETER = /	
07510 1369 0006	TBITNM RMB	6	
07520 136F 20	FCC	/ (HEX)/	
07530 1375 0D0A	FDB	\$0D0A	
07540 1377 43	FCC	/COMPRESSION RATIO = TOTAL /	
07550 1391 44	FCC	/DATA BITS STORED PER ITEM /	
07560 13AA 42	FCC	/BITS AVAILABLE/	
07570 13B8 0D0A	FDB	\$0D0A,\$0D23	
07580 13PC 50	FCC	/PRESS RETURN =/	

PRSTAT-9

07590 13CA 04

FCB 4

07600

*

07610

*

END OF PRSTAT

07620

*

07630

END

```

00030 *****
00040 *
00050 * OVERLAY UNIT : NOCPRS
00060 * AUTHOR : CAPT. MEL TOWNSEND
00070 * VERSION : 1.7
00080 * VERSION DATE : 3 OCT 80
00090 *
00100 *****
00110 *
00120 * OVERLAY DESCRIPTION
00130 *
00140 * THIS OVERLAY SAMPLES THE EKG DATA VIA THE
00150 * A/D CONVERTERS AND STORES THE 8 BIT ROUNDED
00160 * VALUES INTO MEMORY LOCATIONS 3C00-7FFF. THE
00170 * THE X,Y,Z CHANNELS ARE SAMPLED AND STORED
00180 * SEQUENTIALLY STARTING AT 3C00 AND WORKING UP.
00190 *
00200 *****
00210 *****
00220 *
00230 * START OF NOCPRS
00240 *
00250 *****
00260 *
00270 0100          ORG      $0100      OVERLAY START ADDRESS
00280 *
00290          OPT      0          ASSE OPT-COM CRU FILE
00300          OPT      NOG       ASSE OPT-SUPPRESS FULL FCC LI
00310 *
00320 *****
00330 *
00340 * LABEL DECLARATIONS
00350 *
00360 * SUPPORT SUBROUTINE ADDRESSES
00370 *
00380 CA87  OUTPUT EQU      SCAB7      EPGQDOS. ALPH STRING TO CONS
00390 CA8F  OUTCR EQU      SCABF      EPGQDOS. ALPH STRING, NO CR/LF
00400 CA2C  INTD EQU      SCAC2C     EPGQDOS. CONSOLE INPUT ROUTE
00410 CA36  KEND0 EQU      SCAC36     EPGQDOS. CONSOLE INPUT, NO ?
00420 1D00  START EQU      1D00      EPG-EXE C. START EPG-EXEC
00430 *
00440 * DATA BUFFERS
00450 *
00460 3400  HDRSTR EQU      $3400      DATA MEMORY HEADER BUFFER
00470 FFF8  IRQVIC EQU      $FFF8      INTERRUPT VECTOR BUFFER
00480 1C96  STKSAV EQU      $1C96      STACK SAVE BUFFER
00490 1C98  CRTTYP EQU      $1C98      CONNECTION ALPH CTR FILE
00500 3002  ENDBUF EQU      $3002      ADDR OF LAST CRP BUFFER
00510 1C9D  VEC SAV EQU      $1C9D      IRQ VECTOR SAVE BUFFER
00520 1CA1  FILHLC EQU      $1CA1      FILTER ALPH PAGE FILTER
00530 1CA3  SAVHLC EQU      $1CA3      SAVHLC ALPH PAGE FILTER
00540 3490  LOOPCT EQU      $3490      TOTAL NTH TIME TO LOOP
00550 3494  SAVHNO EQU      $3494      NUMBER OF TIMES TO LOOP
00560 3496  LOCAL EQU      $3496      AVG'D MAX LOOP COUNT

```

```

00570 3497 MAXZ EQU S3497 MAX VLU IN CH Z
00580 3498 MAXYLO EQU S3498 MAX VLU LOC IN CH Y
00600 3499 MINZLO EQU S3499 MIN VLU LOC IN CH Z
00610 349D MAXY EQU S349D MAX VLU IN CH Y
00620 349E MAXYLO EQU S349E MAX VLU LOC IN CH Y
00630 34A0 MINY EQU S34A0 MIN VLU IN CH Y
00640 34A1 MINYLO EQU S34A1 MIN VLU LOC IN CH Y
00650 34A3 MAXX EQU S34A3 MAX VLU IN CH X
00660 34A4 MAXXLO EQU S34A4 MAX VLU LOC IN CH X
00670 34A6 MINX EQU S34A6 MIN VLU IN CH X
00680 34A7 MINXLO EQU S34A7 MIN VLU LOC IN CH X
00690 34AC DIBIT EQU S34AC NUM OF BITS USED TO STR DIB
00700 34B0 XBITS EQU S34B0 NUM OF BITS USED TO STR X
00710 34B3 YBITS EQU S34B3 NUM OF BITS USED TO STR Y
00720 34B6 ZBITS EQU S34B6 NUM OF BITS USED TO STR Z
00730 34B9 TBITS EQU S34B9 NUM OF BITS USED TO STR T
00740 34BC ACHLCT EQU S34BC # BITS SET TO VLG IN CHR HCCC
00750 3600 XPDFM1 EQU S3600 0 VLU LOC OF X PDF
00760 3800 YPDFM1 EQU S3800 0 VLU LOC OF Y PDF
00770 3A00 ZPDFM EQU S3A00 0 VLU LOC OF Z PDF
00780 3E00 TPDF EQU S3E00 0 VLU LOC TO TIME PDF
00790 3C00 SECZRO EQU S3C00 SWEEP OF DATA NEW FILE
00800 *
00810 * HARDWARE ADDRESSES
00820 *
00830 E400 ADCZRO EQU SE400 ADC CHANNEL ZERO
00840 E404 ADCYLO EQU SE404 ADC CHANNEL TWO
00850 E500 DACZRO EQU SE500 DAC CHANNEL ZERO
00860 *
00870 *
00880 *FUNCTION :NOCPRS
00890 *INPUTS :STATUS PORTS FROM EKG-EXEC
00900 *OUTPUTS :DATA TO DIB
00910 *CALLS :OUTPUT,PRIME,TIME,QUENCH
00920 * KEYBD0,SAVITL,START
00930 *DESTROYS :ALL REGISTERS
00940 *PURPOSE :TO COLLECT 3 CHANNELS OF EKG
00950 * DATA AND STORE INTO MEMORY.
00960 *
00970 *
00980 0100 0F NOCPRS SET
00990 0101 CE 0233 LDX #STARTG
01000 0104 ED C48F JSR OUTCON "THIS MODULE SAMPLES THE..."
01010 0107 ED C42C JSR KEYBD GET FLAGGED FROM CONSOLE
01020 010A FE 3002 LDX ENDBUT
01030 010D 09 DEX
01040 010F 16 00 LDA B 0,X
01050 0110 C1 59 CMP B #Y IS INPUT YES?
01060 0112 27 03 BLO NOCPRI YFS. KEEP EXECUTING THIS ROUT
01070 0114 7E 1D00 JMP START NO. PIN TO EKG-EXEC
01080 0117 CE 4E43 NOCPRI LDX "S4E43 FLAG COMMISSION TIME (NO)
01090 011A FE 1C98 STX CTIMEP
01100 011D ED 051A JSR FILDER SET UP DATA FILE HEADER

```

```

01110 0120 CE 02B5      LDX      #IROMSG
01120 0121 CE 02B6      STX      #IROMSG      "INITIALIZE SURVIVAL AND REG..."
01130 0122 CE 02B7      STX      #IROMSG
01140 0129 CE 3C00 NOCPRGO LDX      #SECZRO INITIALIZE REG FOR START LOC
01150 012C FF 0230      STX      BUFPTR  STORE IN BUFFER POINTER
01160 012F 86 03      LDA A    #3      PICK UP COUNTER FOR TIME CAL
01170 0131 B7 0227      STA A    CALCNT  STORE IN BUFFER
01180 0134 FE FFF8      LDX      IRQVEC  PICK UP CURRNT IRQ VECTOR
01190 0137 FF 1C9D      STX      VEC SAV  SAVE IN BUFFER
01200 013A CE 051F THICAL LDX      #CALINT GET INTR VECTOR ADDR FOR CAL
01210 013D FF FFF8      STX      IRQVEC  PUT IN VECTOR ADDRESS
01220 0140 85 00      LDA A    #0      INIT COUNTER FOR 256 TEST LOC
01230 0142 B7 0232      STA A    DONTEST STORE IN DONTEST BUFFER
01240 0145 CE 0000      LDX      #0      CLEAR LOOP COUNTER
01250 0148 FF 022C      STX      IRQCNT  IRQCNT
01260 014F FF 022E      STX      IRQCNT+2
01270 014E FF E500      STX      DACZRO  PULSE INT ENABLE CIRCUIT
01280 0151 20 09      BNA      SPOO1  PPA TO COUNTING LOOP
01290
01300      *
01310      *****
01320      *
01330      * BASIC TIMING LOOP FOR EFFICIENCY TEST
01340      *
01350      *****
01360 0153 01      SPOO2  NOP      DELAY TO MATCH TIME
01370 0154 01      NOP      IT TAKES TO EXECUTE ONE
01380 0155 01      NOP      INCREMENT OF DATA 3&4
01390 0156 01      NOP      OF THE LOOP COUNTER
01400 0157 01      NOP      WHEN COUNT CARRIES TO
01410 0158 01      NOP      HIGH 2 BYTES OF 4 BYTE
01420 0159 7E 015F      JMP      SPOO3  JUMP TO CONTINUE LOOP
01430 015C 0E      SPOO1  CLI      PREPARE FOR TIMER INT
01440 015D 3E      WAI      STOP PROCESSOR & WAIT
01450 015E 0E      SPOO4  CLI
01460 015F FE 022E      LDX      IRQCNT+2 RET'ED FROM INT, INC COUNT
01470 0162 08      INX
01480 0163 FF 022E      STX      IRQCNT+2 SAVE IT
01490 0166 26 FB      BNE      SPOO2  COUNT GONE FFFF TO 0000?
01500 0168 FE 022C      LDX      IRQCNT  YES. INC BYTES 3&4
01510 016B 08      INX
01520 016C FF 022C      STX      IRQCNT  SAVE COUNT
01530 016F B6 0232 SPOO3  LDA A    DONTEST IS DONE TEST SATISFIED?
01540 0172 81 00      CMP A    #0
01550 0174 26 E8      BNE      SPOO4  NO. KEEP LOOPING
01560
01570      *
01580      *****
01590      *
01600      * END BASIC TIMING LOOP FOR EFF TEST
01610      *
01620      *****
01630 0176 0F      SEI      PREVENT SER TO INTR LATER
01640 0177 CE 4000      LDX      #S4000  PICK UP INTR COUNTER

```

01650 017A FF B500	STX	DACZRO	DISABLE INTER FLIGH TIMER
01660 017D B7 B400	STA A	ADCZRO	CLR 176800 INT FLIP FLAG
01670 017E 00 0000	LDX	#0	
01680 0181 7A 0000	LDX	CALCNT	IS COUNT CALCULATED
01690 0184 27 38	BEQ	SPXONE	IS COUNT 3? YES, GO SAVE DATA
01700 0185 F6 0227	LDA A	CALCNT	NO. CNT COUNT
01710 0189 81 02	CMF A	#2	IS COUNT 2?
01720 018B 26 22	RNE	PARSAV	NO. SAVE PARAMETER FROM INPUT
01730 018D FE 022E	LDX	IRQCNT+2	YES. SAVE TIME LOOP COUNTER
01740 0190 FF 022A	STX	CALONE	
01750 0193 CE 0316	LDX	#SAMPLE	NOW PUT SAMPLE ADDR IN IRQVEC
01760 0196 FF FFF8	STX	IRQVEC	
01770 0199 86 AA	LDA A	#SAA	SET UP DONTST
01780 019E B7 0232	STA A	DONTST	
01790 019E CE 0000	LDX	#0	CLR LOOP COUNT AND SAMPLE COU
01800 01A1 FF 022C	STX	IRQCNT	
01810 01A4 FF 022E	STX	IRQCNT+2	
01820 01A7 FF 3494	STX	SAMENO	
01830 01AA FF B500	STX	DACZRO	ENABLE TIMER INTERRUPTS
01840 01AD 20 AD	RPA	SPOOL	GO WAIT FOR INTERRUPT
01850 01AF FE 022C	PARSAV LDX	IRQCNT	SAVE LOOP COUNT
01860 01B2 FF 3490	STX	LOOPCT	
01870 01B5 FE 022E	LDX	IRQCNT+2	
01880 01B8 FF 3492	STX	LOOPCT+2	
01890 01BB 7F 013A	JMP	TIFCAL	GO EXECUTE ANOTHER TIME CAL
01900 01BE B6 022E	SPXONE LEA A	CALZRO+1	AVERAGE TWO TIME CAL REUS
01910 01C1 FB 022F	ADD A	IRQCNT+3	ADD 158 BYTES
01920 01C4 B7 0229	STA A	CALZRO+1	STORE IN BUFFER
01930 01C7 F6 022A	LDA A	CALONE	PICK UP 158 OF FIRST
01940 01CA B9 022E	ADC A	IRQCNT+2	ADD WITH CARRY
01950 01CD B7 0228	STA A	CALZRO	STORE IN BUFFER
01960 01D0 77 0228	ASR	CALZRO	DIVIDE BY 2 TO AVERAGE
01970 01D3 76 0229	ROR	CALZRO+1	
01980 01D6 85 08	LEA A	#8	SET UP COUNTER FOR /256
01990 01D8 77 0228	CALSHF ASR	CALZRO	NOW DIVIDE BY 256 FOR
02000 01D3 76 0229	ROR	CALZRO+1	LOOP/INTERSECT
02010 01DE 4A	DEC A		
02020 01DF 26 F7	RNE	CALSHF	8 SHIFTS YET?
02030 01E1 F6 0229	LDA A	CALZRO+1	PICK UP SHIFTED RESULT
02040 01F4 B7 3455	STA A	LOCAL	STORE IN FILE FOR BUFFER
02050 01F7 F6 34AF	LDA A	DATAIT+3	PUT DATA BIT COUNT INTO ACCTC
02060 01FA B7 34FF	STA A	ACCHCT+3	SINCE NO COMPRESSION PERFECT
02070 01FD B6 34AE	LDA A	DATAIT+2	
02080 01F0 B7 34DE	STA A	ACCHCT+2	
02090 01F3 B6 34AD	LDA A	DATAIT+1	
02100 01F6 B7 34DD	STA A	ACCHCT+1	
02110 01F9 B6 34FC	LDA A	DATAIT	
02120 01FC B7 34FC	STA A	ACCHCT	
02130 01FF FE 1C9D	LDX	VECSAV	YES. RESTORE IRQ VECTOR
02140 0202 FF FFF8	STX	IRQVEC	
02150 0205 CE 0312	LDX	#SAVMSG	
02160 0208 FD CA8F	JSR	OUTNCR	"SAMPLING COMPLETE..."
02170 020B FD C4CC	JSR	KEYPD	GET SAVE DECISION
02180 020E FE 3002	LDX	ENDRUF	

NOCTIS-7

```

02190 0211 09          DEX
02200 0212 E6 00      LDA B  0,X
02210 0214 C1 F9      CUP B  #Y      IS DECISION YET?
02220 0216 28 00      FCB      EXJMP  NO. TEST LOG FULL?
02230 0218 ED 0515     JSR      SAVEPL  YES. SAVE THE FILE ON DISK
02240 021F CE 0381     LDY      #SDONE
02250 021E BD CA8F     JSR      OUTCHR
02260 0221 ED CA2C     JSR      KEYDD
02270 0224 7E 1D00     EXJMP  JMP      START  JUMP BACK TO PRG-EXEC
02280                  *
02290 0227 0001      CALCNT RMB  1      TEST CAL LOOP COUNTER
02300 0228 0002      CALZRO RMB  2      AVG'D CAL LOOPS DURING DAT CO
02310 022A 0002      CALONE RMB  2      TEMP RUP FOR PRE OMNISC CAL
02320 022C 0004      IRQCNT RMB  4      TEMP INTR LOOP COUNTER
02330 0230 0002      BUFFPTR RMB  2      BUFFER POINTER (OT NEXT AVAIL
02340 0232 0001      DONTST RMB  1      DONE TEST FLG FOR MEM FULL
02350                  *
02360 0233 1A07      STRMSG FDB  $1A07
02370 0235 54          FCC      /THIS MODULE SAMPLES THE EKG /
02380 0251 49          FCC      /INPUT AND STORES THE DATA WITH/
02390 026F 0D0A      FDB      $0D0A
02400 0271 4E          FCC      /NO DATA COMPRESSION./
02410 0285 0D0A      FDB      $0D0A,$0D0A
02420 0289 44          FCC      /DO YOU WISH TO EXECUTE THIS /
02430 02A5 4D          FCC      /MODULE (Y OR N)/
02440 02B4 04          FCB      4
02450 02B5 1A      IRQMSG FCB  $1A
02460 02B6 49          FCC      /INSURE SUBJECT AND EKG DEVICE /
02470 02D4 52          FCC      /READY!/
02480 02DA 0D0A      FDB      $0D0A,$0D0A,$0D0A
02490 02E0 50          FCC      /PRESS RETURN, THEN CLOSE INTERRUPT/
02500 0302 20          FCC      / ENABLE SWITCH./
02510 0311 04          FCB      4
02520 0312 1A07      SAVMSG FDB  $1A07
02530 0314 53          FCC      /SAMPLING COMPLETE. PLEASE OPEN /
02540 0334 49          FCC      /INTERRUPT ENABLE SWITCH./
02550 034C 0D0A      FDB      $0D0A,$0D0A,$0D0A
02560 0352 44          FCC      /DO YOU WISH TO SAVE THIS DATA ON /
02570 0373 44          FCC      /DISK (Y OR N)/
02580 0380 04          FCB      4
02590 0381 20      SDONE  FCC      / PRESS RETURN/
02600 038E 04          FCB      4
02610                  *
02620                  *****
02630                  END NOCTIS
02640                  *****
02650                  *
02660                  *
02670                  *FUNCTION :SAMPLE
02680                  *INPUT   :STATUS BUFFERS
02690                  *OUTPUTS  :COMPRESSED,ROUNDED DATA IN INT' BUFF
02700                  *CALLS   :NOTHING
02710                  *PURPOSE  :THIS ROUTINE SAMPLES THE EKG LEADS,
02720                  * ROUNDS THE VALUES TO 3 BITS (FIGURE 12)

```



```

02730      * CALCULATES MAX, MIN, # OF BITS, SAMPLES
02740      * ETC., AND SAVES THESE PARAM AND DATA
02750      * TO DATA MEMORY.
02760      *
02770      *
02780 038F 0001 SHFBUF RMB 1      TEMP SHIFT BUFFER
02790 0390 0001 TEMPAL RMB 1      TEMP REG, MSB CH 0
02800 0391 0001 TEMPB1 RMB 1      TEMP REG, LSB CH 0
02810 0392 0001 TEMP2 RMB 1      TEMP REG, MSB CH 1
02820 0393 0001 TEMP3 RMB 1      TEMP REG, LSB CH 1
02830 0394 0001 TEMP3A RMB 1     TEMP REG, MSB CH 2
02840 0395 0001 TEMP3B RMB 1     TEMP REG, LSB CH 2
02850      *
02860 0396 FE 3494 SAMPLE LDX  SAMENO  GET CUR SAMPL COUNT
02870 0399 08      INX
02880 039A FF 3494      STX  SAMENO
02890 039D BF 1C96      STS  STKSAV  SAVE STACK PRT
02900 03A0 8E 3280      LDS  #S3280  INITIALIZE STACK IN UNUSED MEM
02910 03A3 CE E400      LDX  #ADC2R0 NOW PULSE A/D TO START CONV
02920 03A6 A7 00      STA A  0,X      FOR CHANNEL 0
02930 03A8 01      NOP
02940 03A9 A6 00      LDA A  0,X
02950 03AB E6 01      LDA B  1,X
02960 03AD B7 0390      STA A  TEMPAL
02970 03E0 F7 0391      STA B  TEMPB1
02980 03B3 08      INX
02990 03B4 08      INX
03000 03B5 A7 00      STA A  0,X      NOW PULSE A/D FOR CONV
03010 03B7 01      NOP      ON CHANNEL 1
03020 03B8 A6 00      LDA A  0,X
03030 03BA E6 01      LDA B  1,X
03040 03BC B7 0392      STA A  TEMP2
03050 03BE F7 0393      STA B  TEMPB2
03060 03C2 08      INX
03070 03C3 08      INX
03080 03C4 A7 00      STA A  0,X      NOW PULSE A/D FOR CONV
03090 03C6 01      NOP      ON CHANNEL 2
03100 03C7 A6 00      LDA A  0,X
03110 03C9 F6 01      LDA B  1,X
03120 03CB B7 0394      STA A  TEMP3A
03130 03CE F7 0395      STA B  TEMP3B
03140 03D1 CE 0390      LDY  #TEMPAL
03150 03D4 A6 00  SAMPLE1 LDA A  0,X
03160 03D6 E6 01      LDA B  1,X
03170 03D8 47      ASR A      NOW SHIFT 2 BYTE VLD TO ACCUM
03180 03D9 56      ROR B      4 POSITIONS TO RIGHT TO INFO
03190 03DA 47      ASR A      12 BIT TO 8 BIT RIGHT CONV
03200 03DB 56      ROR B
03210 03DC 47      ASR A
03220 03DD 56      ROR B
03230 03DE 47      ASR A
03240 03DF 56      ROR B
03250 03E0 F7 039F      STA B  SHFBUF  SAVE 8 BIT RESULT OF SHIFT
03260 03E3 86 00      LDA A  #0      AND ADD CARRY OUT OF LAST

```

03270 03F5 B9 039F	ADC A	SHRBUF	SHIFT ROUNDING UP OR DOWN
03280 03F8 36	PSH A		SAVE VLU IN STACK TEMPORARILY
03290 03F9 00	LDX		
03300 03FA 00	LDX		
03310 03EB 8C 0396	CPX	#TEMPA3+2	CHANNEL Z ROUNDED TO 8 BITS
03320 03EE 26 E4	BNE	SAMPL1	NO. GO SAMPLE NEXT CHANNEL
03330 03F0 32	PUL A		GET Z DATA
03340 03F1 FE 0230	LDX	BUFPTR	PICK UP CUR MEM FILE PTR
03350 03F4 A7 02	STA A	2,X	SAVE DATA TO MEM FILE
03360 03F6 FE 3494	LDX	SAMPNO	GET CUR SAMPLE COUNT
03370 03F9 B1 3497	CMP A	MAXZ	IS CUR Z MAX OVE SAMPLE SET?
03380 03FC 2F 06	BLE	SPZMIN	NO. GO CHECK FOR MIN
03390 03FE B7 3497	STA A	MAXZ	YES. KEEP CUR VLU
03400 0401 FF 3493	STX	MAXZLO	KEEP CUR SAMPLE NUM
03410 0404 B1 349A	CMP A	MINZ	IS CUR Z MIN OVER SAMPLE SET?
03420 0407 2C 06	BGE	SPYMAX	NO. GO CHECK FOR Y MAX
03430 0409 B7 349A	STA A	MINZ	YES. KEEP CUR VLU
03440 040C FF 349B	STX	MINZLO	KEEP CUR SAMPLE NUM
03450 040F 32	SPYMAX	PUL A	GET CUR Y DATA
03460 0410 FE 0230	LDX	BUFPTR	PICK UP CURRENT MEM FILE PTR
03470 0413 A7 01	STA A	1,X	SAVE DATA TO MEM FILE
03480 0415 FE 3494	LDX	SAMPNO	GET CUR SAMPLE COUNT
03490 0418 B1 349D	CMP A	MAXY	IS CUR Y MAX OVER SAMPLE SET
03500 041B 2F 06	BLE	SPYMIN	NO. GO CHECK FOR MIN
03510 041D B7 349D	STA A	MAXY	YES. KEEP CUR Y VLU
03520 0420 FF 349E	STX	MAXYLO	KEEP CUR SAMPLE NUM
03530 0423 B1 34A0	CMP A	MINY	IS CUR Y MIN OVER SAMPLE SET?
03540 0426 2C 06	BGE	SPXMAX	NO. GO CHECK FOR X MAX
03550 0428 B7 34A0	STA A	MINY	YES. KEEP CUR Y VLU
03560 042B FF 34A1	STX	MINYLO	KEEP CUR SAMPLE NUM
03570 042E 32	SPXMAX	PUL A	GET CUR X DATA
03580 042F FE 0230	LDX	BUFPTR	PICK UP CURRENT MEM FILE PTR
03590 0432 A7 00	STA A	0,X	SAVE DATA TO MEM FILE
03600 0434 FE 3494	LDX	SAMPNO	GET CUR SAMPLE COUNT
03610 0437 B1 34A3	CMP A	MAXX	IS CUR X MAX OVER SAMPLE SET?
03620 043A 2F 06	BLE	SPXMIN	NO. GO CHECK X MIN
03630 043C B7 34A3	STA A	MAXX	YES. KEEP CUR X VLU
03640 043F FF 34A4	STX	MAXXLO	KEEP CUR SAMPLE VLU
03650 0442 B1 34A6	CMP A	MINX	IS CUR X MIN OVER SAMPLE SET?
03660 0445 2C 06	BGE	MNDONE	NO. EXIT MAX,MIN UPDATE
03670 0447 B7 34A6	STA A	MINX	YES. KEEP CUR X VLU
03680 044A FF 34A7	STX	MINXLO	KEEP CUR SAMPLE COUNT
03690 044D 0C	MNDONE	CLC	DONE WITH MAX & MIN'S. NOW
03700 044E 86 18	LDA A	#24	UPDATE ACTUAL DATA BIT STORED
03710 0450 C6 00	LDA B	#0	
03720 0452 BB 34AF	ADD A	DTABIT+3	
03730 0455 B7 34AF	STA A	DTABIT+3	
03740 0458 86 00	LDA A	#0	
03750 045A F9 34AE	ADC B	DTABIT+2	
03760 045D F7 34AE	STA B	DTABIT+2	
03770 0460 C6 00	LDA B	#0	
03780 0462 B9 34AD	ADC A	DTABIT+1	
03790 0465 B7 34AD	STA A	DTABIT+1	
03800 0468 F9 34AC	ADC B	DTABIT	

03810 046B F7 34AC	STA B	DTABIT	
03820 046E 86 08	LDA A	#8	NOW UPDATE XBIT,YBIT,ZBIT COU
03830 0470 C6 0C	LDX B		UPDATE COUPTER TO 0
03840 0472 BB 34B2	ADD A	XBITS+1	X,Y,Z COUNTERS INCR.
03850 0475 B7 34B2	STA A	XBITS+2	
03860 0478 B7 34B5	STA A	YBITS+2	
03870 047B B7 34B8	STA A	ZBITS+2	
03880 047E 86 00	LDA A	#0	
03890 0480 F9 34B1	ADC B	XBITS+1	
03900 0483 F7 34B1	STA B	XBITS+1	
03910 0486 F7 34B4	STA B	YBITS+1	
03920 0489 F7 34B7	STA B	ZBITS+1	
03930 048C B9 34B0	ADC A	XBITS	
03940 048F B7 34B0	STA A	XBITS	
03950 0492 B7 34B3	STA A	YBITS	
03960 0495 B7 34B6	STA A	ZBITS	
03970 0498 7F 34BB	CLR	TBITS+2	
03980 049B 7F 34BA	CLR	TBITS+1	
03990 049E 7F 34B9	CLR	TBITS	
04000 04A1 34	DES		NOW POINT STACK BACK TO DATA
04010 04A2 34	DES		
04020 04A3 34	DES		
04030 04A4 32	PUL A		GET Z DATA
04040 04A5 CE 3A00	LDX	#ZPDFM	
04050 04A8 BD 04D2	JSR	PDFSTR	UPDATE Z PDF BIN COUNTER
04060 04AB 32	PUL A		GET Y DATA
04070 04AC CE 3800	LDX	#YPDFM	
04080 04AF BD 04D2	JSR	PDFSTR	UPDATE Y PDF BIN COUNTER
04090 04B2 32	PUL A		
04100 04B3 CE 3600	LDX	#XPDFM	
04110 04B6 BD 04D2	JSR	PDFSTR	
04120 04B9 FE 0230	LDX	BUFPTR	NOW GET & UPDATE MEM BUFF PTR
04130 04EC 08	INX		
04140 04ED 08	INX		
04150 04EE 08	INX		
04160 04BF FF 0230	STX	BUFPTR	
04170 04C2 BE 1C96	LDS	STKSAV	RETRIEV ENTRY STACK POINTER
04180 04C5 8C 7FFE	CPX	#S7FFE	IS MEM BUF FULL?
04190 04C8 26 03	BNE	SAMRTI	NO. RET FROM INTR & SAMPLE AC
04200 04CA 7F 0232	CLR	DONTST	YES. RESET MEM FULL FLAG
04210 04CD 3B	SAMRTI	RTI	
04220	*		
04230 04CE 0002	TEMPDF	RMB 2	TEMP WORKING BUFFER
04240 04D0 0002	TEMPST	RMB 2	TEMP STACK SAVE BUFFER
04250	*		
04260 04D2 BF 04D0	PDFSTR	STX	TEMPST
04270 04D5 FF 04CE	STX	TEMPDF	SAVE PDF PTR
04280 04D8 16	TAB		SAVE INPUT DATA VLI
04290 04D9 BB 04CF	ADD A	TEMPDF+1	NOW CALCULATE ADDRESS FOR INC
04300 04DC B7 04CF	STA A	TEMPDF+1	IN THE PDF MEM BUFFER
04310 04DF 86 00	LDA A	#0	
04320 04E1 C5 80	BIT B	#S80	IS THE DATA VLI NEG?
04330 04E3 2B 05	EMI	PDFST1	YES. SRC MSB VS ADC
04340 04E5 B9 04CE	ADC A	TEMPDF	NO. ADD WITH CARRY INB

```

04350 04E8 20 03      FRA    PDFST2    AND STORE
04360 04FA B2 04CE PDFST1 SBC A    TEMPDF    VLU IS NEG. SBC
04370 04FC B7 04CE PDFST1 STA A    TEMPDF    AND STORE
04380 04FD 17          TBA          RECOVER VLU AND ADD (CARRY)
04390 04F1 BB 04CF      ADD A    TEMPDF+1 FOR PROPER ADDR CALCULATION
04400 04F4 B7 04CF      STA A    TEMPDF+1
04410 04F7 86 00      LDA A    #0
04420 04F9 C5 80      BIT B    #$80    IS THE DATA NEG?
04430 04FB 2B 05      BMI    PDFST3    YES. SBC MSB VS ADC
04440 04FD B9 04CE      ADC A    TEMPDF    NO. ADD WITH CARRY MSB BYTE
04450 0500 20 03      FRA    PDFST4    AND STORE
04460 0502 B2 04CE PDFST3 SBC A    TEMPDF    VLU IS NEG. SUB WITH CARRY
04470 0505 B7 04CE PDFST4 STA A    TEMPDF    AND STORE
04480 0508 FE 04CE      LDX    TEMPDF    NOW LOAD CALC ADDR FOR INDEX
04490 050B AE 00      LDS    0,X      GRAB VLU IN CALC ADDR
04500 050D 31          INS          INCREMENT IT
04510 050E AF 00      STS    0,X      AND STORE IT BACK IN BUFFER
04520 0510 BE 04D0      LDS    TEMPST    RECOVER STACK POINTER
04530 0513 17          TBA          RECOVER ORIGINAL DATA
04540 0514 39          RTS
04550      *
04560      *FUNC: EXEC JUMPS
04570      *INPUTS: ACCUMULATORS
04580      *OUTPUTS: NONE
04590      *CALLS: SAVEFL, FILHDR (EKG-EXEC) VIA ADDR BUFFERS
04600      *DESTROYS: X,A,B,CC
04610      *PURPOSE: TO JUMP TO DESIRED ROUTINES VIA RELOC ADDR
04620      *
04630 0515 FE 1CA3 SAVEFL LDX    SAVELC    GET ADDR OF SAVEFL
04640 0518 6E 00      JMP    0,X      JUMP TO SAVEFL
04650      *
04660 051A FE 1CA1 FILHDR LDX    FILHLC    GET ADDR OF FILHLC
04670 051D 6E 00      JMP    0,X      JUMP TO SAVEFL
04680      *
04690      *****
04700      * END SAMPLE
04710      *****
04720      *
04730      *
04740      *FUNCTION :CALINT
04750      *INPUTS (REG) :NONE
04760      *OUTPUTS (REG) :NONE
04770      *CALLS :NOTHING
04780      *DESTROYS (REG):NONE (INTERRUPT HANDLER)
04790      *PURPOSE :THIS ROUTINE IS USED FOR CALIBRATING
04800      * THE MAX NUMBER OF LOGS POSSIBLE
04810      * DURING AN INTERRUPT PERIOD. THIS
04820      * ROUTINE UPDATES THE DONTST FLAG AND
04830      * RESETS THE ST6800 INT FLIP FLOP &
04840      * THE RETURNS.
04850      *
04860 051F 7C 0232 CALINT INC    DONTST    INCREMENT DONT TEST FLAG
04870 0522 B7 E400      STA A    ADCZRO    CLR ST6800 INTR FLIP FLAG
04880 0525 01          NOP

```

NOCPRS-7

```
04899 0526 3B          RTI
04900                  *
04901                  *****
04920                  *  END CALINT
04930                  *****
04940                  *
04950                  *
04960                  *****
04970                  *
04980                  *  END OF NOCPRS OVERLAY ROUTINES
04990                  *
05000                  *****
05010                  *
05020                  END
```

TOLAN-8

```

00030 *****
00040 *
00050 * OVERLAY NAME: TOLAN-A
00060 * AUTHOR : CAPT. MEL TOUNGILL
00070 * VERSION : 1.8
00080 * VERSION DATE : 22 OCT 80
00090 *
00100 *****
00110 *
00120 * OVERLAY DESCRIPTION
00130 *
00140 * THIS OVERLAY SAMPLES THE HKG DATA VIA THE
00150 * A/D CONVERTERS, ROUNDS THE DATA TO 8 BITS,
00160 * AND THEN COMPRESSES THE DATA VIA THE TOLAN-A
00170 * ALGORITHM. THE COMPRESSED DATA IS THEN STORED INTO
00180 * MEMORY DATA FILE FROM 3C00-7FFF.
00190 *
00200 *****
00210 *
00220 * START OF TOLAN-A
00230 *
00240 *****
00250 *
00260 0100          ORG      $0100      OVERLAY START ADDRESS
00270 *
00280          OPT      O          ASSB OPT-GLN ORJ FILE
00290          OPT      NOG       ASSB OPT-SUPPRESS FULL FCC LI
00300 *
00310 *****
00320 *
00330 * LABEL DECLARATIONS
00340 *
00350 * SUPPORT SUBROUTINE ADDRESSES
00360 *
00370 CA87  OUTPUT EQU  $CA87  EPRG1DOS. ALPH STRING TO COMS
00380 CA8F  OUTNCR EQU  $CA8F  EPRG1DOS. ALPH STRING, NO CR/LF
00390 CA2C  KEYPD  EQU  $CA2C  EPRG1DOS. CONSOLE INPUT ROUTE
00400 CA36  KEYBD0 EQU  $CA36  EPRG1DOS. CONSOLE INPUT, NO ?
00410 1D00  START  EQU  $1D00  EKG-EXEC. FILE CREATE ROUTINE
00420 *
00430 * DATA BUFFERS
00440 *
00450 3400  HDRSTR EQU  $3400  DATA MEMORY BUFFER HEADER
00460 FFF8  IRQVEC EQU  $FFF8  INTERRUPT VECTOR ADDR
00470 1C96  STKSAV EQU  $1C96  STACK SAVE BUFFER
00480 1C98  CRTTYP EQU  $1C98  COMPRESSION ALGOR FLAG
00490 3002  ENDBUF EQU  $3002  ADDR OF LAST CHAR PRINT
00500 1C9D  VLSAV  EQU  $1C9D  IRQ VECTOR SAVE BUFFER
00510 1CA1  FILHLC EQU  $1CA1  FILTER ADDR FLAG BUFFER
00520 1CA3  SAVHLC EQU  $1CA3  SAVHLC ADDR FLAG BUFFER
00530 3490  LOOPCT EQU  $3490  TOTAL NUM TIMES LOOP EXEC'D
00540 3494  SAMPL0 EQU  $3494  NUMBER OF SAMPLES TAKEN
00550 3496  LRCAL  EQU  $3496  AVG'D MAX LOOP COUNT/ITER
00560 3497  MAXZ   EQU  $3497  MAX VLU IN CH Z

```

00577	3498	MAXZLO EQU	\$3498	MAX VLU LOC IN CH Z
00589	349A	MINZ EQU	\$349A	MIN VLU IN CH Z
00591	349B	MINZLO EQU	\$349B	MIN VLU LOC IN CH Z
00601	349D	MAXY EQU	\$349D	MAX VLU IN CH Y
00610	349E	MAXYLO EQU	\$349E	MAX VLU LOC IN CH Y
00620	34A0	MINY EQU	\$34A0	MIN VLU IN CH Y
00630	34A1	MINYLO EQU	\$34A1	MIN VLU LOC IN CH Y
00640	34A3	MAXX EQU	\$34A3	MAX VLU IN CH X
00650	34A4	MAXXLO EQU	\$34A4	MAX VLU LOC IN CH X
00660	34A6	MINX EQU	\$34A6	MIN VLU IN CH X
00670	34A7	MINXLO EQU	\$34A7	MIN VLU LOC IN CH X
00680	34AC	DPABIT EQU	\$34AC	NUM OF BITS USED TO STR DTA
00690	34B0	XBITS EQU	\$34B0	NUM OF BITS USED TO STR X
00700	34B3	YBITS EQU	\$34B3	NUM OF BITS USED TO STR Y
00710	34B6	ZBITS EQU	\$34B6	NUM OF BITS USED TO STR Z
00720	34B9	TBITS EQU	\$34B9	NUM OF BITS USED TO STR T
00730	34BC	ACELOT EQU	\$34BC	# BITS FED TO VAR LEN CODE
00740	3600	XPDEFM EQU	\$3600	0 VLU LOC OF X PDEF
00750	3800	YPDEFM EQU	\$3800	0 VLU LOC OF Y PDEF
00760	3A00	ZPDEFM EQU	\$3A00	0 VLU LOC OF Z PDEF
00770	3B00	TPDF EQU	\$3B00	0 VAL LOC OF TIME VAR HIST
00780	3C00	SECZRO EQU	\$3C00	DATA STORE ADDR START
00790		*		
00800		* HARDWARE ADDRESSES		
00810		*		
00820	E400	ADCZRO EQU	\$E400	ADC CHANNEL ZERO
00830	E404	ADCTWO EQU	\$E404	ADC CHANNEL TWO
00840	E500	DACZRO EQU	\$E500	DAC CHANNEL ZERO
00850		*		
00860		*		
00870		*FUNCTION :TOLAN-A		
00880		*INPUTS :STATUS BUFFERS FROM EKG-EXEC		
00890		*OUTPUTS :DATA TO DISK		
00900		*CALLS :OUTPUT, FILLDR, KEYPD, OUTNCR		
00910		* KEYBD0, SAVEFL, START		
00920		*DESTROYS :ALL REGISTERS		
00930		*PURPOSE :TO COLLECT 3 CHANNELS OF EKG		
00940		* DATA AND STORE INTO MEMORY.		
00950		*		
00960		*		
00970	0100 0F	TOLANA SEI		
00980	0101 CE 0252	LDX	#STRMSG	
00990	0104 BD CA8F	JSR	OUTNCR	"THIS MODULE SAMPLES THE..."
01000	0107 ED CA2C	JSR	KEYPD	GET RESPONSE FROM CONSOLE
01010	010A FE 3002	LDX	INDRUP	
01020	010D 09	DEX		
01030	010E E6 00	IDA B	0,X	
01040	0110 C1 59	CMP B	#'Y	IS INPUT YES?
01050	0112 27 03	BEQ	TOLANI	YES. KEEP EXECUTING THIS ROUT
01060	0114 7E 1D00	JMP	START	NO. RIN TO EKG-EXEC
01070	0117 CE 5441	TOLANI LDX	#\$5441	FLAG COMPRESSION TYPE (TA)
01080	011A FF 1C98	STX	CRTYP	
01090	011D BD 023B	JSR	FILLDR	SET UP DATA FILE HEADER
01100	0120 CE 02D9	LDX	#IRMSG	

```

01110 0123 BD CA8F      JSR      OUTNCR      "INSURE SUBJECT" AND INC..."
01120 0126 BD CA36      JSR      KEYPDO
01130 0129 C3 3000 TOLAGO INX      $COUNTCO INITIALIZE THE FOR STATE LOG
01140 012C FF 024E      STX      BUFFER      STORE IN BUFFER POINTER
01150 012F 86 04        LDA A      #4          PICK UP COUNTER FOR TIME CAL
01160 0131 B7 0245      STA A      CALCNT      STORE IN BUFFER
01170 0134 86 80        LDA A      #S80       INITIALIZE BIT POINTER TO LEFT
01180 0136 B7 03C7      STA A      BITPTR
01190 0139 4F          CLR A          INIT COMPRESS VAR
01200 013A B7 0251      STA A      FIRST      CLEAR FIRST SAMPLE FLAG
01210 013D B7 03BC      STA A      XSLOPE
01220 0140 B7 03ED      STA A      YSLOPE
01230 0143 B7 03BE      STA A      ZSLOPE
01240 0146 FE FFF8      LDX      IRQVEC      PICK UP CURREN IRQ VECTOR
01250 0149 FF 1C9D      STX      VECSAV      SAVE IN BUFFER
01260 014C CE 06AE TIMCAL LIX      $CALINT     GET INTR VECTOR ADDR FOR CAL
01270 014F FF FFF8      STX      IRQVEC      PUT IN VECTOR ADDRESS
01280 0152 86 00        LDA A      #0          INIT COUNTER FOR 256 TEST LOG
01290 0154 B7 0250      STA A      DONTST      STORE IN DONTST BUFFER
01300 0157 CE 0000      LDX      #0          CLEAR LOOP COUNTER
01310 015A FF 024A      STX      IRQCNT
01320 015D FF 024C      STX      IRQCNT+2
01330 0160 FF E500      STX      DACZRO      PULSE INT ENABLE CIRCUIT
01340 0163 20 09      BRA      SPOO1      BRA TO COUNTING LOOP
01350                      *
01360                      *****
01370                      *
01380                      * BASIC TIMING LOOP FOR EFFICIENCY TEST
01390                      *
01400                      *****
01410                      *
01420 0165 01          SPLOOP NOP          DELAY TO MATCH TIME
01430 0166 01          NOP          IT TAKES TO EXECUTE THE
01440 0167 01          NOP          INCREMENT OF BYTES 3&4
01450 0168 01          NOP          OF THE LOOP COUNTER
01460 0169 01          NOP          WHEN COUNT CARRIES TO
01470 016A 01          NOP          HIGH 2 BYTES OF 4 BYTE
01480 016B 7E 0181      JMP      SPOO3      JUMP TO CONTINUE LOOP
01490 016E 0E          SPOO1 CLI          PREPARE FOR TIMER INT
01500 016F 3E          WAI          STOP PROCESSOR & WAIT
01510 0170 0E          SPOO4 CLI
01520 0171 FE 024C      LDX      IRQCNT+2 RET'ED FROM INT, INC COUNT
01530 0174 08          INX
01540 0175 FF 024C      STX      IRQCNT+2 SAVE IT
01550 0178 26 EB      BNE      SPLOOP      COUNT COME FFFF TO 0000?
01560 017A FE 024A      LDX      IRQCNT      YES, INC BYTES 3&4
01570 017D 08          INX
01580 017E FF 024A      STX      IRQCNT      SAVE COUNT
01590 0181 E6 0250 SPOO3 LDA A      DONTST      IS DONE TEST SATISFIED?
01600 0184 81 00      CMP A      #0
01610 0186 26 E8      BNE      SPOO4      NO, KEEP LOOPING
01620                      *
01630                      *****
01640                      *

```



```

01610                                * END BASIC TIMING LOOP FOR EFT TEST
01620                                *
01630                                *****
01640                                *
01650 0188 0F                          SEI          PREVENT SER TO MORE INTR
01700 0189 CE 4000                  LDX          #S4000    PICK UP INTR OFF WORD
01710 018C FF E500                  STX          DACZRO    DISABLE INTR FROM TIMER
01720 018F B7 E400                  STA A      ADCZRO    CLR ST6800 INT FLIP FLOP
01730 0192 01                      NOP
01740 0193 7A 0245                  DEC          CALCNT    DEC TIME CAL LOOP COUNTER
01750 0196 27 4F                    BRQ          SPDONE    IS COUNT 3? YES, GO SAVE DATA
01760 0198 B6 0245                  LDA A      CALCNT    NO. GET COUNT
01770 019E 81 03                    CMP A      #3          IS COUNT 3?
01780 019D 26 25                    RNE          MAINSP    NO. GO TO MAIN SPL LOOP
01790 019F FE 024C                  LDX          IRQCNT+2 YES. SAVE TIME LOOP COUNTER
01800 01A2 FF 0248                  STX          CALONE
01810 01A5 CE 03C8                  LDX          #SAMPLE    NOW PUT SAMPLE ADDR IN IRQVEC
01820 01A8 FF FFF8                  STX          IRQVEC
01830 01AB 86 00                    LDA A      #80          SET UP DONTST
01840 01AD B7 0250                  STA A      DONTST
01850 01B0 7F 0251                  CLR          FIRST    RESET FIRST SMP FLG
01860 01B3 CE 0000                  LDX          #0          CLR LOOP COUNT AND SAMPLE COU
01870 01B6 FF 024A                  STX          IRQCNT
01880 01B9 FF 024C                  STX          IRQCNT+2
01890 01BC FF 3494                  STX          SAMINO
01900 01BF FF E500                  STX          DACZRO    ENABLE TIMER INTERRUPTS
01910 01C2 20 AA                    BRA          SPOOL     GO WAIT FOR INTERRUPT
01920 01C4 81 02                    MAINSP CMP A      #2          IS COUNT 2?
01930 01C6 26 10                    RNE          PARSAR    NO. SAVE PARAMETER FROM ERG
01940 01C8 86 55                    LDA A      #S55        SET FIRST FLAG
01950 01CA B7 0251                  STA A      FIRST
01960 01CD B7 0250                  STA A      DONTST    SET DONTST FOR SMPL COLLEC
01970 01D0 CE 0000                  LDX          #0          ENABLE INTR CLOCK
01980 01D3 FF E500                  STX          DACZRO
01990 01D6 20 96                    BRA          SPOOL     GO WAIT FOR INTERRUPT
02000 01D8 FE 024A PARSAR LDX          IRQCNT    SAVE LOOP COUNT
02010 01DB FF 3490                  STX          LOOPCP
02020 01DE FE 024C                  LDX          IRQCNT+2
02030 01E1 FF 3492                  STX          LOOPCP+2
02040 01E4 7E 014C                  JMP          TINCAL    GO EXECUTE ANOTHER TIME CAL
02050 01E7 B6 0249 SPDONE LDA A      CALONE+1    AVERAGE TWO TIME CAL RUNS
02060 01EA BB 024D                  ADD A      IRQCNT+3    ADD LSB BYTES
02070 01ED B7 0247                  STA A      CALZRO+1    STORE IN BUFFER
02080 01F0 B6 0248                  LDA A      CALONE    PICK UP MSB OF FIRST
02090 01F3 B9 024C                  ADC A      IRQCNT+2    ADD WITH CARRY
02100 01F6 B7 0246                  STA A      CALZRO    STORE IN BUFFER
02110 01F9 77 0246                  ASR          CALZRO    DIVIDE BY 2 TO AVERAGE
02120 01FC 76 0247                  ROR          CALZRO+1
02130 01FF 86 02          LDA A      #8          SET UP COUNTER FOR 256
02140 0201 77 0246 CALSHF ASR          CALZRO    NOW DIVIDING BY 256 TO GET
02150 0204 76 0247                  ROR          CALZRO+1    LOOPS/INTERUPT
02160 0207 4A                    DEC A
02170 0208 26 F7                    PNE          CALSHF    8 SHIFTS YET?
02180 020A B6 0247                  LDA A      CALZRO+1    PICK UP SHIFTED RESULT

```

TOLAN-8

02190 020D B7 3496	STA A	LICAL	STORE IN FILE FOR PLOTTER
02200 0210 BD 059E	JSR	DTACNT	NO. COUNT # OF DATA BITS SEC
02210 0212 11 1000	LDX	DTACNT	NO. COUNT # OF DATA BITS SEC
02220 0216 FF F1F8	STX	DTACNT	NO. COUNT # OF DATA BITS SEC
02230 0219 CE 0336	LDX	#SAVMSG	
02240 021C BD CA8F	JSR	OUTNCR	"SAMPLING COMPLETE..."
02250 021F BD CA2C	JSR	KEYED	GET SAVE DECISION
02260 0222 FE 3002	LDX	ENDBUF	
02270 0225 09	DEX		
02280 0226 E6 00	LDA B	0,X	
02290 0228 C1 59	CMP B	#'Y	IS DECISION YES?
02300 022A 26 0C	BNE	FXJMP	NO. RTN TO ERG-EXEC
02310 022C BD 0240	JSR	SAVEFL	YES. SAVE THE FILE ON DISK
02320 022F CE 03A4	LDX	#SDONE	
02330 0232 BD CA8F	JSR	OUTNCR	
02340 0235 BD CA2C	JSR	KEYED	
02350 0238 7E 1D00	EXJMP	JMP	START JUMP BACK TO ERG-EXEC
02360	*		
02370	*FUNC: RELOC JUMP		
02380	*INPUTS: A		
02390	*OUTPUTS: NONE		
02400	*CALLS: FILHDR,SAVEFL (ERG-EXEC ROUTINES)		
02410	*DESTROYS: A,R,CC,X		
02420	*PURPOSE : THIS ROUTINE INABLES TOLAN-A TO BE RELOCA		
02430	*WITHOUT WORRY OF CHANGING CALLS TO ERG-EXEC.		
02440	*		
02450 023E FE 1CA1	FILHDR LDX	FILHLC	GET FILHDR ADDRESS
02460 023E 6E 00	JMP	0,X	JUMP TO IT
02470	*		
02480 0240 FE 1CA3	SAVEFL LDX	SAVEFLC	GET SAVEFL ADDRESS
02490 0243 6E 00	JMP	0,X	JUMP TO IT
02500	*		
02510 0245 0001	CALCNT RMB	1	TEST CAL LOOP COUNTER
02520 0246 0002	CALZPO RMB	2	AVG'D CAL LOOPS DURING DAT 00
02530 0248 0002	CALCIE RMB	2	TEMP BUF FOR PRE COLLIC CAL
02540 024A 0004	IRQCNT RMB	4	TEMP INTR LOOP COUNTER
02550 024E 0002	BUTYTR RMB	2	BUFFER POINTER OF NEXT ANAIL
02560 0250 0001	DONTEST RMB	1	DONE TEST FLAG FOR NEW FULL
02570 0251 0001	FIRST RMB	1	FIRST A/D SAMPLE FLAG
02580	*		
02590 0252 1A07	STRMSG FDB	\$1A07	
02600 0254 54	FCC		/THIS MODULE SAMPLES THE ERG /
02610 0270 49	FCC		/INPUT AND STORES THE DATA /
02620 028A 0D0A	FDB		\$0D0A
02630 028C 43	FCC		/COMPRESSED VIA MAJOR TOLAN-A./
02640 02A9 0D0A	FDB		\$0D0A,\$0D0A
02650 02AD 44	FCC		/DO YOU WISH TO INTERRUPT THIS /
02660 02C9 4D	FCC		/MODULE (Y OR N)/
02670 02D8 04	FCB	4	
02680 02D9 1A	IRQMSG FCB	\$1A	
02690 02DA 49	FCC		/INSURE SUBJECT AND ERG DEVICE /
02700 02F8 52	FCC		/READY! /
02710 02FE 0D0A	FDB		\$0D0A,\$0D0A,\$0D0A
02720 0304 50	FCC		/PRESS RETURN, THEN CLOSE INTERRUPT/

TOLAN-8

```

02730 0326 20          FCC      / ENABLE SATEMCH./
02740 0335 04          FCB      4
02750 0336 1A07        SATING FCB      S1A07
02760 0338 53          FCC      /SATING CALIBER. HAVING CALIB./
02770 0357 13          FCC      /INTERMPT ENABLE SWITCH./
02780 036F 0D0A        FDB      S0D0A,S0D0A,S0D0A
02790 0375 44          FCC      /DO YOU WISH TO SAVE THIS DATA ON /
02800 0396 44          FCC      /DISK (Y OR N)/
02810 03A3 04          FCB      4
02820 03A4 20          SDONE FCC      / PRESS RETURN/
02830 03B1 04          FCB      4
02840                  *
02850                  *
02860                  * END TOLAN-A
02870                  *
02880                  *
02890                  *FUNCTION :SAMPLE
02900                  *INPUTS :STATUS BUFFERS
02910                  *OUTPUTS :COMPISSHD,ROUNDED DATA IN MEM BUFF
02920                  *CALLS :NOTHING
02930                  *PURPOSE :THIS ROUTINE SAMPLES THE EKG LEADS,
02940                  * ROUNDS THE VALUES TO 8 BITS (FROM 12)
02950                  * CALCULATES MAX,MIN,# OF BITS, SAMPLES
02960                  * ETC., AND SAVES THESE PARAM AND DATA
02970                  * TO DATA MEM FILE.
02980                  *
02990                  *
03000 03E2 0001        SHFBUF RMB      1          TEMP SHIFT BUFFER
03010 03E3 0001        TEMPA1 RMB      1          TEMP REG, MSB CH 0
03020 03E4 0001        TEMPA1 RMB      1          TEMP REG, LSB CH 0
03030 03E5 0001        TEMPA2 RMB      1          TEMP REG, MSB CH 1
03040 03E6 0001        TEMPA2 RMB      1          TEMP REG, LSB CH 1
03050 03E7 0001        TEMPA3 RMB      1          TEMP REG, MSB CH 2
03060 03E8 0001        TEMPA3 RMB      1          TEMP REG, LSB CH 2
03070 03E9 0001        XDATA RMB      1          TEMP X CH DTA BUF
03080 03FA 0001        YDATA RMB      1          TEMP Y CH DTA BUF
03090 03EB 0001        ZDATA RMB      1          TEMP Z CH DTA BUF
03100 03EC 0001        XSLOPE RMB      1          1ST DIFF X(N)-X(N+1)
03110 03ED 0001        YSLOPE RMB      1          1ST DIFF Y(N)-Y(N+1)
03120 03EE 0001        ZSLOPE RMB      1          1ST DIFF Z(N)-Z(N+1)
03130 03EF 0001        XEXP RMB      1          EXPECTED VAL OF X(N+1)
03140 03F0 0001        YEXP RMB      1          EXPECTED VAL OF Y(N+1)
03150 03F1 0001        ZEXP RMB      1          EXPECTED VAL OF Z(N+1)
03160 03F2 0001        XACCEL RMB      1          DIFF X(N+1)-E(X(N+1))
03170 03F3 0001        YACCEL RMB      1          DIFF Y(N+1)-E(Y(N+1))
03180 03F4 0001        ZACCEL RMB      1          DIFF Z(N+1)-E(Z(N+1))
03190 03F5 0001        SIFLG RMB      1          I HAVE JUST GIVEN 1100
03200 03F6 0001        ICTF RMB      1          TIME CTR VAR (<127)
03210 03F7 0001        BITPTR RMB      1          BIT POINTER FOR SET & RESET
03220                  *
03230 03C8 FF 3494      SAMPLE LDX      SAMPNO    GET CUR SAMPLE COUNT
03240 03C8 08          INX
03250 03CC FF 3494      STX      SAMPNO
03260 03CF DF 1C96      SIS      SIKSAV    SAVE STACK PTR

```

03270 03D2 8E 3280	LDS	#S3280	INITIALIZE STACK IN USED MEM
03280 03D3 15 03F1	LDA A	FIRST	OK FIRST FLAG
03290 03D4 15 03F1	LDX	BUFPTR	IF IT TRUE
03300 03DA 1D 05C8	JSR	SAMPLE0	YES. SAMPLE & IN IN
03310 03DD 1E 024E	LDX	BUFPTR	NOW SET FIRST DATA VALUES
03320 03E0 16 03F9	LDA A	XDATA	
03330 03E3 A7 00	STA A	0,X	
03340 03E5 08	INX		
03350 03E6 16 03FA	LDA A	YDATA	
03360 03E9 A7 00	STA A	0,X	
03370 03EB 08	INX		
03380 03EC 16 03FB	LDA A	ZDATA	
03390 03EF A7 00	STA A	0,X	
03400 03F1 08	INX		
03410 03F2 FF 024E	STX	BUFPTR	UPDATE NEW BUFPTR
03420 03F5 85 55	LDA A	#S55	SET STRFLG & FIRST FLAG
03430 03F7 B7 03C5	STA A	STRFLG	
03440 03FA B7 0251	STA A	FIRST	
03450 03FD 86 08	LDA A	#8	UPDATE X,Y,Z BIT CNTS FOR IN
03460 03FF B7 34E2	STA A	XBITS+2	
03470 0402 B7 34E5	STA A	YBITS+2	
03480 0405 B7 34E8	STA A	ZBITS+2	
03490 0408 7F 0250	CLR	DNTEST	CLEAR DOWNTEST FLAG FOR 1 SML
03500 040B BE 1C96	LDS	STKSAV	
03510 040E 3E	RTI		
03520 040F 7D 03C5	COMPR1	TST STRFLG	HAS DATA JUST BEEN STORED ?
03530 0412 27 05	BEQ	COMPR1	NO. KEEP COUNTING TIME
03540 0414 86 01	LDA A	#1	YES. SET TIME COUNT TO 1
03550 0416 17 03C6	STA A	ICNT	
03560 0419 B6 03F9	COMPR1	LDA A	GET X VLU
03570 041C BB 03FC	ADD A	XSI OPE	ADD {X(N-1)-X(N-2)}
03580 041F B7 03FF	STA A	XEXP	CREATE EXPECTED VAL OF X(N+1)
03590 0422 16 03FA	LDA A	YDATA	GET LAST Y SML VALUE
03600 0425 BB 03FD	ADD A	YSLOPE	ADD {(Y(N-1)-Y(N-2))}
03610 0428 B7 03C0	STA A	XEXP	CREATE EXPECTED VAL OF Y(N+1)
03620 042B 16 03FB	LDA A	ZDATA	GET LAST Z SML VALUE
03630 042E BB 03FE	ADD A	ZSLOPE	ADD {(Z(N-1)-Z(N-2))}
03640 0431 B7 03C1	STA A	XEXP	CREATE EXPECTED VAL OF Z(N+1)
03650 0434 1D 05C8	JSR	SAMPLE0	NOW GO GET NEXT SML. N=N+1
03660 0437 16 03F9	LDA A	XDATA	CALC DIF X(N)-L(X(N))
03670 043A 1D 03FF	SUB A	XEXP	
03680 043D B7 03C2	STA A	XACCEL	SAVE DIF IN VAR
03690 0440 16 03FA	LDA A	YDATA	CALC DIF Y(N)-L(Y(N))
03700 0443 1D 03C0	SUB A	XEXP	
03710 0446 B7 03C3	STA A	YACCEL	SAVE DIF IN VAR
03720 0449 16 03FB	LDA A	ZDATA	CALC DIF Z(N)-L(Z(N))
03730 044C B0 03C1	SUB A	XEXP	
03740 044F B7 03C4	STA A	ZACCEL	SAVE DIF IN VAR
03750 0452 16 03C2	LDA A	XACCEL	
03760 0455 26 21	ENE	COMPR2	ARE ALL INTERPOLATIONS GOOD?
03770 0457 16 03C3	LDA A	YACCEL	IF NO, NEED TO SET ALL DIFF'S
03780 045A 26 1C	ENE	COMPR2	
03790 045C 16 03C2	LDA A	XACCEL	
03800 045F 26 17	ENE	COMPR2	

03810	0461	B6	03C6	LDA A	ICNT	UPDATE RUN TIME COUNTER
03820	0464	AC		INC A		
03830	0467	2F	71	LDZ	CONTR	TEST FOR CONTR
03840	0467	2F	71	LDZ	CONTR	YES. CONTINUE DATA & DATA IS
03850	0469	B7	03C6	STA A	ICNT	NO. INC ICNT & RTI
03860	046C	7F	03C5	CLR	SIRFLG	RESET I HAVE JUST STORED FLG
03870	046F	86	55	LDA A	#555	
03880	0471	B7	0250	STA A	DONEST	KEEP DONE TEST SET
03890	0474	BE	1C96	LDS	STKSAV	RETRIEVE STACK
03900	0477	3B		RTI		
03910						
03920	0478	CE	34B0	COMPR2	LDX	#XBITS
03930	047E	B6	03C2	LDA A	XACCEL	UPDATE X BIT CTR
03940	047E	BD	0568	JSR	BITCNT	
03950	0481	CE	3600	LDX	#XPDEF	UPDATE CH X FREQ OF OCCUR HIS
03960	0484	B6	03C2	LDA A	XACCEL	
03970	0487	BD	066B	JSR	PDFSTR	
03980	048A	B6	03C2	LDA A	XACCEL	
03990	048D	B6	03C2	LDA A	XACCEL	
04000	0490	BB	03FC	ADD A	XSLOPE	CALC & SAVE X(N-1)-X(N-2)
04010	0493	B7	03FC	STA A	XSLOPE	
04020	0496	CE	34B3	LDX	#YBITS	UPDATE Y BIT CTR
04030	0499	B6	03C3	LDA A	YACCEL	
04040	049C	BD	0568	JSR	BITCNT	
04050	049F	CE	3600	LDX	#YPDEF	UPDATE CH Y FREQ OF OCCUR HIS
04060	04A2	B6	03C3	LDA A	YACCEL	
04070	04A5	BD	066B	JSR	PDFSTR	
04080	04A8	B6	03C3	LDA A	YACCEL	
04090	04AB	B6	03C3	LDA A	YACCEL	
04100	04AE	BB	03ED	ADD A	YSLOPE	CALC & SAVE Y(N-1)-Y(N-2)
04110	04B1	B7	03ED	STA A	YSLOPE	
04120	04B4	CE	34B6	LDX	#ZBITS	UPDATE Z BIT CTR
04130	04B7	B6	03C4	LDA A	ZACCEL	
04140	04BA	BD	0568	JSR	BITCNT	
04150	04BD	CE	3600	LDX	#ZPDEF	UPDATE CH Z FREQ OF OCCUR HIS
04160	04C0	B6	03C4	LDA A	ZACCEL	
04170	04C3	BD	066B	JSR	PDFSTR	
04180	04C6	B6	03C4	LDA A	ZACCEL	
04190	04C9	B6	03C4	LDA A	ZACCEL	
04200	04CC	BB	03FE	ADD A	ZSLOPE	CALC & SAVE Z(N-1)-Z(N-2)
04210	04CF	B7	03FE	STA A	ZSLOPE	
04220	04D2	CE	34B9	LDX	#TBITS	NOW UPDATE TIME BIT CTR
04230	04D5	86	06	LDA A	#6	
04240	04D7	BD	0568	JSR	BITCNT	
04250	04DA	B6	03C2	COMPR3	LDA A	XACCEL
04260	04DD	BD	051B	JSR	CONSTR	PICK UP XACCEL VAL
04270	04E0	B6	03C3	LDA A	YACCEL	GO CODE & STORE IT
04280	04E3	BD	051B	JSR	CONSTR	PICK UP XACCEL VAL
04290	04E6	B6	03C4	LDA A	ZACCEL	GO CODE & STORE IT
04300	04E9	BD	051B	JSR	CONSTR	PICK UP XACCEL VAL
04310	04EC	B6	03C6	LDA A	ICNT	GO CODE & STORE IT
04320	04EF	CE	3F00	LDX	#TDEF	NOW GET ICNT FOR TIME CTR
04330	04F2	BD	066B	JSR	PDFSTR	UPDATE TIME HISTOGRAM
04340	04F5	BD	066A	JSR	RESET	PUT IN 0 FOR INITIAL VAL FOR

AD-A100 799 AIR FORCE INST OF TECH WRIGHT-PATTERSON AFB OH SCHOO--ETC F/6 6/5
ANALYSIS AND PERFORMANCE EVALUATION OF ELECTROCARDIOGRAM DATA C--ETC(U)
DEC 80 M D TOWNSEND
UNCLASSIFIED AFIT/GE/EE/80D-46 NL

AIR FORCE INST OF TECH WRIGHT-PATTERSON AFB OH SCHOO--ETC F/G 6/5
ANALYSIS AND PERFORMANCE EVALUATION OF ELECTROCARDIOGRAM DATA C--ETC(U)
DEC 80 M D TOWNSEND
AFIT/GE/EE/80D-46 NL

NL

3 of 3

AG A
100789

END
DATE
FILMED
7-8
DTIC

04350	04F8 79 03C6	ROL	ICNT	ROL 7 BIT ICNT TO MS POSITION
04360	04FB 86 08	LDA A	#8	SET COUNTER FOR TIME STORE LOOP
04370	04FD 4A	COMPR4	DEC A	DECREMENT LOOP COUNTER
04380	04FE 27 0F	BEQ	COMPR6	TIME STORED YET?
04390	0500 79 03C6	ROL	ICNT	NO. ROTATE TIME WORD ONCE LEFT
04400	0503 25 05	BCS	COMPR5	1 IN CARRY?
04410	0505 BD 055A	JSR	RESET	NO. PUT 0 IN MEM FILE
04420	0508 20 F3	BRA	COMPR4	
04430	050A BD 0534	COMPR5	JSR	SET
04440	050D 20 EE	BRA	COMPR4	YES. PUT 1 IN MEM FILE
04450	050F 86 55	COMPR6	LDA A	#55
04460	0511 B7 03C5	STA A	STRFLG	NOW SET "I HAVE JUST STORED"
04470	0514 B7 0250	STA A	DONTST	SET DONTST FLAG
04480	0517 BE 1C96	LDS	STKSAV	
04490	051A 3B	RTI		
04500		*		
04510	051B BD 055A	COMSTR	JSR	RESET
04520	051E 4D	TST A		PUT 0 IN MEM FILE FOR DELIM
04530	051F 2A 06	BPL	COMST1	NOW TEST STATUS OF A
04540	0521 40	NEG A		IS IT +?
04550	0522 BD 0534	JSR	SET	NO. GET 2'S COMP FOR MAG
04560	0525 20 03	BRA	COMST2	NO. STORE 1 FOR SIGN BIT
04570	0527 BD 055A	COMST1	JSR	RESET
04580	052A 4D	COMST2	TST A	YES. STORE 0 FOR SIGN BIT
04590	052B 27 06	BEQ	COMST3	NOW CNT DOWN & STORE DATA
04600	052D BD 0534	JSR	SET	CNT 0?
04610	0530 4A	DEC A		NO. STORE 1 TO MEM FILE
04620	0531 20 F7	BRA	COMST2	DEC COUNTER
04630	0533 39	COMST3	RTS	
04640		*		
04650		*		
04660		*FUNC:	SET	
04670		*INPUTS:	BUFPTR, BITPTR	
04680		*OUTPUTS:	BIT SET POINTED TO BY ABOVE	
04690		*CALLS:	NOTHING	
04700		*DESTROYS:	B, X, CC	
04710		*PURPOSE:	THIS ROUTINE SETS THE BIT POINTED TO	
04720		* BY BITPTR & BUFPTR AND THEN UPDATES THESE COUNTERS		
04730		*		
04740	0534 FE 024E	SET	LDX	BUFPTR
04750	0537 F6 03C7		LDA B	BITPTR
04760	053A EA 00		ORA B	0,X
04770	053C E7 00		STA B	0,X
04780	053E 0C	SET0	CLC	
04790	053F 76 03C7		ROR	BITPTR
04800	0542 24 15		BCC	SETRTS
04810	0544 76 03C7		ROR	BITPTR
04820	0547 08		INX	
04830	0548 FF 024E		STX	BUFPTR
04840	054B F6 024E		LDA B	BUFPTR
04850	054E C1 80		CMP B	#80
04860	0550 26 07		ENE	SETRTS
04870	0552 7F 0250		CLR	DONTST
04880	0555 BE 1C96		LDS	STKSAV

```

04850 0558 3B          RTI
04900 0559 39          SETRTS RTS
04910                  *
04920                  *FUNC: RESET
04930                  *INPUTS: BUFPTR,BITPTR
04940                  *OUTPUTS: BIT RESET IN MEM FILE PT'ED TO BY ABOVE
04950                  *CALLS: NOTHING
04960                  *DESTROYS: B,X,CC
04970                  *PURPOSE: THIS ROUTINE RESETS THE BIT IN THE MEMFILE
04980                  * POINTED TO BY BUFPTR & BITPTR & UPDATES THESE COUN
04990                  *
05000 055A FE 024E RESET LDX      BUFPTR
05010 055D F6 03C7      LDA B      BITPTR
05020 0560 53           COM B
05030 0561 E4 00        AND B      0,X      RESET BIT IN MEM WORD
05040 0563 E7 00        STA B      0,X      SAVE BIT IN MEM FILE
05050 0565 7E 053E      JMP       SET0     UPDATE CNTR'S AND CK IF MEM F
05060                  *
05070                  *FUNC: BITCNT
05080                  *INPUTS: ACCEL VLU IN A, BIT CNTR IN X
05090                  *OUTPUTS: UPDATED BIT COUNT
05100                  *CALLS: NOTHING
05110                  *DESTROYS: A,B,CC
05120                  *PURPOSE: THIS ROUTINE INC 3 BYTE BIT COUNTER WITH
05130                  * NUMBER OF BIT USED TO STORE X,Y,Z, & TIME.
05140                  *
05150 0568 2A 01        BITCNT BPL      BITCN1  IS ACCEL +
05160 056A 40           NEG A          NO. GET ABS VLU
05170 056B 4C          BITCN1 INC A      ADD 1 FOR DELIM 0 BIT
05180 056C 4C          INC A          ADD 1 FOR SIGN BIT
05190 056D C6 00        LDA B      #0     RESET B WITHOUT AFFECTING CAR
05200 056F AB 02        ADD A      2,X     INC LSB BYTE OF COUNT PTED TO
05210 0571 A7 02        STA A      2,X
05220 0573 86 00        LDA A      #0
05230 0575 E9 01        ADC B      1,X
05240 0577 E7 01        STA B      1,X
05250 0579 C6 08        LDA B      #8     NOW UP COUNT OF ACTUAL DATA B
05260 057B A9 00        ADC A      0,X
05270 057D A7 00        STA A      0,X
05280 057F 86 00        LDA A      #0
05290 0581 FB 34BF      ADD B      ACELCT+3 NO. INC ACELCT COUNT
05300 0584 F7 34BF      STA B      ACELCT+3
05310 0587 C6 00        LDA B      #0
05320 0589 B9 34BE      ADC A      ACELCT+2
05330 058C B7 34BE      STA A      ACELCT+2
05340 058F 86 00        LDA A      #0
05350 0591 F9 34BD      ADC B      ACELCT+1
05360 0594 F7 34BD      STA B      ACELCT+1
05370 0597 B9 34BC      ADC A      ACELCT
05380 059A B7 34BC      STA A      ACELCT
05390 059D 39          BITRTS RTS
05400                  *
05410                  *FUNC: DTACNT
05420                  *INPUTS: NUM ON SMPLS COLLECTED

```



```

05430          *OUTPUTS: BITCNT=TO SMPLS*24
05440          *CALLS: NOTHING
05450          *DESTROYS: A,B,X,CC
05460          *PURPOSE: TO CALC NUM OF BITS STORED DURING COLLECTI
05470          *
05480 059E FE 3494 DTACNT LDX      SAMPNO
05490 05A1 86 18   DTACN1 LDA A    #24
05500 05A3 C6 00           LDA B    #0
05510 05A5 BB 34AF          ADD A    DTABIT+3
05520 05A8 B7 34AF          STA A    DTABIT+3
05530 05AB 86 00           LDA A    #0
05540 05AD F9 34AE          ADC B    DTABIT+2
05550 05B0 F7 34AE          STA B    DTABIT+2
05560 05B3 C6 00           LDA B    #0
05570 05B5 B9 34AD          ADC A    DTABIT+1
05580 05B8 B7 34AD          STA A    DTABIT+1
05590 05BB F9 34AC          ADC B    DTABIT
05600 05BE F7 34AC          STA B    DTABIT
05610 05C1 09           DEX
05620 05C2 8C 0000          CPX      #0
05630 05C5 26 DA          BNE      DTACN1
05640 05C7 39           RTS
05650          *
05660          *FUNC: SAMPLO
05670          *INPUTS: DATA VIA A/D CONV
05680          *OUTPUTS: ROUNDED DATA TO XDTA,YDTA,ZDTA
05690          *CALLS:PDFSTR (UPDATE PDF HISTOGRAMS)
05700          *DESTROYS:A,B,X,CC
05710          *PURPOSE: THIS ROUTINES SAMPLES THE A/D AND THEN
05720          * ROUNDS THE DATA TO 8 BITS FROM 12. THE DATA
05730          * IS PLACED IN XDTA,YDTA,ZDTA AND THE COLLEC
05740          * STAT ARE UPDATED.
05750          *
05760 05C8 CE E400 SAMPLO LDX      #ADCZRO  NOW PULSE A/D TO START CONV
05770 05CB A7 00          STA A    0,X      ON CHANNEL 0
05780 05CD 01           NOP
05790 05CE A6 00          LDA A    0,X
05800 05D0 E6 01          LDA B    1,X
05810 05D2 B7 03B3          STA A    TEMPAL
05820 05D5 F7 03B4          STA B    TEMPBL
05830 05D8 08           INX
05840 05D9 08           INX
05850 05DA A7 00          STA A    0,X      NOW PULSE A/D TO CONV
05860 05DC 01           NOP                ON CHANNEL 1
05870 05DD A6 00          LDA A    0,X
05880 05DF E6 01          LDA B    1,X
05890 05E1 B7 03B5          STA A    TEMPAL
05900 05E4 F7 03B6          STA B    TEMPBL
05910 05E7 08           INX
05920 05E8 08           INX
05930 05E9 A7 00          STA A    0,X      NOW PULSE A/D TO CONV
05940 05EB 01           NOP                ON CH 2
05950 05EC A6 00          LDA A    0,X
05960 05EE E6 01          LDA B    1,X

```

TOLAN:-8

05970	05F0	B7	03B7	STA A	TEMPA3	
05980	05F3	F7	03B8	STA B	TEMPB3	
05990	05F6	CE	03F3	IDX	#TEMPA1	
06000	05F9	A6	00	SAMPL1 LDA A	0,X	
06010	05FB	E6	01	LDA B	1,X	
06020	05FD	47		ASR A		12 BIT TO 8 BIT ROUNDED CONV
06030	05FE	56		ROR B		
06040	05FF	47		ASR A		
06050	0600	56		ROR B		
06060	0601	47		ASR A		
06070	0602	56		ROR B		
06080	0603	47		ASR A		
06090	0604	56		ROR B		
06100	0605	F7	03B2	STA B	SHFBUF	SAVE 8 BIT RESULT OF SHIFT
06110	0608	86	00	LDA A	#0	AND ADD CARRY OUT OF LAST
06120	060A	B9	03B2	ADC A	SHFBUF	SHIFT ROUNDING UP OR DOWN
06130	060D	36		PSH A		SAVE TO MEM FILE BUFFER VIA S
06140	060E	08		INX		
06150	060F	08		INX		
06160	0610	8C	03B9	CPX	#TEMPA3+2	CHANNEL Z ROUNDED TO 8 BITS
06170	0613	26	E4	BNE	SAMPL1	NO. GO SAMPLE NEXT CHANNEL
06180	0615	32		PUL A		GET Z DATA
06190	0616	B7	03BB	STA A	ZDATA	
06200	0619	FE	3494	LDX	SAMPNO	GET CUR SAMPLE COUNT
06210	061C	B1	3497	CMP A	MAXZ	IS CUR Z MAX OVE SAMPLE SET?
06220	061F	2F	06	BLE	SPZMIN	NO. GO CHECK FOR MIN
06230	0621	B7	3497	STA A	MAXZ	YES. KEEP CUR VLU
06240	0624	FF	3498	STX	MAXZLO	KEEP CUR SAMPLE NUM
06250	0627	B1	349A	SPZMIN CMP A	MINZ	IS CUR Z MIN OVER SAMPLE SET?
06260	062A	2C	06	BGE	SPYMAX	NO. GO CHECK FOR Y MAX
06270	062C	B7	349A	STA A	MINZ	YES. KEEP CUR VLU
06280	062F	FF	349B	STX	MINZLO	KEEP CUR SAMPLE NUM
06290	0632	32		SPYMAX PUL A		GET CUR Y DATA
06300	0633	B7	03BA	STA A	YDATA	
06310	0636	B1	349D	CMP A	MAXY	IS CUR Y MAX OVER SAMPLE SET
06320	0639	2F	06	BLE	SPYMIN	NO. GO CHEK FOR MIN
06330	063B	B7	349D	STA A	MAXY	YES. KEEP CUR Y VLU
06340	063E	FF	349E	STX	MAXYLO	KEEP CUR SAMPLE NUM
06350	0641	B1	34A0	SPYMIN CMP A	MINY	IS CUR Y MIN OVER SAMPLE SET?
06360	0644	2C	06	BGE	SPXMAX	NO. GO CHEK FOR X MAX
06370	0646	B7	34A0	STA A	MINY	YES. KEEP CUR Y VLU
06380	0649	FF	34A1	STX	MINYLO	KEEP CUR SAMPLE NUM
06390	064C	32		SPXMAX PUL A		GET CUR X DATA
06400	064D	B7	03B9	STA A	XDATA	
06410	0650	B1	34A3	CMP A	MAXX	IS CUR X MAX OVER SAMPLE SET?
06420	0653	2F	06	BLE	SPXMIN	NO. GO CHECK X MIN
06430	0655	B7	34A3	STA A	MAXX	YES. KEEP CUR X VLU
06440	0658	FF	34A4	STX	MAXXLO	KEEP CUR SAMPLE NUM
06450	065B	B1	34A6	SPXMIN CMP A	MINX	IS CUR X MIN OVER SAMPLE SET?
06460	065E	2C	06	BGE	MMDONE	NO. EXIT MAX,MIN UPDATE
06470	0660	B7	34A6	STA A	MINX	YES. KEEP CUR X VLU
06480	0663	FF	34A7	STX	MINXLO	KEEP CUR SAMPLE COUNT
06490	0666	39		MMDONE RTS		
06500				*		

TOLAN-8

```

06510 0667 0002    TEMPDF RMB    2    TEMP WORKING BUFFER
06520 0669 0002    TEMPST RMB    2    TEMP STACK SAVE BUFFER
06530
*
06540 066B BF 0669 PDFSTR STS    TEMPST    SAVE STACK PTR
06550 066E FF 0667    STX    TEMPDF    SAVE PDF PTR
06560 0671 16    TAB
06570 0672 BB 0668    ADD A    TEMPDF+1    NOW CALCULATE ADDRESS IN INC
06580 0675 B7 0668    STA A    TEMPDF+1    IN THE PDF MEM BUFF
06590 0678 86 00    LDA A    #0
06600 067A C5 80    BIT B    #$80    IS THE DATA VLU NEG
06610 067C 2B 05    BMI    PDFST1    YES. SBC MSB VS ADC
06620 067E B9 0667    ADC A    TEMPDF    NO. ADC MSB
06630 0681 20 03    BRA    PDFST2    AND STORE
06640 0683 B2 0667 PDFST1 SBC A    TEMPDF    VLU IS NEG. SBC
06650 0686 B7 0667 PDFST2 STA A    TEMPDF    AND STORE
06660 0689 17    TBA    RECOVER VLU AND ADD (SUB) AGA
06670 068A BB 0668    ADD A    TEMPDF+1    FOR PROPER ADDR CALC
06680 068D B7 0668    STA A    TEMPDF+1
06690 0690 86 00    LDA A    #0
06700 0692 C5 80    BIT B    #$80    IS VLU NEG?
06710 0694 2B 05    BMI    PDFST3    YES. SBC MSB VS ADC MSB
06720 0696 B9 0667    ADC A    TEMPDF    NO. ADC MSB
06730 0699 20 03    BRA    PDFST4
06740 069B B2 0667 PDFST3 SBC A    TEMPDF    VLU IS NEG. SBC
06750 069E B7 0667 PDFST4 STA A    TEMPDF    AND STORE
06760 06A1 FE 0667    LDX    TEMPDF    NO LOAD CALC ADDR FOR INDEX
06770 06A4 AE 00    LDS    0,X    GRAB VLU IN CALC ADDR
06780 06A6 31    INS    INC IT
06790 06A7 AF 00    STS    0,X    AND STORE IT BACK IN BUFFER
06800 06A9 BE 0669    LDS    TEMPST    RECOVER STACK POINTER
06810 06AC 17    TBA
06820 06AD 39    RTS
06830
*
06840    *FUNC: CALINT
06850    *INPUTS: NONE
06860    *OUTPUTS: NONE
06870    *CALLS: NOTHING
06880    *DESTROYS :NO REG (INTR HANDLER)
06890    *PURPOSE: THIS ROUTINE IS USED FOR CALIBRATING
06900    * THE MAX NUMBER OF LOOPS POSSIBLE DURING
06910    * AN INTERRUPT PERIOD. THIS ROUTINE UPDATES THE
06920    * DONTST FLAG AND RESETS THE ST6800 INT FLIP FLOP
06930    * AND THEN RETURNS.
06940
*
06950 06AE 7C 0250 CALINT INC    DONTST    INCREMENT DONE TEST FLAG
06960 06B1 B7 E400    STA A    ADCZRO    CLR ST6800 INTR FLIP FLOP
06970 06B4 01    NOP
06980 06B5 3B    RTI
06990
*
07000    * END OF TOLAN-A OVERLAY ROUTINES
07010
*
07020    END

```

```

00030      *
00040      *      PROGRAM NAME:  DECPRS
00050      *      AUTHOR :  CAPT. MEL TOWNSEND
00060      *      VERSION :  1.3
00070      *      VERSION DATE :  29 NOV 80
00080      *
00090      *****
00100      *
00110      *      PROGRAM DESCRIPTION
00120      *
00130      *      THIS PROGRAM PERFORMS THE DECOMPRESSION OPERATI
00140      *      ON THE DATA COMPRESSED BY THE TOLAN-A,TOLAN-B,DOVE
00150      *      ETC.  THE ROUTINE ASSUMES THAT THE DATA IS IN MEMO
00160      *      AND READS THE COMPRESSION TYPE FROM THE MEMORY FIL
00170      *      BUFFER HEADER.  AFTER PROMPT TO THE TERMINAL, THIS
00180      *      PROGRAM BEGINS DECOMPRESSION AND OUTPUTS THE DECOM
00190      *      DATA TO AN ANALOG DISPLAY DEVICE (IE OSCILLISCOPE)
00200      *      D/A CONVERTER 0.  THE ROUTINE REQUIRES AN INTERRUPT
00210      *      CLOCK VIA THE A/D BOARD (ST6800) AND THE ASSUMPTIO
00220      *      IS MADE THAT THE INTERRUPT CLOCK FREQUENCY IS ADJU
00230      *      TO THE SAMPLE RATE AT WHICH THE DATA WAS TAKEN.
00240      *
00250      *****
00260      *
00270      *      START OF DECPRS
00280      *
00290      *****
00300      *
00310 0500      *      ORG      $0500      PROGRAM ORIGIN
00320      *
00330      *      OPT      O      ASSB OPT. LIST ASSMBLY
00340      *      OPT      NOG      ASSB OPT. SUPRESS FCC LIST
00350      *
00360      *****
00370      *
00380      *      LABLE DESCRIPTIONS
00390      *
00400      *      SUPPORT SUBROUTINES
00410      *
00420      CA8F      OUTINCR EQU      $CA8F      EPROMDOS. OUTPUT STRING
00430      CA2C      KEYBDO EQU      $CA2C      EPROMDOS. INPUT ALPH STRING
00440      0100      DISPLA EQU      $0100      EKG-EXEC. TERMINAL INTFC DRIV
00450      *
00460      *      DATA BUFFERS
00470      *
00480      0020      CHNLBF EQU      $0020      CHNL SELECT FLG FRM DISPLAY
00490      3C00      XINIT EQU      $3C00      INITIAL COND CH X
00500      3C01      YINIT EQU      $3C01      INITIAL COND CH Y
00510      3C02      ZINIT EQU      $3C02      INITIAL COND CH Z
00520      3002      ENDBUF EQU      $3002      END-OF-STRING FRM TERM INPUT
00530      FFF8      IRQVEC EQU      $FFF8      MASKABLE INTR JMP VECTOR ADDR
00540      3400      HDRSTR EQU      $3400      MEM FILE HEADER ADDR
00550      3C00      SECZRO EQU      $3C00      DATA START LOCATION
00560      *

```

```

00570      * HARDWARE ADDRESSES
00580      *
00590      E500  DACZRO EQU    $E500    D/A CH 0
00600      E400  ADCZRO EQU    $E400    A/D CH 0
00610      *
00620      *
00630      *FUNCTION: DECPRS
00640      *INPUTS : CHNL SELECT VIA CHNLBF
00650      *OUTPUTS : DECRPRSED DATA VIA D/A
00660      *CALLS : OUTNCR,KEYBD0,DISPLAY
00670      *DESTROYS (REG) : A,B,X,CC
00680      *PURPOSE : THIS IS THE SET UP ROUTINE WHICH INITIALI
00690      * THE PROGRAM CONSTANTS AND, DETERMINES WHAT TYPE OF
00700      * COMPRESSION WAS USED, AND ENABLES RECONSTRUCTION
00710      * INTERUPPTS.
00720      *
00730 0500 0F      DECPRS SEI
00740 0501 FE 3400      LDX      HDRSTR      GET CMPRS TYPE FROM MEM FILE
00750 0504 FF 0538      STX      NOTDN1      STR CPR TYPE IN ERROR MSG
00760 0507 8C 5441      CPX      #$5441      IS IT TOLAN-A?
00770 050A 27 61      BEQ      TADEC2      YES. GO DECPRS IT.
00780 050C CE 0518      LDX      #NOTDNE     NO. CPRS ROUTINE NOT DONE
00790 050F BD CA8F      JSR      OUTNCR      PRINT ERROR MSG & RETURN
00800 0512 BD CA2C      JSR      KEYBD0
00810 0515 7E 0100      JMP      DISPLA      RETURN TO DISPLAY ROUTINE
00820      *
00830 0518 1A07      NOTDNE FDB      $1A07,$0D0A
00840 051C 44      FCC      /DECOMPRESSION FOR FILE TYPE /
00850 0538 0002      NOTDN1 RMB      2
00860 053A 20      FCC      / NOT YET IMPLEMENTED./
00870 054F 0D0A      FDB      $0D0A,$0D0A
00880 0553 50      FCC      /PRESS ANY KEY TO CONTINUE/
00890 056C 04      FCB      4
00900      *
00910 056D 96 20      TADEC2 LDA A  CHNLBF      GET WHICH CHANNEL TO BE DECPR
00920 056F 81 00      CMP A  #0      IS IT CHNL X?
00930 0571 26 07      TADEC2 BNE TADEC1      NO. CHECK IF Y.
00940 0573 C6 58      LDA B  #'X      PUT ASCII X IN ERROR MSG
00950 0575 F7 0846      STA B  CINASC
00960 0578 20 10      BRA      TADEC3
00970 057A 81 01      TADEC1 CMP A  #1      IS IT CH Y?
00980 057C 26 07      BNE      TADEC2      NO. MUST BE Z.
00990 057E C6 59      LDA B  #'Y      PUT ASCII Y IN ERROR MSG
01000 0580 F7 0846      STA B  CHNASC
01010 0583 20 05      BRA      TADEC3
01020 0585 C6 5A      TADEC2 LDA B  #'Z      PUT ASCII Z IN ERROR MSG
01030 0587 F7 0846      STA B  CINASC
01040 058A CE 0811      TADEC3 LDX      #GOMSG      GET PROMPT MSG LOC
01050 058D BD CA8F      JSR      OUTNCR      NOW PRINT PROMPT
01060 0590 BD CA2C      JSR      KEYBD0
01070 0593 FE 3002      LDX      ENDBUF      GET ANSWER FROM TERMINAL
01080 0596 09      DEX
01090 0597 E6 00      LDA B  0,X
01100 0599 C1 4E      CMP B  #'N      WAS IT NO?

```

DECPRS-3

01110	059B	26	03	BNE	TADEC4	YES. EXIT ROUTINE
01120	059D	7E	0100	JMP	DISPLA	
01130	05A0	CF	3C03	TADEC4 LDX	#SECTPC+3	POSITION POINTER TO DATA STR
01140	05A3	FF	0611	STX	PUPPTR	
01150	05A6	86	80	LDA A	#\$80	INITIALIZE BITPTR TO LEFT BIT
01160	05A8	B7	0613	STA A	BITPTR	
01170	05AB	7F	0617	CLR	XSLOPE	INIT SLOPES TO ZERO (2ND INIT
01180	05AE	7F	0618	CLR	YSLOPE	
01190	05B1	7F	0619	CLR	ZSLOPE	
01200	05B4	B6	3C00	LDA A	XINIT	GET FIRST INIT COND
01210	05B7	B7	0614	STA A	XDTA	& PUT IN DATA BUF
01220	05BA	B6	3C01	LDA A	YINIT	
01230	05BD	B7	0615	STA A	YDTA	
01240	05C0	B6	3C02	LDA A	ZINIT	
01250	05C3	B7	0616	STA A	ZDTA	
01260	05C6	7F	0629	CLR	TARFLG	CLR DECERSION FLG
01270	05C9	7F	061A	CLR	DELT	INIT TIME COUNTER
01280	05CC	FE	FFF8	LDX	IRQVEC	PICK UP CURRENT INT VEC
01290	05CF	FF	061D	STX	VECSAV	SAVE IT IN TEMP BUFFER
01300	05D2	CE	062F	LDX	#TARECN	PUT NEW INT VEC IN
01310	05D5	FF	FFF8	STX	IRQVEC	
01320	05D8	86	AA	LDA A	#\$AA	SET DONTST FLAG
01330	05DA	B7	062A	STA A	DONTST	
01340	05DD	CE	0000	LDX	#0	
01350	05E0	FF	E500	STX	DACZRO	SEND ENABLE PLS TO INT CIRCUIT
01360				*		
01370	05E3	0E		TAINTR CLI		ENABLE CPU FOR RECONS INTR'S
01380	05E4	E6	062A	LDA A	DONTST	CHECK IF DECPRS DONE
01390	05E7	81	00	CMP A	#0	DONE YET?
01400	05E9	26	F8	BNE	TAINTR	NO. KEEP LOOPING
01410	05EB	0F		SEI		INSURE INTR DISABLE
01420	05EC	CE	0878	LDX	#GOAGIN	PRINT GO AGAIN MSG
01430	05EF	BD	CA8F	JSR	OUTNCR	"ENTER CONTINUATION COMMAND"
01440	05F2	BD	CA2C	JSR	KEYBD0	
01450	05F5	FE	3002	LDX	ENDBUF	NOW GET COMMAND
01460	05F8	09		DEX		
01470	05F9	E6	00	LDA B	0,X	
01480	05FB	17		TBA		PUT IT IN A
01490	05FC	90	28	SUB A	\$28	TRANSFORM X TO 0
01500	05FE	97	20	STA A	CHNLBF	PUT IN CHNL FLAG
01510	0600	F7	0846	STA B	CHNASC	PUT ASCII IN ERROR MSG
01520	0603	C1	58	CMP B	#'X	WAS IT .LT. 'X?
01530	0605	2D	07	BLT	TADEC5	YES. EXIT ROUTINE
01540	0607	C1	5A	CMP B	#'Z	WAS CMD .GT. 'Z ?
01550	0609	2E	03	BGT	TADEC5	YES. EXIT ROUTINE
01560	060B	7E	05A0	JMP	TADEC4	NO. GO EXECUTE DECPRS AGAIN
01570	060E	7E	0100	TADEC5 JMP	DISPLA	RETURN TO DISPLAY
01580				*		
01590				*****		
01600				*	END DECPRS INIT ROUTINE	
01610				*****		
01620				*		
01630				*		
01640				*	PARAMETER BUFFER	

01650		*			
01660	0611	0002	RUFPTR	RMB	2
01670	0612	0001	FILEPTR	RMB	1
01680	0614	0001	XDATA	RMB	1
01690	0615	0001	YDATA	RMB	1
01700	0616	0001	ZDATA	RMB	1
01710	0617	0001	XSLOPE	RMB	1
01720	0618	0001	YSLOPE	RMB	1
01730	0619	0001	ZSLOPE	RMB	1
01740	061A	0001	DELT	RMB	1
01750	061B	0001	VALUE	RMB	1
01760	061C	0001	CNT	RMB	1
01770	061D	0002	VECSAV	RMB	2
01780	061F	0001	BITVLU	RMB	1
01790	0620	0001	SIGN	RMB	1
01800	0621	0001	XACCEL	RMB	1
01810	0622	0001	XTEMP	RMB	1
01820	0623	0001	YACCEL	RMB	1
01830	0624	0001	YTEMP	RMB	1
01840	0625	0001	ZACCEL	RMB	1
01850	0626	0001	ZTEMP	RMB	1
01860	0627	0002	WORKBF	RMB	2
01870	0629	0001	TARFLG	RMB	1
01880	062A	0001	DCNTST	RMB	1
01890	062B	0002	DABUF	RMB	2
01900	062D	0002	STKSAV	RMB	2
01910			*		
01920			*FUNCTION :	TARECN	
01930			*INPUTS :	BUFFER VALUES	
01940			*OUTPUTS :	D/A VALUES VIA ST-6800	
01950			*CALLS :	GETVLU,GETTIM,OUTDA	
01960			*DESTROYS	A,B,X,CC	
01970			*PURPOSE :	THIS IS THE INTERRUPT SERVICE ROUTINE WHI	
01980			* PERFORMS	RECONSTRUCTION ON THE TOLAN-A ENCODED DAT	
01990			* FILE IN	MEMORY. THE DATA IS DECODED AND OUTPUT ON	
02000			* D/A (ST6800)	CH 0.	
02010			*		
02020			*****		
02030			*		
02040	062F	B7 E400	TARECN	STA A	ADCZRO
02050	0632	01		NOP	
02060	0633	B7 062D		STS	STKSAV
02070	0636	B6 0846		LDA A	CHNASC
02080	0639	81 58		CMP A	#'X
02090	063B	26 08		BNE	TAREC1
02100	063D	B6 0614		LDA A	XDATA
02110	0640	FD 0796		JSR	OUTDA
02120	0643	20 12		BRA	TARE22
02130	0645	81 59	TAREC1	CMP A	#'Y
02140	0647	26 08		BNE	TAREC2
02150	0649	B6 0615		LDA A	YDATA
02160	064C	BD 0796		JSR	OUTDA
02170	064F	20 06		BRA	TARE22
02180	0651	B6 0616	TAREC2	LDA A	ZDATA
					GET CUR ZDATA VALUE.

02190 0654 BD 0756	JSR	OUTDA	OUTPUT Z VALUE VIA D/A CH 0
02200 0657 B6 0629 TARE22	LDA A	TARFLG	WAS LAST DATA PT EXTROPOLATED
02210 065A 81 00	CLR A		
02220 065C 26 33	BEQ	TAREC3	YES. DON'T GET NEXT ACCEL. M
02230 065E 86 58	LDA A	#'X	PUT ASCII X IN GETVLU ERROR M
02240 0660 B7 0A3D	STA A	OVFLCH	
02250 0663 BD 0700	JSR	GETVLU	NOW GO GET NEXT XACCEL VALUE
02260 0666 B6 061B	LDA A	VALUE	AND PUT XACCEL
02270 0669 B7 0621	STA A	XACCEL	
02280 066C B7 0622	STA A	XTEMP	SAVE XACCEL IN TEMP BUFF
02290 066F 86 59	LDA A	#'Y	PUT 'Y IN GETVLU ERROR MSG
02300 0671 B7 0A3D	STA A	OVFLCH	
02310 0674 BD 0700	JSR	GETVLU	GET YACCEL DTA
02320 0677 B6 061B	LDA A	VALUE	
02330 067A B7 0623	STA A	YACCEL	PUT IN YACCEL VAR
02340 067D B7 0624	STA A	YTEMP	PUT IN YACCEL TEMP VAR
02350 0680 86 5A	LDA A	#'Z	PUT 'Z IN GETVLU ERROR MSG
02360 0682 B7 0A3D	STA A	OVFLCH	
02370 0685 BD 0700	JSR	GETVLU	GET ZACCEL DATA
02380 0688 B7 0625	STA A	ZACCEL	PUT IN ZACCEL VAR
02390 068B B7 0626	STA A	ZTEMP	PUT IN ZACCEL TEMP VAR
02400 068E BD 075E	JSR	GETTIM	NOW GET TIME COPPRESSION RUN L
02410 0691 7F 0629 TAREC3	CLR	TARFLG	RESET FLAGS
02420 0694 7A 061A	DEC	DELT	NOW CK IF RUN LEN ONLY 1
02430 0697 27 09	BEQ	TAREC4	YES. ACCEL VAR NOT ZERO
02440 0699 7F 0621	CLR	XACCEL	ACCEL WAS ZRO SO CLR X-Y-ZACC
02450 069C 7F 0623	CLR	YACCEL	
02460 069F 7F 0625	CLR	ZACCEL	
02470 06A2 B6 0617 TAREC4	LDA A	XSLOPE	NOW GET LAST SLOPE
02480 06A5 BB 0621	ADD A	XACCEL	& CALC NEW SLOPE WITH ACCEL
02490 06A8 B7 0617	STA A	XSLOPE	
02500 06AB B6 0618	LDA A	YSLOPE	
02510 06AE BB 0623	ADD A	YACCEL	
02520 06B1 B7 0618	STA A	YSLOPE	
02530 06B4 B6 0619	LDA A	ZSLOPE	
02540 06B7 BB 0625	ADD A	ZACCEL	
02550 06BA B7 0619	STA A	ZSLOPE	
02560 06BD B6 0614	LDA A	XDTA	NOW CALC DATA VALUE USING SLO
02570 06C0 BB 0617	ADD A	XSLOPE	
02580 06C3 B7 0614	STA A	XDTA	
02590 06C6 B6 0615	LDA A	YDTA	
02600 06C9 BB 0618	ADD A	YSLOPE	
02610 06CC B7 0615	STA A	YDTA	
02620 06CF B6 0616	LDA A	ZDTA	
02650 06D8 B6 061A	LDA A	DELT	WAS RUN LEN .GT. 2 PPS ?
02660 06DB 27 22	BEQ	TARRTI	IF YES, DONT GET NEW RUN FIL
02670 06DD 81 01	CMP A	#1	JUST EXTROPOLATE POINTS UNTIL
02680 06DF 27 07	BEQ	TAREC5	DELTA IS REDUCED TO 1
02690 06E1 86 37	LDA A	#55	SET TARFL1
02700 06E3 B7 0629	STA A	TARFLG	
02710 06E6 20 17	BRA	TARRTI	
02720 06E8 B6 0622 TAREC5	LDA A	XTEMP	RETRIEVE ACCEL VLUS PREVIOUSL
02730 06EB B7 0621	STA A	XACCEL	
02740 06EE B6 0624	LDA A	YTEMP	


```

02750 06F1 B7 0623      STA A  YACCEL
02760 06F4 B6 0626      LDA A  ZTEMP
02770 06F7 B7 0625      STA A  ZACCEL
02780 06FA 9C 55        LDA A  #855      NOW SET TARFL2 FLAG
02790 06FC B7 0629      STA A  TARFLG
02800 06FF 3B          TARRTI RTI
02810                  *
02820                  *****
02830                  *          END OF INTERRUPT HANDLER
02840                  *****
02850                  *
02860                  *          SUBROUTINES
02870                  *
02880                  *
02890                  *FUNCTION : GETVLU
02900                  *INPUTS : BUFPTR,BITPTR,SIGN BUFFERS
02910                  *OUTPUTS :VALUE BUFFER
02920                  *CALLS : GETBIT
02930                  *DESTROYS : A,X,CC
02940                  *PURPOSE : THIS ROUTINE SCANS MEMORY STARTING AT POI
02950                  * DEFINED BY BUFPTR,BITPTR AND GETS ACCEL VALUE. BUF
02960                  * AND BITPTR ARE UPDATED.
02970                  *
02980                  *
02990 07C0 7F 061B GETVLU CLR    VALUE    INSURE VALUE IS CLEAR ON STAR
03000 0703 ED 07ED      JSR      GETBIT   GET DELIMITER BIT
03010 0706 E6 061F      LDA A  BITVLU
03020 0709 81 00        CMP A  #0        IS DELIM BIT =0 ?
03030 070B 27 06        BEQ      GETVLI   YES. GET NEXT BIT
03040 070D CE 097D      LDX      #SYNCHS  NO. SYNC ERROR HAS OCCURED. E
03050 0710 7E 07F5      JMP      ERROR
03060 0713 7F 0620 GETVLI CLR    SIGN     INSURE SIGN CLR
03070 0716 ED 07ED      JSR      GETBIT   GET SIGN BIT
03080 0719 E6 061F      LDA A  BITVLU
03090 071C 81 00        CMP A  #0        IS ACCEL POS?
03100 071E 27 03        BEQ      GETVL2   YES. LEAVE SIGN CLR
03110 0720 73 0620      CCM      SIGN     NO. SET SIGN FLAG
03120 0723 BD 07ED GETVL2 JSR      GETBIT   GET FIRST CODE BIT
03130 0726 F6 061F      LDA A  BITVLU
03140 0729 26 11        RNE      GETVL4
03150 072B FE 0611 GETVL3 LDX      BUFPTR   END OF VALUE DECODE. PACK UP
03160 072E 78 0613      ASL      BITPTR   FOR CORRECT POSITIONING FOR N
03170 0731 24 2A        BCC      GETVL6   VALUE FETCH
03180 0733 79 0613      ROL      BITPTR
03190 0736 09          DEX
03200 0737 FF 0611      STX      BUFPTR
03210 073A 20 21        BRA      GETVL6
03220 073C 7C 051B GETVL4 INC      VALUE   NOW PAST DELIM & SIGN. NOW CO
03230 073F 2A 06        BPL      GETVL5   BITS. IS CODE WORD > 127 ?
03240 0741 CE 09FB      LDX      #OVFLMS YES. OVERFLOW HAS OCCURED. PR
03250 0744 7E 07F5      JMP      ERROR
03260 0747 ED 07ED GETVL5 JSR      GETBIT   NOW GET DATA VALUE CODE BITS
03270 074A B6 061F      LDA A  BITVLU
03280 074D 81 00        CMP A  #0        DELIM DETECTED YET?

```

```

03290 074F 26 EB      BNE      GETVL4    NO. KEEP FETCHING BITS
03300 0751 B6 0620    LDA A  SIGN      YES. GET SIGN AND ADJUST CNT
03310 0753 81 00      CLZ A  #0
03320 0756 27 00      LDR      GETVL3    SIGN 0. VALUE COR. , CO ADJ PT
03330 0756 70 061B    NEG      VALUE     SIGN 1. TAKE 2'S COMP, ADJ PT
03340 075B 20 CE      BRA      GETVL3
03350 075D 39      GETVL6 RTS
03360      *
03370      *****
03380      *          END OF SUBROUTINE GETVL6
03390      *****
03400      *
03410      *
03420      *FUNCTION : GETTIM
03430      *INPUTS : BITVL6 (CURR BIT VALUE PT'ED TO)
03440      *OUTPUTS : DELT (CURRENT TIME COUNT)
03450      *CALLS : GETBIT, ERROR
03460      *DESTROYS : A,X,CC REGISTERS
03470      *PURPOSE : THIS ROUTINE, WHEN CALLED AFTER 3 GETVL6
03480      * RETURNS THE VALUE OF THE RUN LENGTH TIME COUNTER I
03490      * VARIABLE DELT.
03500      *
03510      *
03520 075E BD 07ED GETTIM JSR      GETBIT    GET DELIM BIT
03530 0761 B6 061F      LDA A  BITVL6
03540 0764 81 00      CMP A  #0          IS BIT 0?
03550 0766 27 06      BEQ      GETTI1     YES. PROPER DELIM. CONTINUE
03560 0768 CE 09AA      LDX      #SYNCTI  NO. SYNC ERROR. PRINT & ERR 0
03570 076B 7E 07F5      JMP      ERROR
03580 076E 7F 061A GETTI1 CLR      DELT     INSURE TIME CNT INITIALLY CLR
03590 0771 86 07      LDA A  #7          INIT CNT FOR 7 BIT TIME
03600 0773 B7 061C      STA A  CNT
03610 0776 BD 07ED GETTI2 JSR      GETBIT    GET DELTA T BIT
03620 0779 B6 061F      LDA A  BITVL6
03630 077C BA 061A      ORA A  DELT     SET LSB BIT ACCORD TO T BIT
03640 077F B7 061A      STA A  DELT
03650 0782 78 061A      ASL      DELT     SHIFT DELT FOR NEXT BIT
03660 0785 7A 061C      DEC      CNT     DONE WITH TIME YET?
03670 0788 26 EC      BNE      GETTI2     NO. KEEP FETCHING TIME BITS
03680 078A B6 061A      LDA A  DELT     YES. GET TIME VALUE
03690 078D 26 06      BNE      GETTI3     IS TIME CNT 0 ?
03700 078F CE 0A3F      LDX      #TIMERR  YES. ERROR. PRNT ERR MSG
03710 0792 7E 07F5      JMP      ERROR    AND ERROR OFF
03720 0795 39      GETTI3 RTS
03730      *
03740      *****
03750      *          END GETTIM SUBROUTINE
03760      *****
03770      *
03780      *
03790      *FUNCTION : OUTDA
03800      *INPUTS : DATA VALUE TO BE D/A'ED IN ACC A
03810      *OUTPUTS : VALUE IN ACC A VIA D/A CH 0
03820      *CALLS : NOTHING

```

```

03830      *DESTROYS : B,CC
03840      *PURPOSE : THIS ROUTINE OUTPUTS THE DATA IN
03850      * ACCOUNT FOR A TO THE D/A FOR CONVERSION AND
03860      * TRANSMISSION. THE 8 BIT VALUE IS SHIFTED
03870      * TO FORM THE 12 BIT OPERAND REQUIRED BY THE
03880      * D/A (ST6800).
03890      *
03900      *
03910 0796 36      OUTDA  PSH A          SAVE VLU IN A ACC
03920 0797 B7 062C      STA A  DABUF+1  PUT VALUE IN D/A OUT BUFF
03930 079A 7F 062B      CLR   DABUF     CLR MSB BYTE OF D/A OUT BUFF
03940 079D 78 062C      ASL   DABUF+1   LEFT SHIFT 2 BYTE D/A WORD
03950 07A0 79 062B      ROL   DABUF     FOR CONVERSION OF 8 BIT
03960 07A3 78 062C      ASL   DABUF+1   VALUE TO 12 BIT
03970 07A6 79 062B      ROL   DABUF
03980 07A9 78 062C      ASL   DABUF+1
03990 07AC 79 062B      ROL   DABUF
04000 07AF 78 062C      ASL   DABUF+1
04010 07B2 79 062B      ROL   DABUF
04020 07B5 FE 062B      LDX   DABUF     NOW LD INDEX WITH 12 BIT VLU
04030 07B8 FF E500      STX   DACZRO   AND OUTPUT TO D/A
04040 07BB 32          PUL A          RETRIEV ENTY A ACC VLU
04050 07BC 39          RTS
04060      *
04070      *****
04080      *          END SUBROUTINE OUTDA
04090      *****
04100      *
04110      *
04120      *FUNCTION : GETBIT
04130      *INPUTS : BUFPTR,BITPTR POINTERS
04140      *OUTPUTS : BITVLU, DONIST FLAGS
04150      *CALLS :NOTHING
04160      *DESTROYS :A,CC
04170      *PURPOSE : THIS ROUTINE CHECKS THE STATE OF THE BIT
04180      * POINTED TO BY THE BUFPTR,BITPTR PAIR AND SETS THE
04190      * THE BITVLU FLAG ACCORDINGLY. THIS ROUTINE CHECK
04200      * FOR END OF MEMORY AND FORCES EXIT OF RECONSTRUCTIO
04210      * ROUTINES IF EQ1 IS DETECTED.
04220      *
04230      *
04240 07BD 7F 061F  GETBIT CLR   BITVLU  INSURE BITVLU STARTS RESET
04250 07C0 FE 0611      LDX   BUFPTR  GET MEM WORD POINTER
04260 07C3 F6 0613      LDA B  BITPTR  GET BIT POINTER IN THAT MEM W
04270 07C6 E5 00      BIT B  0,X      CHECK IF BIT SET
04280 07C8 27 03      BEQ   GETBIT1  NO IT WAS RESET. RESET BITVLU
04290 07CA 7C 061F      INC   BITVLU  YES. BIT WAS SET. SET BITVLU
04300 07CD 74 0613  GETBIT1 LSR   BITPTR  NOW UPDATE POINTERS
04310 07D0 24 22      BCC   GETRIS
04320 07D2 76 0613      ROR   BITPTR
04330 07D5 02      INX
04340 07D6 FF 0611      STX   BUFPTR
04350 07D9 F6 0611      LDA B  BUFPTR  CHECK FOR END OF MEMORY
04360 07DC C1 80      CMP B  #$80

```

```

04370 07DE 26 14      BNE    GETRTS    NO. CONTINUE
04380 07E0 0F          SEI                      NOW TURN OFF INTR AND RET
04390 07E1 0F 4000      LDX    #$4000
04400 07E4 FF E500      STX    DACZRO
04410 07E7 FE 061D      LDX    VECSAV
04420 07EA FF FFF8      STX    IRQVEC
04430 07ED 7F 062A      CLR    DONTST    YES. CLR DONTST FLAG AND RTI
04440 07F0 BE 062D      LDS    STKSAV    RETRIEVE INTR STACK PTR
04450 07F3 3B          RTI
04460 07F4 39          GETRTS RTS
04470                  *
04480                  *****
04490                  *                END OF SUBROUTINE GETBIT
04500                  *****
04510                  *
04520                  *
04530                  *FUNCTION : ERROR
04540                  *INPUTS : MESSAGE STRING POINTED TO BY INDX RFG
04550                  *OUTPUTS : ERROR MESSAGE TO CONSOLE DEVICE
04560                  *CALLS : OUTNCR,KEYBD0,DISPLA
04570                  *DESTROYS : A,B,X,CC
04580                  *PURPOSE : THIS ROUTINE PRINTS ERROR MESSAGES UPON
04590                  * ERROR DETECTION BY THE DECOMPRESSION ROUTINES.
04600                  *
04610                  *
04620 07F5 BD CA8F      ERROR JSR    OUTNCR
04630 07F8 CE 09DF      LDX    #ANYKEY    GET ANYKEY MESSAGE
04640 07FB BD CA8F      JSR    OUTNCR    PRINT "TO CONTINUE PRESS ANY
04650 07FE BD CA2C      JSR    KEYBD0
04660 0801 0F          SEI                      NOW CLR INTR CONDITION
04670 0802 CE 4000      LDX    #$4000    AND RETURN TO DISPLAY
04680 0805 FF E500      STX    DACZRO
04690 0808 FE 061D      LDX    VECSAV    RETRIEV IRQ VEC ADDR
04700 080B FF FFF8      STX    IRQVEC
04710 080E 7E 0100      JMP    DISPLA    RETURN TO "DISPLAY"
04720                  *
04730                  *****
04740                  *                END ERROR HANDLING ROUTINE
04750                  *****
04760                  *
04770                  *                MESSAGE STRINGS
04780                  *
04790 0811 1A0D      GOMSG FDB    $1A0D,$0A07
04800 0815 44          FCC    /DECOMPRESSION OVERLAY/
04810 082A 0D0A      FDB    $0D0A,$0D0A
04820 082E 43          FCC    /CURRENT CHANNEL FLAG IS /
04830 0846 0001      CHNASC RMB    1
04840 0847 0D0A      FDB    $0D0A,$0D0A
04850 084B 44          FCC    /DO YOU WISH TO EXECUTE THIS MODULE (
04860 0877 04          FCB    4
04870 0878 0D0A      GOAGIN FDB    $0D0A,$0A07
04880 087C 44          FCC    /DECOMPRESSION AND DISPLAY COMPLETE/
04890 089E 0D0A      FDB    $0D0A,$0D0A
04900 08A2 45          FCC    /ENTER CONTINUATION COMMAND/

```

```

04910 084C 0D0A      FDB      $0D0A
04920 08BE 20        FCC      / X=DISPLAY CHANNEL X ON D-A CHANNE
04930 08C5 0D0A      FDB      $0D0A
04940 08E7 20        FCC      / Y=DISPLAY CHANNEL Y OF D-A CHANNE
04950 090E 0D0A      FDB      $0D0A
04960 0910 20        FCC      / Z=DISPLAY CHANNEL Z ON D-A CHANNE
04970 0937 0D0A      FDB      $0D0A,$0D0A
04980 093B 41        FCC      /ANY OTHER KEY RETURNS CONTROL TO "DI
04990 0966 0D0A      FDB      $0D0A,$0D0A
05000 096A 45        FCC      /ENTER COMMAND NOW=/
05010 097C 04        FCB      4
05020 097D 1A07      SYNCMS FDB      $1A07,$0D0A
05030 0981 53        FCC      /SYNC ERROR DETECTED BY TOLAN-A DECOD
05040 09A7 0D0A      FDB      $0D0A
05050 09A9 04        FCB      4
05060 09AA 1A07      SYNCTI FDB      $1A07,$0D0A
05070 09AE 53        FCC      /SYNC ERROR DETECTED BY /
05080 09C5 54        FCC      /TOLAN-A DURING TIME FETCH/
05090 09DE 04        FCB      4
05100 09DF 0D0A      ANYKEY FDB      $0D0A
05110 09E1 54        FCC      /TO CONTINUE PRESS ANY KEY/
05120 09FA 04        FCB      4
05130 09FB 1A07      OVFLMS FDB      $1A07,$0D0A
05140 09FF 41        FCC      /ACCEL VALUE OVERFLOW DETECTED BY TOL
05150 0A2F 0D0A      FDB      $0D0A
05160 0A31 20        FCC      / IN CHANNEL /
05170 0A3D 0001      OVFLCH FCB      1
05180 0A3E 04        FCB      4
05190 0A3F 1A07      TIMERR FDB      $1A07,$0D0A
05200 0A43 54        FCC      /TIME COUNT ERROR DETECTED BY TOLAN-A
05210 0A6F 0D0A      FDB      $0D0A
05220 0A71 54        FCC      /TIME CNT (DELT)=0/
05230 0A82 04        FCB      4
05240                *
05250                *****
05260                *
05270                *      END OF DECPRS OVERLAY ROUTINES
05280                *
05290                *****
05300                *
05310                END

```

ENTROPY

```

0100 REM *****
0110 REM *
0120 REM *      EKG ENTROPY CALCULATION PROGRAM
0130 REM *
0140 REM *****
0150 REM *
0160 REM *      THIS PROGRAM READS THE FREQUENCY
0170 REM * OF OCCURENCE DATA IN MEMORY FROM AN
0180 REM * EKG DATA COLLECTION AND CALCULATES
0190 REM * THE ENTROPY OF THE X,Y,AND Z DATA
0200 REM * SOURCES.  THE CALCULATED ENTROPY IS
0210 REM * STORED BACK IN THE MEMORY FILE HEADER
0220 REM * IN ASCII. THIS HEADER CAN THEN BE
0230 REM * INSERTED INTO THE DISK FILE USING
0240 REM * MINIDOS.
0250 REM *
0260 REM *****
0270 REM
0280 DIM D(255),N$(8),T(255),M9(6),T$(3)
0290 STRING= 8
0291 DIGITS= 4
0292 T$(1)="X"
0293 T$(2)="Y"
0294 T$(3)="Z"
0295 LINE= 80
0296 L7=6
0297 E1=13507
0313 REM
0325 REM *****
0326 REM * NOW GET FILENAME FROM MEMORY & PRINT IT
0327 REM *****
0328 REM
0330 FOR I=1 TO 8
0340 J=13313+I
0350 N$(I)=CHR$(PEEK(J))
0360 A$=A$+N$(I)
0370 NEXT I
0372 P$=CHR$(12)
0373 PRINT P$
0374 PRINT
0375 PRINT
0376 PRINT
0377 PRINT
0378 PRINT
0380 PRINT
0385 PRINT "EKG ENTROPY CALCULATION"
0386 PRINT
0387 PRINT
0390 PRINT "FILENAME.....";A$
0400 PRINT
0401 PRINT
0402 GOSUB 1696
0403 GOSUB 1900
0404 IF G$ = "S" THEN GO TO 1120

```

ENTROPY

```

0410 REM
0411 REM *****
0412 REM * E1=MEM FILE DESTINATION FOR ASCII RESULTS
0413 REM * 13856=$3620 (HEX)=FLAG MEM USED BY BASIC. DATA REPLACED
0414 REM * 13567=34FF (HEX)=XPDF BUFFER-1
0415 REM * 14079=36FF (HEX)=YPDF BUFFER-1
0416 REM * 14591=38FF (HEX)=ZPDF BUFFER-1
0417 REM * 13500=34BC (HEX)=BASSAV. MEM BUF FOR $3620 DATA SAVE
0419 REM *****
0420 REM
0421 REM
0424 REM
0425 REM *****
0426 REM * NOW GET VLU IN BASSAV & STR TO $3620 (HEX)
0427 REM *****
0428 REM
0429 B=13856
0430 C1=1
0432 N=PEEK(13500)
0433 POKE( B,N)
0435 FOR I=1 TO 255
0436 T(I)=0
0437 NEXT I
0438 T9=0
0439 S2=0
0440 S=0
0441 REM
0442 REM *****
0443 REM * NOW BEGIN LOOP TO CALC X,Y,Z & TOTAL ENTROPY
0444 REM *****
0445 REM
0450 IF C1 <> 1 THEN GO TO 480
0460 K=13567
0480 IF C1 <> 2 THEN GO TO 510
0490 K=14079
0510 IF C1 <> 3 THEN GO TO 540
0520 K=14591
0540 IF C1>= 4 THEN GO TO 800
0542 REM
0543 REM *****
0544 REM * GET 2 BYTE PDF DATA & MERGE INTO FLTING POINT NUMBER
0545 REM *****
0546 REM
0550 FOR I=1 TO 509 STEP 2
0560 M=(I+1)/2
0570 D(M)=256*PEEK(I+K)+PEEK(I+K+1)
0572 T(M)=T(M)+D(M)
0590 S=D(M)+S
0592 S2=D(M)+S2
0600 NEXT I
0605 D9=256*PEEK(K+511)+PEEK(K+512)
0606 S=S+D9
0607 T9=T9+D9
0608 S2=S2+D9

```

ENTROPY

```

0610 E=0
0612 REM
0613 REM *****
0614 REM * NOW CALCULATE PROBABILITY OF OCCURRENCE OF SPECIFIC
0615 REM * VLU'S FROM NUM OF OCCURENCES. THEN CALC ENTROPY SUM.
0616 REM *****
0617 REM
0620 FOR I=1 TO 255
0630 D(I)=D(I)/S
0640 IF D(I)=0 THEN GO TO 660
0650 E=E+(D(I)*(-(LOG(D(I))/.693147)))
0660 NEXT I
0663 D9=D9/S
0664 IF D9=0 THEN GO TO 666
0665 E=E+(D9*(-(LOG(D9)/.693147)))
0666 IF E=0 THEN GOTO 672
0667 C2=8/E
0668 GOTO 683
0672 C2=1
0673 GOTO 680
0674 REM
0675 REM *****
0676 REM * PRINT RESULTS OF CALC TO TERMINAL (PRINTER).
0677 REM *****
0678 REM
0680 PRINT
0681 PRINT "MAXIMUM COMPRESSION FOR CHANNEL ";T$(C1);" IS INFINITE"
0682 GOTO 690
0683 PRINT
0684 PRINT "MAX COMPRESSION RATIO FOR CHANNEL ";T$(C1);" IS ";C2;" ; 1"
0690 PRINT
0700 PRINT "EKG LEAD ";T$(C1);" ENTROPY = ";E;" BITS."
0710 PRINT
0712 REM
0713 REM *****
0714 REM * NOW CONVERT CALCULATED ENTROPY TO ASCII & STORE BACK
0715 REM * IN MEMORY BUFFER
0716 REM *****
0717 REM
0720 LET E$=STR$(E)
0740 GOSUB 1750
0785 C1=C1+1
0790 GOTO 440
0800 REM
0801 REM *****
0802 REM * WITH X,Y,Z ENTROPY CALCULATED, NOW CALCULATE TOTAL,
0803 REM * COMBINED ENTROPY BY ADDING X,Y,Z BIN COUNTS & DIVIDING
0804 REM * BY TOTAL SAMPLE COUNT
0805 REM *****
0806 REM
0807 E5=0
0810 FOR I=1 TO 255
0820 T(I)=T(I)/S2
0822 IF T(I)=0 THEN GO TO 840

```


ENTROPY

```

0830 E5=E5+(T(I)*(-(LOG(T(I))/.693147)))
0840 NEXT I
0850 T9=T9/S2
0852 IF T9=0 THEN GO TO 870
0860 E5=E5+(T9*(-(LOG(T9)/.693147)))
0870 C3=8/E5
0872 REM
0873 REM *****
0874 REM * PRINT COMBINED ENTROPY TO TERM & STR RESULT IN MEM BUFFER
0875 REM *****
0876 REM
0880 PRINT
0890 PRINT "MAX COMPRESSION RATIO FOR 3 LEAD EKG SYSTEM ";C3;" : 1"
0910 PRINT
0920 PRINT "3 LEAD EKG SOURCE ENTROPY = ";E5;" BITS."
0930 PRINT
0940 LET E$=STR$(E5)
0950 GOSUB 1750
1000 LET E$=STR$(C3)
1005 GOSUB 1750
1010 REM
1020 REM *****
1030 REM * NOW GET OTHER STATISTICAL VARIABLES &
1040 REM * CALCULATE CHANNEL MAX,MINS,COMPRESSION RATIO
1050 REM * OBTAINED, AND COMPRESSION TIME EFFICIENCY
1060 REM *****
1070 REM
1080 REM *****
1082 REM * START WITH ACHIEVED COMPRESSION RATIO. THEN
1090 REM * PRINT THE TIME COMPRESSION EFFICIENCY. THIS IS
1091 REM * THE PERCENTAGE OF THE TIME AVAILABLE THAT WAS USED
1092 REM * TO COLLECT,COMPRESS, & CALCULATE STAT VARIABLES.
1100 REM *****
1110 REM
1120 DIGITS= 4
1130 L7=6
1140 M4=PEEK(13483)+PEEK(13482)*256+PEEK(13481)*256*256
1150 D4=PEEK(13487)+PEEK(13486)*256+PEEK(13485)*256*256
1160 D4=PEEK(13484)*256*256*256+D4
1170 C4=D4/M4
1180 LET E$=STR$(C4)
1190 GOSUB 1750
1200 PRINT
1201 PRINT "COMPRESSION RATIO ACHIEVED = ";C4;" : 1"
1210 IF G$="S" THEN GOTO 1260
1211 REM
1212 REM *****
1213 REM * CALCULATE OVERALL COMPRESSION EFFICIENCY
1214 REM *****
1215 REM
1216 DIGITS= 1
1217 C6=((C4-1)/(C3-1))*100
1218 PRINT
1219 PRINT "ACHIEVED COMPRESSION EFFICIENCY = ";C6;" %"

```

ENTROPY

```

1220 LET E$=STR$(C6)
1221 L7=4
1222 GOSUB 1750
1229 REM
1230 REM *****
1240 REM * NOW CALCULATE COMPRESSION TIME EFFECIENCY
1250 REM *****
1252 REM
1260 S4=PEEK(13461)+PEEK(13460)*256
1262 L4=PEEK(13462)
1264 L5=PEEK(13459)+PEEK(13458)*256+PEEK(13457)*256*256
1266 L5=PEEK(13456)*256*256*256+L5
1270 T4=(1-(L5/(L4*S4)))*100
1272 DIGITS= 1
1274 L7=4
1276 LET E$=STR$(T4)
1278 GOSUB 1750
1280 PRINT
1282 PRINT "COLLECTION TIME EFFICIENCY = ";T4;" %"
1284 PRINT
1290 REM
1291 REM *****
1292 REM * CALCULATE COLLECTION DURATION
1293 REM *****
1294 REM
1295 T6=S4/S8
1296 PRINT
1297 PRINT "COLLECTION DURATION = ";T6;" SECONDS"
1298 PRINT "AT A SAMPLE RATE OF ";S8;" HZ"
1299 LET E$=STR$(T6)
1300 GOSUB 1750
1301 REM
1302 REM *****
1310 REM * RETRIEVE DATA MAX & MINS & CALC VOLTS
1320 REM *****
1330 REM
1339 DIGITS= 5
1340 FOR I=6 TO 1 STEP -1
1350 M9(I)=PEEK(13460+I*3)
1360 IF M9(I) > 127 THEN GOTO 1390
1370 M9(I)=M9(I)*.0390625
1380 GOTO 1400
1390 M9(I)=- (10-M9(I)*.0390625)
1400 NEXT I
1401 E1=13408
1402 FOR I4=5 TO 1 STEP -2
1403 LET E$=STR$(M9(I4))
1404 L7=7
1405 GOSUB 1750
1406 I3=I4+1
1407 LET E$=STR$(M9(I3))
1408 GOSUB 1750
1409 NEXT I4
1410 REM

```

ENTROPY

```

1420 REM *****
1430 REM * NOW PRINT MAX & MINS TO TERM (PRINTER)
1440 REM *****
1450 REM
1460 PRINT
1470 FOR I=1 TO 6 STEP 2
1480 PRINT
1490 PRINT "CHANNEL ";T$(I+1)/2;" MINIMUM = ";M9(7-I);" VOLTS."
1500 PRINT
1510 PRINT "CHANNEL ";T$(I+1)/2;" MAXIMUM = ";M9(6-I);" VOLTS."
1520 PRINT
1530 NEXT I
1540 REM
1550 REM *****
1560 REM * FINISH WITH A PRINT OF NUM OF BITS USED TO
1570 REM * STORE X,Y,Z IN MEM WITH THIS COMPRESSION TYPE
1580 REM *****
1590 REM
1600 X3=PEEK(13490)+PEEK(13489)*256+PEEK(13488)*256*256
1610 Y3=PEEK(13493)+PEEK(13492)*256+PEEK(13491)*256*256
1620 Z3=PEEK(13496)+PEEK(13495)*256+PEEK(13494)*256*256
1630 T3=PEEK(13499)+PEEK(13498)*256+PEEK(13497)*256*256
1632 DIGITS= 0
1640 PRINT
1650 PRINT "NUMBER OF BITS USED TO STORE CHANNEL X =";X3
1652 PRINT
1660 PRINT "NUMBER OF BITS USED TO STORE CHANNEL Y =";Y3
1662 PRINT
1670 PRINT "NUMBER OF BITS USED TO STORE CHANNEL Z =";Z3
1672 PRINT
1680 PRINT "NUMBER OF BITS USED TO STORE TIME =";T3
1682 PRINT
1683 PRINT
1690 GOTO 2000
1691 REM
1692 REM *****
1693 REM * SUB TO PRINT OUT COMPRESSION TYPE
1694 REM *****
1695 REM
1696 H=PEEK(13312)*256+PEEK(13313)
1697 IF H=20035 THEN HS(1)="NOT COMP"
1698 IF H=20035 THEN HS(2)="RESSESSED "
1699 IF H=21569 THEN HS(1)="TOLAN A "
1700 IF H=21569 THEN HS(2)=" "
1701 IF H=21570 THEN HS(1)="TOLAN B "
1702 IF H=21570 THEN HS(2)=" "
1703 IF H=17487 THEN HS(1)="DOWER "
1704 IF H=17487 THEN HS(2)=" "
1705 IF H=21584 THEN HS(1)="TURNING "
1706 IF H=21584 THEN HS(2)="POINT "
1707 IF H=18766 THEN HS(1)="2ND ORD "
1708 IF H=18766 THEN HS(2)="INTERPOL"
1709 PRINT "COMPRESSION USED.....";HS(1);HS(2)
1710 PRINT

```

ENTROPY

```

1711 RETURN
1720 REM
1730 *****
1732 REM * SUBROUTINE TO WRITE ASCII DATA IN ES TO MEM
1733 REM * LOCATION POINTED TO BY E1. WILL WRITE 8 CHAR
1734 REM * STRINGS TO MEMORY.
1739 REM *****
1740 REM
1750 FOR I=1 TO L7
1760 F$=MID$(E$,I,L7)
1770 F=ASC(F$)
1780 POKE( E1,F)
1781 E1=E1+1
1790 NEXT I
1800 POKE( E1,4)
1801 E1=E1+1
1810 RETURN
1820 REM
1830 REM *****
1840 REM * SUBROUTINE TO PROMPT COMMAND FOR SHORT RUN OR LONG
1850 REM * RUN. SHORT DOES NOT CALCULATE ENTROPY OR COMPRESSION
1860 REM * RATIOS.
1870 REM *****
1880 REM
1900 PRINT
1910 PRINT "ENTER S FOR SHORT RUN (NO ENTROPY CALCULATED)"
1920 PRINT "ENTER L FOR LONG RUN (WITH ENTROPY CALCULATED)"
1930 INPUT G$
1932 PRINT "ENTER COLLECTION SAMPLE RATE (IE. 500)"
1934 INPUT O$
1935 S8=VAL(O$)
1940 RETURN
1950 REM
1960 REM *****
1970 REM * END OF ENTROPY CALCULATION PROGRAM
1980 REM *****
1990 REM
2000 PRINT
2010 PRINT
2020 PRINT "ENTROPY CALCULATION & STATISTICS PRINTOUT COMPLETE"
2030 END

```

Appendix D

This appendix contains a listing of the data compressed by the TOLAN-A8 module for the thesis experiment.

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1359PA
 SUBJECT PARNELL
 SAMPLING RATE 500 HZ
 DATE OF COLLECTION. . . . 23 OCT 80
 TIME OF COLLECTION. . . . 1359
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 2.4871 BITS
 CHANNEL Y ENTROPY 2.4751 BITS
 CHANNEL Z ENTROPY 2.4217 BITS
 TOTAL SOURCE ENTROPY. . . 2.4796 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 3.2263 : 1
 COMPRESSION RATIO
 ACHIEVED. 2.2595 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 56.5 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED . . . 77.8 % SMP INTERVAL
 COLLECTION DURATION . . . 26.2 SECONDS
 CHANNEL X MAXIMUM 2.50000 VOLTS
 CHANNEL X MINIMUM -0.1171 VOLTS
 CHANNEL Y MAXIMUM 2.50000 VOLTS
 CHANNEL Y MINIMUM -0.1171 VOLTS
 CHANNEL Z MAXIMUM 2.50000 VOLTS
 CHANNEL Z MINIMUM -0.1562 VOLTS
 COMMENTS. X,Y,Z IN COMMON LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 3336 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 0001B9DB (HEX)
 TIME EFFICIENCY = (1-(LOOPCT%(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 40 (HEX) AT SAMPLE NUMBER 0DA7 (HEX)
 XMIN= FD (HEX) AT SAMPLE NUMBER 0171 (HEX)
 YMAX= 40 (HEX) AT SAMPLE NUMBER 0DA7 (HEX)
 YMIN= FD (HEX) AT SAMPLE NUMBER 0171 (HEX)
 ZMAX= 40 (HEX) AT SAMPLE NUMBER 0DA7 (HEX)
 ZMIN= FC (HEX) AT SAMPLE NUMBER 0171 (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FFC (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00040000 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 0004CD10 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 006652 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 0066A8 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006160 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 010000 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABE

FILE TA1359PA
 CHANNEL X AMPLITUDE DISTRIBUTION
 DATA VALUE NUMBER OF OCCURRENCES

80.....	0000 (HEX)
81.....	0000 (HEX)
82.....	0000 (HEX)
83.....	0000 (HEX)
84.....	0000 (HEX)
85.....	0000 (HEX)
86.....	0000 (HEX)
87.....	0000 (HEX)
88.....	0000 (HEX)
89.....	0000 (HEX)
8A.....	0000 (HEX)
8B.....	0000 (HEX)
8C.....	0000 (HEX)
8D.....	0000 (HEX)
8E.....	0000 (HEX)
8F.....	0000 (HEX)
90.....	0000 (HEX)
91.....	0000 (HEX)
92.....	0000 (HEX)
93.....	0000 (HEX)
94.....	0000 (HEX)
95.....	0000 (HEX)
96.....	0000 (HEX)
97.....	0000 (HEX)
98.....	0000 (HEX)
99.....	0000 (HEX)
9A.....	0000 (HEX)
9B.....	0000 (HEX)
9C.....	0000 (HEX)
9D.....	0000 (HEX)
9E.....	0000 (HEX)
9F.....	0000 (HEX)
A0.....	0000 (HEX)
A1.....	0000 (HEX)
A2.....	0000 (HEX)
A3.....	0000 (HEX)
A4.....	0000 (HEX)
A5.....	0000 (HEX)
A6.....	0000 (HEX)
A7.....	0000 (HEX)
A8.....	0000 (HEX)
A9.....	0000 (HEX)
AA.....	0000 (HEX)
AB.....	0000 (HEX)
AC.....	0000 (HEX)
AD.....	0000 (HEX)
AE.....	0000 (HEX)
AF.....	0000 (HEX)
B0.....	0000 (HEX)
B1.....	0000 (HEX)
B2.....	0000 (HEX)

B3.....0000 (HEX)
B4.....0000 (HEX)
B5.....0000 (HEX)
B6.....0000 (HEX)
B7.....0000 (HEX)
B8.....0000 (HEX)
B9.....0000 (HEX)
BA.....0000 (HEX)
BB.....0000 (HEX)
BC.....0000 (HEX)
BD.....0000 (HEX)
BE.....0000 (HEX)
BF.....0000 (HEX)
C0.....0000 (HEX)
C1.....0000 (HEX)
C2.....0000 (HEX)
C3.....0000 (HEX)
C4.....0000 (HEX)
C5.....0000 (HEX)
C6.....0000 (HEX)
C7.....0000 (HEX)
C8.....0000 (HEX)
C9.....0000 (HEX)
CA.....0000 (HEX)
CB.....0000 (HEX)
CC.....0000 (HEX)
CD.....0000 (HEX)
CE.....0000 (HEX)
CF.....0000 (HEX)
D0.....0000 (HEX)
D1.....0000 (HEX)
D2.....0000 (HEX)
D3.....0000 (HEX)
D4.....0000 (HEX)
D5.....0000 (HEX)
D6.....0000 (HEX)
D7.....0000 (HEX)
D8.....0000 (HEX)
D9.....0000 (HEX)
DA.....0000 (HEX)
DB.....0000 (HEX)
DC.....0000 (HEX)
DD.....0000 (HEX)
DE.....0000 (HEX)
DF.....0000 (HEX)
E0.....0000 (HEX)
E1.....0000 (HEX)
E2.....0000 (HEX)
E3.....0000 (HEX)
E4.....0000 (HEX)
E5.....0000 (HEX)
E6.....0000 (HEX)
E7.....0000 (HEX)
E8.....0000 (HEX)

E9.....	0000	(HEX)
FA.....	0000	(HEX)
FB.....	0000	(HEX)
EC.....	0000	(HEX)
ED.....	0000	(HEX)
EE.....	0000	(HEX)
EF.....	0000	(HEX)
F0.....	0000	(HEX)
F1.....	0000	(HEX)
F2.....	0001	(HEX)
F3.....	0001	(HEX)
F4.....	0003	(HEX)
F5.....	000C	(HEX)
F6.....	0012	(HEX)
F7.....	0016	(HEX)
F8.....	000F	(HEX)
F9.....	000E	(HEX)
FA.....	000E	(HEX)
FB.....	000E	(HEX)
FC.....	001B	(HEX)
FD.....	004D	(HEX)
FE.....	0260	(HEX)
FF.....	0961	(HEX)
00.....	06A2	(HEX)
01.....	09A5	(HEX)
02.....	024A	(HEX)
03.....	005F	(HEX)
04.....	0027	(HEX)
05.....	0011	(HEX)
06.....	0003	(HEX)
07.....	0008	(HEX)
08.....	0008	(HEX)
09.....	0004	(HEX)
0A.....	0001	(HEX)
0B.....	0001	(HEX)
0C.....	0006	(HEX)
0D.....	0003	(HEX)
0E.....	0003	(HEX)
0F.....	000A	(HEX)
10.....	0009	(HEX)
11.....	0000	(HEX)
12.....	0003	(HEX)
13.....	0002	(HEX)
14.....	0000	(HEX)
15.....	0000	(HEX)
16.....	0000	(HEX)
17.....	0000	(HEX)
18.....	0000	(HEX)
19.....	0000	(HEX)
1A.....	0000	(HEX)
1B.....	0000	(HEX)
1C.....	0000	(HEX)
1D.....	0000	(HEX)
1E.....	0000	(HEX)

1F.....0000 (HEX)
20.....0000 (HEX)
21.....0000 (HEX)
22.....0000 (HEX)
23.....0000 (HEX)
24.....0000 (HEX)
25.....0000 (HEX)
26.....0000 (HEX)
27.....0000 (HEX)
28.....0000 (HEX)
29.....0000 (HEX)
2A.....0000 (HEX)
2B.....0000 (HEX)
2C.....0000 (HEX)
2D.....0000 (HEX)
2E.....0000 (HEX)
2F.....0000 (HEX)
30.....0000 (HEX)
31.....0000 (HEX)
32.....0000 (HEX)
33.....0000 (HEX)
34.....0000 (HEX)
35.....0000 (HEX)
36.....0000 (HEX)
37.....0000 (HEX)
38.....0000 (HEX)
39.....0000 (HEX)
3A.....0000 (HEX)
3B.....0000 (HEX)
3C.....0000 (HEX)
3D.....0000 (HEX)
3E.....0000 (HEX)
3F.....0000 (HEX)
40.....0000 (HEX)
41.....0000 (HEX)
42.....0000 (HEX)
43.....0000 (HEX)
44.....0000 (HEX)
45.....0000 (HEX)
46.....0000 (HEX)
47.....0000 (HEX)
48.....0000 (HEX)
49.....0000 (HEX)
4A.....0000 (HEX)
4B.....0000 (HEX)
4C.....0000 (HEX)
4D.....0000 (HEX)
4E.....0000 (HEX)
4F.....0000 (HEX)
50.....0000 (HEX)
51.....0000 (HEX)
52.....0000 (HEX)
53.....0000 (HEX)
54.....0000 (HEX)

55.....0000 (HEX)
56.....0000 (HEX)
57.....0000 (HEX)
58.....0000 (HEX)
59.....0000 (HEX)
5A.....0000 (HEX)
5B.....0000 (HEX)
5C.....0000 (HEX)
5D.....0000 (HEX)
5E.....0000 (HEX)
5F.....0000 (HEX)
60.....0000 (HEX)
61.....0000 (HEX)
62.....0000 (HEX)
63.....0000 (HEX)
64.....0000 (HEX)
65.....0000 (HEX)
66.....0000 (HEX)
67.....0000 (HEX)
68.....0000 (HEX)
69.....0000 (HEX)
6A.....0000 (HEX)
6B.....0000 (HEX)
6C.....0000 (HEX)
6D.....0000 (HEX)
6E.....0000 (HEX)
6F.....0000 (HEX)
70.....0000 (HEX)
71.....0000 (HEX)
72.....0000 (HEX)
73.....0000 (HEX)
74.....0000 (HEX)
75.....0000 (HEX)
76.....0000 (HEX)
77.....0000 (HEX)
78.....0000 (HEX)
79.....0000 (HEX)
7A.....0000 (HEX)
7B.....0000 (HEX)
7C.....0000 (HEX)
7D.....0000 (HEX)
7E.....0000 (HEX)
7F.....0000 (HEX)

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1413LU
 SUBJECT LUTZ
 SAMPLING RATE 500 HZ
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1413
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 3.7826 BITS
 CHANNEL Y ENTROPY 3.7950 BITS
 CHANNEL Z ENTROPY 3.8149 BITS
 TOTAL SOURCE ENTROPY. 3.8028 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.1036 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.2519 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 22.8 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED 95.4 % SMP INTERVAL
 COLLECTION DURATION 14.5 SECONDS
 CHANNEL X MAXIMUM 1.79687 VOLTS
 CHANNEL X MINIMUM -1.0156 VOLTS
 CHANNEL Y MAXIMUM 1.75781 VOLTS
 CHANNEL Y MINIMUM -1.0156 VOLTS
 CHANNEL Z MAXIMUM 1.75781 VOLTS
 CHANNEL Z MINIMUM -1.0156 VOLTS
 COMMENTS. X,Y,Z IN COMMON. LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 1C60 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 000031E3 (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 2E (HEX) AT SAMPLE NUMBER 0E7F (HEX)
 XMIN= E6 (HEX) AT SAMPLE NUMBER 11C4 (HEX)
 YMAX= 2D (HEX) AT SAMPLE NUMBER 0E7F (HEX)
 YMIN= E6 (HEX) AT SAMPLE NUMBER 11C4 (HEX)
 ZMAX= 2D (HEX) AT SAMPLE NUMBER 0E7F (HEX)
 ZMIN= E6 (HEX) AT SAMPLE NUMBER 11C4 (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00032040 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 0002A900 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 007786 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 0077F8 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 00764C (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 00CA50 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABE

FILE TA1413LU

CHANNEL X AMPLITUDE DISTRIBUTION

DATA VALUE NUMBER OF OCCURENCES

80.....	0000 (HEX)
81.....	0000 (HEX)
82.....	0000 (HEX)
83.....	0000 (HEX)
84.....	0000 (HEX)
85.....	0000 (HEX)
86.....	0000 (HEX)
87.....	0000 (HEX)
88.....	0000 (HEX)
89.....	0000 (HEX)
8A.....	0000 (HEX)
8B.....	0000 (HEX)
8C.....	0000 (HEX)
8D.....	0000 (HEX)
8E.....	0000 (HEX)
8F.....	0000 (HEX)
90.....	0000 (HEX)
91.....	0000 (HEX)
92.....	0000 (HEX)
93.....	0000 (HEX)
94.....	0000 (HEX)
95.....	0000 (HEX)
96.....	0000 (HEX)
97.....	0000 (HEX)
98.....	0000 (HEX)
99.....	0000 (HEX)
9A.....	0000 (HEX)
9B.....	0000 (HEX)
9C.....	0000 (HEX)
9D.....	0000 (HEX)
9E.....	0000 (HEX)
9F.....	0000 (HEX)
A0.....	0000 (HEX)
A1.....	0000 (HEX)
A2.....	0000 (HEX)
A3.....	0000 (HEX)
A4.....	0000 (HEX)
A5.....	0000 (HEX)
A6.....	0000 (HEX)
A7.....	0000 (HEX)
A8.....	0000 (HEX)
A9.....	0000 (HEX)
AA.....	0000 (HEX)
AB.....	0000 (HEX)
AC.....	0000 (HEX)
AD.....	0000 (HEX)
AE.....	0000 (HEX)
AF.....	0000 (HEX)
B0.....	0000 (HEX)
B1.....	0000 (HEX)
B2.....	0000 (HEX)

B3.....0000 (HEX)
B4.....0000 (HEX)
B5.....0000 (HEX)
B6.....0000 (HEX)
B7.....0000 (HEX)
B8.....0000 (HEX)
B9.....0000 (HEX)
BA.....0000 (HEX)
BB.....0000 (HEX)
BC.....0000 (HEX)
BD.....0000 (HEX)
BE.....0000 (HEX)
BF.....0000 (HEX)
C0.....0000 (HEX)
C1.....0000 (HEX)
C2.....0000 (HEX)
C3.....0000 (HEX)
C4.....0000 (HEX)
C5.....0000 (HEX)
C6.....0000 (HEX)
C7.....0000 (HEX)
C8.....0000 (HEX)
C9.....0000 (HEX)
CA.....0000 (HEX)
CB.....0000 (HEX)
CC.....0000 (HEX)
CD.....0000 (HEX)
CE.....0000 (HEX)
CF.....0000 (HEX)
D0.....0000 (HEX)
D1.....0000 (HEX)
D2.....0000 (HEX)
D3.....0000 (HEX)
D4.....0000 (HEX)
D5.....0000 (HEX)
D6.....0000 (HEX)
D7.....0000 (HEX)
D8.....0000 (HEX)
D9.....0000 (HEX)
DA.....0000 (HEX)
DB.....0000 (HEX)
DC.....0000 (HEX)
DD.....0000 (HEX)
DE.....0000 (HEX)
DF.....0000 (HEX)
E0.....0000 (HEX)
E1.....0000 (HEX)
E2.....0000 (HEX)
E3.....0000 (HEX)
E4.....0000 (HEX)
E5.....0000 (HEX)
E6.....0000 (HEX)
E7.....0000 (HEX)
E8.....0000 (HEX)

E9.....	0000	(HEX)
EA.....	0001	(HEX)
EB.....	0000	(HEX)
EC.....	0000	(HEX)
ED.....	0000	(HEX)
EE.....	0000	(HEX)
EF.....	0001	(HEX)
F0.....	0000	(HEX)
F1.....	0000	(HEX)
F2.....	0001	(HEX)
F3.....	0006	(HEX)
F4.....	0007	(HEX)
F5.....	0012	(HEX)
F6.....	0012	(HEX)
F7.....	0022	(HEX)
F8.....	0035	(HEX)
F9.....	0050	(HEX)
FA.....	0082	(HEX)
FB.....	00CC	(HEX)
FC.....	0136	(HEX)
FD.....	01D9	(HEX)
FE.....	02EF	(HEX)
FF.....	03D0	(HEX)
00.....	0164	(HEX)
01.....	040A	(HEX)
02.....	02BB	(HEX)
03.....	01E4	(HEX)
04.....	0132	(HEX)
05.....	00AD	(HEX)
06.....	006B	(HEX)
07.....	0058	(HEX)
08.....	0041	(HEX)
09.....	001F	(HEX)
0A.....	0015	(HEX)
0B.....	0013	(HEX)
0C.....	000C	(HEX)
0D.....	000A	(HEX)
0E.....	0002	(HEX)
0F.....	0002	(HEX)
10.....	0000	(HEX)
11.....	0000	(HEX)
12.....	0001	(HEX)
13.....	0000	(HEX)
14.....	0000	(HEX)
15.....	0000	(HEX)
16.....	0000	(HEX)
17.....	0001	(HEX)
18.....	0000	(HEX)
19.....	0000	(HEX)
1A.....	0000	(HEX)
1B.....	0000	(HEX)
1C.....	0000	(HEX)
1D.....	0000	(HEX)
1E.....	0000	(HEX)

1F.....0000 (HEX)
20.....0000 (HEX)
21.....0000 (HEX)
22.....0000 (HEX)
23.....0000 (HEX)
24.....0000 (HEX)
25.....0000 (HEX)
26.....0000 (HEX)
27.....0000 (HEX)
28.....0000 (HEX)
29.....0000 (HEX)
2A.....0000 (HEX)
2B.....0000 (HEX)
2C.....0000 (HEX)
2D.....0000 (HEX)
2E.....0000 (HEX)
2F.....0000 (HEX)
30.....0000 (HEX)
31.....0000 (HEX)
32.....0000 (HEX)
33.....0000 (HEX)
34.....0000 (HEX)
35.....0000 (HEX)
36.....0000 (HEX)
37.....0000 (HEX)
38.....0000 (HEX)
39.....0000 (HEX)
3A.....0000 (HEX)
3B.....0000 (HEX)
3C.....0000 (HEX)
3D.....0000 (HEX)
3E.....0000 (HEX)
3F.....0000 (HEX)
40.....0000 (HEX)
41.....0000 (HEX)
42.....0000 (HEX)
43.....0000 (HEX)
44.....0000 (HEX)
45.....0000 (HEX)
46.....0000 (HEX)
47.....0000 (HEX)
48.....0000 (HEX)
49.....0000 (HEX)
4A.....0000 (HEX)
4B.....0000 (HEX)
4C.....0000 (HEX)
4D.....0000 (HEX)
4E.....0000 (HEX)
4F.....0000 (HEX)
50.....0000 (HEX)
51.....0000 (HEX)
52.....0000 (HEX)
53.....0000 (HEX)
54.....0000 (HEX)

55.....0000 (HEX)
56.....0000 (HEX)
57.....0000 (HEX)
58.....0000 (HEX)
59.....0000 (HEX)
5A.....0000 (HEX)
5B.....0000 (HEX)
5C.....0000 (HEX)
5D.....0000 (HEX)
5E.....0000 (HEX)
5F.....0000 (HEX)
60.....0000 (HEX)
61.....0000 (HEX)
62.....0000 (HEX)
63.....0000 (HEX)
64.....0000 (HEX)
65.....0000 (HEX)
66.....0000 (HEX)
67.....0000 (HEX)
68.....0000 (HEX)
69.....0000 (HEX)
6A.....0000 (HEX)
6B.....0000 (HEX)
6C.....0000 (HEX)
6D.....0000 (HEX)
6E.....0000 (HEX)
6F.....0000 (HEX)
70.....0000 (HEX)
71.....0000 (HEX)
72.....0000 (HEX)
73.....0000 (HEX)
74.....0000 (HEX)
75.....0000 (HEX)
76.....0000 (HEX)
77.....0000 (HEX)
78.....0000 (HEX)
79.....0000 (HEX)
7A.....0000 (HEX)
7B.....0000 (HEX)
7C.....0000 (HEX)
7D.....0000 (HEX)
7E.....0000 (HEX)
7F.....0000 (HEX)

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1448L
 COLLECTION. 10Y14 LN
 SAMPLING RATE. 500 HZ
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1449
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY. 2.9299 BITS
 CHANNEL Y ENTROPY. 2.9243 BITS
 CHANNEL Z ENTROPY. 2.9564 BITS
 TOTAL SOURCE ENTROPY. 2.9436 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.7176 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.6008 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 34.9 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED. 90.6 % SMP INTERVAL
 COLLECTION DURATION. 18.5 SECONDS
 CHANNEL X MAXIMUM. 1.79687 VOLTS
 CHANNEL X MINIMUM. -0.7421 VOLTS
 CHANNEL Y MAXIMUM. 1.79687 VOLTS
 CHANNEL Y MINIMUM. -0.7421 VOLTS
 CHANNEL Z MAXIMUM. 1.79687 VOLTS
 CHANNEL Z MINIMUM. -0.7812 VOLTS
 COMMENTS. X,Y,Z IN COMMON. LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 2448 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 0000345C (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 2E (HEX) AT SAMPLE NUMBER 0D0A (HEX)
 XMIN= ED (HEX) AT SAMPLE NUMBER 1D0A (HEX)
 YMAX= 2E (HEX) AT SAMPLE NUMBER 0D0A (HEX)
 YMIN= ED (HEX) AT SAMPLE NUMBER 1D0A (HEX)
 ZMAX= 2E (HEX) AT SAMPLE NUMBER 0D0A (HEX)
 ZMIN= EC (HEX) AT SAMPLE NUMBER 1D0A (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 0003A500 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 000366C0 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 0069AA (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 00693E (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006739 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 00EB98 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABLE

FILE TA1448LT
CHANNEL X AMPLITUDE DISTRIBUTION

DATA VALUE	NUMBER OF OCCURENCES
80.....	0000 (HEX)
81.....	0000 (HEX)
82.....	0000 (HEX)
83.....	0000 (HEX)
84.....	0000 (HEX)
85.....	0000 (HEX)
86.....	0000 (HEX)
87.....	0000 (HEX)
88.....	0000 (HEX)
89.....	0000 (HEX)
8A.....	0000 (HEX)
8B.....	0000 (HEX)
8C.....	0000 (HEX)
8D.....	0000 (HEX)
8E.....	0000 (HEX)
8F.....	0000 (HEX)
90.....	0000 (HEX)
91.....	0000 (HEX)
92.....	0000 (HEX)
93.....	0000 (HEX)
94.....	0000 (HEX)
95.....	0000 (HEX)
96.....	0000 (HEX)
97.....	0000 (HEX)
98.....	0000 (HEX)
99.....	0000 (HEX)
9A.....	0000 (HEX)
9B.....	0000 (HEX)
9C.....	0000 (HEX)
9D.....	0000 (HEX)
9E.....	0000 (HEX)
9F.....	0000 (HEX)
A0.....	0000 (HEX)
A1.....	0000 (HEX)
A2.....	0000 (HEX)
A3.....	0000 (HEX)
A4.....	0000 (HEX)
A5.....	0000 (HEX)
A6.....	0000 (HEX)
A7.....	0000 (HEX)
A8.....	0000 (HEX)
A9.....	0000 (HEX)
AA.....	0000 (HEX)
AB.....	0000 (HEX)
AC.....	0000 (HEX)
AD.....	0000 (HEX)
AE.....	0000 (HEX)
AF.....	0000 (HEX)
B0.....	0000 (HEX)
B1.....	0000 (HEX)
B2.....	0000 (HEX)

B3.....	0000	(HEX)
B4.....	0000	(HEX)
B5.....	0000	(HEX)
B6.....	0000	(HEX)
B7.....	0000	(HEX)
B8.....	0000	(HEX)
B9.....	0000	(HEX)
BA.....	0000	(HEX)
BB.....	0000	(HEX)
BC.....	0000	(HEX)
BD.....	0000	(HEX)
BE.....	0000	(HEX)
BF.....	0000	(HEX)
C0.....	0000	(HEX)
C1.....	0000	(HEX)
C2.....	0000	(HEX)
C3.....	0000	(HEX)
C4.....	0000	(HEX)
C5.....	0000	(HEX)
C6.....	0000	(HEX)
C7.....	0000	(HEX)
C8.....	0000	(HEX)
C9.....	0000	(HEX)
CA.....	0000	(HEX)
CB.....	0000	(HEX)
CC.....	0000	(HEX)
CD.....	0000	(HEX)
CE.....	0000	(HEX)
CF.....	0000	(HEX)
D0.....	0000	(HEX)
D1.....	0000	(HEX)
D2.....	0000	(HEX)
D3.....	0000	(HEX)
D4.....	0000	(HEX)
D5.....	0000	(HEX)
D6.....	0000	(HEX)
D7.....	0000	(HEX)
D8.....	0000	(HEX)
D9.....	0000	(HEX)
DA.....	0000	(HEX)
DB.....	0000	(HEX)
DC.....	0000	(HEX)
DD.....	0000	(HEX)
DE.....	0000	(HEX)
DF.....	0000	(HEX)
E0.....	0000	(HEX)
E1.....	0000	(HEX)
E2.....	0000	(HEX)
E3.....	0000	(HEX)
E4.....	0000	(HEX)
E5.....	0000	(HEX)
E6.....	0000	(HEX)
E7.....	0000	(HEX)
E8.....	0000	(HEX)

E9.....	0000	(HEX)
EA.....	0000	(HEX)
EB.....	0000	(HEX)
EC.....	0000	(HEX)
ED.....	0000	(HEX)
EE.....	0000	(HEX)
EF.....	0000	(HEX)
F0.....	0000	(HEX)
F1.....	0000	(HEX)
F2.....	0000	(HEX)
F3.....	0000	(HEX)
F4.....	0000	(HEX)
F5.....	0000	(HEX)
F6.....	0002	(HEX)
F7.....	0001	(HEX)
F8.....	0003	(HEX)
F9.....	0007	(HEX)
FA.....	001C	(HEX)
FB.....	0034	(HEX)
FC.....	008F	(HEX)
FD.....	0180	(HEX)
FE.....	03F7	(HEX)
FF.....	06A2	(HEX)
00.....	030D	(HEX)
01.....	0714	(HEX)
02.....	03F0	(HEX)
03.....	018A	(HEX)
04.....	008A	(HEX)
05.....	002A	(HEX)
06.....	0013	(HEX)
07.....	0006	(HEX)
08.....	0002	(HEX)
09.....	0002	(HEX)
0A.....	0002	(HEX)
0B.....	0000	(HEX)
0C.....	0000	(HEX)
0D.....	0000	(HEX)
0E.....	0000	(HEX)
0F.....	0000	(HEX)
10.....	0000	(HEX)
11.....	0000	(HEX)
12.....	0000	(HEX)
13.....	0000	(HEX)
14.....	0000	(HEX)
15.....	0000	(HEX)
16.....	0000	(HEX)
17.....	0000	(HEX)
18.....	0000	(HEX)
19.....	0000	(HEX)
1A.....	0000	(HEX)
1B.....	0000	(HEX)
1C.....	0000	(HEX)
1D.....	0000	(HEX)
1E.....	0000	(HEX)

1F.....	0000	(HEX)
20.....	0000	(HEX)
21.....	0000	(HEX)
22.....	0000	(HEX)
23.....	0000	(HEX)
24.....	0000	(HEX)
25.....	0000	(HEX)
26.....	0000	(HEX)
27.....	0000	(HEX)
28.....	0000	(HEX)
29.....	0000	(HEX)
2A.....	0000	(HEX)
2B.....	0000	(HEX)
2C.....	0000	(HEX)
2D.....	0000	(HEX)
2E.....	0000	(HEX)
2F.....	0000	(HEX)
30.....	0000	(HEX)
31.....	0000	(HEX)
32.....	0000	(HEX)
33.....	0000	(HEX)
34.....	0000	(HEX)
35.....	0000	(HEX)
36.....	0000	(HEX)
37.....	0000	(HEX)
38.....	0000	(HEX)
39.....	0000	(HEX)
3A.....	0000	(HEX)
3B.....	0000	(HEX)
3C.....	0000	(HEX)
3D.....	0000	(HEX)
3E.....	0000	(HEX)
3F.....	0000	(HEX)
40.....	0000	(HEX)
41.....	0000	(HEX)
42.....	0000	(HEX)
43.....	0000	(HEX)
44.....	0000	(HEX)
45.....	0000	(HEX)
46.....	0000	(HEX)
47.....	0000	(HEX)
48.....	0000	(HEX)
49.....	0000	(HEX)
4A.....	0000	(HEX)
4B.....	0000	(HEX)
4C.....	0000	(HEX)
4D.....	0000	(HEX)
4E.....	0000	(HEX)
4F.....	0000	(HEX)
50.....	0000	(HEX)
51.....	0000	(HEX)
52.....	0000	(HEX)
53.....	0000	(HEX)
54.....	0000	(HEX)

55.....0000 (HEX)
56.....0000 (HEX)
57.....0000 (HEX)
58.....0000 (HEX)
59.....0000 (HEX)
5A.....0000 (HEX)
5B.....0000 (HEX)
5C.....0000 (HEX)
5D.....0000 (HEX)
5E.....0000 (HEX)
5F.....0000 (HEX)
60.....0000 (HEX)
61.....0000 (HEX)
62.....0000 (HEX)
63.....0000 (HEX)
64.....0000 (HEX)
65.....0000 (HEX)
66.....0000 (HEX)
67.....0000 (HEX)
68.....0000 (HEX)
69.....0000 (HEX)
6A.....0000 (HEX)
6B.....0000 (HEX)
6C.....0000 (HEX)
6D.....0000 (HEX)
6E.....0000 (HEX)
6F.....0000 (HEX)
70.....0000 (HEX)
71.....0000 (HEX)
72.....0000 (HEX)
73.....0000 (HEX)
74.....0000 (HEX)
75.....0000 (HEX)
76.....0000 (HEX)
77.....0000 (HEX)
78.....0000 (HEX)
79.....0000 (HEX)
7A.....0000 (HEX)
7B.....0000 (HEX)
7C.....0000 (HEX)
7D.....0000 (HEX)
7E.....0000 (HEX)
7F.....0000 (HEX)

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1545T
 SERVICE POLYMAP
 SAMPLING RATE 500 Hz
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1546
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 3.2801 BITS
 CHANNEL Y ENTROPY 3.2716 BITS
 CHANNEL Z ENTROPY 3.3188 BITS
 TOTAL SOURCE ENTROPY. 3.3026 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.4223 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.3617 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 25.4 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED 93.9 % SMP INTERVAL
 COLLECTION DURATION 15.8 SECONDS
 CHANNEL X MAXIMUM 3.59375 VOLTS
 CHANNEL X MINIMUM -1.4453 VOLTS
 CHANNEL Y MAXIMUM 3.59375 VOLTS
 CHANNEL Y MINIMUM -1.4453 VOLTS
 CHANNEL Z MAXIMUM 3.55468 VOLTS
 CHANNEL Z MINIMUM -1.4453 VOLTS
 COMMENTS. X,Y,Z, IN COMMON, LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 1EDD (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 26 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 0000470E (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 5C (HEX) AT SAMPLE NUMBER 03DD (HEX)
 XMIN= DB (HEX) AT SAMPLE NUMBER 145D (HEX)
 YMAX= 5C (HEX) AT SAMPLE NUMBER 03DD (HEX)
 YMIN= DB (HEX) AT SAMPLE NUMBER 145D (HEX)
 ZMAX= 5B (HEX) AT SAMPLE NUMBER 03DD (HEX)
 ZMIN= DB (HEX) AT SAMPLE NUMBER 145D (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00034820 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 0002E4B8 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 0071DF (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 007199 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 007013 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 00D220 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABLE

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1548T
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1548
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 3.0171 BITS
 CHANNEL Y ENTROPY 3.0114 BITS
 CHANNEL Z ENTROPY 3.0373 BITS
 TOTAL SOURCE ENTROPY. 3.0355 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.6354 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.4998 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 30.5 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED 91.4 % SMP INTERVAL
 COLLECTION DURATION 17.4 SECONDS
 CHANNEL X MAXIMUM 2.77343 VOLTS
 CHANNEL X MINIMUM -1.3671 VOLTS
 CHANNEL Y MAXIMUM 2.73437 VOLTS
 CHANNEL Y MINIMUM -1.3671 VOLTS
 CHANNEL Z MAXIMUM 2.73437 VOLTS
 CHANNEL Z MINIMUM -1.3671 VOLTS
 COMMENTS. X,Y,Z IN COMMON, LEAD AVL
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 21FE (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 00007134 (HEX)
 TIME EFFICIENCY = (1-(LOOPCT%(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 47 (HEX) AT SAMPLE NUMBER 0A9E (HEX)
 XMIN= DD (HEX) AT SAMPLE NUMBER 1447 (HEX)
 YMAX= 46 (HEX) AT SAMPLE NUMBER 095A (HEX)
 YMIN= DD (HEX) AT SAMPLE NUMBER 1447 (HEX)
 ZMAX= 46 (HEX) AT SAMPLE NUMBER 095A (HEX)
 ZMIN= DD (HEX) AT SAMPLE NUMBER 1447 (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00036000 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 00032FD0 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 006EE8 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 006E90 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006F68 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 00DB48 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABLE

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME TA1559B
 BLANK
 SAMPLING RATE 500 HZ
 DATE OF COLLECTION 23 OCT 80
 TIME OF COLLECTION 1559
 COMPRESSION USED TOLAN-A
 CHANNEL X ENTROPY 3.0165 BITS
 CHANNEL Y ENTROPY 3.0038 BITS
 CHANNEL Z ENTROPY 3.0261 BITS
 TOTAL SOURCE ENTROPY 3.0233 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE 2.6460 : 1
 COMPRESSION RATIO
 ACHIEVED 1.5280 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY 32.0 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED 93.0 % SMP INTERVAL
 COLLECTION DURATION 17.7 SECONDS
 CHANNEL X MAXIMUM 2.92968 VOLTS
 CHANNEL X MINIMUM 0.19531 VOLTS
 CHANNEL Y MAXIMUM 2.89062 VOLTS
 CHANNEL Y MINIMUM 0.19531 VOLTS
 CHANNEL Z MAXIMUM 2.89062 VOLTS
 CHANNEL Z MINIMUM 0.15625 VOLTS
 COMMENTS X,Y,Z IN COMMON, LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 22A2 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 26 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 00005B1E (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 4B (HEX) AT SAMPLE NUMBER 1FA2 (HEX)
 XMIN= 05 (HEX) AT SAMPLE NUMBER 00BD (HEX)
 YMAX= 4A (HEX) AT SAMPLE NUMBER 1FA2 (HEX)
 YMIN= 05 (HEX) AT SAMPLE NUMBER 005C (HEX)
 ZMAX= 4A (HEX) AT SAMPLE NUMBER 1FA2 (HEX)
 ZMIN= 04 (HEX) AT SAMPLE NUMBER 01D5 (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN COMP = 0003A780 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 00033F30 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 006BE5 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 006AC9 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006857 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETER = 00E9E0 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABE

:EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1511S
 SAMPLING RATE 500 HZ
 DATE OF COLLECTION. . . . 1511
 TIME OF COLLECTION. . . . 1511
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 3.3255 BITS
 CHANNEL Y ENTROPY 3.3341 BITS
 CHANNEL Z ENTROPY 3.3560 BITS
 TOTAL SOURCE ENTROPY. . . 3.3457 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.3911 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.4319 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 31.0 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED . . . 94.0 % SMP INTERVAL
 COLLECTION DURATION . . . 16.6 SECONDS
 CHANNEL X MAXIMUM 1.32812 VOLTS
 CHANNEL X MINIMUM -2.0703 VOLTS
 CHANNEL Y MAXIMUM 1.32812 VOLTS
 CHANNEL Y MINIMUM -2.0703 VOLTS
 CHANNEL Z MAXIMUM 1.32812 VOLTS
 CHANNEL Z MINIMUM -2.0703 VOLTS
 COMMENTS. X,Y,Z IN COMMON. LEAD AVL
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 2074 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 00004ADC (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 22 (HEX) AT SAMPLE NUMBER 1171 (HEX)
 XMIN= CB (HEX) AT SAMPLE NUMBER 021D (HEX)
 YMAX= 22 (HEX) AT SAMPLE NUMBER 1423 (HEX)
 YMIN= CB (HEX) AT SAMPLE NUMBER 021D (HEX)
 ZMAX= 22 (HEX) AT SAMPLE NUMBER 1423 (HEX)
 ZMIN= CB (HEX) AT SAMPLE NUMBER 00E8 (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00039000 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 00030AE0 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 00717A (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 00713D (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006127 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETERS = 00F000 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM BITS AVAILABLE

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1520B
 SUBJECT. BALSALO
 SAMPLING RATE. 500 HZ
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1520
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY. 2.6012 BITS
 CHANNEL Y ENTROPY. 2.5808 BITS
 CHANNEL Z ENTROPY. 2.6229 BITS
 TOTAL SOURCE ENTROPY. 2.6104 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 3.0645 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.7247 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 35.1 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED. 90.9 % SMP INTERVAL
 COLLECTION DURATION. 20.0 SECONDS
 CHANNEL X MAXIMUM. 2.07031 VOLTS
 CHANNEL X MINIMUM. -0.5468 VOLTS
 CHANNEL Y MAXIMUM. 2.07031 VOLTS
 CHANNEL Y MINIMUM. -0.5468 VOLTS
 CHANNEL Z MAXIMUM. 2.07031 VOLTS
 CHANNEL Z MINIMUM. -0.5468 VOLTS
 COMMENTS. X,Y,Z IN COMMON. LEA D 2
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 2717 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 00008A6D (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 35 (HEX) AT SAMPLE NUMBER 1F9E (HEX)
 XMIN= F2 (HEX) AT SAMPLE NUMBER 0A63 (HEX)
 YMAX= 35 (HEX) AT SAMPLE NUMBER 1F9E (HEX)
 YMIN= F2 (HEX) AT SAMPLE NUMBER 054B (HEX)
 ZMAX= 35 (HEX) AT SAMPLE NUMBER 1F9E (HEX)
 ZMIN= F2 (HEX) AT SAMPLE NUMBER 054B (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00000000 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 0002AA28 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 00066E (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 000666 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 0006CC (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER INFORMATION = 00F630 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER 1024 BITS AVAILABLE

EKG SAMPLE COLLECTION STATISTICS : PAGE 1

FILENAME. TA1439ST
 SUBJECT STROUP
 SAMPLING RATE 500 HZ
 DATE OF COLLECTION. 23 OCT 80
 TIME OF COLLECTION. 1439
 COMPRESSION USED. TOLAN-A
 CHANNEL X ENTROPY 3.2669 BITS
 CHANNEL Y ENTROPY 3.2756 BITS
 CHANNEL Z ENTROPY 3.3094 BITS
 TOTAL SOURCE ENTROPY. 3.2912 BITS
 PRESS RETURN FOR PAGE 2 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS : PAGE 2

APPROX MAX COMPRESSION
 RATIO POSSIBLE. 2.4307 : 1
 COMPRESSION RATIO
 ACHIEVED. 1.4314 : 1
 ACHIEVED COMPRESSION
 EFFICIENCY. 30.1 % OF MAXIMUM
 COMPRESSION TIME
 EFFICIENCY OBTAINED 93.4 % SMP INTERVAL
 COLLECTION DURATION 16.6 SECONDS
 CHANNEL X MAXIMUM 3.35937 VOLTS
 CHANNEL X MINIMUM -0.5468 VOLTS
 CHANNEL Y MAXIMUM 3.35937 VOLTS
 CHANNEL Y MINIMUM -0.5468 VOLTS
 CHANNEL Z MAXIMUM 3.35937 VOLTS
 CHANNEL Z MINIMUM -0.5859 VOLTS
 COMMENTS. X,Y,Z IN COMMON. LEAD 1
 PRESS RETURN FOR PAGE 3 OF STATISTICS=

EKG SAMPLE COLLECTION STATISTICS: PAGE 3

NUMBER OF SAMPLES TAKEN (SAMPNO)= 2071 (HEX)
 MAXIMUM LOOP COUNT PER INTERRUPT (LPCAL)= 27 (HEX)
 TOTAL WAITING LOOP COUNTS DURING COLLECTION (LOOPCT)= 0000526C (HEX)
 TIME EFFICIENCY = (1-(LOOPCT*(SAMPNO*LPCAL)))*100
 CHANNEL MAXIMUMS AND MINIMUMS
 XMAX= 56 (HEX) AT SAMPLE NUMBER 1472 (HEX)
 XMIN= F2 (HEX) AT SAMPLE NUMBER 1E5A (HEX)
 YMAX= 56 (HEX) AT SAMPLE NUMBER 1472 (HEX)
 YMIN= F2 (HEX) AT SAMPLE NUMBER 0FD8 (HEX)
 ZMAX= 56 (HEX) AT SAMPLE NUMBER 1472 (HEX)
 ZMIN= F1 (HEX) AT SAMPLE NUMBER 1E5A (HEX)
 COMPRESSION STATISTICS:
 NUMBER OF MEMORY BITS AVAILABLE = 021FF0 (HEX)
 NUMBER OF BITS AVAILABLE TO VAR LEN CODER = 00037700 (HEX)
 TOTAL NUMBER OF DATA BITS STORED = 00030A98 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL X = 006E88 (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Y = 006EDC (HEX)
 NUMBER OF BITS USED TO STORE CHANNEL Z = 006CF0 (HEX)
 NUMBER OF BITS USED TO STORE TIME OR OTHER PARAMETERS = 00DDC0 (HEX)
 COMPRESSION RATIO = TOTAL DATA BITS STORED PER MEM. BITS AVAILABLE

Appendix E

This appendix contains a photocopy of the specifications for the pertinent equipment used in this thesis.


DK-VL

Copied from Vn. VanSleet's Plan. (1/22)



Bld 79

550150

RESEARCH RECORDER



OPERATING INSTRUCTIONS



ELECTRONICS FOR MEDICINE INC

30 VIRGINIA ROAD, WHITE PLAINS, N.Y. 10613 TEL. (914) 941-6112


DR-8 RESEARCH RECORDER MONITOR

The Model DR-8 Research Recorder is designed to operate from a 100-130 volt, 60 Hz line. Power consumption is about 700 watts. The circuit is protected by a 10 ampere circuit breaker which is incorporated in the power ON-OFF switch on the camera panel.

It is essential for safety and to minimize 60 Hz interference in the tracings that the Research Recorder be connected to a good "ground". A third wire is brought out at the power plug for this purpose. A good "ground" can usually be obtained at a cold water pipe.

The instrument combines as many as eight amplifiers with a cathode ray tube recording camera and two cathode ray tube monitors. The multitrace monitor displays any number of traces up to eight. The single monitor displays any one of the traces enlarged or reduced in size. It is also used for the presentation of loops, such as the vectorcardiogram or the lung compliance and lung resistance loops, and for balancing strain gauge amplifiers. All phenomena can be monitored continuously. Recording is done at will, merely by turning a switch. The display on the cathode ray tube in the camera can be identical to that on either the multitrace tube or the single trace tube. It is therefore possible to photograph either scalar tracings or loops. In actual operation, size adjustments and preliminary measurements are made by observing the monitor screens without recording. Thus, the camera runs for only short intervals.

The multitrace tube has only one beam, but because the amplifiers are connected to it through an electronic switch operating at a rate of 45,000 Hz/sec., it is possible to display as many as eight traces simultaneously. The beam switches rapidly from one trace to the next, with blanking of the cathode ray tube taking place during the switching period. Thus, each beam is on screen only 1/8th of the time, although giving the impression of being continuous. A switch is provided to make possible more frequent sampling for high speed phenomena, by reducing the number of channels displayed.



POWER switch, located on the right side of camera panel, turns the Research Recorder on. A minute or two is required before the traces on the monitor screens are stabilized. Adjustments may now be made.

INTENSITY and FOCUS control adjust the cathode ray tube beams. The upper two controls are for the multitrace tube, the lower two for the single trace tube. The intensity is set according to the brightness of the room in which the instrument is being used. The traces will be sharper at lower intensities. Focus is adjusted for the sharpest spot and should be done when the beams are approximately 1/3rd of the distance from either side of the tube to allow best overall focusing.

FUNCTION selector switch determines the presentation on the single trace tube; the multitrace tube will display scalar traces continuously. The single trace tube will display the trace from any amplifier when the FUNCTION selector switch is in the SWEEP position. In the BALANCE position, a pattern is displayed on the single trace screen which permits rapid balancing of the pressure transducers. In LOOPS, any channel can be presented on the horizontal axis of the single trace tube, and at the same time any other channel can be presented on the vertical axis. This position is used for the VCG (vertical) loops (X-Y plot).

SWEEP switch determines the speed at which the traces move across the screen. The proper setting depends on the repetition rate of the phenomena displayed. Nine monitoring speeds are available: 2.5, 5, 10, 25, 50, 75, 100, 150 and 200 mm/sec.

The MULTITRACE MONITOR permits the display of up to eight channels at one time. The positions of the individual traces on the multitrace tube are controlled by the CENTER knob on the individual channels. Each of the traces can be positioned anywhere on the screen, superimposed on other traces, or can be moved off screen in either direction if it is not being used. The amplitude of the traces on the multitrace screen will depend upon the setting of the AMPLITUDE controls on the individual amplifiers.

The CHANNELS switch determines the number of amplifiers that will be displayed. The multitrace presentation is the result of an 8 channel electron switch, whose basic frequency is 48 K Hz. In the 1-8 position, each amplifier will be sampled equally 6,000 times/second, or every 166 microseconds. In the 1-5 position, channels 1, 2 and 3 are sampled at 12,000 times/second or every 82 microseconds. Channels 4 and 5 are sampled at the normal 6,000 times/second and channels 6, 7 and 8 are not displayed at all. Placing the multitrace switch at 1-6 will sample channels 1 and 2 at the 12,000 rate and 3, 4, 5 and 6 at the normal 6,000 rate. Channels 7 and 8 will not appear on the multitrace screen. Amplifiers providing fast rise time phenomena, such as heart sounds, action and nerve potentials, should be in positions 1, 2 and 3 to take advantage of the faster sampling rate and to minimize the dotting effect.

When the SUBTRACT SELECTOR switch is in the OFF position, the built-in baseline, gradient and subtract functions are inoperative. Turning switch to BASELINE will allow built in baseline to be positioned to any point on screen with POSITION control. An event MARK button and jack (for marking from a remote position) are provided. To display either baseline or gradient 1-5, 1-6, 1-8 switch must be in 1-8, since baseline or gradient trace replace channel 8. When BASELINE or GRADIENT position is selected channel 8 amplifier will not be seen.

A PRESSURE GRADIENT can be displayed on the multitrace screen together with the two pressures. The adjustment procedure is as follows:

- 1) Attach pressure transducers to pressure amplifiers; allow them to warm up, balance the amplifiers. Position the traces to baseline from which measurements are to be made; equalize the amplitudes. (See SGA-SGM instructions)
- 2) Turn the SUBTRACT SELECTOR switch to "GRADIENT". The gradient trace will now be substituted for Channel 8 and its position controlled by the baseline POSITION control.
- 3) Turn HORIZONTAL SIZE and SUBTRACT SIZE to zero (counterclockwise)

- 4) Select the channel to be subtracted from with the Horizontal switch; the channel to be subtracted with the Subtract switch.
- 5) Position the gradient trace to the baseline position of the two pressure channels. There are now 3 or 4 lines superimposed i.e., 2 pressure amplifiers, the difference trace (gradient) and optionally a baseline, which can be a baseline, marker trace, or any unused pressure amplifier.
- 6) Turn the +/- switch to +.
- 7) Turn HORIZONTAL SIZE slowly clockwise; note in which direction trace moves. If it moves upwards, turn HORIZONTAL CENTER counter-clockwise to return trace to its original position. Advance HORIZONTAL SIZE fully clockwise and readjust HORIZONTAL CENTER, if necessary. When the center is optimally adjusted, turning HORIZONTAL SIZE should not affect the position of trace. Return HORIZONTAL SIZE to zero. Repeat procedure with SUBTRACT SIZE and CENTER controls, returning SIZE to zero after obtaining balance.
- 8) Turn CALIBRATE control on pressure amplifier which was selected by HORIZONTAL SWITCH so that its trace moves approximately full screen; adjust HORIZONTAL SIZE until gradient trace coincides with the pressure trace. Remove the CALIBRATE signal from the pressure amplifier and both lines will return to their original positions.
- 9) Turn the CALIBRATE on that amplifier selected by SUBTRACT switch until the pressure trace moves approximately full screen. Adjust SUBTRACT SIZE until lines are superimposed. Return CALIBRATE switch to zero and both lines will return to their original positions.
- 10) The "Gradient" trace will now represent the sum of the two pressures. To observe the gradient, turn the - - switch to -. Readjust gradient POSITION, if necessary.

For the H and V sections of SUBTRACT SELECTOR switch, see loops.

SINGLE TRACE MONITOR presentation is determined by the FUNCTION switch. Any individual amplifier can be displayed when the FUNCTION switch is in the SWEEP position. Turn the VERTICAL select or switch to the channel number representing the desired amplifier. Make certain this amplifier is visible on multitrace screen. The position of the selected trace can be controlled by the CENTER control on the channel, or by the VERTICAL CENTER control. The amplitude on the screen is adjusted by the VERTICAL SIZE control and can be 2 1/2 times the amplitude of the corresponding trace on the multitrace tube, or it

can be reduced as low as zero amplitude. When the SIZE control is at zero (completely counterclockwise), the CENTER control will have no effect. If trace is more than 1" from mid-screen, (SIZE control fully counterclockwise) turn screwdriver adjust below VERTICAL CENTER to return baseline to mid-screen. Turning the SIZE control clockwise may move trace off screen. Note in which direction it moves (up or down) and return trace to screen with VERTICAL CENTER control. If trace moves upward, turn CENTER counterclockwise; if trace moves downward, turn CENTER clockwise. "

The BALANCE position of the FUNCTION switch will display a lissajous figure, which is used as an aid in balancing a transducer in a pressure amplifier.

Turn VERTICAL selector switch to pressure amplifier (with transducer connected) to be balanced. The horizontal deflection represents the transducer excitation voltage, the vertical represents the output voltage from the transducer through a pressure amplifier.

The vertical SIZE and CENTER controls do not have any effect in the BALANCE position. When the transducer is brought into balance, there will not be any output voltage, and the lissajous presentation will be a straight horizontal line.

In LOOPS position, the VERTICAL selector switch selects the channel to be presented on the vertical axis, and the vertical SIZE and CENTER controls are used to adjust vertical amplitude and position.

The HORIZONTAL switch selects the channel to be presented on the horizontal axis, and the horizontal SIZE and CENTER controls are used to adjust horizontal amplitude and position. When the HORIZONTAL (or VERTICAL) SIZE control is turned to minimum counterclockwise, the CENTER controls will not have any effect. A dot should appear at mid-screen. If it is more than 1" from mid-screen, adjust screwdriver controls below VERTICAL CENTER or HORIZONTAL SIZE until dot is at mid-screen.

Advance HORIZONTAL SIZE clockwise and adjust HORIZONTAL CENTER to return dot to its original position near mid-screen. Proper adjustment of center is achieved when it is possible to turn horizontal SIZE from minimum counterclockwise to maximum clockwise without moving dot when there isn't any applied signal. Do the same with the VERTICAL and SUBTRACT switches: return the subtract SIZE to minimum counterclockwise.

The H or V positions of SUBTRACT SELECTOR make it possible to add or subtract the signals from two channels on either the horizontal or vertical axis when measuring loops. This is especially useful when measuring Lung Compliance or Lung Resistance. For example: Let us suppose it is desired to subtract from the horizontal axis. The HORIZONTAL switch selects the amplifier to be subtracted from. Then, the SUBTRACT selector switch is used to select the channel to be subtracted. - switch is turned to -. The SUBTRACT SELECTOR is turned from C (Flow) to Horizontal. The SUBTRACT, SIZE and CENTER controls will now also be able to control horizontal amplitude and position but in a direction opposite to that of the HORIZONTAL controls. See instructions for Lung Resistance and Lung Compliance loops.

Timing and direction marking can be introduced to the loop patterns by turning the camera selector to LOOPS and choosing the appropriate TIMING rate. The camera ON-OFF switch need not be ON to display timing on loop pattern, even when photographing loop. Camera timing should be in .004 or .02 to properly display LC/LR loops; or .004 for VCG.


PMB
9/69

ELECTRONICS for MEDICINE, Inc., 30 Virginia Rd., White Plains, N.Y.

MONITOR
MODEL DR-12

The DR-12 recorder is a modified version of the DR-8. The DR-8 instruction plus the following information covers the operation of the DR-12.

The amplifiers in the DR-12 are numbered from 1 through 12. Ten of these, 1-10, starting from the top, are located on the left side of the recorder; numbers 11 and 12 are located beneath the monitor panel. The eight displayable channels of the recorder are lettered from A through H with a dial type switch supplied for each channel which permits any of the 12 amplifiers to be selected on any of the 8 channels. Since it is possible to select the same amplifier on more than one channel, the user should be certain that each of the dials is set to a different number or the same trace will appear on the screen more than once.



The A-E, A-F, A-H switch on the monitor panel can be used to sample two or three channels more frequently than the others. The electronic switching rate is 48,000 Hz/second; each of the 8 channels is sampled 6,000 times per second when the switch is in the A-H position. In the A-E position, channels A, B, and C are sampled 12,000 times/second; the channels D and E are sampled 6,000 times/second. Only

5 channels are operative in this position. In the A-F position, channels A and B are sampled 12,000 times/second while channels C, D, E and F are sampled 6,000 times/second. Operation is limited to 6 channels in this position. Any channel can be sampled at these higher rates by selecting them on channels A or B.

SUBTRACT SELECTOR switch marked BASELINE, GRADIENT, OFF, H and V, when positioned on BASELINE or GRADIENT, will remove channel H (A-E, A-F, A-H switch must be on A-H) from display on multitrace screen and replace it with a baseline or gradient trace. With switch in BASELINE, the POSITION control, below the switch, can be used to set the baseline to any desired point. A MARK push button and remote MARK jack can be used as an event marker. A hand/foot switch can be patched into MARK jack, to mark events from a remote point. Placing switch on GRADIENT position will permit the display of gradients using built in gradient amplifier (POSITION control is the same as for baseline). The subtraction is done horizontally. See section on "subtraction" in DR-8 instructions. The H and V position of the SUBTRACT SELECTOR switch, permit subtraction on vector screen, either vertically or horizontally in the OFF position, subtract or GRADIENTS and BASELINE are inoperative. See DR-8 instructions.


PMB
9/69

ELECTRONICS for MEDICINE, Inc., 30 Virginia Road, White Plains, N.Y.

SINETRAC SERIES
MODEL ST-6800
A/D-D/A PERIPHERAL SYSTEMS
INSTRUCTION MANUAL

Part No. 58-12140-25

JANUARY 1979



SPECIFICATIONS

DATA ACQUISITION SECTION

Specifications

Typical @ +25°C, dynamic conditions, unless otherwise specified.

Analog Inputs

Number of Channels	32 Single-ended or 16 differential
Channel Expansion	Up to 128 single ended or 112 differential using ADC-Slave Expander Cards (Model ST-6800 ADX)
Full Scale Input Ranges	0 to +5 Volts 0 to +10 Volts -5 to +5 Volts -10 to +10 Volts
Current Input Channel Range (8 Channels)	4-20 mA type
Common Mode Range	+10 Volts
Input Overvoltage	+35 Volts Max. continuous
Input Impedance	100 Megohms differential or to ground
Input Bias Current	3nA typ., 10nA max.
Input Capacitance	5pF, OFF CHANNEL, 100 pF ON CHANNEL to ground

Performance

Accuracy @ +25C	Within $\pm 0.025\%$ of input range
Resolution	12 Binary bits (1 part in 4096)
Nonlinearity	$\pm 1/2$ LSB maximum
Differential non-linearity	$\pm 1/2$ LSB maximum
Gain Error	Adjustable to zero
Offset or Zero Error	Adjustable to zero
Gain Temperature Drift (Bipolar)	Within ± 10 ppm of FSR/°C
Zero Temperature Drift (Unipolar)	Within ± 5 ppm of FSR/°C max.
Common Mode Rejection	70 dB min, DC to 1 kHz with 1k Ω unbalance
Power Supply Rejection	100 dB to +5V bus

Dynamic Characteristics

Typical Data Transfer	
I/O Period (Total)	36 microseconds
Throughput Period	20 microseconds
Acquisition Time	8 microseconds
A/D Conversion Time	12 microseconds
Aperture Time	100 nanoseconds
Sample/Hold Switch	.01% Max.
Feedthrough	
MUX Crosstalk from OFF Channels	.007% @ 1 kHz, $R_s = 1K$

Digital Outputs

Output Coding	Straight Binary (Unipolar) } Offset Binary (Bipolar) } Jumper 2's Complement (Bipolar) } Selected
Output Format	2-Byte group electrically compatible to Motorola's EXORciser bidirectional bus. Sign extension is jumper selected on bits 12 thru 15 for 2's complement units. Bits 12 thru 15 are logic zero for all other units.
Channel Addressing	Random channel addressing may be started by external interrupt input for event operation or by internal program control.
Base Address	Prewired by PC Board jumpers for one of 128 base addresses.

Data Distribution Section (D/A Analog Outputs)

Number of Channels	2 (Expandable only by stand-alone ST-6800DA Boards.)
Resolution	12 Bits
Full Scale Output	0 to +5 Volts
Voltage Ranges	0 to +10 Volts -5 to +5 Volts -10 to +10 Volts
Input Coding	Straight Binary (Unipolar) Offset Binary (Bipolar) 2's Complement (Bipolar)
Output Impedance	.05 ohm
Output Current	± 5 mA min

Performance

Non linearity	$\pm 1/2$ LSB, maximum
Differential Nonlinearity	$\pm 1/2$ LSB, maximum
Gain Error	Adjustable to zero using pot. for each channel
Offset or Zero Error	Adjustable to zero using pot for each channel
Gain Temperature Drift	± 20 ppm of output/ $^{\circ}\text{C}$
Zero Temperature Drift	(Unipolar output) ± 5 ppm of FSR/ $^{\circ}\text{C}$
Offset Temperature Drift	(Bipolar output) ± 10 ppm of FSR/ $^{\circ}\text{C}$
Settling Time (20V change)	4 microseconds to $\pm 1/2$ LSB
Slew Rate	20V/usec
Power Supply Rejection	$\pm 0.02\%$ of FSR per 1% variation

Power Consumption

1.2 ampstypical @ +5 vdc supplied from MPU bus connector. On-board DC to DC Converter supplies ± 15 vdc to linear circuits.

Physical

Operating Temperature Range	0° to $+70^{\circ}\text{C}$
Storage Temperature Range	-25°C to $+85^{\circ}\text{C}$
Card Size	9.75"W x 5.75"H x .062"D

10/1/57

10/1/57

Room 243

Mr. Carlson

NYT/55318

BRUSH

OPERATING INSTRUCTIONS

Brush Recorder Mark II
Model RD 2522 20

General Information

Specifications

Sensitivity	10 millivolts per chart line (mm). Full scale deflection from chart center : 200 millivolts.
Sensitivity steps01, .02, .05, .1, .2, .5, 1, 2, 5 and 10 volts per chart line (mm). Maximum attenuator error 1% with balanced input.
Measurement range	
Single-ended input (5 megohm)010 volt to 400 volts
Balanced input (10 megohm)010 volt to 400 volts, side to side, allowable voltage off-ground at any attenuator step is 1000 x volts per chart line switch setting up to 500 volts off-ground.
Common mode rejection	Better than 1000 to 1, attenuator set in .01 volt per chart line position.
Zero line stability	Less than 1/4 chart line (mm) per hour. Total drift over eight hour period not more than 1/2 chart line (mm).
Noise	Not noticeable on chart with shorted input.
Frequency response	The recorded peak-to-peak amplitude of a constant voltage sine wave will be within ± 1 chart line (mm) of a nominal 10 lines from 0.01 to 100 cps.
Maximum amplitude	40 lines (mm) peak-to-peak, D.C. to 40 cps. 20 lines (mm) peak-to-peak, D.C. to 70 cps. 10 lines (mm) peak-to-peak, D.C. to 100 cps.

General Information

Specifications

Pen Bias	Permits positioning of pen on chart, ± 20 chart lines (mm). Effective for either single-ended or balanced input.
Trace linearity	D. C. within 2% full chart width. A. C. within 3% full chart width, any frequency within limits of maximum amplitude for electric writing.
Trace width	0.006 with Model RA 2522-31 pen.
Writing method	Electric stylus.
Number of recording channels	2
Number of event channels	1 actuated by external switch. 1 actuated by panel switch.
Channel width	40 mm (40 divisions).
Chart Supply	150 feet
Chart speeds	1, 5, 25, and 125 mm per second
Chart speed regulation	Synchronous motor, direct drive.

General Information

Specifications

Operating Temperature range, ambient 0°C to 55°C.

Power requirements 105-125 volts, 60 cps, 135 watts at 115 volts.

Transistors 1-5651 and 2-5687

Tubes 1-6BW4, 2-12B4A, 2-12AT7, 2-12AX7

Input terminals

Front (signal) Binding posts.

Rear (event Marker) Binding posts.

Supplies

DESCRIPTION	BRUSH PART NO.
Chart Paper, 2 channel	RA 2922-22
Electric Styli (4)	RA 2822-31
Gram Gage Assembly	227416-910
Pen Mounting Tool	126717

CHAPTER 1 GENERAL DESCRIPTION

1-1 INTRODUCTION

The M6800 EXORciser (M68SDT) is a system development tool used in the design and development of M6800 Microcomputer Systems. The EXORciser Debug and the user's system are built around the M6800 Microcomputer Family of Parts. The M6800 Microcomputer Family of Parts are discussed in the following documents:

- M6800 Microprocessor Programming Manual
- M6800 Microprocessor Applications Manual
- M6800 Microprocessor Family of Parts Data Sheets

The EXORciser may be configured in a variety of applications and with various EXORciser options. This manual, rather than discussing every possible configuration and option, discusses only the basic EXORciser with each option except the Wirewrap and Extender Modules discussed in a supplement to this manual. The basic EXORciser is discussed in Paragraph 1-5 and the options are identified in Paragraph 1-6.

This manual provides general information, installation instructions, applications information, operating

procedures, and theory of operation, for Motorola's M6800 EXORciser. The M6800 EXORciser is illustrated in Figure 1-1.

1-2 EXORciser FEATURES

The features of the basic EXORciser include:

- Flexible, adaptable, and expandable design development tool
- Easy to use
- Provides the Microprocessing Unit capability for both the EXORciser and the user's system
- Saves system design and development time
- Decreases system design and development costs
- Evaluates and debugs the user's program
- Evaluates and debugs the user's system hardware

1-3 EXORciser SPECIFICATIONS

Table 1-1 identifies the basic EXORciser specifications.

TABLE 1-1. Basic EXORciser Specifications

CHARACTERISTICS	SPECIFICATIONS
Power Requirements	95-135 205-250 V AC 47-420 Hz, 250W
Word Size	
Data	8 bits
Address	16 bits
Instructions	8, 16, or 24 bits
Memory Size	64k bytes
Instruction Set	72 variable length instructions
Clock cycle time	Selectable, 1 μ s or an external clock between 1 and 10 μ s
Interrupt	Maskable real time interrupt
Physical Characteristics	
Table top	
Length	19.25 in
Depth	17.50 in
Height	7.00 in
Rack Mountable	
Length	19.00 in
Depth	17.00 in
Height	7.00 in
Baud Rates (Switch Selectable)	110, 150, 300, 600, 1200, 2400, 4800, and 9600

1.4 FUNCTIONAL DESCRIPTION

The M6800 EXORciser, illustrated in Figure 1-2, is a general purpose emulator for the shelf systems (single board computers) of the M6800 Microcomputer Family of Parts. The EXORciser may be easily tailored to meet the user's needs in the design and development of his system. Its modular design reduces the time required to develop a system and at the same time, provides great flexibility in configuring an emulation (functional representation) of the user's system. The EXORciser's EXbug Firmware, through its debug and program control features, minimizes the time required to develop user's systems. The EXbug firmware provides the EXORciser with the capability to:

- Display the contents of the MPU registers
- Step through the program
- Trace through user's programs to locate problem areas
- Stop the program on a selected program step
- Provide an oscilloscope trigger signal on a selected program step
- Abort from the user's program and return to the EXbug control program on command
- Reinitialize the EXORciser on command
- Change the contents of memory

The user communicates with the EXORciser in one of two ways:

- Through a RS 232C or TTY terminal
- Through the EXORciser front panel controls and indicators

The terminal device permits the user to communicate directly with the EXORciser's EXbug Firmware. The EXORciser's front panel permits the user to apply power to the EXORciser, to abort (exit) the EXORciser from a routine, and to initialize and restart the EXORciser. The EXORciser's unique front panel was designed to incorporate future EXORciser options such as data keys and displays.

1-5 BASIC EXORciser SYSTEM (FIGURES 1-2 AND 1-3)

The basic EXORciser (M68SDF) consists of the MPU Module, the Debug Module, the Baud Rate Module, the Power Supply, and the chassis. These modules are built around the M6800 Microcomputer Family of Parts (MC6800 Microprocessing Unit, MC6820 Peripheral Interface Adapter, MC6850 Asynchronous Communications Interface Adapter, MCM6810 Random Access Memory, and MCM6830 Read Only Memory devices).




FIGURE 1-2 EXORciser Simplified Block Diagram

The MPU Module (Figure 1-4) incorporates the MC6800 Microcomputer Unit (MPU) and the system clock. This module, as illustrated in Figure 1-2, serves a dual function in that the module provides the MPU and clock to the EXORciser. It is part of the basic system. The MPU Module also initiates an EXORciser restart operation and initializes the EXORciser. The MC6800 Microprocessing Unit is an 8-bit parallel device capable of addressing 64k bytes of memory. In addition,

the MPU addresses its input and output as memory. The MPU also provides the EXORciser with 72 variable length instructions and the capability of reading and storing data to memory. When the EXORciser enters system start-up, the MPU (after a delay) the MPU Module appears exactly like a MC6800 Microprocessing Unit with unlimited TTL bus drive capability.




FIGURE 1-3. Basic EXORciser System




FIGURE 1-4. MPU Module

Vita

Melvin D. Townsend was born 5 April 1952 in Colorado Springs, Colorado. Upon graduation from Gen. William Mitchell High School in 1970, he enrolled in undergraduate engineering study at Colorado State University. In December, 1974, he was awarded the BSEE degree and commissioned as a second lieutenant, USAF, via AFROTC. Entering active duty immediately, he was assigned to the Foreign Technology Division (AFSC), WPAFB, Ohio. After three and one half years working technical intelligence engineering, he was transferred PCA to the AF Avionics Laboratory (AFSC) where he performed R&D engineering on electro-optical target acquisition systems. In June, 1979, he entered the Air Force Institute of Technology. Captain Townsend is a member of Eta Kappa Nu and IEEE.

Permanent Address

621 Prairie Rd.

Colorado Springs, CO

80909

UNCLASSIFIED

2. REPORT DOCUMENTATION PAGE

REPORT DOCUMENTATION PAGE		READ IN 15 MIN BRIEFING CENTER FORM
1. REPORT NUMBER AF/ED/10D-46	2. GOVT ACCESSION NO. AF-1100-111	3. REPORTING DATA 100-111
4. TITLE (and Subtitle) Analysis and Performance Evaluation of Electrocardiogram Data Compression Techniques		5. TYPE OF REPORT MS Thesis
7. AUTHOR MELVIN D. TOWNSEND Captain USAF		8. CONTRACT OR GRANT NUMBER
9. PERFORMING ORGANIZATION NAME (and Address) Air Force Institute of Technology (AFIT/EN) Wright-Patterson AFB, Ohio 45433		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
11. CONTROLLING OFFICE NAME AND ADDRESS Clinical Sciences Division USAF School of Aerospace Medicine/NSF (UAAFSAM/NSF), Brooks AFB, Texas 78225		12. REPORT DATE December 1980
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		13. NUMBER OF PAGES 252
		15. SECURITY CLASS. of this report Unclassified
16. DISTRIBUTION STATEMENT of this Report Approved for public release; distribution unlimited		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
17. DISTRIBUTION STATEMENT (in the abstract enter 1 or Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES Approved for public release; LAW AER 100-17 <i>Fred C. Lynch</i> FREDIC C. LYNCH, Major, USAF Director of Public Affairs 16 JUN 1981		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) Electrocardiogram Data Reduction EKG Data Compression Source Encoding		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) EKG data compression techniques were investigated for potential real time implementation on an 8 bit Motorola 6800 microprocessor. Research indicated entropy reduction transform technique such as the Fast Fourier Transform and the discrete Karhunen-Loeve Transform were not feasible for implementation on the 6800. Two redundancy reduction (RR) techniques (TOLAN and DOWNS) utilizing 2nd order difference		

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

operations in conjunction with variable length encoding were studied in detail. One such technique (CLIR) can be implemented and tested with in vivo ECG data. Analysis revealed compression ratios ranging from 1.25:1 to 2.26:1. Investigation of the poor performance of the compression algorithm showed significant degradation of the 2nd order difference "decorrelator" due to a noisy collection environment. It was concluded that real time EKG data compression is feasible on the 6800 but that time compression techniques which store a zero value sequence counter versus the value of zero are not efficient in a high noise environment.

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

DATE
FILMED
-18-